

Introduction

Selenium Isotope Dilution Analysis using Innovative Reaction cell ICP-MS

Selenium (Se) is one of the major elements required to be quantified in various types of samples due to its nutritional essentiality in animals and humans, and its toxicity when ingested in excess. However, for several reasons, it is a difficult element to accurately quantify. The first is that the required Limit of Quantification (LOQ) for Se is commonly as low as ppb to sub ppb levels. The second is that Se has high ionization potential (IP=9.75e V), causing low sensitivity due to poor ionization in plasma, susceptibility to matrix suppression, and difficulty in finding suitable Internal Standard (ISTD) elements. The third is that all Se isotopes suffer from severe interferences as summarized in Table 1. Isotope Dilution (ID) methodology is well known for its absolute accuracy due to its ability to minimize the effects of signal drift and matrix suppression, which are common problems in ICP-MS. However the application of isotope dilution to Se quantification has been limited because of the third reason stated above: interferences.

Introduction

The innovative reaction cell of the new Agilent triple quadrupole ICP-MS, enables measurement of more than one Se isotope free from interferences, making it possible to quantify Se using isotope dilution. Recently, Giuseppe et al.¹ applied the ID method for multi element analysis focusing on less susceptibility to matrix suppression, rather than on absolute accuracy. They proposed Online Isotope Dilution Analysis (OIDA) which removes the time-consuming step of spiking enriched-isotope standards into the samples.

Experimental




Table 1. Interference on major Se isotopes

m/z	Abundance %	Isobar	argide	oxide	hydride	chloride	doubly charged	dimer
77	7.63		³⁷ Cl ⁴⁰ Ar+, ³⁹ K ³⁸ Ar+	⁶¹ Ni ¹⁶ O+, ⁶⁵ Co ¹⁸ O+	⁷⁶ GeH+, ⁷⁶ SeH+	³⁷ Cl ⁴⁰ Ca+	¹⁵⁴ Sm ⁺⁺ , ¹⁵⁴ Gd ⁺⁺	
78	23.77	⁷⁸ Kr+	³⁸ Ar ⁴⁰ Ar+, ³⁸ Ar ⁴⁰ Ca+	⁶² Ni ¹⁶ O+	⁷⁷ SeH+	⁴¹ K ³⁷ Cl+	¹⁵⁶ Gd ⁺⁺ , ¹⁵⁶ Dy ⁺⁺	³⁸ K ³⁹ K+
80	49.61	⁸⁰ Kr+	⁴⁰ Ar ⁴⁰ Ar+, ⁴⁰ Ca ⁴⁰ Ar+	⁶⁴ Ni ¹⁶ O+, ⁶⁴ Zn ¹⁶ O+	⁷⁹ BrH+	⁴⁵ Sc ³⁵ Cl+	¹⁶⁰ Gd ⁺⁺ , ¹⁶⁰ Dy ⁺⁺	³⁸ K ⁴¹ K+
82	8.73	⁸² Kr+	⁴² Ca ⁴⁰ Ar+	⁶⁶ Zn ¹⁶ O+	⁸¹ BrH+	⁴⁶ Sc ³⁷ Cl+	¹⁶² Dy ⁺⁺ , ¹⁶² Er ⁺⁺	

Figure 1. Configuration of Agilent ICP-QQQ

Experimental

Agilent Technologies has developed a new triple quadrupole ICP-MS (ICP-QQQ). As shown in Figure 1, it has two quadrupoles, one before and one after the collision/reaction cell. The first quadrupole selects ions entering the cell, providing **consistent reaction conditions** during changing sample composition. It solves the problems associated with current reaction cell technologies, allowing the analyst to use the reaction cell for more elements and applications, more effectively. The Agilent ICP-QQQ has the following features:

- Superior interference removal:
 - Advanced reaction cell technology overcomes the weaknesses of existing reaction cell instruments due to the unique QQQ configuration.
 - Versatile He collision mode with performance better than the Agilent 7700 Series ICP-MS due to MS/MS scanning ability.
- Sensitivity is much higher than the single quadrupole 7700 ICP-MS. A lab-prototype achieved 810 Mcps/ppm for Yttrium with oxide ratio (CeO+/Ce+) of 1.47% when using the high matrix x-lens and 1.2 Gcps/ppm with oxide ratio of 2.47% when using the high sensitivity s-lens.
- Random background noise is as low as <0.2cps
- Abundance sensitivity is immeasurably low: <10-9
- Maintains Agilent 7700 compatible robustness and durability, being applicable to developed Agilent 7700 applications

Results and Discussion

Cell gas mode investigation

The interference removal capability was studied to determine the optimum cell gas conditions and the best Se isotopes for the ID method. Experimental results are given in Figure 2. The figure shows the background equivalent concentrations (BECs) using four reaction gas modes in various synthetic matrix samples. All modes used MS/MS scanning that operates the first Quad as a 1 amu-wide bandpass filter, selecting ions entering reaction cell. While Se isotopes were detected at the natural mass in the three cell gas modes, He, H₂ and NH₃+H₂ modes, they were detected as SeO⁺ in O₂+H₂ mode. In this mode, Se⁺ reacts with O₂, converting Se⁺ to SeO⁺. We refer to this as the mass shift method. Figure 3 illustrates the mass shift method. The reaction between Se⁺ and O₂ is endothermic ($\Delta H_f = -0.769$ eV), however, the reaction was promoted through the application of an octopole bias of -18 V, effectively forming SeO⁺. Addition of small amount of H₂ into the cell was found to be effective at further reducing the BEC.

From Figure 2, we can see O₂+H₂ mass shift method gives the best interference removal capability, allowing for the measurement of three Se isotopes, ⁷⁸Se, ⁸⁰Se and ⁸²Se nearly free from interference. O₂ cell gas can be used with existing single-quadropole cell ICP-MS, though the method does not work well when the sample contains either Zr, Mo or Ru, as their isotopes overlap SeO⁺. With ICP-QQQ, due to ion selection by Q1, the method can be applied to samples regardless of the composition. This is a clear advantage of ICP-QQQ over existing cell ICP-MS, with the developed method being applicable to a variety of sample types.




Figure 3. Illustration of O₂ mass shift method of ICP-QQQ

Selenium Online Isotope Dilution Analysis (OIDA)

A ⁸²Se enriched standard (⁸²Se 97.43%, ⁸⁰Se 1.65% and ⁷⁸Se 0.51%) was purchased from Oak Ridge National Laboratory. The standard was dissolved in ultrapure HNO₃, and diluted to the appropriate concentration. It was spiked to all tested samples online using the Agilent online-ISTD mixing kit. Isopropyl alcohol (IPA) was added to the spike solution to enhance Se⁺ signal. Refer to Reference 1 for details of OIDA. The Se concentration of the sample was calculated using the measured isotope ratios of two Se isotopes according to the formula presented below.

$$C_x = C_s \left(\frac{R_m - R_n}{R_m - R_s} \right) \left(\frac{R_s - R_n}{R_s - R_m} \right) \quad (1)$$

R_m : measured isotope ratio of mixed reference standard and spike solution
 R_n : measured isotope ratio of mixed unknown sample and spike solution
 C_x : concentration of unknown sample
 C_s : concentration of natural reference standard
 R_s : isotope ratio of spike




Figure 2. Result of interference removal capability study

