Application Note Clinical Research



Analytical Development of a Four-Stream Multiplexed LC/MS Method for the Simultaneous Determination of SDMA and ADMA

Authors

Haiqing Ding and Cory Bystrom, Ph.D Cleveland HeartLab Cleveland, Ohio

Abstract

A profound technical challenge to laboratories today is achieving the highest productivity while meeting the demands for precise, accurate, and cost-effective test results. The Agilent StreamSelect LC/MS System can help deliver all this, by seamlessly multiplexing up to four concurrent HPLC separations with one triple quadrupole mass spectrometer. The system can accommodate a single method mirrored across all HPLC systems, or can have an unique method assigned to each HPLC. With MassHunter-based software-driven scheduling and acquisition, a four-fold increase in throughput can be realized without any loss of analytical fidelity. This Application Note demonstrates the development and confirmation of an analytical method for asymmetric and symmetric dimethyl arginine (ADMA and SDMA respectively).

Introduction

ADMA and SDMA are two modified amino acids that are emerging biomarkers. A method for the simultaneous determination of ADMA and SDMA by LC/MS/MS was developed and rigorously confirmed on a four-stream Agilent StreamSelect LC/MS system. The data demonstrate excellent and equivalent analytical performance across all four streams, and show that the system can meet strict analytical performance targets.

Experimental

Stock solutions of ADMA and SDMA were prepared at 1 mg/mL in water. Calibrators for ADMA and SDMA were prepared in charcoal-stripped human serum at 10, 20, 50, 100, 250, and 500 ng/mL by appropriate dilution of stock solutions, and stored frozen. A 50-µL volume of human serum or plasma was mixed with 150 µL of methanol containing 60 ng/mL of heavy-labeled ADMA to precipitate proteins and provide an internal standard. The precipitate was removed by centrifugation, and the supernatant transferred to a clean vial, and subsequently diluted with 200 µL of 0.5 % heptafluorobutyric acid in water prior to injection.

Instrumentation

An Agilent 6490 triple quadrupole LC/MS system with Agilent StreamSelect software was configured with an HTC-PAL autosampler with four injection ports, four Agilent 1260 Infinity binary pumps, and a 9-port/8-position stream selection valve.

Instrument conditions

LC conditions							
Analytical column	Phenomenex Kinetix 2.6 μ M 2.1 × 50 mm, alternatively, an Agilent Poroshell						
Column temperature	Ambient						
Injection volume	3 µL						
Mobile phase	A) H ₂ O + 0.005 % HFBA B) Methanol + 0.005 % HFBA						
Flow rate	0.4 mL minutes						
Stop time	5 minutes						
Isocratic separation	boratic separation 5 % mobile phase B for 2.5 minutes with a wash at 95 %B for 2.5 minutes						
MS/MS conditions							
Instrument	Agilent 6490 triple quadrupole LC/MS						
lon mode	Positive						
Drying gas temperature	250 °C						
Gas flow	14 L/min						
Nebulizer	20 psi						
Sheath gas temperature	300 °C						
Shealth gas flow	12						
Capillary voltage	3,000						
EMV	+100 V						

MRM transitions

	Analyte	Precursor	Product	Fragmentor	Dwell	Collision energy
	ADMA-d6-quant	209.2	69.9	380	100	27
	ADMA-d6-qual	209.2	52	380	100	22
	SDMA-quant	203.2	172.1	380	100	10
	SDMA-qual	203.2	69.9	380	100	10
	ADMA-quant	203.2	69.9	380	100	27
	ADMA-qual	203.2	46	380	100	22

Results and discussion

Imprecision studies were carried out running eight replicate injections from each of five unique batches over a span of 14 days. In the absence of reference material, bias was determined against assigned values established for synthetic and sample pools from calibrators and reagents prepared separately from the ones used in this study. Each individual stream was evaluated separately. For both analytes, the mean bias was below 4.0% for a high QC pool, and below 2.5% for a low QC pool. Interbatch imprecision was below 3.0% CV, and intrabatch imprecision was below 2.4% CV, as shown in Table 1.

Table 1. Assay performance.

	ADMA	SDMA		
Intra-assay imprecision	<1.7%	<2.4%		
Inter-assay imprecision	<3.0%	<2.8%		
Mean Bias at 30ng/mL	0.7%	2.5%		
Mean Bias at 90 ng/mL	1.2%	3.9%		

Calibration stability

Calibration samples were run with each batch. Linear regression with 1/x weighting was used to generate calibration curves. To illustrate stability, calibration statistics for slope, intercept, R², and maximal back-calculated deviation at day 1 and day 14 are found along with mean internal standard intensity in Table 2. The individual HPLC streams generate calibration curves that are essentially indistinguishable from each other.

Table 2. ADMA calibration statistics. Slope Stream Day 1 Day 14 1 1.0000 0.9996 2 1.0000 1.0010 3 1.0020 1.0010 4 0.9997 0.9950 Intercept 1 0.0180 -0.0184 2 -0.0158 0.0093 3 -0.0493 -0.0207 4 -0.0312 -0.0160 R² 1 0.9997 0.9999 2 0.9998 1.0000 3 0.9998 0.9998 4 0.9997 0.9999 Maximum % backcalculated deviation 1 3.0 -4.0 2 2.0 2.0 3 -3.0 -4.0 4 -3.0 3.0 Mean IS area for calibrator 1 33433 38601 2 40489 43467 3 34397 40618 4 48449 37296

Stream-by-stream performance

Figures 1A and 1B show the individual data points for each observation for sample pool with eight observations per day per stream carried out over 14 days, and emphasizes the quality

of the analytical data generated in multistream mode. The data revealed no dependence on injection order and very good performance, with only 12 of 320 analyses having a bias of more than 5 %.



Figure 1. Low QC pool precision and accuracy for ADMA (A) and SDMA (B).

The suitablity for the four-stream system to be fully multiplexed and thought of as one system with no distinction between streams can conveniently be evaluated using total error, which incorporates estimates of both bias and imprecision. To ensure that each stream performed acceptably across a wide dynamic range, we examined the total error for each stream as a function of concentration. Figure 2 shows that the total error for each of four streams remains below the threshold for ADMA and SDMA: between 7.5 and 400 ng/mL for ADMA, and 5.0-400 ng/mL for SDMA. This performance indicates that full multiplex operation where calibrators, unknowns, and QC specimens are distributed across all four streams is acceptable, providing the highest throughput.

Conclusions

Deploying analytical methods that retain optimal performance as throughput increases is a challenging task. The StreamSelect LC/MS System provides an effective multiplexing platform. We have developed a four-stream method that meets the stringent analytical performance criteria for the simultaneous analysis of SDMA and ADMA. The data support the view that the four-stream HPLC system can be considered as one system, allowing distribution of calibrators, QC samples, and unknowns across all available streams for highest throughput.





www.agilent.com/chem

For Research Use Only. Not for use in diagnostic procedures.

This information is subject to change without notice.

© Agilent Technologies, Inc. 2016, 2018 Printed in the USA, October 24, 2018 5991-6764EN

