MACHEREY-NAGEL

CHROMAFIL® Xtra





Contents

Introduction	3
Technical information	5
Specifications	11
Quality control	12
Packages and membranes	13
Chemical compatibility	15
Optimal use of CHROMAFIL®	16
How to select the optimal CHROMAFIL® syringe filter	17
Ordering information	18
CHROMAFIL® – Disposable filters from MACHEREY-NAGEL	31

CHROMAFIL Finder online

The cross reference for syringe filters www.mn-net.com/chromafilfinder







Introduction

CHROMAFIL® – Ideal for GC, HPLC and UHPLC sample filtration

By using CHROMAFIL® syringe filters unwanted particles are removed which can otherwise cause contamination and clogging of sensitive instrumentation. So, filtering your samples leads to an increase of lifetime and less downtime of chromatographic columns and equipment, gives more consistent and reproducible results and thereby saves time and money.

Features

- Different membrane types to meet multiple filtration needs
- Low content of extractable compounds
- Fast flow geometry for easy filtration
- Low hold-up volume for maximum filtrate recovery
- HPLC certified
- Designed to be compatible for the use on all common automated filtration systems, e.g., SOTAX® dissolution systems

Your benefits when using CHROMAFIL® syringe filters

- Protection of sensitive instrumentation
- Less system downtime and increased column lifetime
- Reproducible and reliable analytical results

CHROMAFIL® Xtra

Labeled for method validation and certification

- Xtra imprint for the direct identification of the membrane type, diameter and pore size
- Xtra low bleeding PP housing
- Xtra color-free plain polypropylene





 Low content of extractable substances due to a high density polypropylene housing combined with ultrapure filtration membranes.

In comparison to filters made of polycarbonate, polyacrylate or polystyrene, all CHROMAFIL® housings are resistant to nearly all organic solvents. (see list of chemical compatibility on page 15)




HPLC test


Conditions	2 mL of the solvent (specified on top of the chromatograms) were applied to the filter; 100 μL of the filtrate were injected into the HPLC
Eluent A	water
Eluent B	acetonitrile
Gradient	10 % → 95 % B in 25 minutes
Flow rate	0.5 mL/min
Sensitivity	-5 to 10 mAU at 254 nm
Column	125/4 NUCLEODUR® C ₁₈ Gravity 5 μm (REF 760100.40)

Acetonitrile:





Blank:




Low Bleeding PP housing

Even a treatment with very aggressive solvents/reagents does not lead to significant blind values or extractables.

For proof a filter was rinsed with 0.5 mL MSTFA (N-methyl-N-trimethylsilyl-trifluoroacetamide), a very powerful silylation reagent. The result is shown in the GC/MS chromatogram.







Lowest content of extractable substances

- The housing of every CHROMAFIL® filter is ultrasonically sealed.
- The filters are welded, not glued, because glue may have extractable ingredients.
- The welding leads to a tight connection between both parts, thus the filter can be used in both directions. No fluid can leak from the filter housing.

The special thick rim of the housing is ideal for use in laboratory robots (e.g., $SOTAX^{\otimes}$, $Benchmate^{TM}$, ...).



Luer lock on inlet

 For a safe connection on the "high pressure" side every CHROMAFIL® filter provides a Luer lock on inlet.

Luer outlet

 For the 3, 13 and 25 mm diameter filters: standard Luer outlet. This Luer configuration offers low hold-up volume and easy filtration into autosampler vials and NMR tubes

Filter inlet and filter outlet can be fitted to all CHROMABOND® columns and accessories for selective sample preparation with the aid of a special adaptor.









No breakage of the membrane

- Stabilizing "crash" plate
- The sample fluid is deviated in four lanes by the "crash" plate and does not directly hit the membrane. The resulting pressure distribution protects the membrane against breakage.

Optimal flow geometry

- By star-shaped distribution plate
- The fluid penetrates the membrane on the whole surface, not only on a small area; the filter will not clog rapidly, which guarantees a high flow efficiency.

Different pore sizes

- For multiple filtration application
- Available pore sizes 0.2 and 0.45 μm (additional: PET filters with 1.2 μm, glass fiber filters with 1 μm, PES filters with 5 μm). Filters with 0.45 μm pore size remove fine particles which can clog chromatography columns. 0.2 μm pore size filters are recommended for the filtration of UHPLC samples.

Combi filters

For solutions with a high load of particulate matter







Filter Sizes

3, 13 and 25 mm are the effective membrane diameters. The small diameter filters are especially recommended for very small samples, which require extremely low dead volumes.


Sample volume	Recommended membrane diameter	Dead volume	Filtration area
≤ 1 mL	3 mm	5 μL	0.07 cm ²
1–5 mL	13 mm	30 µL	1.33 cm ²
5–100 mL	25 mm	80 µL	4.91 cm ²

All filters can be autoclaved at 121 °C and 1.1 bar for 30 min.

25 mm CHROMAFIL® filters are designed to be 100 % compatible and reliable for use with the SOTAX® AT70 smart fully automated dissolution testing systems.



Outer dimensions of 25 mm and 13 mm CHROMAFIL® filters





Specifications


Enhanced quality control for better results


MN certifies CHROMAFIL® syringe filters to be low in UV absorbing extractables.

All filters and membrane types have been HPLC tested for compatibility with the most common HPLC solvents (methanol, water, acetonitrile, see test chromatograms).

HPLC-test certificates are available for every membrane type.

Please visit: www.mn-net.com





A specification data sheet is available for all membranes and filter diameters.

Quality control

Pressure stability of CHROMAFIL®

The syringe filter housing is stable up to 12 bar. The "blue" test - membrane, pressure and filtration batch test - with blue colored silica particles in matching particle sizes provides an excellent method to find leaks or membrane deviations.







Packages and membranes

Package sizes

- Packs of 100 or 400 (BIGboxes) for 25 mm Ø filter
- Packs of 100 for 13 mm Ø filter
- Packs of 100 for 3 mm Ø filter
- Packs of 50 for sterile filter



Different membrane materials

Depending on your filtration task you can choose filter membranes made from different materials:

Polyester	(PET) with or without
•	glass fiber prefilter*
	<u> </u>
Regenerated cellulose	(RC) with or without
	glass fiber prefilter*
Polytetrafluoroethylene	(PTFE) with or without
	glass fiber prefilter*
Hydrophilized Polytetrafluoroethylene	(H-PTFE)
Cellulose mixed esters	(MV)
Cellulose acetate	(CA) sterile and non-sterile
Polyamide / Nylon	(PA) with or without
	glass fiber prefilter*
Polyethersulfone	(PES)
Polyvinylidene difluoride	(PVDF) with or without
	glass fiber prefilter*
Glass fiber	(GF)

^{*} Filters with (nom. 1 µm) GF prefilter provide a 2-4 times greater throughput than filters without prefilter for extremely viscous and most difficult-to-filter samples

Packages and membranes

Combi syringe filters

Benefits

- A coarse glass fiber (GF) prefilter and a small-pore membrane as main filter
- For solutions with a high load of particulate matter: lower back pressure, easy filtration
- For high yields of filtrate: more mL of pure filtrate per filter


The technology

The glass fiber membrane (1 μ m) removes coarse particles, before they can block the fine main membrane. This results in a better filtration efficiency, particularly for highly contaminated samples.

Housing	solvent-resistant, ultra low bleed polypropylene
Inlet	Luer lock
Outlet	Luer
Pore diameter	1.0/0.20 μm or 1.0/0.45 μm
Filter diameter	25 mm
Void volume	< 80 µL
Packing unit	100 filters/BIGboxes with 400 filters

Available membranes with GF-prefilter

- Polyester (PET)
- Regenerated cellulose (RC)
- Polyvinylidene difluoride (PVDF)
- Polytetrafluoroethylene (PTFE)
- Polyamide / Nylon (PA)



Chemical compatibility

Solvent	Materia MV	al CA	RC	PA	PTFE	H-PTFE	PVDF	PES	PET	GF	IC	PP	
Acetaldehyde	-	-	+	0	+	+	+		+	+		0	1
Acetic acid, 100 %	_	-	_	_	+	+	+	+	+	+		+	1
Acetone	-	-	+	+	+	+	-	-	+	+		+	1
Acetonitrile	-	_	+	+	+	+	+	+	+	+		+	1
Ammonia, 25 %	-	-	0	-	+	+	+	+	0	+	-	+	1
Benzene	+	+	+	+	+	+	0	+	+	+		0	1
n-Butanol	+	+	+	0	+	+	+	+	+	+		+	1
Cyclohexane	+	+	+	0	+	+	+	+	+	+		+	1
Dichloromethane	+	-	+	-	+	+	+	-	+	+		_	1
Diethyl ether	0	0	+	+	+	+	+	+	+	+		0	1
Dimethylformamide	-	-	0	+	+	+	-	-	+	+		+	1
1,4-Dioxane	-	-	+	+	+	+	0	-	+	+		0	1
Ethanol	-	+	+	+	+	+	+	+	+	+		+	1
Ethyl acetate	-	-	+	+	+	+	+	+	+	+		0	1
Ethylene glycol	0	0	+	+	+	+	+	+	+	+		+	Data not guaranteed
Formic acid, 100 %	+	-	0	_	+	+	+	+	0	+		+	+ resistant
Hydrochloric acid, 30 %	-	_	-	-	+	+	+	+	-	+	-	+	- not resistant
Methanol	-	-	+	+	+	+	+	+	+	+		+	o limited resistance
Nitric acid, 65 %	_	-	_	_	0	0	0	_	0	+	_	_	Material
Oxalic acid, 10 % aqueous	+	-	+	_	+	+	+	0	+	+		+	Membranes
Petroleum ether	+	+	+	+	+	+	+	+	+	+		+	MV = cellulose mixed
Phosphoric acid, 80 %	-	-	0	-	+	+	0	-	+	+	-	+	CA = cellulose aceta RC = regenerated ce
Potassium hydroxide, 1 mol/L	-	-	0	+	+	+	0	0	0	0	+	+	PA = polyamide
2-Propanol	+	+	+	+	+	+	+	+	+	+		+	PTFE = polytetrafluo
Sodium hydroxide, 1 mol/L	-	_	0	+	+	+	0	0	0	0	+	+	H-PTFE = hydrophili
Tetrachloromethane	+	-	+	+	+	+	0	-	+	+		0	ethylene
Tetrahydrofuran	-	-	+	0	+	+	+	-	+	+		0	PVDF = polyvinylidene PES = polyethersulfon PET = polyester GF = glass fiber IC = special filter for io Housing material
Toluene	+	-	+	+	+	+	+	+	+	+		0	
Trichloroethene	+	+	+	0	+	+	+	0	+	+		0	
Trichloromethane (chloroform)	+	-	+	-	+	+	+	-	+	+		-	
Urea	+	+	+	+	+	+	+	+	+	+		+	
Water	+	+	+	+	+	+	+	+	+	+	+	+	
Xylene	+	+	+	+	+	+	0	0	+	+		0	PP = polypropylene

esters

ellulose

roethylene

zed polytetrafluoro-

ne difluoride

ion chromatography

Optimal use of CHROMAFIL®

For achieving the full benefits of filtration we recommend the following instructions.











Draw up the sample into the syringe. Then draw approximately 1 mL of air into the syringe. The air helps to minimize the remaining fluid in the filter.

Plug the CHROMAFIL® syringe filter onto the syringe with the Luer connection. Ensure a tight connection by turning gently.

Start with gentle pressure to filter your sample into a vial*. This helps to assure maximum throughput.

Tips/additional information

We recommend either discarding the first 1 mL or rinsing the filter unit with 1 mL of primary solvent before sample filtration.

In order to avoid the breakage of the membrane only syringes with volumes of 10 mL or higher should be used.

Do not reuse syringe filters


Do not use at temperatures above 55 °C (131°F)

Warning: CHROMAFIL® syringe filters are intended for laboratory use only. Do not use CHROMAFIL® syringe filters for direct patient care applications.

* MACHEREY-NAGEL offers a wide range of vials and caps. More information at www.mn-net.com/vials

How to select the optimal CHROMAFIL® syringe filter

1. Filter size



2. Pore size of filter membrane

For general purpose HPLC columns packed with particles ≥ 3 μm, GC, SFC, ...

0.45 µm

Recommended for

UHPLC-, core-shell and HPLC columns, packed with particles ≤ 3 µm, GC, SFC, ...



3. Membrane type

Properties of sample	Recommended	Alternatives
Aqueous, polar hydrophilic		
low particle-load	PET	H-PTFE, MV, RC
high particle-load, prefiltration required	GF/PET	GF/RC, GF/PVDF
Mid-polar e.g. HPLC eluents		
low particle-load	PET	PA, RC
high particle-load, prefiltration required	GF/PET	GF/PA
Proteins		
low binding capacity of proteins	CA	PVDF, PES
high binding capacity of proteins	GF	GF/PVDF, GF/PET
Strong acids and bases		
low particle-load	H-PTFE	PTFE
high particle-load, prefiltration required	GF/PTFE	GF
Organic, nonpolar, hydrophob		
low particle-load	PTFE	PET
high particle-load, prefiltration required	GF/PTFE	GF/PET, GF/PVDF
Aqueous, for ion chromatography determinations	IC	

Polyester (PET)

- Hydrophilic multipurpose membrane
- For polar as well as nonpolar solvents
- The HPLC filter, especially suited for mixtures of water and organic solvents for TOC/DOC determination, not cytotoxic, does not inhibit the growth of microorganisms and higher cells
- PET filters with integrated glass fiber prefilter* (GF/PET) are recommended for solutions with a high load of particulate matter or for highly viscous solutions



	Туре	Pore size [µm]	Membrane diameter [mm]	Color code top	bottom	Standard pack filters / pack	REF	BIGbox filters/pack	REF
	Xtra PET-20/25	0.20	25	labeled		100	729221	400	729221.400
	Xtra PET-45/25	0.45	25	labeled		100	729220	400	729220.400
Ш	Xtra PET-120/25	1.2	25	labeled		100	729229	400	729229.400
	Xtra PET-20/13	0.20	13	labeled		100	729222		
	Xtra PET-45/13	0.45	13	labeled		100	729223		
	Combi Filters								
	GF/PET-20/25	1.0/0.20	25	blue	orange	100	729032	400	729032.400
	GF/PET-45/25	1.0/0.45	25	black	orange	100	729033	400	729033.400

^{*} Glass fiber exhibits a high protein-binding capacity.

Regenerated cellulose (RC)

- Hydrophilic membrane with very low adsorption
- For aqueous and organic/aqueous liquids, i.e. polar and medium polar sample solutions
- Binding capacity for proteins 84 µg per 25 mm filter
- RC filters with integrated glass fiber prefilter* (GF/RC) are recommended for solutions with a high load of particulate matter or for highly viscous solutions



Туре	Pore size [µm]	Membrane diameter [mm]	Color code top	bottom	Standard pack filters / pack	REF	BIGbox filters/pack	REF
Xtra RC-20/25	0.20	25	labeled		100	729230	400	729230.400
Xtra RC-45/25	0.45	25	labeled		100	729231	400	729231.400
Xtra RC-20/13	0.20	13	labeled		100	729236		
Xtra RC-45/13	0.45	13	labeled		100	729237		
Combi filters								
GF/RC-20/25	1.0/0.20	25	blue	blue	100	729050	400	729050.400
GF/RC-45/25	1.0/0.45	25	black	blue	100	729051	400	729051.400

^{*} Glass fiber exhibits a high protein-binding capacity.

Polytetrafluoroethylene (PTFE)

- Hydrophobic membrane
- For nonpolar liquids and gases
- Very resistant to all kinds of solvents as well as acids and bases; flushing with alcohol, followed by water, makes the originally hydrophobic membrane more hydrophilic
- PTFE filters with integrated glass fiber prefilter* (GF / PTFE) are recommended for solutions with a high load of particulate matter or for highly viscous solutions



Туре	Pore size [µm]	Membrane diameter [mm]		Standard pack filters / pack	REF	BIGbox filters/pack	REF
Xtra PTFE-20/25	0.20	25	labeled	100	729207	400	729207.400
Xtra PTFE-45/25	0.45	25	labeled	100	729205	400	729205.400
Xtra PTFE-100/25	1.00	25	labeled	100	729247		
Xtra PTFE-20/13	0.20	13	labeled	100	729208		
Xtra PTFE-45/13	0.45	13	labeled	100	729209		
PTFE-20/3	0.20	3		100	729014		
PTFE-45/3	0.45	3		100	729015		
Combi Filters							
Xtra GF/PTFE-20/25	1.0/0.20	25	labeled	100	729270		
Xtra GF/PTFE-45/25	1.0/0.45	25	labeled	100	729271		

^{*} Glass fiber exhibits a high protein-binding capacity.

Hydrophilized polytetrafluoroethylene (H-PTFE)

- Hydrophobic membrane with additional hydrophilic properties
- For polar and nonpolar sample solutions
- Resistant to all kinds of solvents as well as acids and bases



Type	Pore size [µm]	Membrane diameter [mm]		Standard pack filters/pack	REF	BIGbox filters/pack	REF
Xtra H-PTFE-20/25	0.20	25	labeled	100	729245		
Xtra H-PTFE-45/25	0.45	25	labeled	100	729246	400	729246.400
Xtra H-PTFE-20/13	0.20	13	labeled	100	729256		
Xtra H-PTFE-20/13	0.45	13	labeled	100	729257		

Cellulose mixed ester (MV)

- Hydrophilic membrane with very low adsorption
- For aqueous or polar solutions





	Туре	Pore size [µm]	Membrane diameter [mm]		Standard pack filters/pack	REF	BIGbox filters/pack	REF
	Xtra MV-20/25	0.20	25	labeled	100	729206		
	Xtra MV-45/25	0.45	25	labeled	100	729204	400	729204.400

Cellulose acetate (CA)

- Hydrophilic membrane
- For the filtration of water-soluble oligomers and polymers, especially suited for biological macromolecules
- Very high stability in aqueous solutions
- Binding capacity for proteins 21 μg per 25 mm filter
- Also available in a sterile package (S) for filtration under sterile conditions (each filter individually sealed)



	Туре	Pore size [µm]	Membrane diameter [mm]	Color code top	bottom	Standard pack filters / pack	REF	BIGbox filters/pack	REF
	Xtra CA-20/25	0.20	25	labeled		100	729226	400	729226.400
	Xtra CA-45/25	0.45	25	labeled		100	729227	400	729227.400
	Xtra CA-20/13	0.20	13	labeled		100	729254		
	Xtra CA-45/13	0.45	13	labeled		100	729255		
	Sterile filters								
	CA-20/25 (S)	0.20	25	yellow	red	50	729024		
	CA-45/25 (S)	0.45	25	colorless	red	50	729025		

Polyamide (PA) = Nylon

- Moderately hydrophilic membrane
- For aqueous and organic/aqueous medium polar liquids
- PA filters with integrated glass fiber prefilter* (GF/PA) are recommended for solutions with a high load of particulate matter or for highly viscous solutions



	Туре	Pore size [µm]	Membrane diameter [mm]		Standard pack filters / pack	REF	BIGbox filters/pack	REF
	Xtra PA-20/25	0.20	25	labeled	100	729212	400	729212.400
	Xtra PA-45/25	0.45	25	labeled	100	729213	400	729213.400
	Xtra PA-20/13	0.20	13	labeled	100	729248		
	Xtra PA-45/13	0.45	13	labeled	100	729249		
	PA-20/3	0.20	3		100	729010		
	PA-45/3	0.45	3		100	729011		
7	Combi Filters							
	Xtra GF/PA-20/25	1.0/0.20	25	labeled	100	729260		
	Xtra GF/PA-45/25	1.0/0.45	25	labeled	100	729261		

Polyethersulfone (PES)

- Hydrophilic membrane
- For aqueous and slightly organic liquids with higher flow rates
- Very low adsorption for pharmaceuticals and proteins
- Good stability against organic acids and bases
- Binding capacity for proteins 29 µg per 25 mm filter



	Туре	Pore size [µm]	Membrane diameter [mm]		Standard pack filters / pack	REF	BIGbox filters/pack	REF
	Xtra PES-20/25	0.20	25	labeled	100	729240		
	Xtra PES-45/25	0.45	25	labeled	100	729241	400	729241.400
Ш	Xtra PES-500/25	5.0	25	labeled	100	729242		

Polyvinylidene difluoride (PVDF)

- Hydrophilic membrane
- For 100 % aqueous solutions, water-soluble oligomers and polymers like proteins
- Low binding capacity for proteins 20 µg per 25 mm filter
- PVDF filters with integrated glass fiber prefilter* (GF/PVDF) are recommended for the filtration of biological samples with high particle loads




	Туре	Pore size [µm]	Membrane diameter [mm]	Color code top	bottom	Standard pack filters / pack	REF	BIGbox filters/pack	REF
	Xtra PVDF-20/25	0.20	25	labeled		100	729218	400	729218.400
	Xtra PVDF-45/25	0.45	25	labeled		100	729219	400	729219.400
-	Xtra PVDF-20/13	0.20	13	labeled		100	729243		
	Xtra PVDF-45/13	0.45	13	labeled		100	729244		
	Combi filters								
	GF/PVDF-45/25	1.0/0.45	25	black	white	100	729039	400	729039.400

^{*} Glass fiber exhibits a high protein-binding capacity.

Glass fiber (GF)

- Inert filter, nominal pore size 1 µm, allows higher flow rates than small pore filters
- For solutions with high loads of particulate matter or for highly viscous solutions (e.g. soil samples, fermentation broths)
- As prefilters for other CHROMAFIL[®] filters, they prevent the clogging of the membrane
- High binding capacity for proteins





Ion chromatography (IC)

- Special filter for ion chromatography
- For the filtration of aqueous liquids
- For optimal results with blind values < 5 ppb we recommend to pre-wash the filter with de-ionized water



Туре	Pore size [µm]	Membrane diameter [mm]		Standard pack filters / pack	REF
Xtra IC-45/25	0.45	25	labeled	100	729258

Polypropylene vials for ion chromatography With various interior shapes and volumes

Ideal protection for your samples

- 1.5 mL PP screw neck vials N 9 (702500) with closure 702288.1 for Dionex® instruments
- 0.3 mL transparent/amber PP vials, inner cone (screw N 9: 702009/702172; snap ring/crimp neck N 11: 702809/702173)
- 0.7 mL transparent PP vials, round bottom insert (702010, screw N 9/702174 snap ring/crimp neck vials N 11)



Filtration cartridges

- Filtration cartridges for sample clarification under vacuum (e.g., using the CHROMABOND® vacuum manifold or SPE automation systems like Gilson Aspec™, Rapidtrace) or by gravity flow.
- Cartridge sizes 3 mL and 6 mL
- Different membranes (PET, RC, PTFE, PVDF, GF) and pore sizes (0.20, 0.45 and 1.0 µm). The membrane materials correspond to the respective CHROMAFIL® syringe filters.



	Туре	Pore size [µm]	Pack of	Column volume REF 3 mL	REF 6 mL
	PET (polyester)	0.20	100	730578.320	730578.620
	PET (polyester)	0.45	100	730578.345	730578.645
	RC (regenerated cellulose)	0.20	100	730068.320	730068.620
	RC (regenerated cellulose)	0.45	100	730068.345	730068.645
U	PTFE (polytetrafluoroethylene)	0.20	100	730570.320	730570.620
	PTFE (polytetrafluoroethylene)	0.45	100	730570.345	730570.645
	PVDF (polyvinylidene difluoride)	0.20	100	730579.320	730579.620
	PVDF (polyvinylidene difluoride)	0.45	100	730579.345	730579.645
	GF (glass fiber)	nom. 1.0	100	730517.3100	730517.6100

MULTI 96 filter plates

- 96-well polypropylene plates for the simultaneous filtration of 96 samples
- Advantages of this high-throughput system:
 Economical by saving time and solvent
 Use of multi-channel pipettors facilitates liquid transfer steps
 Readily adaptable to all common automated/robotic handling systems
 Minimized dead volume (= 40 µL)
- Membrane materials correspond to the respective CHROMAFIL® syringe filters





Description	Pack of	REF
Filter plates with cellulose mixed ester filter elements (0.20 µm)	1	738770.M
Filter plates with cellulose mixed ester filter elements (0.45 µm)	1	738771.M
Filter plates with RC filter elements (regenerated cellulose, 0.20 µm)	1	738656.M
Filter plates with RC filter elements (regenerated cellulose, 0.45 µm)	1	738657.M
Filter plates with PTFE filter elements (0.20 µm)	1	738660.M
Filter plates with PTFE filter elements (0.45 µm)	1	738661.M
Filter plates with PTFE filter elements (1.0 µm)	1	738662.M
Filter plates with PTFE filter elements (3.0 µm)	1	738663.M
Filter plates with PE filter elements (40–100 μm)	1	738659.M
Filter plates with glass fiber filter elements (nominal 1 µm)	1	738655.2M
Filter plates with glass fiber filter elements (nominal 3 µm)	1	738658.M

CHROMAFIL® - Disposable filters from MACHEREY-NAGEL















SPE and Flash



Vials and caps

Germany

MACHEREY-NAGEL GmbH & Co. KG Valencienner Str. 11

52355 Düren · Germany Tel.: +49 24 21 969-0 info@mn-net.com

France

MACHEREY-NAGEL SAS 1, rue Gutenberg - BP135 67722 Hoerdt · France Tél.: +33 388 68 22 68

sales-fr@mn-net.com

MACHEREY-NAGEL SAS (Société par Actions Simplifi ée) au capital de 186600 € Siret 379 859 531 00020 · RCS Strasbourg B379859531

Nº intracommunautaire FR04 379 859 531

Switzerland

MACHEREY-NAGEL AG Hirsackerstr. 7 4702 Oensingen · Switzerland Tel.: +41 62 388 55 00

sales-ch@mn-net.com

USA

MACHEREY-NAGEL Inc. 924 Marcon Blvd. Suite 102 Allentown PA. 18109 · USA Toll-free: 888 321 6224 (MACH) sales-us@mn-net.com

Registered and common law trademarks

ASPEC® Gilson Inc. (USA)

CHROMAFII® MACHEREY-NAGEL GmbH & Co. KG (Germany) CHROMABOND® MACHEREY-NAGEL GmbH & Co. KG (Germany)

Dionex® Dionex Corporation SOTAX® Sotax AG (Switzerland)

Benchmate™ Caliper Life Sciences Inc. (USA)

www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG DE Tel.: +49 24 21 969-0 info@mn-net.com Valencienner Str. 11 52355 Düren · Germany

CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com

FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com US Tel.: +1 888 321 62 24 sales-us@mn-net.com

