# MACHEREY-NAGEL

# CHROMABOND® HLB



- Hydrophilic-lipophilic balanced phase
- Enhanced retention for polar compounds
- Consistent recoveries



#### Introduction

The enrichment of hydrophilic analytes from polar matrices is a challenge because they often cannot be retained by  $C_{18}/C_8$ -modified silica or hydrophobic polymer adsorbents. To overcome this issue MACHEREY-NAGEL developed CHROMABOND® HLB, a hydrophilic-lipophilic balanced N-vinylpyrrolidone-divinylbenzene copolymer. Its lipophilic backbone interacts with nonpolar hydrocarbon residues of the analytes while the linked hydrophilic groups interact with polar functional groups to provide enhanced retention.


#### Typical applications

Polar organic molecules from polar matrices e.g.,

- Sulfonamides
- Pesticides
- Chloramphenicol
- lodinated contrast media

## Advantages of CHROMABOND® HLB

- Applicable for a wide range of analyte polarities
- Enhanced retention for polar compounds
- High loadability and outstanding performance
- Water wettable even if bed runs dry, SPE can be continued
- The alternative to Oasis® HLB



#### Technical data

Hydrophilic-lipophilic balanced N-vinylpyrrolidone-divinylbenzene copolymer

Particle shape: spherical pH stability: 1–14

Particle size: 60 µm and 30 µm

Pore size: 65 Å Specific surface: 750 m²/g

# Standard SPE protocol (subsequent HPLC analysis) MN Appl. No. 306300

Column: CHROMABOND® HLB, 3 mL, 200 mg

MN REF: 730924

Column conditioning: 5 mL methanol, 5 mL dest. water

Sample application: slowly aspirate sample through column

Washing: 5 mL dest. water

Drying: 10 min with applied vacuum

Elution: 8 mL methanol Evaporation: under nitrogen

Reconstitution: in 1 mL dest. water + 0.1 % formic acid

# Standard SPE protocol (subsequent GC analysis) MN Appl. No. 306310

CHROMABOND® HLB, 3 mL, 200 mg

....

MN REF: 730924

Column conditioning: 5 mL solvent (e.g., ethyl acetate), 5 mL

methanol, 5 mL dest. water

Sample application: slowly aspirate sample through column

Washing: 5 mL dest. water

Drying: 10 min with applied vacuum

Elution: solvent 1) (typical solvents: ethyl acetate,

MTBE, methylene chloride)

Evaporation: under nitrogen, dry with sodium sulfate 2),

adjust to final volume

1) usually nonpolar, therefore often 10 % methanol are added

2) e.g., with CHROMAFIX® Dry



#### Pharmaceuticals from serum

MN Appl. No. 306510

Columns\*: CHROMABOND® HLB, 1 mL, 30 mg

Oasis® HLB, 1 mL, 30 mg

MN REF: 730921

Column conditioning: 1 mL methanol, 1 mL dest. water

Sample application: 1 mL serum (spiked with 50 ng of each

analyte)

Washing: 1 mL dest. water

Drying: 10 min with applied vacuum

Elution: 2 mL methanol Evaporation: under nitrogen, 40 °C

Reconstitution: in 1 mL dest. water – acetonitrile (95:5, v/v)

#### Recovery rates ± RSD [%], n = 5

Compound	CHROMABOND® HLB	Oasis® HLB
Amitriptyline	77.7 ± 1.2	26.9 ± 0.8
Atenolol	71.2 ± 1.3	71.1 ± 1.2
Atropine	$84.8 \pm 0.9$	80.5 ± 1.4
Carbamazepine	97.7 ± 0.3	57.7 ± 4.4
Chlorpheniramine	85.9 ± 2.2	91.7 ± 1.2
Clomipramine	73.6 ± 6.7	48.1 ± 0.8
Diphenhydramine	88.3 ± 2.1	94.7 ± 1.3
Indapamide	87.7 ± 3.3	$49.0 \pm 2.0$
Ketamine	90.8 ± 1.8	$88.4 \pm 2.0$
Ketoprofen	84.1 ± 3.1	48.9 ± 1.6

#### Further analysis: LC-MS/MS, according to MN Appl. No. 128200

Column: EC 50/2 NUCLEOSHELL® PFP, 2.7 μm

MN REF: 763532.20

Eluent: A: dest. water + 0.1 % formic acid

B: acetonitrile + 0.1 % formic acid 5–95 % B in 7.5 min, 95 % B for 1 min, 95–5 % B in 0.5 min, 5 % B for 5 min

Flow rate: 0.3 mL/min Temperature:  $30 \,^{\circ}\text{C}$ 

Detection: MS, Selected Reaction Monitoring (SRM)

Injection: 5 µL

Compound	CHROMABOND® HLB	Oasis <sup>®</sup> HLB
Nortriptyline	$76.6 \pm 2.1$	$14.9 \pm 3.5$
Propanolol	107.7 ± 1.4	108.3 ± 1.7
Sulfachloropyridazine	85.8 ± 1.6	84.0 ± 1.4
Sulfadoxine	99.8 ± 2.0	91.2 ± 2.0
Sulfamethoxazole	94.3 ± 1.6	81.2 ± 1.6
Sulfapyridine	64.6 ± 1.8	61.6 ± 3.9
Sulfaquinoxaline	127.1 ± 3.4	104.8 ± 2.8
Sulfamerazine	67.3 ± 0.8	63.4 ± 3.7
Trimipramine	81.5 ± 2.3	37.3 ± 1.2
Verapamil	107.5 ± 1.7	48.9 ± 0.9

# Drugs from tab water

MN Appl. No. 306330

Columns\*: CHROMABOND® HLB, 3 mL, 200 mg

Oasis® HLB, 3 mL, 200 mg

MN REF: 730924

Column conditioning: 5 mL methanol, 5 mL dest. water

Sample application: 1000 mL tap water (spiked with 5  $\mu g/L$  of

each analyte), ~10 mL/min

Washing: 5 mL dest. water

Drying: 10 min with applied vacuum

Elution: 8 mL methanol
Evaporation: under nitrogen, 40 °C

Reconstitution: in 1 mL dest. water + 0.1 % formic acid

# Further analysis: HPLC, according to MN Appl. No. 128110 see Drugs from serum, page 7

#### Recovery rates $\pm$ RSD [%], n = 5

Compound	CHROMABOND® HLB	Oasis® HLB
Azidothymidine	$98.2 \pm 0.7$	$96.9 \pm 0.8$
Caffeine	84.2 ± 0.9	67.9 ± 0.9
trans-doxepin	$78.5 \pm 0.9$	64.5 ± 0.9
cis-doxepin	81.3 ± 1.3	62.7 ± 0.6
Propanolol	93.7 ± 0.9	83.7 ± 1.2
Protriptyline	78.2 ± 2.9	60.2 ± 2.3

#### CHROMABOND® HLB for polar analytes

CHROMABOND® HLB provides high recovery rates for polar drugs from water.



#### Sulfa drugs from serum

MN Appl. No. 306340

CHROMABOND® HLB, 1 mL, 30 mg Columns\*:

Oasis® HLB, 1 mL, 30 mg

MN RFF: 730921

Column conditioning: 1 mL methanol, 1 mL dest. water

Sample application: 1 mL serum (spiked with 10 µg/mL of each

analyte)

Washing: 1 mL dest. water

Drying: 10 min with applied vacuum

Elution: 2 mL methanol under nitrogen, 40 °C Evaporation:

Reconstitution: in 1 mL dest, water + 0.1 % formic acid

#### Equivalence to Oasis® HLB

CHROMABOND® HLB shows equivalent recovery rates to

Oasis® HLB for the three tested sulfa drugs.

#### Further analysis: HPLC, according to MN Appl. No. 128130

Column: EC 150/2 NUCLEODUR® C<sub>18</sub> Pyramid,

3 um

MN REF: 760261.20

Eluent: dest. water + 0.1 % formic acid - methanol

+ 0.1 % formic acid (85:15, v/v), 5 min

Flow rate: 0.6 mL/min Temperature: 25 °C Detection UV. 254 nm Injection: 5 μL

Recovery rates  $\pm$  RSD [%], n = 5

Compound	CHROMABOND® HLB	Oasis® HLB	
Compound	CHROMABOND® HLB	Oasis <sup>®</sup> HLB	
Sulfadiazine	$97.3 \pm 2.9$	$92.0 \pm 3.8$	
Sulfamerazine	94.4 ± 1.8	92.8 ± 1.6	
Sulfathiazole	90.3 ± 2.9	89.6 ± 1.5	

#### Chloramphenicol from honey

MN Appl. No. 306350

Columns\*: CHROMABOND® HLB, 3 mL, 200 mg

Oasis® HLB, 3 mL, 200 mg

MN REF: 730924

Sample preparation: Weigh out 5 g of honey. Add 4 mL water

> and shake rigorously for 30 sec. Spike with 1 mL standard solution (c = 5 ng/mL in methanol) and shake rigorously for 30 sec. Add 15 mL ethyl acetate and shake rigorously for 30 sec. Centrifuge at 3000 rpm for 10 min. Take 12 mL of supernantant for eluent exchange. Evaporate extracts to dryness at 40 °C under a stream of nitrogen.

Redissolve residue in 10 mL water.

3 mL methanol (dispensing speed 1 mL/min), Column conditioning:

5 mL dest. water (disp. speed 1 mL/min)

9 mL water sample (disp. speed 3 mL/min Sample application:

over sample loop)

Washing: 10 mL dest. water (disp. speed 3 mL/min) Drying: 100 mL air (disp. speed 100 mL/min) Elution: 5 mL ethyl acetate - methanol (80:20, v/v) Drying: 100 mL air (disp. speed 100 mL/min)

under nitrogen, 40 °C Evaporation:

Reconstitution: in 1 mL dest. water - acetonitrile (95:5, v/v)

The SPE application was performed with a FREESTYLE® SPE

automation system.

#### Further analysis: LC-MS/MS, according to MN Appl. No. 128140

Column: EC 150/2 NUCLEODUR® π<sup>2</sup>. 5 μm

MN REF: 760624.20 A: dest. water Eluent: B: acetonitrile

> 5-95 % B in 7.5 min, 95 % B for 1 min, 95-5 % B in 1 min, 5 % B for 5 min

Flow rate: 0.3 mL/min Temperature: 35 °C

MS, Selected Reaction Monitoring (SRM) Detection:

Injection: 5 μL

#### Recovery rates $\pm$ RSD [%], n = 5

Compound	CHROMABOND® HLB	Oasis® HLB	
Chloramphenicol-d5	$90.9 \pm 5.4$	$90.0 \pm 9.3$	

#### Good to know

Antibiotics and pesticides contamination of agricultural products such as honey has been an issue in the recent years and resulted in stricter guidelines in food safety control.



<sup>\*</sup>Same conditions for all used columns. Due to a better comparability CHROMABOND® HLB and Oasis® HLB adsorbents (60 µm) were packed into equal column hardware. The shown chromatograms may not be representative of other applications

#### Pesticides from tap water

MN Appl. No. 306360

Columns\*: CHROMABOND® HLB, 3 mL, 200 mg

Oasis® HLB, 3 mL, 200 mg

MN REF: 730924

Column conditioning: 5 mL methanol, 5 mL dest. water

Sample application: 1000 mL tap water (spiked with 50 ng of

each analyte)

Washing: 10 mL dest. water

Drying: 5 min with applied vacuum ( -15 psi)

Elution: 6 mL acetonitrile
Evaporation: under nitrogen, 40 °C

Reconstitution: in 1 mL dest. water - acetonitrile (95:5, v/v)

#### Recovery rates $\pm$ RSD [%], n = 5

Compound	CHROMABOND® HLB	Oasis <sup>®</sup> HLB
Acetamiprid	73.3 ± 5.0	112.1 ± 9.9
Atrazine	110.3 ± 17.8	114.0 ± 11.6
Azoxystrobin	74.7 ± 5.4	98.1 ± 10.8
Carbaryl	65.7 ± 5.4	69.1 ± 7.1
Chlorotoluron	82.7 ± 5.7	101.2 ± 3.8
Chlorpyrifos	50.3 ± 5.4	47.0 ± 3.7
Clofentezine	27.8 ± 2.7	21.4 ± 3.7
Clothianidin	69.4 ± 6.5	$52.9 \pm 2.9$
Coumaphos	69.8 ± 4.8	82.3 ± 5.2
Cyanazine	99.8 ± 9.3	85.1 ± 7.2
Desethylatrazine	94.8 ± 15.1	87.4 ± 11.4
Desisopropylatrazine	92.5 ± 7.6	N/A
Diazinon	71.5 ± 7.9	73.3 ± 4.7
Difenoconazole	83.9 ± 6.5	28.8 ± 5.0
Diuron	70.0 ± 4.8	80.1 ± 8.4
Ethoprophos	72.4 ± 9.3	85.4 ± 7.2
Hexazinone	88.4 ± 7.7	104.3 ± 7.4

#### Further analysis: LC-MS/MS, according to MN Appl. No. 128150

Column: EC 50/2 NUCLEOSHELL® PFP, 2.7 μm

MN REF: 763532.20

Eluent: A: dest. water + 0.1 % formic acid

B: acetonitrile + 0.1 % formic acid 5–95 % B in 15 min, 95 % B for 5 min, 95–5 % B in 1 min, 5 % B for 9 min

Flow rate: 0.3 mL/min Temperature: 40  $^{\circ}$ C

Detection: MS, Selected Reaction Monitoring (SRM)

Injection: 5 µL

Compound	CHROMABOND® HLB	Oasis® HLB
Imazalil	27.3 ± 15.7	N/A
Imidacloprid	93.4 ± 5.1	$40.3 \pm 5.2$
Isoproturon	100.2 ± 4.2	$102.8 \pm 13.0$
Linuron	84.5 ± 7.6	88.3 ± 9.5
Methabenzthiazuron	72.5 ± 5.3	$48.0 \pm 3.7$
Methomyl	78.8 ± 5.4	83.6 ± 5.6
Metobromuron	73.8 ± 5.6	85.6 ± 9.3
Metolachlor	79.0 ± 5.2	89.2 ± 5.0
Monolinuron	75.4 ± 6.2	97.9 ± 7.2
Myclobutanil	101.8 ± 11.4	88.7 ± 14.5
Phosalone	63.8 ± 7.7	74.0 ± 4.0
Piperonylbutoxide	101.4 ± 8.6	$99.7 \pm 7.9$
Propazine	102.1 ± 13.6	90.9 ± 9.4
Propyzamide	84.8 ± 7.1	86.4 ± 10.6
Terbuthylazine	107.9 ± 13.3	100.0 ± 13.6
Thiacloprid	74.1 ± 6.3	86.5 ± 10.8



<sup>\*</sup>Same conditions for all used columns. Due to a better comparability CHROMABOND® HLB and Oasis® HLB adsorbents (60 µm) were packed into equal column hardware. The shown chromatograms may not be representative of other applications.

#### lodinated contrast media from serum

#### MN Appl. No. 306370

Columns\*: CHROMABOND® HLB, 1 mL, 30 mg

Oasis® HLB, 1 mL, 30 mg

MN REF: 730921

Column conditioning: 1 mL methanol, 1 mL dest. water

Sample application: 1 mL serum (spiked with 10  $\mu$ g/mL of each

analyte)

Washing: 1 mL dest. water


Drying: 10 min with applied vacuum

Elution: 2 mL methanol Evaporation: under nitrogen, 40 °C

Reconstitution: in 1 mL dest. water + 5 mM ammonium

formate

#### Recovery rates $\pm$ RSD [%], n = 5



#### Further analysis: HPLC, according to MN Appl. No. 128160

Column: EC 150/2 NUCLEODUR® C<sub>18</sub> Gravity-SB,

3 µm

MN REF: 760608.20

Eluent: A: dest. water + 5 mM ammonium formate

B: methanol/acetonitrile, (1:2, v/v) + 5 mM

ammonium formate

7-10 % B in 5 min, 10-100 % B in 5 min,

100 % B for 10 min

Flow rate: 0.2 mL/min
Temperature: 40 °C
Detection: UV, 254 nm
Injection: 5 µL

#### Superior to Oasis® HLB

CHROMABOND® HLB provides higher recovery rates for the five given analytes from serum in comparison to Oasis® HLB.

■ CHROMABOND® HLB

■ Oasis® HLB

#### Tetracyclines and alkaloids from serum at pH 5

#### MN Appl. No. 306380

Columns\*: CHROMABOND® HLB, 1 mL, 30 mg

Oasis® HLB, 1 mL, 30 mg

MN REF: 730921

Column conditioning: 1 mL methanol, 1 mL dest. water

Sample application: 1 mL serum pH 5, adjusted with formic acid

(spiked with 20 µg/mL of each analyte)

Washing: 1 mL dest. water

Drying: 10 min with applied vacuum

Elution: 2 mL methanol Evaporation: under nitrogen, 40 °C

Reconstitution: in 1 mL dest. water + 0.1 % formic acid

#### Further analysis: HPLC, according to MN Appl. No. 128170

Column: EC 50/2 NUCLEOSHELL® RP 18plus,

2.7 µm

MN REF: 763232.20

Eluent: A: dest. water + 0.1 % formic acid

B: acetonitrile + 0.1 % formic acid 2–60 % B in 4 min, 60 % B for 1 min, 60–2 % B in 0.5 min, 2 % B for 3 min

Flow rate: 0.75 mL/min Temperature: 22 °C Detection: UV, 330 nm Injection:  $5 \mu L$ 

#### Recovery rates ± RSD [%], n = 4

Compound	CHROMABOND® HLB	Oasis <sup>®</sup> HLB
Berberine	$85.4 \pm 0.3$	$82.5 \pm 0.6$
Chlortetracycline	72.1 ± 1.4	66.3 ± 2.8
Hydrastine	88.9 ± 2.6	99.3 ± 5.7
Oxytetracycline	82.3 ± 1.4	78.7 ± 1.4
Tetracycline	78.1 ± 1.4	70.7 ± 2.6

<sup>\*</sup>Same conditions for all used columns. Due to a better comparability CHROMABOND® HLB and Oasis® HLB adsorbents (60 µm) were packed into equal column hardware. The shown chromatograms may not be representative of other applications.

#### Drugs from serum MN Appl. No. 306320

Columns\*: CHROMABOND® HLB, 1 mL, 30 mg Hydrophobic polymer, 1 mL, 30 mg

MN REF: 730921

Column conditioning: 1 mL methanol, 1 mL dest. water

Sample application: 1 mL serum (spiked with 10 µg/mL of each

analyte)

Washing: 1 mL dest. water

Drying: 10 min with applied vacuum

Elution: 2 mL methanol Evaporation: under nitrogen, 40 °C

Reconstitution: in 1 mL dest. water + 0.1 % formic acid

### Further analysis: HPLC, according to MN Appl. No. 128110

Column: EC 150/2 NUCLEODUR® C<sub>18</sub> Pyramid,

3 µm

MN REF: 760261.20


Eluent: A: dest. water + 0.1 % formic acid

B: methanol + 0.1 % formic acid

30-45 % B in 15 min

Flow rate: 0.3 mL/min
Temperature: 30 °C
Detection: UV, 254 nm
Injection: 5 µL

#### Recovery rates ± RSD [%], n = 5



#### Good to know

No conditioning step is needed when using CHROMABOND® HLB due to its excellent water wettability. Remark: Conditioning often improves analyte recovery. Therefore we recommend comparing results.

- CHROMABOND® HLB (conditioned)
- CHROMABOND® HLB (dry)
- Hydrophobic polymer (conditioned)
- Hydrophobic polymer (dry)



<sup>\*</sup>Same conditions for all used columns. Due to a better comparability CHROMABOND® HLB and Oasis® HLB adsorbents (60 µm) were packed into equal column hardware. The shown chromatograms may not be representative of other applications.

## Ordering information

Volume	Adsorbent weight →							Pack of
	30 mg	60 mg	100 mg	150 mg	200 mg	500 mg	1 g	
CHROMABO	ND® HLB polyprop	ylene columns (6	60 μm)					
1 mL	730921		730922					30
3 mL		730923			730924	730925		30
6 mL				730944	730926	730927		30
15 mL						730928	730929	20
CHROMABO	ND® HLB polyprop	ylene columns (6	60 μm) · BIGpacks					'
3 mL		730923.250			730924.250			250
6 mL					730926.250	730927.250		250
CHROMABO	ND® HLB polyprop	ylene columns (3	80 μm)					'
1 mL	730921P30		730922P30					30
3 mL		730923P30			730924P30			30
6 mL				730944P30				30
CHROMABO	ND® LV-HLB (30 µr	n)						
15 mL	732140	732141						30

Size →	S	М	L	Pack of
Minimum adsorbent weight →	50 mg	120 mg	350 mg	
CHROMAFIX® HLB cartridges (60 μm)				
	731921	731922	731923	50
Adsorbent weight →	96 x 10 mg	96 x 30 mg	96 x 60 mg	
CHROMABOND® MULTI 96 HLB (60 µm)				
			738920.060M	1
CHROMABOND® MULTI 96 HLB (30 µm)				
	738921.010M	738921.030M		1

#### Registered trademarks

Oasis® Waters Corp. (USA)

CHROMABOND® MACHEREY-NAGEL GmbH & Co. KG (Germany) CHROMAFIX® MACHEREY-NAGEL GmbH & Co. KG (Germany)

FREESTYLE® LCTech GmbH (Germany)

# Your local distributor

www.mn-net.com

## **MACHEREY-NAGEL**



MACHEREY-NAGEL GmbH & Co. KG 52355 Düren · Deutschland

DE Tel.: +49 24 21 969-0

CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com

info@mn-net.com US Tel.: +1 888 321 62 24 sales-us@mn-net.com

