# Spe-ed<sup>TM</sup> Helix





www.appliedseparations.com

## Helix



## Basic Requirements for All Extractions

All extractions require at least these Applied Separations, Inc. components:

#### Base unit

CO<sub>2</sub> Pump Pressure vessel assembly Recirculating bath (chiller) and these utilities:

A source of air delivered at 7 BAR Electrical power: 240v or 120v Source of liquid CO<sub>2</sub>

The Helix is made up of several "base" components. The basic components are put together in a variety of standard or custom configurations to make a unit to perform a specific function.

With the base system you will be able to use the same components to do separations and extractions as well as make nanoparticles, but not at the same time. This means on one day you can do extractions and on another day you can make nanoparticles.

#### **Base Unit**

The compact Base Unit, measuring 10"w X 16"d X 34"h is the starting platform for operations. Pressure vessel assemblies up to 1 liter are placed on the base unit's shelf. Input, output and vent lines are controlled by shutoff valves located on the front of the unit. Digital temperature and pressure indicators also on the front of the unit show pressure and temperature.

Pressure vessels are heated by specially designed band heaters which are plugged into the front of the base unit making for easy access. An additional  $CO_2$  preheater is employed to ensure that the  $CO_2$  is at the designed temperature before entering the pressure vessel. A back pressure regulator controls the flow of gaseous  $CO_2$  if exiting to ambient collection or regulating the pressure in a downstream pressure vessel (e.g. cyclone separator).

A base unit with its vessel assembly can be linked to other base units for additional processing capabilities: cyclone separators, precipitation vessel, expansion vessel, etc.



## Laboratory System

## Helix

## **Helix Configuration Options**

menk connigeranon opnom	
Basic Configuration	#7409
Helix SCF Base Unit 240v	#7305
Touchpad Controller and	
Standard CO <sub>2</sub> Pump Module 240v	#7401
1 Liter Vessel Assembly with basket	#7322
500mL Vessel Assembly	#7323
300mL Vessel Assembly	#7324
100mL Vessel Assembly	#7329
50mL Vessel Assembly	#6414
32mL Vessel Assembly	#6413
24mL Vessel Assembly	#6412
Standard Flow Meter	#7927
Standard Collector Vessel	#7735

## **Basic Configuration Options**

Modifier/Liquid Pump, Helix 240v	#7172
Recirculating Bath, Helix 240v	#7027
800mL CO, Pump Module 240v	#7316
Stirring Assembly - stirrer, controller	#7320

## Variety of Flow Meters

CO <sub>2</sub> Recycle Module	
Chiller, Level Indicator, etc.	#7399
	#7027



Vessel with Stirrer Assembly



Collection



Pressure Relief Device



Touchpad System Control



**RESS** collector



Basic Helix system with the separator module.



## **Laboratory System**

Helix


 $CO_2$ In

Light Port

Liquid Reactant In

Emergency Pressure Relief Device\*

Endoscope UV/IR Particle Sizing Turbidity



Stirrer

Pressure Probe

Gas Reactant In

- $H_2$
- O,
- CO<sub>2</sub>
- $NH_3^2$ , etc.

Reactor Vessel with Agitator

Electric In/Out e.g. Ultrasound

Out

Temperature Probe



Endoscope



UV/IR



Monitor Meters

• Spectrophotometer



Stirrer Shaft/Impeller

- Variable length
- Interchangable impellers



Close-up of the top and bottom vessel ports.

\*Patent Pending



## **Laboratory System**

## Helix

### Cryo Cooler

Low pressure/ambient aluminum containers are available to collect a variety of extracts. Should the extract be highly viscous heat can be applied, or if volatile, chilling can be applied.



### **Pump Module**

There are two standard CO<sub>2</sub> pumps (#7321 and #7316). These are air-driven, liquid pumps capable of delivering CO<sub>2</sub> from a cylinder (nominally 60 BAR at room temperature) to a pressure of 690 BAR. Because they are pneumatic, they are inherently more compact, safer, cleaner, quieter and requiring less maintenance than either electric or hydraulic pumps. In both pumps, the pressure is set and shown by a digital readout. There is another gauge to show the air pressure.

### **Reciruclating Bath**



The *Spe-ed* RCB for Helix (#7029) (820 BTU/hr, 15LPM, -10 C to +40C) is a specifically designed recirculating cooling bath that chills the CO<sub>2</sub> to liquidity. It is microprocessor controlled, with a small footprint and nearly noise free.

#### Stirrer

Applied Separations now offers new stirrers to go into their extraction/ reaction vessels via the 5-Port cap that allows



access to the inside of the vessel during your process.

The stirrers are rated to 10,000 PSI, 650°F, and go up to 3,000 rpm. The instrument control panel includes a digital display to monitor rpm. Stirrers are available with different shaft lengths, with a wide variety of impellers available to attach to the bottom of the stirrer.

### Modifier/Liquid Pump

Liquid pumps may be necessary in several operations when using the Helix: adding polar modifiers, introducing solvents during PCA and for the operation of the countercurrent column.

Spe-ed MAX (#7298) is a microprocessor controlled pump delivering 690 BAR at adjustable flow rates up to 24 ml/min.



### **Pressure Vessel Assemblies**

316 Stainless steel pressure vessels for the Helix are hand-tightenet range in size from 24ml to 1,000ml. The assembly is comprised of the pressure vessel, heating elements, electrical input cable, and insulation.

#### Standard sizes are

Standa	ard sizes are	
7972	5mL Vessel	1" O.D. x 5.125" O.L390" I.D. x 2.24" I.L.
7972	10mL Vessel	1" O.D. x 5.125" O.L560" I.D. x 2.24" I.L.
7973	24mL Vessel	1" O.D. x 8.875" O.L560" I.D. x 5.9" I.L.
7974	32mL Vessel	1" O.D. x 10.5" O.L560" I.D. x 8" I.L.
7975	50mL Vessel	1" O.D. x 15.25" O.L560" I.D. x 12.72" I.L.
7329	100mL Vessel	2.25" O.D. x 9.57" O.L. 1.25 " I.D. x 4.97" I.L.
7324	300mL Vessel	3.5" O.D. x 11.42" O.L. 2" I.D. x 5.87" I.L.
7323	500mL Vessel	4.75" O.D. x 9.49" O.L. 3" I.D. x 4.49" I.L.
7322	1000mL Vessel	4.75" O.D. x 13.63" O.L. 3" I.D. x 8.62" I.L.



