


# An introduction to the Centri portfolio

A breakthrough in automated sample extraction and enrichment




# What is Centri®?

A breakthrough in automated sample extraction and enrichment for GC(–MS)



# The Centri family



Centri 90




Centri 180



Centri 360


# Analyte preconcentration

## Step 1: Desorption of sampled analytes and focusing of volatiles



# GC injection

## Step 2: Rapid backflush desorption and transfer/injection to the GC



# Types of sampling techniques

Compatible on Centri

## Automated High Capacity Sorptive Extraction

Robust immersive or headspace sampling from liquids and solids using HiSorb probes.



## Static Headspace

Versatile sampling from solids and liquids contained in regular headspace vials



## Tube-based TD


Adaptable sampling onto sorbent packed TD tubes or direct desorption of samples from empty tubes



Centri 380

## SPME

Fast and sensitive sample extraction, with a range of fiber/arrow types offering analyte selectivity



## Liquid Injection


Direct liquid injection to GC compatibility – rail and configuration dependent.

Note: Does not use Centri Injector



# HiSorb, high-capacity sorptive extraction

- Similar extraction principles as solid-phase microextraction (SPME)
- Analytes are extracted via an **absorptive or adsorptive process** depending on the phase type used.
- Intrinsically **robust** in design due to inert metal-core probes supporting the sorptive phase (30-65  $\mu\text{L}$  phase vs.  $\sim 0.5 \mu\text{L}$  on SPME fiber)
- Can be used in **headspace** mode for both solid and liquid samples and robust **immersive** sampling of liquids






# Multi-phase HiSorb

## Phase combinations

- Polar compounds or compounds with  $\log K_{o/w}$  values  $< 3$  are not extracted so well and typically require a different phase type such as:
  - PDMS/CWR (Carbon Wide Range)
  - PDMS/DVB (Divinylbenzene)



For **non-targeted** analysis or to cover the **widest range** possible, a combination of **DVB/CWR/PDMS** is used






# Thermal desorption injection on Centri


- Fully compliant with numerous standard TD methods.
- Cryogen-free, multi-sorbent, back-flushing trap operation facilitates wide analyte volatility range.
- Automated re-collection of HS, SPME, HiSorb or TD sample split flows onto sorbent-packed TD tubes.
- Re-collected samples can be stored for re-analysis without using new sample aliquots, thus preserving valuable or 'unique' samples.
- Complete automated sample traceability ensured at all times.



Sample tubes are placed in the 50-tube autosampler, and leak-tests are performed on every tube to confirm sample integrity.



Analytes are released from the tubes in a flow of heated gas, and the vapours collected on the focusing trap.



The trap is thermally desorbed and the sample injected into the GC-MS.

# Speed up analysis and reducing cost per sample


## Multi-gas enabled

- **Future-proof** thanks to independent certification of Centri 90 and Centri 180 for safe use with hydrogen, **safeguarding against helium shortages** and rising costs by switching now or in the future.
- **Remove dependency on gas bottles**, together with their associated costs, using gas generators to provide a consistent gas supply, eliminate any risk of downtime while **lowering cost-of-ownership**.
- **Speed-up sample analysis** to maximise your return on investment with **faster GC separations** enabling higher sample throughput and so increasing earning potential.

48 samples per system, per day with **helium** carrier


~50% (25) more samples per day with **hydrogen** = 73 samples in total

\$100/sample, therefore, **additional earnings of \$2500 per system, per day**




# Food safety: Detecting fumigants in spices and matrices

- Used as a fumigant to eliminate insects in seasonings, seeds and spices
- In the EU, EtO use is banned due to its highly toxic carcinogenic, mutagenic and reprotoxic properties
- Produces 2-chloroethanol (2-CE) as a by-product<sup>1</sup>, which is also a toxic chemical restricted by the EU
- The EU imposes strict maximum residue limits (MRLs) for EtO and 2-CE in imported products



Ethylene oxide




2-Chloroethanol

# Expanding the compound range for residual fumigants testing

- Other fumigant compounds are reported and may be used instead to circumvent tests for ethylene oxide use.
- 5 fumigants, spanning a volatility range from 3 °C to 131 °C, were analysed using the same methodology as ETO & 2-CE.

Compound	Formula	B.p. (°C)	Quant ion (m/z)	Hazardous properties	EU MRL (mg/kg)
Bromomethane	CH <sub>3</sub> Br	3.6	94	Highly toxic	0.01 <sup>a</sup>
Bromoethane	C <sub>2</sub> H <sub>5</sub> Br	38	108	Toxic	0.01 <sup>a</sup>
Carbon tetrachloride	CCl <sub>4</sub>	77	117	Highly toxic	0.01 <sup>b</sup>
1,2-Dichloroethane	C <sub>2</sub> H <sub>4</sub> Cl <sub>2</sub>	83	62	Carcinogen	0.01 <sup>c</sup>
1,2-Dibromoethane	C <sub>2</sub> H <sub>4</sub> Br <sub>2</sub>	131	107	Highly toxic, carcinogen	0.01 <sup>c</sup>

**Table 1:** Target analyte information with EU MRLs shown, and regulation number denoted with superscript. <sup>a</sup> Regulation (EC) No. 1005/2009. <sup>b</sup> Regulation (EU) 2021/155. <sup>c</sup> Regulation (EU) 2012/649.



- A Bromomethane
- B Bromoethane
- C Carbon tetrachloride
- D 1,2-Dichloroethane
- E 1,2-Dibromoethane

**Figure 2:** Extracted-ion chromatogram (merged EIC of m/z 94, 108, 31, 117, 62 and 107) for all target compounds, spiked at 0.05 mg/kg on sesame seeds.



# One platform, many applications...







Thanks for listening!

Any questions?