



Agilent Technologies

AA Furnace Performance Figures

Introduction

This table provides data on Detection Limits and Characteristic Masses available from Agilent graphite furnace instruments.

These figures have been determined using an AA220 Zeeman with UltrAA lamps (where available for the elements). The sample volume used was 20 μL . All Agilent Zeeman systems, including the AA240Z, AA280Z and AA880Z should at least match these performance figures.

Although based on Zeeman performance, these figures are also representative of the performance achieved when using a deuterium furnace system. The deuterium furnace system should at least match these figures, and should provide better performance for those elements that have MSR's < 100 %.

Element	Wavelength λ (nm)	Characteristic Mass (μg)	Absolute Detection Limit (μg)	Detection Limit ($\mu\text{g/L}$)
Ag	328.1	1.5	0.40	0.020
Al	309.3	7.5, 8.0	1.2	0.060
Al	396.2	12, 14	2.0	0.10
As	193.7	13, 20, 22	4.4	0.22
Au	242.8	9, 11	6.1	0.30
B	249.8			

Ba	553.6	11		
Be	234.9	0.8	0.20	0.010
Bi	306.8	18	5.8	0.29
Ca	422.7	0.51		
Cd	228.8	0.5, 0.7	0.18	0.009
Ce				
Co	240.7	8	4.0	0.20
Cr	357.9	2.5, 3.5	0.75	0.038
Cs	852.1			
Cu	324.8	5, 6		
Cu	327.4	10, 11	2.9	0.15
Dy	421.2			
Er	400.8			
Eu	459.4			
Fe	248.3	4, 5.5	3.1	0.16
Ga	294.4			
Gd	368.4			
Ge	265.2			
Hf	307.3			
Hg	253.7			

Ho	410.4			
In	303.9			
Ir	208.9			
K	766.5	1.0		
La	550.1			
Li	670.8			
Lu	336.0			
Mg	285.2	0.4	0.08	0.004
Mn	279.5	1.7, 2.0	0.75	0.038
Mo	313.3	7.5, 9.0		
Na	589.0	0.67		
Nb	334.9			
Nd	492.5			
Ni	232.0	9, 11	7.5	0.38
Os	290.9			
P	213.6	5600	1800	90
Pb	217.0	6.5, 8	2.2	0.11
Pb	283.3	12, 16, 18	4.0	0.20
Pd	247.6	16.8	7.5	0.38
Pr	495.1			

Pt	265.9	140		
Rb	780.0			
Re	346.1			
Rh	343.5			
Ru	349.9			
Sb	217.6	25, 26, 30	10	0.50
Sc	391.2			
Se	196.0	23, 30	4.6	0.23
Si	251.6	38	20	1.00
Sm	429.7			
Sn	286.3	23	10	0.50
Sr	460.7	1.3		
Ta	271.5			
Tb	432.7			
Te	214.3	16, 18, 20	12	0.60
Th				
Ti	364.3			
Tl	276.8	20, 22	5.0	0.25
Tm	371.8			
V	318.5	22	5.2	0.26

W	255.1			
Y	410.2			
Yb	398.8			
Zn	213.9	0.61, 0.68		
Zr	360.1			

Notes:

- These detection limits are 3 sigma values - that is, the detection limit is defined as the concentration of a solution which gives an absorbance equal to 3 times the standard deviation of the blank. This value establishes the minimum determinable quantity for an element.
- The absolute detection limit is simply the detection limit (in concentration terms) converted to a mass basis. This conversion can be made by multiplying the detection limit (in $\mu\text{g/L}$) by the sample volume (in μL).
- The characteristic mass is simply the characteristic concentration converted to a mass basis. This conversion can be made by multiplying the characteristic concentration (in $\mu\text{g/L}$) by the sample volume (in μL). The characteristic concentration is the slope of the calibration curve, i.e., the concentration required to give 0.0044 absorbance.
- Blank and standard solutions were prepared in a 1 % (v/v) nitric acid and diluted using de-ionized water.

Acknowledgements:

- This data has been gratefully provided by Jean-Pierre Lener of Agilent Technologies, France.



Agilent Technologies