



Agilent 2100 Bioanalyzer System


APPLICATIONS FOR DNA, RNA, PROTEIN AND CELL ANALYSIS

Application Compendium

The Measure of Confidence



Agilent Technologies



Agilent 2100 Bioanalyzer System

APPLICATIONS FOR DNA, RNA, PROTEIN AND CELL ANALYSIS

Application Compendium

Agilent Technologies is a leading global provider of life science and chemical analysis solutions, including instrumentation, supplies, software and services. In Life Sciences, Agilent helps investigators advance basic research, unlock the genetic causes of disease and accelerate the discovery and development of new drugs. Agilent's life sciences business addresses the needs of academic, institutional and pharmaceutical scientists in five key areas:

- Genomics
- Proteomics
- Metabolomics
- Bioinformatics
- Pharmaceutical analysis



Let us provide you with the right solution for your success.

Agilent 2100 Bioanalyzer system

A wide range of applications

The Agilent 2100 Bioanalyzer system is an easy-to-use benchtop platform with ready-to-run kits for a wide range of applications.

- **On-chip flow cytometry**
The easy way to acquire dual-color, cell-based fluorescence data
- **DNA size and quantity**
High resolution separation and quantitation of DNA down to pg/ μ L sensitivity
- **RNA quality check with RIN**
The industry standard for RNA analysis offering total RNA, mRNA and Small RNA data including RIN algorithm (RNA Integrity Number)
- **SDS-PAGE replacement for protein analysis**
The fast and reliable way to determine the quantity and purity of proteins from Coomassie down to silver stain sensitivity

When combined with any one of our kits, you will discover how lab-on-a-chip technology can revolutionize your laboratory.

Microfluidics utilizes interconnected networks of micro-channels and wells for the analysis of various sample types. The technology reduces overall space and volume requirements, and allows online integration of many workflow steps, such as sample enrichment, separation, staining, de-staining, and detection. Advantages of microfluidics include dramatically reduced sample and reagent consumption, significantly faster analysis time and less hands-on activities during sample preparation and data analysis. The Agilent 2100 Bioanalyzer system represents the first microfluidic lab-on-a-chip platform for the electrophoretic analysis of DNA, RNA, and proteins and the flow-cytometric analysis of cells. This versatility makes the 2100 Bioanalyzer system an indispensable tool for the molecular biologist and biochemist.



Lab-on-a-chip technology increases quality and efficiency of your analysis

Lab-on-a-chip technology has many advantages over conventional techniques. These advantages include improved data precision and reproducibility, short analysis times, minimal sample consumption, improved automation and integration of complex workflows.

One of the many benefits the Agilent 2100 Bioanalyzer system has over conventional bioanalytical methods is the elimination of time consuming procedures – you enjoy standardized handling and interpretation of data. It simplifies the process of data gathering and analysis down to three quick and easy steps:

Load sample, run analysis, view data

- The Agilent 2100 Bioanalyzer system and the various **DNA kits** are the tools of choice for automated sizing and quantitation of products generated by RT-PCR and any type of multiplex PCR with unprecedented accuracy and reproducibility. It not only provides the detection of the presence or absence of a PCR product, but also offers quantitation of this product and detection of unspecific amplification. The Agilent 2100 Bioanalyzer system therefore helps to optimize PCR reactions for gene expression, sequencing, cloning and typing. The DNA 12000 kit also offers the accurate analysis of restriction digests.
- The **High Sensitivity DNA kit** provides sizing and quantitation of DNA fragments and DNA smears in the 50 to 7000 bp size range down to pg/ μ L sensitivity. This is especially useful for sample quality control and the monitoring of critical steps in next-generation sequencing (NGS) workflows, including DNA fragmentation, target enrichment, and DNA library amplification.
- Due to the omnipresence of RNases, and the instability of RNA, integrity checks and sample quantitation are essential steps before any RNA-dependent experiment. The 2100 Expert software generates the unambiguous RNA Integrity Number (**RIN**), provides a quantitation estimate, calculates ribosomal ratios of total RNA samples and automatically detects ribosomal RNA contamination in mRNA. The **RNA 6000 Nano kit** is a well-established standard for RNA sample QC. The **RNA 6000 Pico kit** allows detection of RNA degradation with sample amounts as low as 200 pg of total RNA.
- The high resolution **Small RNA kit** allows separation, verification and optimization of miRNA after extraction procedures. By consuming only 1 μ L sample, even pg amounts of purified small RNA can be measured reproducibly and comparably within 30 minutes. By effectively staining single- and double-stranded oligonucleotides at the same time, the assay is a versatile tool.

- The **Protein 80 and Protein 230 kits** provide a fast and easy way to analyze a wide range of samples, whether expressing recombinant proteins, purifying proteins, performing stability studies or checking antibody quality. The on-chip electrophoresis provides size, purity and concentration information for ten protein samples in less than 30 minutes. The lab-on-a-chip approach eliminates handling SDS-PAGE gels, staining or imaging steps.
- With the **High Sensitivity Protein 250 kit** it is possible to analyze proteins down to 1 pg/μL (on chip), which is equivalent or superior to silver stain SDS-PAGE. It provides quantitation over a dynamic range of up to 4 orders of magnitude with the reproducibility and ease-of-use only associated with the Agilent 2100 Bioanalyzer system.
- The **Flow Cytometry set** and the **Cell kit** for the Agilent 2100 Bioanalyzer system allows scientists to perform simple on-chip flow cytometry assays, thereby making the 2100 Bioanalyzer system the industry's only platform able to run RNA, DNA, proteins and cell analysis. Six cell samples each 10 μL with 20,000 prestained cells are loaded onto the chip and the fluorescence intensities in two channels for about 750 single cells per sample measured within 25 minutes.
- The optional **Agilent 2100 Security Pack software** ensures full **21 CFR Part 11 compliance** of your 2100 Bioanalyzer system for regulated environments such as pharmaceutical QA/QC labs or manufacturing and addressing requirements such as electronic signatures, audit trails and user authentication. Along with **IQ and OQ/PV** support services and Declarations of Conformity for all components offered for all assays and kits, your Agilent 2100 Bioanalyzer system will be compliant in no time.



Contents

Page

Agilent 2100 Bioanalyzer system

Cell fluorescence analysis	9
Protein expression monitoring	
Cell surface antibody staining – CD4 in CCRF CEM T-cells	10
Cell surface antibody staining – CD3 in T-cell leukemia	11
CD3 expression in T-cell leukemia via on-chip staining	12
Intracellular glucocorticoid receptor (GR) antibody staining in H4 hepatocytes	13
Analyzing a limited number of cells	14
Baculovirus titre determination	15
Upregulated gene expression in primary cells	16
Transfection efficiency monitoring	
Green fluorescent protein in CHO cells	17
On-chip staining of GFP expression for optimizing transfection conditions with different DNA:lipid ratios	18
Verification of stable transfected cell clones by on-chip antibody staining	19
Transfection of primary cells	20
Apoptosis detection	
Detection of phosphatidylserine on the cell surface via Annexin V binding	21
Intracellular Caspase-3 antibody staining assay	22
Apoptosis detection in primary cells	23
Fast Annexin protocol for time course of apoptosis induction via anti-FAS antibody	24-25
Gene silencing in cell culture	
siRNA transfection optimization	26
Monitoring of gene silencing experiments	27
Cells – others	
Identification and counting of blood cells in whole blood samples	28
Bacteria counting	29
DNA analysis	30
Restriction digest analysis	
Sizing range exemplified by the separation of Adenovirus 2/Dra I	31
Detection of single base mutations	32-33
PCR product analysis	
Separation of 3 different mixtures of PCR products	34
Optimize QPCR assay design	35
Determination of PCR product impurity	36-37
DNA sample quality from automated PCR purification	38
Multiplex PCR analysis of bacteria in chicken	39
Multiplex PCR with 19 products	40
Gene expression analysis	
mRNA expression study by comparative multiplex PCR	41
Standardized end-point RT-PCR	42

Co-amplification of GAPDH and hsp72	43
Co-amplification of GAPDH and hsp72 – response curves	44
Competitive PCR	45
Food analysis	
Estimation of non-basmati rice amounts in basmati rice products	46
Strawberry and raspberry fruit differentiation	47
Development of meat specific assays	48
Detection of non-fish species in fish samples with PCR-RFLP	49
Discrimination of sturgeon and related species by PCR-RFLP	50
GMO detection	
Development of a multiplex assay for soya	51
DNA stability during food processing	52
GMO detection by nested multiplex PCR	53
Oncology	
Tumor cell detection from carcinoma patient blood	54
SNP analysis in cancer related P16 gene	55
K-ras gene SNP detection	56
METH-2 downregulation in lung carcinomas	57
Label-free analysis of microsatellite instability in carcinoma	58
Clinical research	
Detection of a point mutation in the prothrombin gene with PCR-RFLP	59
Genotyping of <i>H. pylori</i>	60
Duplications and deletions in genomic DNA	61
Forensic testing	
Optimization of PCR on mtDNA	62
Pitfalls in mtDNA sequencing	63
Next-generation sequencing	
DNA library quantity and quality	64
DNA library QC in target enrichment and next-gen sequencing workflows	65
Sizing and quantitation of DNA libraries and fragmented DNA	66
Quality control of FFPE DNA samples	67
Analysis of limited DNA material on the Pippin Prep system	68
RNA analysis	69
Analysis of total RNA	
RNA integrity	70
Reproducibility of quantitation	71
Standardization of RNA quality control	72-73
QC of automated RNA isolation from HeLa cells	74
RNA quality from freshly frozen or formalin-fixed paraffin-embedded tissue	75
Assessing integrity of plant RNA	76
Genomic DNA contamination	77


Low amounts of total RNA	
Detection of low levels of RNA	78
RNA integrity with the RNA 6000 Pico kit	79
RNA quality after staining and microdissection	80
Analysis of minimum RNA amounts	81
Genomic DNA in RNA extracts of low abundance	82
Low RNA amounts from kidney sections	83
Analysis of mRNA	
RNA integrity	84
Ribosomal RNA contamination in mRNA samples	85
Analysis of Cy5-labeled samples	
Analysis of cRNA with and without dye in gel matrix	86
Optimization of labeling reactions	87
cRNA fragmentation	88
Analysis of T7-RNA transcripts	
Size estimation	89
Analysis of Small RNA	
Analysis of miRNA content in total RNA samples	90
Analysis of small RNAs from <i>Drosophila</i> Schneider	91
Effect of total RNA quality on mRNA and miRNA expression profiles	92
Protein analysis	93
Protein expression	
Analysis of cell lysates – protein induction	94
Protein purification	
Comparison between lysate and flow through	95
Analysis of protein purification	96
GFP Streptag fusion protein purification	97
Analysis of column capacity	98
Analysis of column fractions to optimize conditions	99
His-tag protein purification with Ni ⁺⁺ ZipTips®	100
Determination of protein expression level and recovery	101
Optimization of on-column cleavage	102
Enzymatic removal of His-tags from recombinant proteins	103
Complementing RP-HPLC protein purification	104
High sensitivity protein detection	
Detection of low protein amounts	105
Highly specific and sensitive alternative to Western blotting	106
Pico labeling protocol for samples below 1 ng/μL protein	107
Antibody analysis	
Analysis of antibodies under reducing and non-reducing conditions	108
High sensitivity analysis of antibodies	109
QA/QC of IgG under reducing conditions	110

Contents

Page

Quantitation of the half-antibody content in IgG4 preparations	111
Comparison of SDS-PAGE, CGE and Agilent 2100 Bioanalyzer for humanized monoclonal antibody analysis	112
Absolute quantitation of IgG	113
Quality control of stressed antibodies	114
Analysis of IgG2 under non-reducing conditions	115
IgG2 analysis under reducing conditions	116
Monitoring antibody charge variants	117
Separation of bispecific antibody chains	118
Protein quantitation	
Quantitation of limited protein with an internal standard	119
Absolute protein quantitation	120
Food analysis	
Rapid wheat varietal identification	121
Protein analysis in milk	122
Protein pattern of different transgenic seedlines	123
Protein – others	
Glycoprotein sizing	124
OFFGEL electrophoresis combined with high-sensitivity on-chip protein detection	125
Protein quality control prior to MS-analysis	126
Depletion of high abundant proteins from blood samples	127
Increased sensitivity by desalting protein samples	128
Analysis of PEGylated proteins	129
Agilent 2200 TapeStation system	130
DNA analysis	
PCR product analysis	
Automated DNA sample analysis from 96-well sample plates	131
Next-generation sequencing	
QC of DNA from FFPE tissues for next-generation sequencing	132
Genomic DNA analysis	
Analysis of high molecular weight genomic DNA	133
RNA analysis	
Analysis of total RNA	
RNA quality control with the RIN [®] quality metric	134
Protein analysis	
Antibody analysis	
Reproducible integrity and purity testing of antibodies	135
Literature	136

Cell fluorescence analysis




Protein expression monitoring

Cell surface antibody staining – CD4 in CCRF CEM T-cells

Flow cytometer (10,000 events)

2100 Bioanalyzer system (500 events)



Kit: Cell kit

Assay: Antibody staining assay


Application: CCRF-CEM cells were stained with hCD4-APC labeled antibodies and calcein live dye. 65 % of all CCRF-CEM live cells (yellow curve) are expressing CD4 protein which is good in comparison to conventional flow cytometer results.

Application note: 5988-4322EN

Protein expression monitoring

Cell surface antibody staining – CD3 in T-cell leukemia

Averaged data per instrument



Mean % CD3+ cells				
2100-1	2100-2	2100-3	2100-4	Flow cyt.
60.9	67.8	66.6	65.0	60.9
34.4	36.7	36.7	34.3	29.8
17.3	17.6	18.7	17.2	13.8
8.9	9.4	9.9	8.3	6.5
5.1	4.4	5.3	4.9	3.2
0.8	0.6	0.3	0.3	0.0

Kit: Cell kit

Assay: Antibody staining assay

Application: Jurkat (T-cell leukemia) cells were stained with calcein alone or with calcein and APC-labeled anti-CD3 antibody. To mimic different subpopulation sizes, mixtures of both populations were prepared at various ratios.


Samples were analyzed with 4 2100 Bioanalyzer instruments on 5 chips and compared to a flow cytometer reference instrument. Interestingly, small subpopulations (like 10 – 20 %) could be analyzed with good accuracy and reproducibility.

Application note: 5988-4322EN

Protein expression monitoring


CD3 expression in T-cell leukemia via on-chip staining

A. On-chip flow cytometry



Sample	% of Gated	Sample	% of Gated
1	-	4	64.6
2	66.9	5	66.7
3	67.2	6	72.0

B. Conventional flow cytometry



Kit: Cell kit

Assay: Antibody staining assay

Application: Jurkat cells were stained on-chip with anti hCD3-APC prediluted 1:5.5 in cell buffer and Calcein (1:50 in cell buffer). After an incubation time of 25 minutes in the chip, samples were measured in the 2100 Bioanalyzer instrument. The faster and easier on-chip staining procedure has the advantage here of reducing cell consumption 17 fold and antibody reagent costs 80 fold.


- A) Overlay of representative histograms of calcein and antibody treated cells.
- B) Comparison between on-chip staining data and data obtained by measuring cells stained by conventional staining on a flow cytometer.

Application note: 5988-7111EN


Protein expression monitoring

Intracellular glucocorticoid receptor (GR) antibody staining in H4 hepatocytes

Chip histogram overlay from 700 cells/sample



Correlation of chip vs. flow cytometer results



Kit: Cell kit

Assay: Generic assay

Application: H4 hepatocytes cells were stained with SYTO16 DNA dye alone or with SYTO16 and GR primary antibody. After washing, both cell preparations were stained with APC-labeled secondary antibody. Mixtures of both populations were prepared at various ratios.

The insert in the left picture shows the overlay of all six cell samples in the blue reference color. The black histogram represents data from the control sample, no GR detected.


All other 5 samples have significant staining above marked fluorescence intensity in the red.

Good chip to chip reproducibility and comparison to flow cytometer is demonstrated.

Application note: 5988-4322EN

Protein expression monitoring

Analyzing a limited number of cells



Cells	Live-CD3+	STD(n=4)
20,000	83.7 %	3.5 %
10,000	85.6 %	4.1 %
5,000	87.7 %	4.2 %
2,500	84.0 %	3.0 %
1,250	89.8 %	6.5 %
625	90.0 %	9.3 %

Kit: Cell kit

Assay: On-chip antibody staining assay


Application: The direct comparison of different input numbers of cells (down to 625 cells in 10 µL) for the on-chip staining protocol reveals that even with a much lower number than the recommended 20000 cells/10 µL for the standard protocol reliable and meaningful results can be achieved with good reproducibility. The data shown were generated with CD3-positive Jurkat cells stained with an anti-CD3 antibody for the CD3 protein and counterstained with the live cell stain Calcein AM. Similar results were obtained with primary human dermal fibroblasts (PHDF) indicating the usefulness of this method for scarce specimen. The lack of sensitivity, automation and convenient quantitation found with other methods can be circumvented easily by using the 2100 Bioanalyzer system.

Application note: 5989-0746EN



Protein expression monitoring

Baculovirus titre determination

Fluorescent Light



Transmitted Light



Kit: Cell kit


Assay: GFP Assay

Application: A fast and convenient method exists for the calculation of baculovirus titre for expression systems facilitating insect cells. Using GFP-linked co-expression plasmids, the 2100 Bioanalyzer instrument and the flow cytometry set allows the calculation of the viral titre for six samples in approximately 90 minutes. It is superior to traditional plaque assays in terms of labor time, automation and user-to-user variability.

Application note: 5989-1644EN

Protein expression monitoring

Upregulated gene expression in primary cells



Kit: Cell kit

Assay: On-chip antibody staining assay

Application: Flow cytometric analysis of primary cells can present a challenge for researchers due to limited availability and life span of primary cells. A dose-responder upregulation of protein expression in primary cells using only a minimum number of cells in a fast on-chip-staining approach is shown here. Activation of peripheral blood lymphocytes by phorbol-12-myristate-13-acetate (PMA) leads to increased expression of the T cell receptor CD3 (Figure A, mean from 3 experiments).



For HUVECs (human umbilical vein endothelia cells) the induction of E-selectin (CD62E) expression upon IL-1 β treatment is shown (Figure B, white bars) in comparison to results from a conventional flow cytometer (white bars).

Application note: 5989-2718EN




Transfection efficiency monitoring

Green fluorescent protein in CHO cells

Mock transfected cells



GFP transfected cells




Kit: Cell kit
Assay: GFP assay

Application: Chinese hamster ovary (CHO-K1) cells were transfected with EGFP DNA by a lipofection method. The upper panel shows the control mock transfection; here cells don't express GFP. Examples for data evaluation in dotplot view and histogram view are shown in comparison to the microscopy view. For analysis on the 2100 Bioanalyzer instrument, cells were stained with a red dye for live cells (reference stain). The transfection efficiency of 56 % can be easily determined with the 2100 Bioanalyzer system.

Application note: 5988-4320EN

Transfection efficiency monitoring

On-chip staining of GFP expression for optimizing transfection conditions with different DNA:lipid ratios



Kit: Cell kit

Assay: On-chip GFP assay


Application: Chinese hamster ovary (CHO-K1) cells were transfected with EGFP DNA by alipofection method. Optimization of transfection conditions were done on one chip. Several DNA:lipofectamine ratios were tried. A ratio of 1:8 gave the best transfection efficiency. All cells were reference stained with a red live dye. On-chip staining was applied, minimizing the staining time, reagent usage and cell consumption.

Application note: 5988-7296EN


Transfection efficiency monitoring

Verification of stable transfected cell clones by on-chip antibody staining

Hek 293 control cells



CD 80 stable clone



Kit: Cell kit


Assay: On-chip antibody staining assay

Application: Verification of CD80 protein expression in stable transfected Hek 293 cells with the 2100 Bioanalyzer system. Control (left dot plot) and CD80 transfected cells (right) are stained on-chip with blue calcein live dye and anti-CD80-CyChrome antibody. Red region marks CD80 protein expressing 293 cells within live cell population – confirming expression in the CD80 stable clone Hek 293 cells.

Application note: 5988-7111EN

Transfection efficiency monitoring

Transfection of primary cells



Kit: Cell kit


Assay: GFP assay

Application: Monitoring the transfection efficiency in primary cells requires low cell consumption, high reproducibility of results, a fast on-chip staining procedure and ease-of-use all provided by the 2100 Bioanalyzer system. The transfection efficiency using a GFP-coding plasmid (pEGFP-C2) at varying plasmid:lipofectamine ratios (DNA:LP ratio) obtained with human umbilical vein endothelial cells (HUVEC) is measured in this optimization series. Images from a fluorescence microscope (A) and dot plots (B), as well as histograms (C) of control- and GFP-transfected cells are shown. Using increasing ratios, better transfection efficiency was achieved, whereas the toxicity of LP caused decreased recovery of living cells. Such data facilitates optimizing transfection conditions.

Application note: 5988-8154EN

Apoptosis detection

Detection of phosphatidylserine on the cell surface via Annexin V binding



Kit: Cell kit

Assay: Apoptosis assay

Application: Apoptosis (programmed cell death) in Jurkat cells was induced with camptothecin. Cells treated for 16 hours and untreated cells were stained with calcein and Annexin-Cy5.


Annexin-V binds to phosphatidylserine – a membrane lipid which is kept to the inner leaflet of the cell membrane of intact cells. Exposure of phosphatidylserine on the outer leaflet is an early indicator of apoptotic processes. Annexin-V binding is made detectable by Cy5 staining of the Annexin-V via a biotin-streptavidin interaction. Calcein staining of cells is used as a live control to distinguish living and apoptotic cells from dead cells. Calcein enters the cell via the membrane as a non-fluorescent ester. The ester is cleaved inside the cell which results in fluorescence.

The histograms on the left show the number and intensity value of all events which generated a signal in the blue channel, corresponding to calcein-stained cells. The histogram on the right shows all events which generated a signal in the red channel, corresponding to Annexin-V binding to apoptotic cells. While the control shows only low intensity values (background noise), the treated sample shows high intensity values (within the red markers) corresponding to apoptotic cells. The dot plot of the treated sample nicely shows the subpopulation of all live cells which are apoptotic.

Application note: 5988-4319EN

Apoptosis detection

Intracellular Caspase-3 antibody staining assay



Kit: Cell kit


Assay: Generic assay

Application: Induction of apoptosis in Jurkat cells was done with anti-FAS antibody treatment. Intracellular staining with specific antibodies against 'active' Caspase-3 were performed. Reference staining was done with SYTO16 DNA dye. Good chip to chip reproducibility and good comparison to conventional flow cytometer results were obtained.

Application note: 5988-4319EN

Apoptosis detection

Apoptosis detection in primary cells



Kit: Cell kit


Assay: On-chip antibody staining assay

Application: The 2100 Bioanalyzer system has been used to study induced apoptosis by monitoring annexin V-binding in primary human endothelial cells (HUVEC, not shown) and human dermal fibroblasts (NHDF, shown). A simple and fast assay protocol was used on cells left untreated or treated for 5 hours with different concentrations of staurosporine, which induces apoptosis. See row A for dot blots and B for histograms at different concentrations. Evaluation of the same samples on a conventional flow cytometer (row C) yielded similar results.

Application note: 5989-2934EN

Apoptosis detection


Fast Annexin protocol for time course of apoptosis induction via anti-FAS antibody



Kit: Cell kit

Assay: Apoptosis assay

Application: Apoptosis (programmed cell death) in Jurkat cells was induced with anti-FAS antibody. Cells treated for 0,1,2,3,4 and 6 hours were stained with calcein and Annexin-Cy5. Annexin-V binds to phosphatidylserine – a membrane lipid which is kept to the inner leaflet of the cell membrane of intact cells. Exposure of phosphatidylserine on the outer leaflet is an early indicator of apoptotic processes. Annexin V binding is detectable by Cy5 staining of the Annexin-V via a biotin-streptavidin interaction. Calcein staining of cells is used as a live control to distinguish living and apoptotic cells from dead cells. Calcein enters the cell via the membrane as non-fluorescent ester. The ester is cleaved inside the cell which results in fluorescence and indicates apoptosis.




The histograms on page 24 (A) show the number and intensity value of all events which generated a signal in the blue channel, corresponding to calcein-stained cells. The histograms on the right show all events which generated a signal in the red channel, corresponding to Annexin-V binding to apoptotic cells. While the control shows only low intensity values (background noise), the treated sample shows high intensity values (within the red markers) corresponding to apoptotic cells.

(B) Time course of the induction of apoptosis by anti-FAS antibody in Jurkat cells. Apoptosis is detectable in a significant amount of cells after 2 hours. Following a treatment of 4 hours, approximately 95 % of the cells are apoptotic.

Application note: 5988-7297EN

Gene silencing in cell culture

siRNA transfection optimization



Kit: Cell kit


Assay: On-chip antibody staining assay

Application: In gene silencing experiments (HeLa cells) we found that increasing amounts of transfection reagent (TransMessenger™) to a constant amount of siRNA leads to a plateau of transfection viability (panel A). Transfection viability reflects the product of the viability of the transfected cells and the transfection efficiency. With a constant siRNA/transfection reagent ratio of 1:4 and increasing total amounts of introduced siRNA (panel B) the viability of transfected cells decreases at a certain point although the transfection efficiency increases. Thus, there are experimental conditions where the number of living and transfected cells are at a maximum. The 2100 Bioanalyzer system features on-chip staining and leads to excellent results with a minimal consumption of cells and reagents.

Application note: 5988-9872EN

Gene silencing in cell culture

Monitoring of gene silencing experiments



Kit: Cell kit


Assay: GFP Assay

Application: After co-transfection of a GFP plasmid and Cy5-labeled siRNA (GFP-specific), GFP expression and viability of cells were detected. The course of GFP expression in control (GFP only) and siRNA/GFP transfected cells was measured on the 2100 Bioanalyzer system. Accurate results were obtained fast and in an automated manner. They easily allow the efficiency and reliability of a given protocol and transfection reagents to be judged. Thus, such an experiment provides efficient monitoring and optimization of any gene silencing experiment.

Application note: 5989-0103EN

Cells – others

Identification and counting of blood cells in whole blood samples



Kit: Cell kit

Assay: Generic assay


Application: Nucleic acid dye and antibody staining was used in accordance with the dual-fluorescence detection system of the 2100 Bioanalyzer system to identify several blood cell populations and count blood cells from whole blood samples without a washing step. The red blood cells (RBC) are stained with an anti-CD235a monoclonal antibody conjugated with the fluorochrome.

- A) On-chip antibody staining results as dot plots for 6 samples including a region with spiked-in beads for quantitation.
- B) Quantitative determination of RBC counts with the 2100 Bioanalyzer system and an automatic blood cell counter. The microvolumes of samples and reagents, the low number of cells required for the analysis as well as the easy use of the 2100 Bioanalyzer system are the specific advantages of this microfluidic chip-based technology in comparison with the cell counter.

Application note: 5989-7171EN

Cells – others

Bacteria counting



Kit: Cell kit

Assay: Generic Series II assay

Application: Bacterial cell counting is a widely performed task. The hemocytometer counting method and the spectrophotometric analysis are error prone and show high variation. Therefore, on-chip flow cytometry with the Agilent 2100 Bioanalyzer system was evaluated as alternative process for bacterial cell counting. Three different kinds of prokaryotes: *Escherichia coli* (gram-negative, rod-shaped), *Pseudomonas aeruginosa* (gram-negative, rod-shaped), and *Staphylococcus aureus* (gram-positive, coccal) were tested. As a reference for bacterial concentration, optical density of the bacterial suspension was measured, and hemocytometer counts were obtained. The figure shows the dot blot for *E. coli* (3×10^6 cells/mL) obtained after on-chip staining with SYTO 16. Furthermore, the number of events was plotted against the cell concentration determined with the hemocytometer to obtain a correlation curve. Interestingly, bacteria with different morphology show different standard curves. *E. coli*, as a short rod, shows a behavior between the other two bacteria.

The Agilent 2100 Bioanalyzer system provides a very efficient and highly sensitive method for bacterial counting. All types of bacteria can be counted, whether gram-positive or gram-negative.

Application note: 5991-2582EN

DNA analysis

Forensics/Homeland Security
 Proteomics
 Pharmaceuticals
 Genomics
 Drug Manufacturing QA/QC
 Drug Discovery/Drug Development
 Agriculture/Food

Restriction digest analysis

- Sizing range exemplified by the separation of Adenovirus 2/Dra I
- Detection of single base mutations (I)
- Detection of single base mutations (II)

PCR product analysis

- Separation of 3 different mixtures of PCR products
- Optimize QPCR assay design
- Determination of PCR product impurity
- DNA sample quality from automated PCR purification
- Multiplex PCR analysis of bacteria in chicken
- Multiplex PCR with 19 products

Gene expression analysis

- mRNA expression study by comparative multiplex PCR
- Standardized end-point RT-PCR
- Co-amplification of GAPDH and hsp72
- Co-amplification of GAPDH and hsp72 – response curves
- Competitive PCR

Food analysis

- Estimation of non-basmati rice amounts in basmati rice products
- Strawberry and raspberry fruit differentiation
- Development of meat specific assays
- Detection of non-fish species in fish samples with PCR-RFLP
- Discrimination of sturgeon and related species by PCR-RFLP

GMO detection

- Development of a multiplex assay for soya
- DNA stability during food processing
- GMO detection by nested multiplex PCR

Oncology

- Tumor cell detection from carcinoma patient blood
- SNP analysis in cancer related P16 gene
- K-ras gene SNP detection
- METH-2 downregulation in lung carcinomas
- Label-free analysis of microsatellite instability in carcinoma

Clinical research


- Detection of a point mutation in the prothrombin gene with PCR-RFLP
- Genotyping of *H. pylori*
- Duplications and deletions in genomic DNA

Forensic testing

- Optimization of PCR on mtDNA
- Pitfalls in mtDNA sequencing


Next-generation sequencing

- DNA library quantity and quality
- DNA library QC in target enrichment and next-gen sequencing workflows
- Sizing and quantitation of DNA libraries and fragmented DNA
- Quality control of FFPE DNA samples
- Analysis of limited DNA material on the Pippin Prep system



Restriction digest analysis

Sizing range exemplified by the separation of Adenovirus 2/Dra I




Kit: DNA 12000 kit
Assay: DNA 12000 assay

Application: Restriction digest analysis of Adenovirus 2/Dra I. For restriction fragment analysis the large linear dynamic range of the lab-on-a-chip approach is very advantageous. Analyzing samples with large and short fragments on slab gels can be difficult because of bands running off the gel and insufficient staining (or over-staining) of bands.

Application note: 5968-7501EN

Restriction digest analysis

Detection of single base mutations (I)



Kit: DNA 7500 kit

Assay: DNA 7500 assay


Application: Mutation detection by RFLP highlights the use of the 2100 Bioanalyzer system.

Two different regions of the p53 gene were amplified with specific primers and digested with Hpa II, which cuts in a location that is prone to mutations. In the presence of a point mutation, the enzyme Hpa II does not cleave the DNA, leaving larger fragments that can be revealed by gel electrophoresis or by analysis with the DNA 7500 kit (see next page).

Data not published

Restriction digest analysis

Detection of single base mutations (II)



Kit: DNA 7500 kit


Assay: DNA 7500 assay

Application: Analysis on the chip showed an identical pattern of digest fragments as seen on the slab gel for the wildtype and Exon 7 & 8 PCR products. Comparison of the calculated sizes of the bands shows 1-2 % variance with the assay, which allows fast and accurate detection of point mutations.

Application note: 5968-7496EN

PCR product analysis

Separation of 3 different mixtures of PCR products



Kit: DNA 500 kit*

Assay: DNA 500 assay*


Application: Overlay of three different electropherograms, which are mixtures of PCR samples ranging from 25 to 500 base pairs in size. The two closest eluting bands (50 bp and 53 bp) are partially separated and identified by the software as two separate peaks. The DNA 500* assay achieves a resolution of five base pairs from 25 to 100 base pairs and a 5 % resolution from 100 to 500 base pairs where the sizing error is less than 10% over the entire size range.

Application note: 5988-3041EN

* replaced with DNA 1000 kit and assay

PCR product analysis

Optimize QPCR assay design




Kit: DNA 1000 kit
Assay: DNA 1000 assay

Application: To determine the best possible RT priming strategy, QPCR and 2100 Bioanalyzer system tests were conducted for three genes (GAPDH, HPRT1 and YWHAZ) by using high quality RNA (RIN = 10) as template. The RNA was reverse transcribed with either oligo-dT or random priming (A: 5'-assay, B: 3'-assay for GAPDH). To assess size and purity of QPCR amplicons, 1 μ L of the QPCR reactions was analyzed with the DNA 1000 assay. In addition, No-template-controls (NTC) were used to assess contamination and potential primer dimer formation, since already a small amount of contaminating template can lead to amplification. For the HPRT1 5' oligo-dT assay one of the NTC was positive, showing a peak in the melt curve at a similar melting temperature compared to the positive control (A). To verify that no contamination of the well has occurred, the NTC was analyzed using the DNA 1000 assay on the 2100 Bioanalyzer instrument (B). Two minor peaks could be detected at 21 and 51 bp which are most probably related to primer and primer dimers. This highlights the high information content obtained by the 2100 Bioanalyzer system and the poor discrimination capabilities of a SYBR Green based melt curve.

Application note: 5989-7730EN

PCR product analysis

Determination of PCR product impurity



Kit: DNA 7500 kit


Assay: DNA 7500 assay

Application: Comparison between the analysis of two PCR reactions (300 and 3000 bp products) using the DNA 7500 kit vs. an agarose gel. Two different concentrations are shown side by side for each PCR reaction (undiluted and 1:4 dilution). The 2100 Bioanalyzer system shows superior performance in locating impurities over a broader concentration range than the gel. The 300 bp fragment appears to be uncontaminated in both the gel and on the 2100 Bioanalyzer system. The 3000 bp fragment shows few impurities on the gel, which become invisible at the 1:4 dilution. These impurities can easily be detected with the 2100 Bioanalyzer system.


Application note: 5968-7496EN

PCR product analysis

Determination of PCR product impurity



Sample	c (DNA) All peaks	Main peak (300 bp)
Pure PCR	41.4 ng/μL	40.7 ng/μL
1:4 dilution	9.6 ng/μL	9.6 ng/μL



Sample	c (DNA) All peaks	Main peak (3000 bp)
Pure PCR	61.9 ng/μL	40.7 ng/μL
1:4 dilution	14.8 ng/μL	9.8 ng/μL

Kit: DNA 7500 kit


Assay: DNA 7500 assay

Application: The quantitative data generated by the 2100 Bioanalyzer system indicate the amount of impurity or non-specific products in the PCR reactions from the previous page. Even in the 300 bp fragment a small impurity can be detected, while the 3000 bp fragment shows more than 50 % impurities.

Application note: 5968-7496EN

PCR product analysis

DNA sample quality from automated PCR purification



Kit: DNA 1000 kit


Assay: DNA 1000 assay

Application: A protocol for automated PCR purification in a 96-well format employing the Bravo Automated Liquid Handling Platform, Automated Centrifuge and StrataPrep 96 PCR Purification kit from Agilent was compared to a manual method. The DNA recovery quality from four wells from different portions of the plate was assessed with the 2100 Bioanalyzer system and the DNA 1000 kit. The figure shows overlaid samples from both manually-purified and robot-purified samples, demonstrating identical sizing and purity. The single peak at just over 600 bp is correctly sized within 1 % of the predicted PCR product size. There are no contaminating peaks or primers formed during PCR of Lambda DNA. The smaller peaks at 15 and 1500 bp are DNA markers.

Article: 5990-3948EN

PCR product analysis

Multiplex PCR analysis of bacteria in chicken



Data kindly provided by GenPoint, NL

Kit: DNA 500 kit*

Assay: DNA 500 assay*


Application: Multiplex PCR with four primer pairs, each one specific for a certain DNA sequence from one of the 4 bacteria to be tested for. Total DNA was extracted from chicken and subjected to PCR. The gel-like image shows traces from different chicken samples with bands showing up when an amplicon could be detected. The electropherogram is one example where bacterial DNA from two species of the *Campylobacter* genus could be detected.

Data not published

* replaced with DNA 1000 kit and assay

PCR product analysis

Multiplex PCR with 19 products



Data kindly provided by Qiagen, Germany

Kit: DNA 7500 kit


Assay: DNA 7500 assay

Application: Many molecular applications include PCR multiplexing as shown above with a PCR that yields 19 products. Applications are genotyping of transgenic organisms, detection of pathogens or GMs and microsatellite genotyping (e.g. short tandem repeat (STR) and variable number tandem repeat (VNTR) analyses). The sample shows optimization of PCR conditions (Mg^{2+} concentration) performed to ensure annealing of the multiple primers under identical conditions. Visualization and evaluation of the results can be performed efficiently with the 2100 Bioanalyzer system because of the high resolution, the accurate sizing, quantitation and extended linear range.

Application note: 5988-9342EN

Gene expression analysis

mRNA expression study by comparative multiplex PCR



Data kindly provided by the Roy Castle Centre

Kit: DNA 1000 kit


Assay: DNA 1000 assay

Application: Two genes were co-amplified in this study. A tumor specific gene (upper band) along with a housekeeping gene (lower band). The upregulation of the tumor gene is visualized via analysis on the 2100 Bioanalyzer system. Building the ratio of the concentration values obtained from the 2100 Bioanalyzer system, numerical values are obtained that are normalized with regard to the RT-PCR amplification efficiency. This way tumor tissue can be distinguished from normal tissue more unambiguously.

Data not published

Gene expression analysis

Standardized end-point RT-PCR



Data kindly provided by the Medical College of Ohio

Kit: DNA 7500 kit


Assay: DNA 7500 assay

Application: Complementary DNA from bronchial epithelial cells (BEC) was analyzed by a Standardized RT-PCR (StaRT) for the expression of 15 different genes. This analysis can be performed at the end-point of PCR without the need for real-time measurement at each cycle of PCR. Three methods for evaluation of representative results were compared (see above). The coefficient of variance (CV) from at least 3 measurements was calculated. The direct comparison of the reproducibility for agarose gel analysis (A, CV = 0.50) and the ABI Prism310 Genetic Analyzer (C, CV = 0.39) with the 2100 Bioanalyzer system (B, CV = 0.29) reveals that the 2100 Bioanalyzer system is superior. It is a reliable and valuable tool in quantitative gene expression analysis.

Application note: 5988-3674 EN

Gene expression analysis

Co-amplification of GAPDH and hsp72



Data kindly provided by Dr. Eric Gottwald, Forschungszentrum Karlsruhe, Germany

Kit: DNA 1000 kit


Assay: DNA 1000 assay

Application: Gel-like image and electropherograms showing the results of separate amplifications and co-amplifications of GAPDH and hsp72 in unstimulated HepG2 cells. Primers for GAPDH yield a PCR product of 443 bp (lane 1), primers for hsp72 yield PCR products of 384 and 650 bp (lane 2 and 3). Lane 4 and 5 show the results of the co-amplification reactions. Due to the competitiveness of the reaction, very little hsp72 products could be detected in lane 4 (insert) and no product was detected in lane 5 (lane 6 = negative control). The broad linear dynamic range of the analysis allows detection of weak bands next to strong bands and helped in the determination of gene expression in this case.

Article: 5988-4556EN

Gene expression analysis

Co-amplification of GAPDH and hsp72 – response curves



Data kindly provided by Dr. Eric Gottwald, Forschungszentrum Karlsruhe, Germany

Kit: DNA 1000 kit



Assay: DNA 1000 assay

Application: The optimized PCR conditions were used to monitor the response of a stimulus to hsp. Gene expression was monitored by comparing the RT-PCR amplification of a housekeeping gene with the co-amplification of hsp. In the current case, the highest gene expression was measured after about 10 minutes. As a comparison, the same set of samples was analyzed using the DNA 500 kit*. Virtually identical results are obtained with both kits, demonstrating that lab-on-a-chip technology can serve as a standardized approach to gel electrophoresis.

Article: 5988-4556EN

Gene expression analysis

Competitive PCR



Data kindly provided by Dr. Eric Gottwald, Forschungszentrum Karlsruhe, Germany

Kit: DNA 1000 kit


Assay: DNA 1000 assay

Application: Two genes were reverse transcribed and co-amplified in one reaction tube. The PCR products were analyzed using the DNA 1000 kit. Primers for hsp72 were present from the beginning of the PCR reactions, while primers for GAPDH were added after various cycle numbers ranging from 20 to 40 cycles (primer dropping method). This allowed optimization of this competitive PCR reaction. The left graph displays the dynamic range (arrow) in the gel like view, whereas the right graph indicates conditions with greatest sensitivity (red line).

Article: 5988-4556EN

Food analysis

Estimation of non-basmati rice amounts in basmati rice products



Kit: DNA 1000 kit


Assay: DNA 1000 assay

Application: Ensuring integrity of raw food materials, ingredients, and products is both a product quality and regulatory compliance concern. Food ingredient suppliers, manufacturers, and retailers can suffer economic and legal damages if proven to be supplying materials or products that are incorrectly labeled due to substitution or contamination. For example, EU Commission 1549/04 grants lower import tax on nine basmati rice varieties. Furthermore, the level of non-basmati rice in a basmati rice product must not exceed 7 %. Therefore it is required to identify basmati rice varieties and determine the level of non-basmati rice in basmati rice products. The 2100 Bioanalyzer system and the DNA 1000 kit can be used as a quick and cost-effective analytical method to differentiate approved and non-approved rice varieties using three primer sets and to estimate the level of non-basmati rice using reference rice admixtures. Four reference samples with 0 to 40 % non-basmati rice content were used to estimate the non-basmati amount in an unknown sample to be 10-20 %.

Application note: 5989-6836EN

Food analysis

Strawberry and raspberry fruit differentiation



Kit: DNA 1000 kit

Assay: DNA 1000 assay


Application: Food authentication is important to help detect fraudulent replacement of expensive food ingredients, and to ensure correct ingredient content in pre-packed foods. A PCR method followed by on-chip electrophoresis with the Agilent 2100 Bioanalyzer system and the DNA 1000 kit was used to distinguish strawberry and raspberry samples.

The figure shows the characteristic profiles obtained by the analysis of strawberry sauce through PCR with the microsatellite primer pairs Fvi11 and Fvi20 for PCR. This method allows clear differentiation between samples that contain strawberry (*Fragaria*) and raspberry (*Rubus*) DNA and may prove useful in food authentication.

Application note: 5990-3327EN

Food analysis

Development of meat specific assays



Data kindly provided by CCFRA, UK

Kit: DNA 500 kit*

Assay: DNA 500 assay*


Application: For detection of individual species in processed food, PCR assays with specific sets of primers can be developed. Example: turkey specific primers do not amplify any other meat species, including beef, chicken, lamb, or pork (see lane 5 and respective electropherogram).

Application note: 5988-4069EN

* replaced with DNA 1000 kit and assay

Food analysis

Detection of non-fish species in fish samples with PCR-RFLP



Kit: DNA 1000 kit

Assay: DNA 1000 assay

Application: The Agilent Fish ID solution was developed to identify fish species from processed foods. It is a PCR-RFLP method using a Cytb PCR target sequence with analysis of the restriction pattern fragments on the Agilent 2100 Bioanalyzer system with DNA 1000 assay. Feasibility for detection of non-fish species, especially mammalian or avian DNA originating from dairy products or meat, was tested. The Fish ID solution was successfully used to differentiate pork, beef and lamb. Pork and wild boar resulted in the same pattern. The figure shows the restriction digestion analysis on the Agilent 2100 Bioanalyzer system and the DNA 1000 assay. The gel-like images show the typical pattern obtained after digestion with three restriction enzymes (Dde I, Hae III, Nla III) for two meat samples per species. Of avian samples, only turkey DNA resulted in easily identifiable patterns.

Application note: 5990-8452EN

Food analysis

Discrimination of sturgeon and related species by PCR-RFLP

Group	Species	Csp6 I		Dde I		Hae III		Nla III		
A	<i>Acipenser gueldenstaedtii</i>	325	122	462		180	148	291	193	
	<i>Acipenser medirostris</i>	325	122	458		180	150	293	192	
	<i>Huso dauricus</i>	326	122	463		180	151	294	193	
B ₁	<i>Acipenser transmontanus</i>	446		455		179	149	293	192	
	<i>Acipenser schrenckii</i>	445		457		180	149	292	192	
B ₂	<i>Acipenser stellatus</i>	489		458		180	149	291	193	
C	<i>Huso huso</i>	323	120	455		179	150	77	291	191
D	<i>Polyodon spathula</i>	481		480		323	150	291	165	
E	<i>Oncorhynchus keta</i>	389	86	359	352	121	433	280	192	
F ₁	<i>Acipenser nudiventris</i>	325	93	462		180	149	293	133	52
F ₂	<i>Acipenser ruthenus</i>	349	94	455		178	149	292	134	62
G	<i>Salmo salar</i>	397	85	362	356	121	330	107	46	460

Kit: DNA 1000 kit

Assay: DNA 1000 assay


Application: Species determination of sturgeon is an important issue to enforce regulations, monitor fraud in caviar shipments and assign import fees. The Agilent Fish ID Ensemble is a PCR-RFLP method using a Cytb PCR target sequence and analysis of restriction pattern fragments on the Agilent 2100 Bioanalyzer system and the DNA 1000 assay. This solution was tested for identification of species in fish eggs from caviar shipments. A modified PCR-RFLP protocol provided reproducible and accurate sizing with high resolution, which is critical to assign the right species to an unknown sample.

The table shows the grouping of the samples based on the number of fragments observed with each restriction enzyme. Although not all species in the study could clearly be differentiated from each other, the achieved grouping is sufficient for a first screening in routine analysis, allowing the analysis from sample to result within one working day.

Application note: 5990-8454EN

GMO detection

Development of a multiplex assay for soya



Data kindly provided by CCFRA, UK

Kit: DNA 500 kit*

Assay: DNA 500 assay*

Application: Multiplex assay for genetically modified (GM) soya. The aim was to develop a model assay that could be used to assess the quality of DNA extracted from heat-processed soya flour samples, in particular, to investigate differences in PCR amplification between small DNA targets. A single multiplex PCR assay was developed that enabled three GM soya targets and one control to be analyzed in a single reaction mix. Primer concentration was optimized in order to obtain four PCR products resolved by gel electrophoresis which corresponded in size to the soya lectin gene target of 80 bp, and the EPSPS (5-enolpyruvyl-shikamate-3-phosphate synthase) gene targets of 117 bp, 150 bp and 202 bp respectively. These latter targets are only found in Roundup Ready GM soya. Figure A: Peaks produced by the four PCR targets when analyzed with the 2100 Bioanalyzer system and DNA 500 kit*. Figure B: Analysis of certified reference materials containing known amounts of GM soya.

Application note: 5988-4070EN

* replaced with DNA 1000 kit and assay

GMO detection

DNA stability during food processing

Time at 100°C and pH 3.3 (min)	Amount of PCR product*			
	80 bp	118 bp	150 bp	202 bp
0	100	100	100	100
3	74	77	73	67
6	57	58	21	6
9	36	23	24	15
12	67	33	47	21
15	48	27	16	0
18	0	0	0	0
21	0	0	0	0

* % product determined relative to the amount at 0 minutes

Data kindly provided by CCFRA, UK

Kit: DNA 500 kit*

Assay: DNA 500 assay*


Application: The multiplex PCR assay was applied to soya flour samples containing approx. 1.3 % GM soya and boiled at either pH 3.3, 4.3 or 6.7 for up to 21 minutes. For accurate determination of the quantity of each PCR product, the samples were analyzed with the DNA 500 kit*. The concentration of each PCR product was calculated using the 2100 Bioanalyzer software. At pH 3.3 where an effect of heating time was observed, the amount of each PCR product at each time point was compared to the amount of each product at 0 minutes (Table 2). At pH 3.3, the relative amount of the 80 bp product was reduced to 48 % after 15 minutes and no product was detected at 18 or 21 minutes. After 15 minutes, the relative amounts of products of 118 bp and 150 bp were reduced to 27 % and 16 % respectively and the 202 bp product was not detected. None of the products were detected after 18 or 21 minutes.

Application note: 5988-4070EN

* replaced with DNA 1000 kit and assay

GMO detection

GMO detection by nested multiplex PCR




Kit: DNA 1000 kit
Assay: DNA 1000 assay

Application: GMO detection by multiplex PCR is widely used for soy and corn. Often sequences from the transgene and species specific controls or internal standard are co-amplified by endpoint PCR in a screening procedure. Multiple products can be analyzed with the 2100 Bioanalyzer system at high resolution and sensitivity. Quantitation and comparison of product amounts may already lead to qualification of a positive screening result prior to analysis by expensive quantitative real time PCR.

Application note: 5989-0124EN

Oncology

Tumor cell detection from carcinoma patient blood



Data kindly provided by AdnaGen

Kit: DNA 500 kit*

Assay: DNA 500 assay*

Application: A combined method of specific tumor cell enrichment and a high sensitivity tumor cell detection by multiplex PCR allows analysis of several tumor marker genes. The method is so sensitive that it allows the detection of only a few tumor cells per 5 mL EDTA-blood. The 2100 Bioanalyzer system provides the performance to detect the PCR products with high sensitivity and automated result flagging. This method offers new possibilities for monitoring and prognosis in routine diagnosis, and may facilitate an appropriate selection of patients for adjuvant therapy.

Application note: 5988-9341EN

* replaced with DNA 1000 kit and assay

Oncology

SNP analysis in cancer related P16 gene

Sample	1	2	3	4	5	6	7	8	9	10	11
Genotype316	C	C	CG	CG	CG	CG	CG	C	C	C	G
Genotype356	C	C	CT	CT	CT	C	C	CT	CT	T	C
gel-like view Agilent 2100 bioanalyzer											
Main Band	198	198	197	198	198	198	198	198	199	198	195
Extra Band 1			204	204	204						
Extra Band 2			216	215	215	214	214				

Data kindly provided by SAIC-Frederick

Kit: DNA 1000 kit


Assay: DNA 1000 assay

Application: Mutations in the exon 3 region of P16 gene are closely related to human cancer. A PCR yields 198 bp fragments with single, expected bands or additional, multiple bands in the 2100 Bioanalyzer system analysis. These observations correspond perfectly to genotyping sequencing data of normal and mutant tissues. The pattern of bands is visible due to slower mobility of the heteroduplex formed by heterozygote mutant of the samples. The method provides fast and reliable acquisition of genetic diagnostic data from cancer patients, also on single nucleotide polymorphisms (SNP).

Application note: 5989-0487EN

Oncology

K-ras gene SNP detection



Data kindly provided by SAIC-Frederick

✂ – Restriction
○ – PCR product

Kit: DNA 1000 kit


Assay: DNA 1000 assay

Application: Mutations in the K-ras gene coding 12 region can lead to cancer in different human tissues. A dedicated combination of PCR and specific restrictions (BstNI digest) reveals the underlying single nucleotide polymorphisms (SNPs). The integral element within this test is the rapid and precise analysis of short amplicons (135 bp, see PCR-product lanes above) and fragments (106 bp, visible in lanes labeled with restriction) with the lab-on-a-chip technique. The test was used to ultimately determine a cancer patient's eligibility for a clinical trial for a peptide vaccine.

Application note: 5989-0487EN

Oncology

METH-2 downregulation in lung carcinomas



Data kindly provided by Roy Castle Lung Cancer Research Programme, University of Liverpool, UK

Kit: DNA 1000 kit


Assay: DNA 1000 assay

Application: Microarray analysis reveals under- or over-representations of transcripts. Screening of several cell lines for independent validation of such observations can be done with different techniques such as comparative multiplex PCR. This application shows the downregulation of a characteristic antiangiogenetic factor (METH-2) for a series of patient samples. Expression in normal tissue and tissue from the non small lung carcinomas is compared. Results from the array experiments were confirmed on a broad basis. Fast and convenient analysis with the 2100 Bioanalyzer system with given quantitation capability fit perfectly in such analytical workflow.

Application note: 5989-3514EN

Oncology

Label-free analysis of microsatellite instability in carcinoma



Kit: DNA 1000 kit


Assay: DNA 1000 assay

Application: Microsatellite instability (MSI) is caused by a failure of the DNA mismatch repair system and occurs frequently in various types of cancer. Given that conventional techniques used for MSI detection, for example, polyacrylamide gel electrophoresis (PAGE) or capillary electrophoresis, turned out to be laborious or expensive, this study aimed to develop a simple and efficient procedure of MSI detection. Detection of MSI could be demonstrated by microsatellite loci-associated, well defined deviations in the electropherogram profiles of tumor and non-tumor material and confirmed the classification of the MSI cases performed by conventional technology (95 % concordance rate). Whereas the results of the MSI detection were comparable to conventional techniques, the on-chip electrophoresis on the 2100 Bioanalyzer system was superior in terms of speed, usability and data management.

Application note: 5989-2626EN

Clinical research

Detection of a point mutation in the prothrombin gene with PCR-RFLP



Kit: DNA 1000 kit
Assay: DNA 1000 assay


Application: A single nucleotide polymorphism (SNP), i.e. a point mutation in the prothrombin gene results in a common hereditary predisposition to venous thrombosis. PCR-RFLP was used for reliable detection of a specific SNP in the prothrombin gene. This SNP allows the introduction of a new Hind III restriction site during PCR. If the Hind III restriction site is present this is evidence for prothrombin G20210A mutation. Heterozygous (two bands) and homozygous (one band) genotypes can be specifically identified. The two fragments potentially generated in this PCR-RFLP protocol are only slightly different in size and are difficult to separate with traditional slab gel electrophoresis, even if optimized separation time of 2 hours is applied.

In contrast, the 2100 Bioanalyzer system and the DNA 1000 kit allow precise baseline resolved separation of the expected fragments of 345 and 322 bp. The visualization of the 23 bp fragment cleaved off during the PCR-RFLP protocol is also possible, which further confirms a successful reaction. This is not possible with agarose gel electrophoresis. The results obtained with the 2100 Bioanalyzer system are superior in terms of sizing accuracy, quantitation capability, reproducibility and resolution compared to slab gel analysis.

Application note: 5989-4313EN

Clinical research

Genotyping of *H. pylori*



Data kindly provided by Institute for Pathology, Cologne

Kit: DNA 1000 kit


Assay: DNA 1000 assay

Application: Different allelic variants are associated with different stages of *H. pylori* virulence. Multiplex PCR on five alleles with products in the range of 102 to 301 bp were used to analyze DNA from paraffin embedded tissues. Agarose gel (A) yields only limited distinctiveness, whereas gel-like images (B) and electropherograms (C) show good resolution and superior reproducibility allowing convenient analysis of all desired products in parallel (D). An extended spectrum of prognostic or therapeutic relevant information is now routinely accessible for simultaneous analysis.

Application note: 5989-0078EN

Clinical research

Duplications and deletions in genomic DNA



Data kindly provided by Center for Human and Clinical Genetics, Leiden

Kit: DNA 500 kit*

Assay: DNA 500 assay*


Application: Multiplex amplifiable probe hybridization (MAPH) and multiplex ligation-dependent amplification (MLPA) are high throughput techniques for the detection of reordered genomic segments. These methods include hybridization of amplifiable probes with either stringent washing or ligation events prior to amplification. Exact and reproducible sizing and quantitation of multiple products are important prerequisites which are delivered by the 2100 Bioanalyzer system and lead to quick and simple analysis of genetically related diseases.

Application note: 5989-0192EN


* replaced with DNA 1000 kit and assay

Forensic testing

Optimization of PCR on mtDNA



Amplified areas
in human mtDNA



Kit: DNA 500 kit*

Assay: DNA 500 assay*



Application: Human mitochondrial DNA (mtDNA) is amplifiable even from small or badly degraded samples, even if genomic DNA is not available. Lanes B and C show homogenous PCR products which can subsequently be sequenced for identification. However, careful optimization of PCR parameters, like pH, Mg^{2+} concentration or polymerase amount is necessary and shown in detail in this application note. For example, a high Taq concentration increased the yield but also increased the level of byproducts. PCR for samples in lane D (impurities) and lane A (C-heteroplasmy) need to be improved. The 2100 Bioanalyzer system provides a rapid quantitative analysis over the broad size and concentration range needed for optimization and QC. It has proven to be an indispensable tool for forensic labs.

Application note: 5989-3107EN

* replaced with DNA 1000 kit and assay

Forensic testing

Pitfalls in mtDNA sequencing



Kit: DNA 500 kit*

Assay: DNA 500 assay*


Application: Analysis of the non-coding sequence of human mitochondrial DNA (mtDNA) is performed for the purpose of identification in forensics. PCR amplification of limited or degraded mtDNA is done prior to sequencing. Quantitation and quality control of these PCR products (10-100 ng/mL, homogenous fragment in the range of 200-500 bp) was performed. Difficult PCR templates may cause G-stutters or other unintended byproducts of higher or lower mass (left). This may lead to indistinct sequence readings (right). Therefore, e.g. FBI guidelines enforce a 10 % impurity level at the most. Fulfillment of this prerequisite can be satisfactorily verified with the 2100 Bioanalyzer system.

Application note: 5989-0985EN

* replaced with DNA 1000 kit and assay

Next-generation sequencing

DNA library quantity and quality



Kit: DNA 1000 kit


Assay: DNA 1000 assay

Application: An automated library construction protocol for the Illumina sequencing platform based on the Agilent Bravo Automated Liquid Handling Platform and the Agencourt Ampure PCR purification kit was used for next-generation sequencing and compared to other methods. The 2100 Bioanalyzer system was used to quantify the amount of sample which was retained, as well as to analyze the purity of the sample. The overlaid electrophoretic traces are showing similar size distributions and sample purity using three different purification methodologies.

Application note: 5990-4942EN

Next-generation sequencing

DNA library QC in target enrichment and next-gen sequencing workflows



From [bp]	To [bp]	Corrected area	% of total	Average size [bp]	Size distribution in CV [%]	Concentration [pg/μL]	Molarity [pmol/L]
100	2,000	395.8	51	254	12.4	283.33	1,713.8

Quantitation after 4 PCR cycles

Kit: High Sensitivity DNA kit


Assay: High Sensitivity DNA assay

Application: Here, the High Sensitivity DNA kit was used for quality control of amplified and purified DNA samples from the post-hybridization PCR amplification step prior to sequencing during the SureSelect Target Enrichment workflow. The electropherograms of typical PCR amplified DNA libraries show a typical smear from 150 to 350 nucleotides. The key observation clearly shown in figure 1B, is that the quality of the PCR product depended on the number of PCR cycles performed. After 14 PCR cycles, an additional DNA smear at approximately 500 bp was detected in the electropherogram. The excellent sensitivity of the High Sensitivity DNA kit allowed the amplified DNA to be detected and reliably quantified, even after only four PCR cycles. Thus, the numbers of library PCR cycles can be reduced, removing amplification bias and significantly improving the data quality with increased accuracy.

Application note: 5990-5008EN

Next-generation sequencing

Sizing and quantitation of DNA libraries and fragmented DNA



Kit: High Sensitivity DNA kit


Assay: High Sensitivity DNA assay

Application: The High Sensitivity DNA kit provides sizing and quantitation of DNA fragments and DNA smears in the 50 to 7000 bp size range down to pg/ μ L sensitivity. This is especially useful for sample quality control and the monitoring of critical steps in next-generation sequencing (NGS) workflows, including DNA fragmentation, target enrichment, and DNA library amplification. The analysis of a dilution series from two typical NGS samples, (A) Illumina DNA library and (B) fragmented DNA was performed. For both DNA sample types, the double logarithmic plot demonstrates an excellent linearity with $r^2 = 0.993$. This linear dynamic range depends on the library type and fragment distribution. The broad linear dynamic range of the High Sensitivity DNA kit enables the detection of less abundant products, such as PCR artifacts and impurities.

Application note: 5990-4417EN

Next-generation sequencing

Quality control of FFPE DNA samples



Sample	Average size [bp]	Peak height [bp]	Concentration [ng/ μ L]
FFPE DNA Q	307	264	34.8
FFPE DNA A	309	254	22.6
Frozen tissue DNA	331	264	89.5
Control DNA	348	271	23.2

Average size, peak height, and quantification of precaptured amplified libraries.

Kit: DNA 1000 kit, High Sensitivity DNA kit

Assay: DNA 1000 assay, High Sensitivity DNA assay


Application: The Agilent 2100 Bioanalyzer system was used for quality control of DNA samples from formalin-fixed paraffin-embedded (FFPE) and fresh-frozen tissues prior to and during the SureSelect target enrichment workflow. The figure shows the electropherogram overlay of pre-capture amplified samples after five PCR cycles run on the Agilent 2100 Bioanalyzer system with the Agilent DNA 1000 kit. Similar profiles were observed for all DNA samples. No amplification artifacts or primer dimers were seen. FFPE DNA samples gave comparable results to DNA from fresh-frozen tissue and control cell line DNA, appropriate for downstream sequencing on the Illumina platform.

Reliable DNA on-chip electrophoresis with the Agilent 2100 Bioanalyzer system provided smear profiles and details for library statistics, such as peak heights, average smear size, size distribution, and DNA concentration.

Application note: 5991-0483EN

Next-generation sequencing

Analysis of limited DNA material on the Pippin Prep system



Kit: DNA 1000 kit, High Sensitivity DNA kit

Assay: DNA 1000 assay, High Sensitivity DNA assay


Application: The Agilent 2100 Bioanalyzer system and the High Sensitivity DNA kit complement the Pippin Prep automated size selection workflow (Sage Science, Inc.).

The figure demonstrates the DNA analysis of a restriction digest of *E. coli* genomic DNA to simulate a sheared Pippin Prep input sample. Using the High Sensitivity DNA kit allows the use of only 1.2 ng DNA to achieve a roughly equivalent signal to 50 ng on the DNA 1000 kit.

Even lower DNA amounts (0.375 ng) yield reasonable electropherograms showing input size distribution. This allows tailoring fractionation settings for the Pippin Prep system, and maximizing chances for successful library construction. Afterwards, the 2100 Bioanalyzer system can be used to confirm size ranges, quality, and purity of the Pippin Prep process. The 2100 Bioanalyzer and Pippin Prep systems work well in concert to enable fine control and increase the efficiency of library generation prior to next generation sequencing.

Application note: 5990-8382EN



RNA analysis



Analysis of total RNA

Reproducibility of quantitation

Reproducibility for 12 consecutive runs



Kit: RNA 6000 Nano kit


Assay: Eukaryote Total RNA Nano assay

Application: Alongside the quality control of RNA samples, measurement of RNA concentration is important for (bio-)chemical reactions, such as labeling reactions in the context of microarray experiments. With the RNA 6000 Nano kit good reproducibility can be achieved (here 6 % CV), which is little affected by sample contaminants, such as phenol.

Application note: 5988-7650 EN

Analysis of total RNA

Standardization of RNA quality control (I)



Kit: RNA 6000 Nano kit
Assay: Eukaryote Total RNA Nano assay



Application: Total RNA samples were degraded at 70 °C for varying times and were analyzed with the 2100 Bioanalyzer system to determine the RNA integrity number (RIN). RIN is a reliable software tool to determine the integrity of RNA samples automatically.

The relative quantities of amplicons positioned at the 5' or 3' end of three target genes (GAPDH, HPRT1 and YWHAZ) differ when using differentially degraded RNA templates or highly intact RNA (RIN = 8.9). The data indicates that the integrity of RNA templates can significantly influence the outcome of a real-time quantitative PCR (QPCR) experiment. Furthermore, depending on the position of the amplicon major differences of the calculated fold-change in a comparative quantitation can occur.

Application note: 5989-7730EN

Analysis of total RNA

Standardization of RNA quality control (II)



Kit: RNA 6000 Nano kit



Assay: Eukaryote Total RNA Nano assay

Application: The RNA integrity number (RIN) is calculated by a dedicated software algorithm to assess the quality of RNA preparations. The RIN tool is a major step in the standardization of user-independent RNA evaluation and delivers more meaningful information than simple ratio calculations for ribosomal RNA peaks. It is not influenced by instrument, sample integration and most important, concentration variability, thereby facilitating the comparison of samples and avoiding cost-intensive experiments with low quality RNA preparations. The RIN algorithm is based on a large collection of RNA data of various tissues and qualities. Furthermore, anomalies like genomic DNA contaminations are indicated with weighted error messages (critical/non-critical) to achieve a maximum of reliability.

Application note: 5989-1165EN

Analysis of total RNA

QC of automated RNA isolation from HeLa cells



Kit: RNA 6000 Nano kit


Assay: Eukaryote Total RNA Nano assay

Application: An automated RNA isolation protocol in a 96-well format employing the Bravo Automated Liquid Handling platform and Absolutely RNA 96 Microprep kit from Agilent was used. The total RNA quality of 12 samples isolated from HeLa cell culture was assessed on the 2100 Bioanalyzer system with the RNA 6000 Nano kit. Representative electropherograms of samples from the four corners (top) and the middle row (bottom) from the 96-well tissue culture dish are overlaid. The RNA is fully intact with sharp peaks for 18S and 28S ribosomal RNA. Furthermore, the presence of small RNA, and the absence of genomic DNA contamination can be observed. The RNA integrity numbers ranged from 9.2 to 9.6.

Application note: 5990-3558EN

Analysis of total RNA

RNA quality from freshly frozen or formalin-fixed paraffin-embedded tissue



Kit: RNA 6000 Nano kit

Assay: Eukaryote Total RNA Nano assay


Application: RNA extracted from freshly frozen tissues is optimal for microarray analysis; however, in many cases, formalin-fixed paraffin-embedded (FFPE) tissues are the only samples available. Here, the effectiveness of an FFPE sample protocol optimized for use with Agilent gene expression microarrays was verified. For this purpose, total RNA was extracted from FFPE or fresh frozen tissues from colon tumor (adenocarcinoma) and normal colon and subjected to microarray analysis. The RNA quality was determined with the 2100 Bioanalyzer system. The majority of RNA fragments isolated from the FFPE tissues were between 100–4,000 bp. The low RNA integrity number (RIN ~2.0), indicating degraded RNA, was typical for FFPE extractions as shown in the top two traces.

However, the quality of the gene expression information generated from FFPE samples and from high quality fresh frozen RNA was comparable. Despite the expected loss of sensitivity when analyzing degraded RNA, the differential expression results generated with FFPE and fresh frozen samples are highly concordant.

Application note: 5990-3917EN

Analysis of total RNA

Assessing integrity of plant RNA



Kit: RNA 6000 Nano kit

Assay: Plant RNA 6000 Nano assay


Application: High quality RNA, free of genomic DNA, is a critical determinant for the success of many downstream techniques used in functional genomics, such as RT-PCR and microarray-based experiments. The dedicated Agilent 2100 Bioanalyzer Plant RNA assay included in the 2100 Expert software (version B.02.07 or higher) allows rapid assessment of plant RNA integrity from multiple plant sources and differing degradation states with excellent precision. The figure shows the electrophoretic separation of different plant total RNA using the RNA 6000 Nano kit with the Plant RNA assay. In all samples, the abundant 25S and 18S ribosomal RNA peaks are well resolved and automatically identified by the software. Compared to the root samples, the leaf and lettuce extracts exhibit additional fast migrating peaks corresponding to smaller chloroplast ribosomal RNAs showing that total RNA profiles can vary depending on species and tissue types.

The Plant RNA assay and the RIN algorithm provide a convenient, user-independent assessment of total plant RNA integrity.

Application note: 5990-8850EN

Analysis of total RNA

Genomic DNA contamination




Kit: RNA 6000 Nano kit
Assay: Eukaryote Total RNA Nano assay

Application: Gel representation of a chip run with total RNA samples (mouse brain) spiked with varying amounts of herring sperm genomic DNA before and after treatment with RNase. The left panel shows the intact RNA with broad bands in the low MW region stemming from the genomic DNA. After the RNase digest (right panel) only the DNA bands remain, ranging in intensity according to the amount of DNA spiked into the sample.

Data not published

Low amounts of total RNA

Detection of low levels of RNA



Analysis of mouse brain RNA at two different concentrations

Kit: RNA 6000 Pico kit



Assay: Eukaryote Total RNA Pico assay

Application: The RNA 6000 Pico kit is complementary to the RNA 6000 Nano kit and is suitable for all applications where the amount of RNA (or cDNA) is limited, e.g. for biopsy samples, samples from microdissection experiments, QC of cDNA made from total RNA, microarray samples, etc. Here 2100 Bioanalyzer system results obtained from mouse brain RNA (Ambion) at 200 and 1000 pg/ μ L are shown. By analysis in repetitions the reproducibility of quality control is demonstrated. Detection of 200 pg total RNA could be achieved without problems.

Data not published

Low amounts of total RNA

RNA integrity with the RNA 6000 Pico kit



Kit: RNA 6000 Pico kit


Assay: Eukaryote Total RNA Pico assay

Application: Detection of RNA degradation with the RNA 6000 Pico kit. Sample: mouse liver total RNA (Ambion) concentration: 1 ng. Degradation was accomplished by adding a low amount of RNase. In Figure A the upper electropherogram and gel-like image show the analysis of high quality total RNA with the 18S and 28S subunit as two distinct bands. Figure B shows the analysis of a partially degraded total RNA sample. Many degradation products appear between the two ribosomal bands and below the 18S band.




Data not published

Low amounts of total RNA

RNA quality after staining and microdissection



Laser micro dissection



Check and optimize RNA quality and yield

Kit: RNA 6000 Pico kit


Assay: Eukaryote Total RNA Pico assay

Application: RNA derived from laser-microdissected tissue isolated by the PALM[®]MicroBeam system was shown to be of high quality by convenient analysis with the RNA 6000 Pico assay. RNA-purification kits from different manufacturers and various common staining procedures have been tested and yielded 130-700 pg/ μ L RNA from 1000 cells with different quality (see above). The RNA 6000 Pico kit was well suited to show differences in RNA quality and yield and, therefore, is an ideal tool to optimize and adapt experimental conditions to individual tissue. The experiments were accompanied by a more laborious real time PCR that revealed similar results. Due to its unprecedented sensitivity, the RNA 6000 Pico assay is an indispensable tool for quality control in the context of microdissection experiments, ensuring successful gene expression profiling experiments.

Application note: 5988-9128EN

Low amounts of total RNA

Analysis of minimum RNA amounts



Kit: RNA 6000 Pico kit



Assay: Eukaryote Total RNA Pico assay

Application: The challenge of analysis of minimal amounts of RNA from e.g. laser micro dissections calls for detailed knowledge of extraction conditions. Some commonly used RNA isolation kits and buffer components were assessed in detail. The majority of the kits had no negative effect on the performance of the analysis, whereas, some kits include buffers which lead to shifted, missing and diminished RNA-peaks. In figure A, RNA isolated after microdissection shows lack of the 28S-peak due to high salt concentration introduced during the isolation process. In figure B, a standard RNA was diluted in water and subsequently extracted with a commercially available RNA extraction kit. The original samples (red) and the eluates after extraction are shown. These data show the importance of evaluating the individual method used for RNA extraction to exclude misleading conclusions.

Application note: 5989-0712EN

Low amounts of total RNA

Genomic DNA in RNA extracts of low abundance



Kit: RNA 6000 Pico kit


Assay: Eukaryote Total RNA Pico assay

Application: Laser capture microdissection enables collection of cells from small tissue areas. A low RNA yield is in the nature of the extraction method from such a specimen that usually complicates quality assessment – a fact that can be circumvented by taking advantage of the 2100 Bioanalyzer system capabilities. A comparative study using mouse kidney cryosections showed that on-column DNase digestion is indispensable to obtain a reasonable result for integrity and yield (figure A). Experiments with omitted on-column DNA digestion confirmed that the peak visible in the inter-region consists of genomic DNA which caused overestimation of extracted RNA amounts (figure B).


Application note: 5989-0991EN

Low amounts of total RNA

Low RNA amounts from kidney sections



Renal medulla RNA



Kit: RNA 6000 Pico kit


Assay: Eukaryote Total RNA Pico assay

Application: High sensitivity quality control of RNA samples using the RNA 6000 Pico kit are demonstrated for microdissected samples (0.1 mm³). DNase I digestion revealed that DNA contamination was present in the sample. Removal of DNA revealed total RNA with a low degree of degradation. Under ideal conditions, the RNA Pico assay has a linear response curve and, therefore, allows estimation of RNA concentrations.


Application note: 5988-8554EN

Analysis of mRNA

RNA integrity



Highly enriched
Poly (A)+ RNA



Progressive degradation
of Poly (A)+ RNA

Kit: RNA 6000 Nano kit


Assay: mRNA Nano assay

Application: Progressive degradation of Poly (A)+ RNA. Poly (A)+ RNA (60 ng/ μ L) from cultured Jurkat cells was incubated for 15 minutes at room temperature with very dilute RNase A (1×10^{-6} and 2×10^{-6} mg/mL, respectively). A progressive shift towards shorter fragment sizes can be observed. Even with a mild degradation, the absence of very long transcripts can be noticed.

Application note: 5968-7495EN

Analysis of mRNA

Ribosomal RNA contamination in mRNA samples



Kit: RNA 6000 Nano kit



Assay: mRNA Nano assay

Application: During the isolation of mRNA, varying amounts of ribosomal RNA can remain in a sample. Since the purity of mRNA is of importance for a number of downstream applications, samples should be checked on the 2100 Bioanalyzer system. This slide shows the analysis of 6 commercially available RNA samples from different suppliers. Analysis on the 2100 Bioanalyzer system reveals large differences in the purity of the mRNA samples.

Application note: 5968-7495EN

Analysis of Cy5 labeled samples

Analysis of cRNA with and without dye in gel matrix



Kit: RNA 6000 Nano kit



Assay: mRNA Nano and Cy5 Labeled Nucleic Acids Nano assay

Application: Analysis of Cy5 labeled and non-labeled cRNA samples. Cy5-labeled samples show the combined signals of the fluorescent label and the RNA signal created by the fluorescence of the RNA 6000 dye. If the RNA 6000 dye is omitted from the gel matrix, only the signal created by Cy5 is detected, allowing the determination of dye incorporation after a labeling reaction. Please note that for Cy3 labeled samples the intactness of the sample can be verified, but the dye incorporation can not be checked.

Application note: 5980-0321EN

Analysis of Cy5 labeled samples

Optimization of labeling reactions



Kit: RNA 6000 Nano kit


Assay: Cy5 Labeled Nucleic Acids Nano assay

Application: An experiment was designed to check the influence of Cy5 dCTP concentration on labeling efficiency. Lane 2 represents the negative control (primer omitted from the reaction mixture), while lane 3 shows the analysis of a reaction with a 6-fold increased Cy5 dCTP concentration. A look at the electropherograms reveals that the high Cy5 dCTP concentration not only gave a high peak of unincorporated Cy5, but also the labeling efficiency for longer fragments was very low. This approach allows the optimization of labeling reactions.

Application note: 5980-0321EN

Analysis of Cy5 labeled samples

cRNA fragmentation



Kit: RNA 6000 Nano kit


Assay: Eukaryote Total RNA Nano assay

Application: The RNA 6000 Nano kit can be used to monitor completion of a cRNA fragmentation reaction. In this example, the profile of a Cy5 labeled cRNA sample was monitored at different time points during a fragmentation reaction. It can be seen that after 10 minutes most of the fragments are in the desired size range. After 20 minutes, no further shift of fragmentation can be observed indicating completion of the fragmentation reaction.

Application note: 5988-3119EN

Analysis of T7 RNA transcripts

Size estimation



Kit: RNA 6000 Nano kit


Assay: Eukaryote Total RNA Nano assay

Application: A number of RNA transcripts, ranging from 350 to 1400 nt in size, were analyzed on the RNA 6000 Nano kit. Although the assay runs under native conditions and the transcripts exhibit a certain degree of secondary structure, a good size estimation can be achieved.

Data not published

Analysis of small RNAs

Analysis of miRNA content in total RNA samples



Kit: Small RNA and RNA 6000 Nano kit

Assay: Small RNA and Eukaryote Total RNA Nano assay

Application: Several total RNA samples containing small RNAs (including miRNAs, siRNA, and snRNA) were analyzed using the RNA 6000 Nano assay to determine concentration and quality, including RIN. All total RNA samples were then analyzed with the Small RNA assay to measure miRNA concentration. The relative amount of miRNA was manually calculated as a ratio of the concentration of miRNA in total RNA.

A) total RNA analyzed with the RNA 6000 Nano assay


B) small RNA analyzed with the Small RNA assay

C) miRNA analysis workflow

Application note: 5989-7870EN

Analysis of small RNAs

Analysis of small RNAs from *Drosophila* Schneider



Kit: Small RNA kit


Assay: Small RNA assay

Application: Small RNA samples were prepared from *Drosophila* Schneider cells according to different protocols and were analyzed using polyacrylamide gel electrophoresis, northern blot, or the 2100 Bioanalyzer system. Only 10 ng of total RNA was run on the Small RNA kit. Two sharp peaks, corresponding to 5.8S and 5S RNA, were clearly resolved. Transfer RNA (tRNA) formed a third large peak. A small well-defined peak at 30 nucleotides indicated that 2S RNA could be detected with only 10 ng of total RNA. Therefore, the Small RNA kit was significantly more sensitive than a denaturing 15 % polyacrylamide gel stained with ethidium bromide or SYBR Green II, respectively. Hence, the Small RNA kit appeared to be the method of choice to assess small RNA species contained in a total RNA preparation.

Application note: 5989-8539EN

Analysis of small RNAs

Effect of total RNA quality on mRNA and miRNA expression profiles




Kit: Small RNA kit

Assay: Small RNA assay

Application: Total RNA quality has a significant impact on quantitative PCR results. The MIQE guidelines also emphasize the importance of RNA quality assessment. The importance of RNA quality should also be considered for the investigation of miRNA expression profiles. Eight different miRNAs were quantified using reverse transcription quantitative PCR (RT-qPCR). For all tissues and in all eight genes, a highly significant, negative correlation ($p < 0.001$) between the RNA integrity number (RIN) and the quantitation cycle (Cq) occurred (see figure). With ongoing RNA degradation the Cq increased, showing an impairing influence of the RNA integrity on the performance of the qPCR. Analog experiments with mRNA showed that the expression analysis of miRNAs is influenced by RNA integrity to a lesser degree than mRNA expression. The RIN=5 threshold level for reliable PCR results for mRNA and miRNA was determined.

Application note: 5990-5557EN

Protein analysis



Protein expression

Analysis of cell lysates – protein induction

Protein purification

Comparison between lysate and flow through

Analysis of protein purification

GFP Streptag fusion protein purification

Analysis of column capacity

Analysis of column fractions to optimize conditions

His-tag protein purification with Ni⁺⁺ ZipTips®

Determination of protein expression level and recovery

Optimization of on-column cleavage

Enzymatic removal of His-tags from recombinant proteins

Complementing RP-HPLC protein purification

High sensitivity protein detection

Detection of low protein amounts

Highly specific and sensitive alternative to Western blotting

Pico labeling protocol for samples below 1 ng/μL protein

Antibody analysis

Analysis of antibodies under reducing and non-reducing conditions

High sensitivity analysis of antibodies

QA/QC of IgG under reducing conditions

Quantitation of the half-antibody content in IgG4 preparations

Comparison of SDS-PAGE, CGE and 2100 Bioanalyzer system for humanized monoclonal antibody analysis

Absolute quantitation of IgG

Quality control of stressed antibodies

Analysis of IgG2 under non-reducing conditions

IgG2 analysis under reducing conditions

Monitoring antibody charge variants

Separation of bispecific antibody chains

Protein quantitation

Quantitation of limited protein with an internal standard

Absolute protein quantitation

Food analysis

Rapid wheat varietal identification

Protein analysis in milk

Protein pattern of different transgenic seedlines

Protein – others

Glycoprotein sizing

OFFGEL electrophoresis combined with high-sensitivity on-chip protein detection

Protein quality control prior to MS-analysis


Depletion of high abundant proteins from blood samples

Increased sensitivity by desalting protein samples

Analysis of PEGylated proteins

Protein expression

Analysis of cell lysates – protein induction



Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*


Application: Two cell lysates, induced and non-induced were compared to verify the induction of protein expression. The overlay feature of the 2100 Expert software allows quick sample comparison. The blue electropherogram trace shows the cell lysate highly expressing β-galactosidase (128 kDa).

Data not published

* replaced with Protein 230 kit and assay

Protein purification

Comparison between lysate and flow through



Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*


Application: Cells were lysed using the Pierce B-Per kit and then loaded onto an affinity column. The protein of interest, a 38 kDa protein, should bind to the column and not show up in the flow through. By overlaying the 2 electropherograms from both samples, the lysate and the flow through, it is visible that the protein of interest has bound to the column as expected. In addition, a 77 kDa protein has bound to the column, which could be attributed to unspecific binding or the binding of a dimer.

Data not published

* replaced with Protein 230 kit and assay

Protein purification

Analysis of protein purification



Courtesy of P. Sebastian and S.R. Schmidt GPC-Biotech AG, Martinsried, Germany

Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*


Application: A 18 kDa protein was purified using affinity chromatography. The starting material and the column fraction were analyzed with the protein assay. The protein of interest was determined to be 99 % pure and the concentration in the column fraction was 167 ng/ μ L. The protein assay allows protein purity and concentration to be determined in one step, in addition it calculates protein size for reconfirmation.

Data not published

* replaced with Protein 230 kit and assay

Protein purification

GFP Streptag fusion protein purification



Courtesy of P. Sebastian and S.R. Schmidt GPC-Biotech AG, Martinsried, Germany

Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*


Application: This example shows the analysis of various steps during the purification workflow of a GFP Streptag fusion protein (28 kDa). The protein was expressed in *E. coli* and purified via affinity chromatography with Strep Tactin Poros as the column matrix. The protein assay allows each purification step from the cell lysis to the elution of the purified protein to be monitored and optimized.

Application note: 5988-5025EN

* replaced with Protein 230 kit and assay

Protein purification

Analysis of column capacity



Kit: Protein 200 Plus kit*
Assay: Protein 200 Plus assay*


Application: The binding of a recombinant antibody Fab fragment to a Sepharose column with immobilized Protein G was analyzed to determine the column capacity and prevent column overloading. The protein assay allows this purification step to be monitored and quickly optimized.

Application note: 5988-4022EN

* replaced with Protein 230 kit and assay

Protein purification

Analysis of column fractions to optimize conditions



Courtesy of P. Sebastian and S.R. Schmidt GPC-Biotech AG, Martinsried, Germany

Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*


Application: Different column conditions were tested to optimize the purification conditions for a 30 kDa protein. The column fractions were analyzed for protein purity and concentration to identify the optimal conditions providing a highly purified protein in a good yield. Using the protein assay it was possible to determine the optimum purification conditions in a short time frame.

Data not published

* replaced with Protein 230 kit and assay

Protein purification

His-tag protein purification using Ni⁺⁺ZipTips®



Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*


Application: ZipTips loaded with a Ni²⁺-resin (in development by Millipore) were used to purify a His-tagged protein expressed in *E. coli*. Both the cell lysate and the purified protein were analyzed with the 2100 Bioanalyzer system to demonstrate the performance of the tips. The purification with the tips takes approximately 5 minutes, usually followed by the analysis of the samples with SDS-PAGE analysis which takes a further 2 hours. The SDS-PAGE analysis was substituted by the much faster Protein 200 Plus* assay run on the 2100 Bioanalyzer system.

Data not published

* replaced with Protein 230 kit and assay

Protein purification

Determination of protein expression level and recovery



Sample	Rel. Conc. [ng/ μ L]	% Total	Volume [mL]	POI total [mg]
Total lysate	1,368.7	48.2	20	274
Lysate supernatant	551.2	38.2	20	110
Flow through	325.1	29.9	20	65
Wash 1	435.7	48.3	30	13
Wash 2	405.9	84.1	30	12
SenP2 cleavage	4,784.3	67.4	2	10
Imidazole eluate	9,745.4	55.2	1	10

Kit: Protein 230 kit


Assay: Protein 230 assay

Application: Monitoring the protein of interest (POI) during purification in a fast and reliable manner is fundamental. Here, protein purification of a His6-Sumo3-tagged DNA-binding protein was performed with Immobilized Metal Affinity Chromatography (IMAC). Protein expression and recovery was monitored with the Protein 230 assay. Based on the relative quantitation with the 2100 Bioanalyzer system, the overall yield for the POI was estimated to 30 mg/g biomass. The expression level was 51 % of the total with a relative concentration of 390 ng/ μ L. The table illustrates the protein recovery as determined with the 2100 Bioanalyzer system. Less than half of the POI in the total lysate was soluble and recovered in the lysate supernatant. As anticipated from the bead capacity, 45 mg were depleted from the lysate supernatant, and subsequently released by washes, SenP2 on-column cleavage and final imidazole elution. Together with the flow through, this adds up to 110 mg POI within the lysate supernatant.

Application note: 5990-6153EN

Protein purification

Optimization of on-column cleavage



Sample	Size [kDa]	Rel. Conc. [ng/ μ L]	% Total
Unoptimized on-column cleavage	71.9	310.2	38.3
	81.9	414.7	51.3
Optimized on-column cleavage	70.9	993.0	78.0
	83.2	11.4	0.9

Kit: Protein 230 kit




Assay: Protein 230 assay

Application: Monitoring the protein of interest (POI) during purification in a fast and reliable manner is fundamental. Here, protein purification of a His6-Sumo3-tagged DNA-binding protein was performed with Immobilized Metal Affinity Chromatography (IMAC). The POI was N-terminally fused to His6-Sumo3-tag, bound to Ni beads, and released from the tag and beads by treatment with His-tagged SenP2 protease. Although the protease was used in excess for on-column cleavage, most of the protein was still found to be fused to its tag. Therefore, the on-column cleavage step was optimized by determining the optimal ratio of protein-to-beads-to-protease by monitoring the percentage of cleaved protein versus its tagged precursor with the Protein 230 assay. The molecular weight difference of about 10 kDa is sufficient to achieve a baseline separation between the tagged and tagless protein. As illustrated in the blue electropherogram, almost all tagged protein can be cleaved under the optimized conditions.

Application note: 5990-6153EN

Protein purification

Enzymatic removal of His-tags from recombinant proteins



Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*

Application: For some applications, it might be necessary to remove the His-tag after the protein purification because of its effects on enzymatic activity or protein structure. Here the TAGZyme system (Qiagen) was used to remove the N-terminal His-tag from two different proteins, a GFP variant and a recombinant Interleukin 1 β . Samples were taken at different time points to study the kinetics of the enzymatic cleavage. The dipeptide cleavage can be detected by a size shift on the gel-like images and the electropherograms. The fast analysis with the bioanalyzer allows multiple kinetic studies to be done in one day instead of waiting until the next day for the results from SDS-PAGE analysis.


Poster presented at ABRF Conference, March 2002 by
F. Schäfer, K. Steinert, C. Feckler, J. Drees, and J. Ribbe, QIAGEN GmbH, Hilden, Germany

Application note: 5988-8144 EN

* replaced with Protein 230 kit and assay

Protein purification

Complementing RP-HPLC protein purification



Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*


Application: Protein purification and characterization was carried out facilitating an 1100 Series purification system for reverse phase HPLC assisted by the 2100 Bioanalyzer system. The final polishing of a 56 kDa protein by RP HPLC from a pre-purified sample (starting material, right: red electropherogram and gel) and the analysis of three HPLC-fractions containing the major components are shown (fractions 1-3). No impurity is visible by RP HPLC reanalysis (left chromatogram, fraction 2) of the fraction containing the target protein. However, because the 2100 Bioanalyzer system is an orthogonal technique compared to reverse phase HPLC a 20 kDa protein could be found as an impurity (see insert). The reverse phase HPLC purification leads to a purity of only 76 % for the protein of interest and the 2100 Bioanalyzer system reveals the necessity of further purification.

Application note: 5988-8630EN

* replaced with Protein 230 kit and assay

High sensitivity protein detection

Detection of low protein amounts



Kit: High Sensitivity Protein 250 kit


Assay: High Sensitivity Protein 250 assay

Application: For highest sensitivity, SDS-PAGE gels are commonly silver stained using a laborious procedure with low reproducibility and insufficient quantitation capabilities. The High Sensitivity Protein 250 kit is a superior alternative to silver staining protocols due to reproducible and fast staining, automated separation and data analysis. It overcomes the critical limitations of traditional silver staining of SDS-PAGE gels by providing high sensitivity and a linear dynamic quantitation range of four orders of magnitude combined with excellent reproducibility. The kit analyzes proteins from 10 to 250 kDa down to an on chip concentration of 1 pg/µL. It is based on the detection of fluorescently labeled proteins that are separated by on-chip electrophoresis.

Technical note: 5989-8940EN

High sensitivity protein detection

Highly specific and sensitive alternative to Western blotting



Kit: High Sensitivity Protein 250 kit

Assay: High Sensitivity Protein 250 assay


Application: A targeted protein analysis employing immunoprecipitation in combination with the High Sensitivity Protein 250 assay (IP/HSP250 method) was evaluated. Sensitivity and specificity of the IP/HSP250 method were investigated using *E. coli* lysate spiked with GST-tagged phosphatase and tensin homolog (PTEN). Samples were first fluorescently labeled, then immunoprecipitated, and finally the complexes formed were directly eluted from magnetic Protein A beads and analyzed with the High Sensitivity Protein 250 assay. Electropherograms of samples with 0.1 % to 0.0001 % PTEN and a negative control are shown. The zoom shows the main peak of 0.001 % PTEN. The limit of detection was determined to be 0.001 % or 100 pg PTEN in 10 μ g *E. coli* lysate. For comparison, Western blots were performed. The PTEN blots showed a high non-specific background in all lanes, due to the secondary antibody. A specific band was observed only at 1 % or 100 ng PTEN in 10 μ g *E. coli* lysate.

Thus, the IP/HSP250 method showed both higher sensitivity and specificity than the Western blot, resulting in a 1,000 fold lower limit of detection for PTEN.

Application note: 5990-4097EN

High sensitivity protein detection

Pico labeling protocol for samples below 1 ng/ μ L protein



Kit: High Sensitivity Protein 250 kit


Assay: High Sensitivity Protein 250 assay

Application: The High Sensitivity Protein 250 assay for the 2100 Bioanalyzer system allows detection of fluorescently-labeled proteins down to 1 pg/ μ L on-chip thus offering a sensitivity superior to silver-stained SDS-PAGE. However, the standard fluorescent labeling (Nano) method requires a minimum total protein concentration of 1 ng/ μ L in the initial sample. The alternative Pico labeling protocol extends the applicability to initial protein samples with concentrations from 1 ng/ μ L down to 10 pg/ μ L. The two labeling protocols were compared analyzing 1 ng/ μ L carbonic anhydrase (CA) labeled according to the Nano (red) and the Pico (blue) labeling protocol. When CA labeled with the Nano protocol but without the 200-fold dilution (green) was analyzed, the lower marker peak exceeded the measuring range of the 2100 Bioanalyzer system. If the labeling reaction is performed with the Pico labeling protocol using 100-fold less fluorescent dye, no dilution of the labeled sample is necessary. This allows analyzing highly diluted proteins, such as secreted proteins in culture media, or the analysis of particularly small amounts of protein derived from laser micro dissections.

Application note: 5990-3703EN

Antibody analysis

Analysis of antibodies under reducing and non-reducing conditions



Kit: Protein 230 kit


Assay: Protein 230 assay

Application: The Protein 230 kit allows simultaneous analysis of antibodies under both conditions reduced and non-reduced on the same chip. Analyzing both in one experiment is not possible using SDS-PAGE, as the reducing agent will diffuse within the gel and will also reduce other samples. The Protein 230 assay sensitivity is in the same range as non-colloidal Coomassie staining. A systempeak (SP) as well as upper (UM) and lower (LM) marker are indicated. Under non-reducing conditions the intact antibody (size 149 kDa) represents 83 % of the total detectable protein. Besides other contaminants the light (lc) and heavy chains (hc) as well as half-antibodies are detected and characterized in size and quantity. Under reducing conditions light (28.2 kDa; 31 % of total) and heavy chain (58.2 kDa; 62 % of total) detection is dominant while higher aggregates are also detected.

Demo data file from P230 assay (Expert revision B.02.06 software)

Antibody analysis

High sensitivity analysis of antibodies



Kit: High Sensitivity Protein 250 kit

Assay: High Sensitivity Protein 250 assay


Application: The Protein 250 kit provides analysis of antibodies at highest sensitivity which is equal or better than silver staining sensitivity in SDS-PAGE. The high dynamic range of the Protein 250 assay facilitates impurity detection down to the $\mu\text{g}/\mu\text{L}$ range and low percentage according to regulatory requirements. Both conditions, reduced (blue) and non-reduced (red), can be analyzed parallel on the same chip.

For this preparation the intact antibody (size 147.5 kDa) represents 92.62 % of the total detectable protein under non-reducing conditions. Besides other contaminants the light (lc, 0.74 %) and heavy chains (hc, 0.16 %) as well as half-antibodies (89 kDa; 1.05 %) are precisely quantified in relation of the total protein. Under reducing conditions light (28.9 kDa; 35.14 % of total) and heavy chain (61.2 kDa; 61.51 % of total) detection is dominant while higher aggregates can be identified. For this example antibody non-glycosylated heavy chain fragments (ngHC, 1.05 % of total) could be separated.

Demo data file from P250 assay (Expert revision B.02.06 software)

Antibody analysis

QA/QC of IgG under reducing conditions



Kit: Protein 80 kit



Assay: Protein 80 assay

Application: The Protein 80 kit and the 2100 Bioanalyzer system can be used for on-chip protein electrophoresis of IgG under reducing conditions. Due to the reduction of the disulfide bonds, in this case with DTT, the IgG is separated in its light and heavy chains. The electropherogram and the gel-like image show three main peaks/bands; the light chain, and two heavy chains, with and without glycosylation. No degradation was detected. This type of analysis is required for protein QA/QC to monitor protein stability, e.g. degradation and integrity, and to determine protein contaminations.

Data not published

Antibody analysis

Quantitation of the half-antibody content in IgG4 preparations



Kit: Protein 200 Plus kit*
Assay: Protein 200 Plus assay*

Application: In the given host cell line for antibody production usually up to 30 % of IgG4 is secreted as half molecule (half antibody). The half-molecule has only a single disulfide bond between the heavy and light chains, the inter-heavy chain disulfide bonds are absent. The protein assay allows the half-antibody content in IgG4 preparations to be determined automatically. In addition, the sizing provided by the 2100 Bioanalyzer system compares very well to the theoretical size and is superior to SDS-PAGE in terms of accuracy and reproducibility.




Poster presented at WCBP Conference, January 27-30, 2002 by E. Vasilyeva, H. Fajardo, P. Bove, F. Brown and M. Kretschmer, BIOGEN, Cambridge, MA, USA

Data not published

* replaced with Protein 230 kit and assay

Antibody analysis

Comparison of SDS-PAGE, CGE and 2100 Bioanalyzer system for humanized monoclonal antibody analysis



Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*

Application: The analysis of a humanized monoclonal antibody under reducing condition was compared using 3 different techniques, the 2100 Bioanalyzer system, 4-20 % SDS-PAGE, stained with Coomassie, and capillary gel electrophoresis. All 3 techniques result in a similar separation pattern showing the light and the heavy chain of the antibody. In addition, the determined sizes of the light and heavy chain were comparable for all 3 techniques and compared well to the molecular weights determined by MALDI-TOF (light chain: 23762 Da, heavy chain: 51003 Da). However, the 2100 Bioanalyzer system provides significant time saving compared to the other techniques.



Poster presented at WCBP Conference, January 2002 by
S.H. Bowen, M. Chan, P. McGeehan, J. Smith, L. Inderdass, R. Strouse, M. Schenerman
MedImmune Inc., Gaithersburg, MD, USA

Data not published

* replaced with Protein 230 kit and assay

Antibody analysis

Absolute quantitation of IgG



Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*


Application: The calibration feature of the software allows determination of the absolute antibody concentration in comparison to user defined standards with known concentration, accurate determination of IgG concentrations and carrying out batch comparison during antibody QA/QC.

Application notes: 5988-4021EN and 5988-6576EN

* replaced with Protein 230 kit and assay

Antibody analysis

Quality control of stressed antibodies



Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*



Application: A quality control step in pharmaceutical QA/QC departments is to trigger typical degradation and aggregation patterns for a specific antibody. The given samples from heat stress stability studies show expected protein byproducts after aging at elevated temperatures. The content of heavy and light chain, representing the intact antibody, is reduced by 5 % or 13 % within 1 month or respectively 12 weeks. Excellent reproducibility in the range from 0.6 to 1.7 % CV for this quantization was achieved in a validation study with three different users and two 2100 Bioanalyzer instruments over several days.

Application note: 5988-9648EN

* replaced with Protein 230 kit and assay

Antibody analysis

Analysis of IgG2 under non-reducing conditions



Kit: Protein 230 and High Sensitivity Protein 250 kit




Assay: Protein 230 and High Sensitivity Protein 250 assay

Application: Human myeloma IgG2 was analyzed with the 2100 Bioanalyzer system and the Protein 230 (P230) and High Sensitivity Protein 250 (HSP-250) assay under non-reducing conditions. One representative electropherogram per assay is shown. The intact IgG2 antibody is detected at 156.6 kDa, which is in close agreement with its theoretical molecular mass of about 150 kDa. A unique feature of the HSP-250 assay is the size and concentration measurement beyond the size range of the ladder, that is, 250 kDa. Therefore, high molecular weight aggregates or impurities above 250 kDa are sized and quantified as well. The P230 and HSP-250 assays clearly resolve light chain (LC), heavy chain (HC) and a mixture of LC and HC peaks including the non-glycosylated form of IgG2 (ngAb).

Application note: 5990-5283EN

Antibody analysis

IgG2 analysis under reducing conditions



Kit: Protein 80, 230 and High Sensitivity Protein 250 kit


Assay: Protein 80, 230 and High Sensitivity Protein 250 assay

Application: Human myeloma IgG2 was analyzed with the 2100 Bioanalyzer system and the Protein 80 (P80), the Protein 230 (P230) and the High Sensitivity Protein 250 (HSP-250) assay in the presence of dithiothreitol (DTT) as reducing agent. One representative electropherogram per assay is shown. Under reducing conditions, the IgG2 light chain (LC) and heavy chain (HC) are well resolved with all three available Protein assays. Aggregates of higher molecular weight are observed with the P230 and HSP-250 assays whereas the P80 assay resolves low molecular weight (MW) impurities associated with the IgG2 sample.

Application note: 5990-5283EN

Antibody analysis

Monitoring antibody charge variants



Kit: High Sensitivity Protein 250 kit


Assay: High Sensitivity Protein 250 assay

Application: During production and purification, antibodies can undergo modifications leading to charge heterogeneity, impacting stability, activity and causing immunologically adverse reactions. The investigation of monoclonal antibody (mAb) heterogeneity is critical for quality control, preferably combining charge and size analysis. In the first dimension, 3100 OFFGEL fractionation, an isoelectric focusing technique delivering fractions in solution, was performed and in the second dimension, the fractions were size separated using the High Sensitivity Protein 250 assay. Under native conditions various mAb charge variants exhibiting acidic pI values were detected (A). The electropherogram overlays show the acidic charge variants, the main product and the initial sample load (B). The charge variants differ from the main mAb by a maximum of 4 kDa and represent only 10 % of the total mAb preparation. Therefore, only the combination of both methods allows the detection of the mAb main product and its charge variants.

Application note: 5990-6521EN

Antibody analysis

Separation of bispecific antibody chains



Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*


Application: In general, antibodies are biopharmaceuticals of great interest. Especially bispecific antibodies often require high resolution to allow analysis of both sets of chains (2100 Bioanalyzer system: A, gel like view, resolved heavy chains; C electropherogram). A labor intensive SDS-PAGE could not resolve the heavy chains (B, marked by an arrow) in the given sample. In contrast, the 2100 Bioanalyzer system is a superior tool for antibody quality control since it is a convenient, fast and easy to standardize method which additionally enables quantitative analysis.

Application note: 5988-9651EN

* replaced with Protein 230 kit and assay

Protein quantitation

Quantitation of low protein amounts with an internal standard



Kit: High Sensitivity Protein 250 kit


Assay: High Sensitivity Protein 250 assay

Application: The High Sensitivity Protein 250 assay is based on the detection of fluorescently labeled proteins that are separated by on-chip electrophoresis. The quantitation and sizing of sample peaks is done relative to a ladder as an external standard on the same chip, in contrast to the Agilent Protein 80 and Protein 230 assays where the quantitation is based on internal standards. With the High Sensitivity Protein 250 kit it is also possible to do quantitation based on the addition of a suitable internal standard. In this example a small standard protein, β-Lactoglobulin (βLG, 18.4 kDa), was used. βLG was spiked into a BSA solution prior to the labeling reaction and the analysis was performed according to the standard protocol. The electropherogram shows an overlay of 10 runs and a zoom on the βLG peak in the insert. The quantitation data for the βLG and BSA monomer is summarized in the table. Quantitation reproducibility is improved using an internal standard, because well-to-well variations for sample injection and matrix influences during staining are excluded.

Application note: 5989-8941EN

Protein quantitation

Absolute protein quantitation



Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*


Application: A comparative analysis of different techniques used for absolute protein quantitation was performed analyzing 3 different proteins (CA, BSA, OV) in 4 different concentrations (40 – 1250 µg/mL). The same samples were quantitated using the 2100 Bioanalyzer system, two commonly used total protein quantitation assays, Lowry and Bradford, and SDS-PAGE, stained with Coomassie. The relative standard deviation (CV) and the error compared to the target concentration were determined. A comparison shows that the CV and error for the 2100 Bioanalyzer system are better than for the SDS-PAGE by a factor of 2. This data demonstrates that the 2100 Bioanalyzer system is a viable alternative for protein quantitation. It allows the quantitation of individual proteins and simultaneous determination of protein purity and size.

Application notes: 5988-4021EN and 5988-6576EN

* replaced with Protein 230 kit and assay

Food analysis

Rapid wheat varietal identification



Kit: Protein 230 kit


Assay: Protein 230 assay

Application: Total wheat proteins (including high molecular weight glutenins, HMW) were extracted from individual grains. The extract samples were separated on the 2100 Bioanalyzer system. The Protein 230 assay produced well-resolved protein profiles, suitable for varietal discrimination (A). Electropherogram profiles were processed using the Phoretix 1D Advanced and 1D Database (Nonlinear Dynamics) software for pattern-matching purposes. Replicates of three different wheat varieties can be correctly grouped in a dendrogram (B). The study has demonstrated that using the 2100 Bioanalyzer system with the Phoretix system offers a standardized, objective method for rapid varietal discrimination. The ease of use and short analysis times of less than one hour from sampling to 2100 Bioanalyzer system result makes it most suitable for mill intake use.

Application note: 5989-7735EN

Food analysis

Protein analysis in milk




Kit: Protein 80 kit
Assay: Protein 80 assay

Application: Qualitative and quantitative information about the main protein species in milk is an important quality criterion. Cow, goat, and sheep milk were analyzed by on-chip electrophoresis with the Agilent Protein 80 kit and the Agilent 2100 Bioanalyzer system. The figure shows the analysis of commercially available cow milk. The five major milk proteins were analyzed separately and overlaid with the milk electropherogram to identify the milk proteins. The electropherograms allow an easy and reliable comparison of different samples. Furthermore, milk from different species could be distinguished based on their protein pattern which facilitates incoming inspection in routine labs. The Agilent 2100 Bioanalyzer system and the Protein 80 kit provide a fast, reproducible and robust qualitative and quantitative method for milk analysis.

Application note: 5990-8125EN

Food analysis

Protein pattern of different transgenic seedlines



SeedLine	Extracted protein level µg/mL	7S/11S Ratio
Control	14,000	0.39±0.004(n=5)
7S	5,200	3.4
11S	14,000	0.04
Unknowns	13,000	0.72±0.1(n=20)

Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*




Application: Determination of protein size and concentration with sufficient accuracy and precision allows the highly efficient characterization of transgenic seed lines. Expressed protein was available after grinding, extraction of seeds and dilution with buffer. Electropherograms were evaluated by integration of regions specific for 7S or 11S seed storage proteins. The elevated ratio of 7S/11S for the analyzed unknown line shows significant changes in the expression profile in comparison with the control.

Application note: 5988-9441EN

* replaced with Protein 230 kit and assay

Protein – others

Glycoprotein sizing



Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*


Application: Due to large carbohydrate moieties glycosylated proteins can differ in amount of incorporated SDS and shape of the protein/SDS-complex from non-glycosylated proteins. This may lead to different migration times in SDS-PAGE, as well as in the Protein 200 Plus assay* run on the 2100 Bioanalyzer system. The data compare deglycosylation of a mixture of three proteins (electropherogram on the left) with a commercial N-Glycosidase F Deglycosylation kit. Sizing experiments comparing glycosylated and non glycosylated states for additional proteins are compared and summarized on the right. Such an approach avoids misinterpretation of sizing due to glycosylation and allows detection of a posttranslational modification of unknown proteins.

Application note: 5989-0332EN

* replaced with Protein 230 kit and assay

Protein – others

OFFGEL electrophoresis combined with high-sensitivity on-chip protein detection



Kit: Protein 250 kit



Assay: Protein 250 assay

Application: Two-dimensional gel electrophoresis (2D-GE) is a tedious and time-consuming procedure. The combination of fractionation with the Agilent 3100 OFFGEL fractionator based on isoelectric point (pI), and high-sensitivity on-chip electrophoresis with the 2100 Bioanalyzer system performs an analytical 2D-GE-type analysis with excellent sensitivity and reproducibility. The ease of use is significantly improved compared with the traditional 2D-GE method. The OFFGEL fractionation method takes advantage of the impressive resolving power of immobilized pH gradient-based isoelectric focusing (IPG-IEF) and delivers the sample in liquid phase. The High Sensitivity Protein 250 assay is based on the detection of fluorescently labeled proteins that are electrophoretically separated, allowing separation of proteins from 10 to 250 kDa, with a sensitivity equivalent to or better than silver staining, and delivers a linear dynamic range across four orders of magnitude.

Application note: 5989-8419EN

Protein – others

Protein quality control prior to MS-analysis



Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*


Application: By applying soft ionization methods like MALDI or ESI mass spectroscopy (MS), mass information from proteins up to 300 kDa can be obtained. However, proper sample preparation is an important precondition. Concentration, purity and assumed size are valuable ex ante information usually given by biochemists to MS-analysis services. Two different examples for proteins analyzed by an LC/MS-method (right panel) with good results for sample A and discrepancies for sample B are shown. An analysis of the samples with the Protein 200 Plus assay* (left panel) showed an impure protein preparation for sample B. Here, two major peaks at higher masses (66 and 132 kDa) potentially caused by aggregates were encountered. The protein of interest (30 kDa) yielded high noise background in the MS. A quality check of the sample with the Protein 200 Plus assay, therefore, may avoid an unproductive MS analysis or data evaluation.

Application note: 5989-0771EN

* replaced with Protein 230 kit and assay

Protein – others

Depletion of high abundant proteins from blood samples



Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*


Application: Depletion of high abundant proteins in human blood plasma was facilitated by the Agilent Multiple Affinity Removal System. Unprocessed serum, the flow through (i.e. the immunodepleted serum) and the bound proteins after specific elution were analyzed on the 2100 Bioanalyzer instrument. Equivalent amounts were analyzed and resulting electropherograms (left) and gel like images (right) show in comparison the completeness of separation. The 2100 Bioanalyzer system, in combination with the Protein 200 Plus kit*, proves to be an excellent method for evaluation of serum processed with albumin removal kits. The system offers a rapid and accurate method to detect proteins both quantitatively and qualitatively.

Data not published

* replaced with Protein 230 kit and assay

Protein – others

Increased sensitivity by desalting protein samples



Kit: Protein 200 Plus kit*

Assay: Protein 200 Plus assay*


Application: Protein purification steps, such as ion exchange chromatography, often implicate high salt concentrations. Nevertheless, quantitation of these samples is effective since the upper marker serves as internal protein standard and is subjected to the same conditions. However, under high salt conditions lower amounts of protein are injected into the microfluidic channels for analysis. Therefore, the sensitivity can be increased by usage of convenient desalting spin columns. Comparably higher and homogeneous peak intensities are obtained while the recovery after such treatment is good.

Application note: 5989-0228EN

* replaced with Protein 230 kit and assay

Protein – others

Analysis of PEGylated proteins



Kit: Protein 230 kit

Assay: Protein 230 assay

Application: Monitoring the degree of PEGylation is important to achieve a balance between retention of bioactivity, stability, and immunogenicity of the PEGylated species. SDS-PAGE is routinely used for PEGylation analysis, however, it cannot precisely identify the extent of the PEGylation. In addition, it can be tedious and time consuming. The Agilent 2100 Bioanalyzer system was evaluated as an alternative to SDS-PAGE. A chimeric protein was subjected to PEGylation at 30 °C for 3 hours using different concentrations of two PEGylating reagents, mPEG-pNP (MW 5000) and mPEG-tresylate (MW 5000). Levels of PEGylation was monitored by SDS-PAGE (A) and the 2100 Bioanalyzer system (B). The profiles of PEGylated protein obtained with the Agilent 2100 Bioanalyzer system were superior and allowed differentiation between levels of PEGylation. The Agilent 2100 Bioanalyzer with Protein 230 assay is an easy-to-use tool that provides high resolution of PEGylated protein species, allowing efficient optimization of reaction conditions as well as fast and quantitative monitoring of production batches.

Application note: 5990-9593EN

AGILENT 2200 TAPESTATION SYSTEM

Building on the success of the 2100 Bioanalyzer system, the Agilent 2200 TapeStation system offers scalable throughput and rapid results, making it an ideal solution for quality control of biological samples in next-generation sequencing, microarray and qPCR workflows, as well as in protein purification and antibody production.

The 2200 TapeStation system uses credit card-sized, ScreenTape consumables that are available for DNA, RNA and protein applications. Sample analysis has never been so easy – simply load the 2200 TapeStation instrument with ScreenTape, loading tips, and your samples in tube strips or 96-well microtiter plates – and you will be reviewing results in as little as 1 minute per sample!

The D1K ScreenTape assay facilitates the separation and analysis of DNA fragments, fragmented DNA and NGS libraries from 35 bp to 1000 bp. Choose between the D1K ScreenTape and the High Sensitivity D1K ScreenTape depending on the sensitivity requirements of your application.


The Genomic DNA ScreenTape assay was developed for the separation and analysis of genomic DNA from 200 to > 60,000 bp. It provides accurate quantification data and allows a quality assessment of genomic DNA starting material, which is ideal for NGS and array comparative genomic hybridization (aCGH) workflows.

The R6K ScreenTape assay provides efficient and reliable RNA analysis for RNA characterization and quality assessment. The RNA integrity number equivalent (RIN[®]) delivers an instant and objective evaluation of eukaryotic total RNA degradation. Choose between the R6K ScreenTape and the High Sensitivity R6K ScreenTape depending on the sensitivity requirements of your application.



DNA – PCR product analysis

Automated DNA sample analysis from 96-well sample plates



Kit: D1K ScreenTape and Reagents


Application: Techniques for DNA sample QC, such as slab gel and capillary electrophoresis, in NGS and microarray studies are often slow and involve more manual steps causing a bottleneck. The 2200 TapeStation system with the D1K ScreenTape assay meets this demand by providing electrophoretic analysis of DNA samples from a 96-well sample plate in less than 100 minutes and automates result interpretation with easy data reporting. The analysis with the D1K ScreenTape assay for quantification and sizing of DNA samples from 96-well sample plates is accurate and highly reproducible. The figure shows how the 2200 TapeStation system was used to optimize the annealing temperature for eliminating amplification of non-target sequence. A PCR amplification using LPIN3 primer was carried out over a temperature gradient from 50 to 60 °C. Non-specific amplification was observed at various temperatures. The amplicon with highest purity and quantity is obtained at 55.7 °C.

This consistency of the D1K ScreenTape assay, even over extended periods when using a 96-well plate, is ideal for successful optimization of PCR amplification reactions on the Agilent SureCycler 8800 instrument.

Application note: 5991-2350EN

DNA – Next-generation sequencing

QC of DNA from FFPE tissues for next-generation sequencing



	Avg size (bp)	Concentration (ng/ μ L)
1. Patient FFPE	305	31.7
2. Patient FFPE	301	20.0
3. Patient fresh frozen tumor	334	89.3
4. Control Cell Line	329	12.4

Kit: D1K ScreenTape and Reagents; High Sensitivity D1K ScreenTape and Reagents


Application: A huge number of formalin-fixed paraffin-embedded (FFPE) tissue samples is archived in biobanks worldwide and are valuable resources for molecular genetic studies. However, the challenges of DNA extraction from FFPE tissues result in low amounts of usable high quality material for downstream assays. Given this, assessing quality of DNA samples is essential. Analysis for sample size distribution, purity, and concentration throughout the sample preparation workflow is critical for successful next-generation sequencing (NGS). The sample preparation workflow for SureSelect target-enrichment includes 3 steps requiring quality control (QC): post-shearing, pre-capture amplification, and post-capture amplification. The figure shows the electropherogram profiles of pre-capture amplified libraries run on the 2200 TapeStation system with the D1K ScreenTape assay. The table summarizes the obtained average size and quantitation of libraries.

The 2200 TapeStation system is an ideal tool for NGS sample quality control of FFPE samples, providing peak profiles, the calculation and determination of library statistics, specifically for peak size, distribution and concentration. With its high degree of automation the 2200 TapeStation system offers fast results, full scalability, and a simplified QC workflow.

Application note: Data not published

DNA – Genomic DNA analysis

Analysis of high molecular weight genomic DNA



Kit: Genomic DNA ScreenTape and Reagents

Application: Genomic DNA (gDNA) is very often used as experimental starting material. Information on gDNA quality and concentration and therefore the efficiency of the used extraction method is important for successful downstream applications, particularly array comparative genome hybridization and next generation sequencing.


The Genomic DNA ScreenTape assay in conjunction with the 2200 TapeStation system was used for sizing, quantification and quality assessment of gDNA samples. The figure shows the gel image representing triplicate analyses of gDNA samples obtained from four extraction kits. The bubble chart shows the relationship between the average molecular weight size of the overall sample trace versus the molecular weight size of the peak maxima, with the bubble size indicating the concentration of the sample for each of the four extraction kits. This information illustrates the overall extraction efficiency for each kit, with a large bubble located in the top right hand quadrant of the chart the most desirable outcome.

The 2200 TapeStation system in conjunction with the Genomic DNA ScreenTape assay provides an excellent solution for assessing the quantity, integrity and overall quality of genomic DNA starting material from a single run.

Application note: 5991-1797EN

RNA – Analysis of total RNA

RNA quality control with the RIN^e quality metric



Kit: R6K ScreenTape and Reagents

Application: The Agilent 2100 Bioanalyzer system and its RNA Integrity Number (RIN) is a well-established tool for the assessment of RNA quality. The Agilent 2200 TapeStation system builds on this success with faster analysis times, more flexibility of throughput at a constant cost per sample, and 96-well plate capabilities.


The figure shows the analysis of different total RNA samples on the 2200 TapeStation system and the R6K ScreenTape (A). The gel image shows the separation profile of each of the individual samples showing 28S, 18S, small rRNAs and lower marker. The RNA quality is presented as RIN^e value for each individual sample below the gel image. In addition a representative electropherogram of total RNA from HeLa cells treated with PMA is shown (B).

The Agilent 2200 TapeStation system offers an easy to use system for analyzing RNA samples with minimum manual intervention. RIN^e generated by the 2200 TapeStation system highly correlates with the RIN quality assessment of the 2100 Bioanalyzer system. Combined with the added advantages of high throughput compatibility, ease of use, minimal manual intervention, and constant cost per sample regardless of throughput, the 2200 TapeStation system is ideal for the assessment of RNA quality.


Application note: 5991-0023EN

Protein – Antibody analysis

Reproducible integrity and purity testing of antibodies



P200 ScreenTape assay



SDS-PAGE/Coomassie

Kit: P200 ScreenTape and Reagents

Application: Monoclonal antibodies (mAbs) are an essential tool in the life sciences and are used in many cutting-edge applications such as antibody therapy, diagnostic testing, and protein microarray analysis. Stringent testing of mAbs throughout production is required. The figure shows the separation of three different mAb preparations with the P200 ScreenTape, or 4-12% SDS-PAGE gel stained with Coomassie blue under reducing and non-reducing conditions. The 2200 TapeStation software automatically presented analyzed protein results for 16 samples in less than one minute per sample. The mAb samples that were run in reducing conditions (lanes 2 to 4) showed two major peaks corresponding to mAb heavy and light chains. The analysis under non-reducing conditions (lanes 5 to 7) showed one major peak (whole Ig) smaller bands corresponding to various heavy/light chain combinations. Similar separations were achieved with SDS-PAGE after 3 hours, including buffer preparation, manual gel loading, and staining.

The P200 ScreenTape assay provides the necessary reliability and speed of analysis for quality measurements. It offers a more accurate, precise, and convenient method for mAb analysis than traditional SDS-PAGE methods, as both reducing and non-reducing analysis can be run side-by-side and as few as two samples can be analyzed without sacrificing unused lanes. The 2200 TapeStation system makes manufacturing quality control procedures such as in-process QC and stability testing for mAbs more efficient, more cost effective and less labor intensive.

Application note: 5990-9052EN

Literature

Application notes and other publications

To download an application note or to find other literature on the 2100 Bioanalyzer and 2200 TapeStation systems visit our websites:

www.agilent.com/genomics/bioanalyzer and

www.agilent.com/genomics/tapestation

Cell Applications

Publication number

Protein expression monitoring

Detection of antibody-stained cell surface and intracellular protein targets with the Agilent 2100 Bioanalyzer 5988-4322EN

Detecting cell surface proteins with the Agilent 2100 Bioanalyzer by on-chip antibody staining 5988-7111EN

A new tool for routine testing of cellular protein expression: integration of cell staining and analysis of protein expression on a microfluidic chip-based system 5989-0021EN

Flow cytometric analysis of a limited number of cells using the Agilent 2100 Bioanalyzer 5989-0746EN

A new method for the calculation of baculovirus titre using the Agilent 2100 Bioanalyzer and the flow cytometry kit 5989-1644EN

Cytometric analysis of upregulated functional gene expression in primary cells by on-chip staining 5989-2718EN

Transfection efficiency monitoring

Monitoring transfection efficiency by green fluorescence protein (GFP) detection 5988-4320EN

Detecting cell surface proteins with the Agilent 2100 Bioanalyzer by on-chip antibody staining 5988-7111EN

Monitoring transfection efficiency in cells using an on-chip staining protocol 5988-7296EN

Flow cytometric analysis of human primary cells using the Agilent 2100 Bioanalyzer and on-chip staining 5988-8154EN

Apoptosis detection

Apoptosis detection by annexin V and active caspase-3 with the Agilent 2100 Bioanalyzer 5988-4319EN

A fast protocol for apoptosis detection by annexin V with the Agilent 2100 Bioanalyzer 5988-7297EN

Measuring multiple apoptosis parameters with the Agilent 2100 Bioanalyzer 5988-8028EN

Detection of apoptosis in primary cells by annexin v binding using the Agilent 2100 Bioanalyzer 5989-2934EN

Gene silencing in cell culture

siRNA transfection optimization with the Agilent 2100 Bioanalyzer 5988-9872EN

Confirming gene silencing mechanism by pGFP/GFP22 – siRNA co-transfection 5989-0103EN

Cells – others

Cell fluorescence assays on the Agilent 2100 Bioanalyzer – general use	5988-4323EN
Identification of red and white blood cells from whole blood samples using the Agilent 2100 Bioanalyzer	5989-7171EN
Bacteria counts with on-chip flow cytometry on the Agilent 2100 Bioanalyzer system	5991-2582EN

DNA Applications

Publication number

Restriction digest analysis

Quantitative analysis of PCR fragments with DNA 7500 kit	5968-7496EN
High precision restriction fragment sizing with DNA 12000 kit	5968-7501EN

PCR product analysis

Quantitative analysis of PCR fragments with DNA 7500 kit	5968-7496EN
High resolution DNA analysis with the DNA 500 and DNA 1000 kits	5988-3041EN
Highly efficient multiplex PCR using novel reaction chemistry	5988-9342EN
Optimizing real-time quantitative PCR experiments with the Agilent 2100 Bioanalyzer	5989-7730EN
Interference of SybrGreen in detecting PCR amplicons with the Agilent 2100 Bioanalyzer	5989-4458EN
Complete automation of the Stratagene StrataPrep 96 PCR Purification kit with the Agilent Bravo Automated Liquid Handling platform and Agilent Automated Centrifuge	5990-3948EN

Gene expression analysis

Quantitative end-point RT-PCR gene expression using DNA 7500 kit	5988-3674EN
Semiquantitative reverse transcription-polymerase chain reaction with the Agilent 2100 Bioanalyzer	5988-4556EN
DNA QC for oligonucleotide array CGH (aCGH) with the Agilent 2100 Bioanalyzer	5989-2487EN

Food Food analysis

Development of meat speciation assays using the Agilent 2100 bioanalyzer	5988-4069EN
Analysis of genetically modified soya using the Agilent 2100 Bioanalyzer	5988-4070EN
Nested multiplex PCR for the determination of DNA from genetically modified corn and soy beans using the Agilent 2100 Bioanalyzer	5989-0124EN
Use of the Agilent 2100 Bioanalyzer for basmati rice authenticity testing	5989-6836EN
Strawberry and raspberry fruit differentiation using the Agilent 2100 Bioanalyzer	5990-3327EN
Identification of different meat species by the Agilent Fish ID solution on the Agilent 2100 Bioanalyzer	5990-8452EN
Discrimination of sturgeon and related species by PCR-RFLP using the Agilent 2100 Bioanalyzer system	5990-8454EN

GMO detection

Analysis of genetically modified soya using the Agilent 2100 Bioanalyzer	5988-4070EN
Detecting genetically modified organisms with the Agilent 2100 Bioanalyzer	5988-4847EN
Nested multiplex PCR for the determination of DNA from genetically modified corn and soy beans using the Agilent 2100 Bioanalyzer	5989-0124EN

Oncology

Agilent 2100 Bioanalyzer replaces gel electrophoresis in prostate cancer research	5988-1086EN
Sensitive detection of tumor cells in peripheral blood of carcinoma patients by a reverse transcription PCR method	5988-9341EN
Mutation detection for the K-ras and P16 genes	5989-0487EN
Label-free analysis of microsatellite instability in colorectal carcinoma by on-chip electrophoresis	5989-2626EN
Measuring the METH-2 promoter hypermethylation and transcript downregulation in non-small cell lung carcinomas with the Agilent 2100 Bioanalyzer	5989-3514EN

Clinical research

Microfluidic analysis of multiplex PCR products for the genotyping of Helicobacter pylori	5989-0078EN
Rapid detection of genomic duplications and deletions using the Agilent 2100 Bioanalyzer	5989-0192EN
Detection of a point mutation in the prothrombin gene with the Agilent 2100 Bioanalyzer	5989-4313EN

Forensic testing

Use of the Agilent 2100 Bioanalyzer and the DNA 500 in the analysis of PCR amplified mitochondrial DNA	5989-0985EN
Using the Agilent 2100 Bioanalyzer to optimize the PCR amplification of mitochondrial DNA sequences	5989-3107EN

Next-generation sequencing

Automation of Agencourt AMPure Purification kit for the purification of Next-Generation Sequencing sample preparation reactions on the Agilent Bravo Automated Liquid Handling platform	5990-4942EN
Improving sample quality for target enrichment and next-gen sequencing with the Agilent High Sensitivity DNA Kit and the Agilent SureSelect Target Enrichment platform	5990-5008EN
Analysis of limited DNA material on the Pippin Prep system	5990-8382EN
Quality control of FFPE DNA samples	5991-0483EN

Technical

Comparing the Agilent 2100 Bioanalyzer performance to traditional DNA analysis	5980-0549EN
New Series II DNA assays for the Agilent 2100 Bioanalyzer	5989-4652EN
Performance characteristics of the High Sensitivity DNA kit for the Agilent 2100 Bioanalyzer	5990-4417EN

RNA Applications

Publication number

Analysis of total RNA

Analysis of total RNA using the RNA 6000 kit	5968-7493EN
Characterization of RNA quality using the RNA 6000 kit	5980-0472EN
The total RNA story	5988-2281EN
Quantitation comparison of total RNA using the Agilent 2100 Bioanalyzer, ribogreen analysis and UV spectrometry	5988-7650EN
Advancing the quality control methodology to assess isolated total RNA and generated fragmented cRNA	5988-9861EN
Stringent RNA quality control using the Agilent 2100 Bioanalyzer	5989-1086EN
RNA Integrity Number (RIN) - Standardization of RNA quality control	5989-1165EN
High-Purity RNA isolation from a wide range of plant species and tissue types	5989-2271EN
Isolation of high purity total cellular RNA from muscle tissues	5989-2312EN
Optimizing real-time quantitative PCR experiments with the Agilent 2100 Bioanalyzer	5989-7730EN
Automation of Stratagene Absolutely RNA 96 Microprep Kit with the Bravo Automated Liquid Handling Platform	5990-3558EN
Gene expression microarray analysis of archival FFPE samples	5990-3917EN
Assessing integrity of plant RNA with the Agilent 2100 Bioanalyzer	5990-8850EN

Low amounts of total RNA

Gene expression profiling of esophageal cancer using laser capture microdissected samples	5989-1088EN
High sensitivity quality control of RNA samples using the RNA 6000 Pico kit	5988-8554EN
Quality assurance of RNA derived from laser microdissected tissue samples obtained by the PALM MicroBeam system using the RNA 6000 Pico kit	5988-9128EN
Successful analysis of low RNA concentrations with the Agilent 2100 Bioanalyzer and the RNA 6000 Pico kit	5989-0712EN
Assessing genomic DNA contaminations of total RNA isolated from kidney cells obtained by laser capture microdissection using the Agilent RNA 6000 Pico assay	5989-0991EN

Analysis of mRNA

Analysis of messenger RNA using the RNA 6000 kit	5968-7495EN
Characterization of RNA quality using the RNA 6000 kit	5980-0472EN
Interpreting mRNA electropherograms	5988-3001EN

Analysis of Cy5 labeled samples

Analysis of Cy5-labeled cRNAs and cDNAs using the RNA 6000 kit	5980-0321EN
Optimizing cRNA fragmentation for microarray experiments using Agilent 2100 Bioanalyzer	5988-3119EN
Monitoring extraction efficiency of small RNAs with the Agilent 2100 Bioanalyzer and the Small RNA kit	5990-6935EN

Analysis of Small RNA

Analysis of miRNA content in total RNA preparations using the Agilent 2100 Bioanalyzer	5989-7870EN
Analysis of small RNAs from <i>Drosophila Schneider</i> cells using the Small RNA assay on the Agilent 2100 Bioanalyzer	5989-8539EN
RNA quality control in miRNA expression analysis	5990-5557EN

Technical

Comparing performance of the Agilent 2100 Bioanalyzer to traditional RNA analysis	5980-2206EN
New Series II RNA 6000 assays for the Agilent 2100 Bioanalyzer	5989-4571EN
Performance of the Agilent Small RNA assay	5989-7002EN

Protein Applications

Publication number

Protein purification

Optimization of protein purification using the Agilent 2100 Bioanalyzer (analysis of column capacity)	5988-4022EN
Comparison of different methods for purification analysis of a green fluorescent strep-tag fusion protein	5988-5025EN
Using the Agilent 2100 Bioanalyzer for analysis of His-tag removal from recombinant proteins	5988-8144EN
Protein purification and characterization using the Agilent 1100 Series purification system and the Agilent 2100 Bioanalyzer	5988-8630EN
Evaluation of albumin removal using the Agilent 2100 Bioanalyzer	5988-9911EN
Monitoring protein fate during purification with the Agilent 2100 Bioanalyzer	5990-6153EN
Monitoring antibody charge variants using a combination of 3100 OFFGEL fractionation by isoelectric point and high sensitivity protein detection with the 2100 Bioanalyzer	5990-6521EN
High throughput purification of human IgG using the Agilent Bravo for protein purification and AssayMAP protein A cartridges	5990-7203EN

High sensitivity protein detection

Performance characteristics of the High Sensitivity Protein 250 assay for the Agilent 2100 Bioanalyzer	5989-8940EN
Additional Pico protocol for the High Sensitivity Protein 250 assay with the Agilent 2100 Bioanalyzer	5990-3703EN
Immunoprecipitation and the High Sensitivity Protein 250 Assay	5990-4097EN

Antibody analysis

Absolute quantitation of IgG with the Protein 200 kit	5988-4021EN
Comparison of different protein quantitation methods	5988-6576EN

Quality control of antibodies using the Agilent 2100 Bioanalyzer and the Protein 200 Plus assay	5988-9648EN
Analysis of bispecific antibodies using the Agilent 2100 Bioanalyzer and the Protein 200 Plus assay	5988-9651EN
Protein analysis with the Agilent 2100 Bioanalyzer – An overview of the protein kit portfolio	5990-5283EN
Monitoring antibody charge variants using a combination of Agilent 3100 OFFGEL Fractionation by isoelectric point and high sensitivity protein detection with the Agilent 2100 Bioanalyzer	5990-6521EN
Protein quantitation	
Absolute quantitation with the Protein 200 kit	5988-4021EN
Comparison of different protein quantitation methods	5988-6576EN
Quantification strategies with the High Sensitivity Protein 250 assay	5989-8941EN
Protein sizing and quantitation with the Agilent Protein 80 and Protein 230 kits on the Agilent 2100 Bioanalyzer	5990-5721EN
Immunoprecipitation combined with microchip capillary gel electrophoresis for detection and quantification of b-galactosidase from crude <i>E. coli</i> cell lysate in bioprocess technology	5990-7008EN
Food analysis	
Characterization of transgenic soybean seedlines by protein expression with the Agilent 2100 Bioanalyzer	5988-9441EN
Rapid wheat varietal identification using the Agilent 2100 Bioanalyzer and automated pattern-matching	5989-7735EN
Milk protein analysis with the Agilent 2100 Bioanalyzer and the Agilent Protein 80 kit	5990-8125EN
Protein – others	
Increased sensitivity by desalting protein samples prior to analysis on the Agilent 2100 Bioanalyzer	5989-0228EN
Glycoprotein sizing on the Agilent 2100 Bioanalyzer	5989-0332EN
Alternative to 2D gel electrophoresis – OFFGEL electrophoresis combined with high-sensitivity on-chip protein detection	5989-8419EN
Analysis of PEGylated proteins using the Agilent 2100 Bioanalyzer	5990-9593EN
Using the Agilent 2100 Bioanalyzer for quality control of protein samples prior to MS-analysis	5989-0771EN
Technical	
Differences and similarities between Protein 200 assay and SDS-PAGE	5988-3160EN
Use of lab-on-a-chip technology for protein sizing and quantitation	5988-8604EN
Series II Protein 80 and 230 Kits	5989-5175EN
Performance characteristics of the High Sensitivity Protein 250 assay	5989-8940EN
Protein sizing and quantitation with the Agilent Protein 80 and Protein 230 kits	5990-5721EN

2200 TapeStation system

DNA Applications

	Publication number
D1K ScreenTape simplifies Mycobacterium tuberculosis genotyping	5990-8760EN
D1K ScreenTape allows Mycoplasma PCR analysis for cell line screening	5990-8761EN
Fast and efficient HLA typing analysis using D1K ScreenTape	5990-9051EN
Analysis of high molecular weight genomic DNA using the Agilent 2200 TapeStation and Genomic DNA ScreenTape	5991-1797EN
Automated DNA sample analysis using a D1K ScreenTape assay and an Agilent 2200 TapeStation system	5991-2350EN

RNA Applications

Comparison of RIN and RIN ^e algorithms for the Agilent 2100 Bioanalyzer and the Agilent 2200 TapeStation systems	5990-9613EN
RNA quality control with the RIN ^e quality metric	5991-0023EN

Protein Applications

Suppression of antibody reduction artifacts with NEM on P200 ScreenTape	5990-8762EN
Optimization and control of protein purification	5990-9053EN
Reproducible integrity and purity testing of antibodies	5990-9052EN
Understanding the effects of proteins and buffers on staining, denaturation, and electrophoresis when analyzing proteins with Agilent P200 ScreenTape	5990-9603EN

For more information

Learn more

www.agilent.com/genomics/bioanalyzer

Buy online

www.agilent.com/genomics/contact

Find an Agilent customer center

www.agilent.com/genomics

© 2009-2014 Agilent Technologies Inc.

Printed in Germany, January 1, 2014

Publication Number 5991-3322EN



5991-3322EN



Agilent Technologies

