



AGILENT 8800 ICP-QQQ APPLICATION HANDBOOK

Primer

The Measure of Confidence



Agilent Technologies

Foreword

Quadrupole-based inductively coupled plasma–mass spectrometry (ICP-QMS) was introduced commercially in 1983. The technique immediately aroused considerable interest in the analytical community as it provided multi-element capabilities and better limits of detection compared to atomic absorption spectrometry (AAS) and ICP optical emission spectroscopy (ICP-OES). In the three decades since it was introduced, ICP-MS has been widely adopted as the technique of choice for trace element analysis in samples ranging from ultrapure reagents used in the semiconductor industry, through food commodities, pharmaceutical products, body fluids and tissues, to waters, soils and rocks measured in the environmental monitoring and mining industries. The ever increasing demand to quantitatively determine (ultra-trace) elements in various sectors and in a myriad of matrices ensures that ICP-MS is still gaining popularity.

As is the case with every other technique, ICP-MS has its drawbacks, and spectral interferences, especially those that jeopardize the accurate trace element determination at mass-to-charge ratios ≤ 80 amu, were rapidly recognized as the technique's proverbial Achilles' heel. Early ICP-MS users learnt how to overcome this obstacle in some instances, but some important elements simply remained inaccessible. With the introduction of sector-field ICP-MS with higher mass resolution, the majority of spectral interferences encountered can be resolved. Sector-field instrumentation is, however, substantially more expensive, and higher mass resolution comes at the cost of reduced ion transmission and thus, lower signal intensity. At that time, ICP-QMS had no "general approach" to offer as an alternative to overcome spectral overlaps. Cool plasma conditions formed a first step in the right direction, then with the introduction of multipole collision/reaction cells, the gap between sector-field and quadrupole-based ICP-MS instrumentation began to close. These collision/reaction cells already provide a versatile means for tackling spectral overlaps, as interferences can be removed either through differences in reaction chemistry (with a reactive cell gas), or by kinetic energy discrimination using a bias voltage (when an inert cell gas is used).

With the introduction of the Triple Quadrupole ICP-MS (ICP-QQQ) instrument, equipped with an octopole collision cell located between two quadrupole mass filters, this approach has been further perfected. Specifically, the capabilities of chemical resolution based on selective reactions between ions extracted from the ICP and the molecules of a reactive gas can now be fully exploited. The tandem mass spectrometer configuration provides two separate mass selection steps, which provides full control over the ion/molecule chemistry that occurs in the cell. The first quadrupole allows only the ions of a given mass-to-charge ratio into the gas-pressurized octopole cell, rejecting all the ions at all other mass-to-charge ratios. The second quadrupole then selects only the ion of interest emerging from the cell, and rejects the ions at all other mass-to-charge ratios. This approach not only provides a clear insight into the cell processes – via precursor ion and/or product ion scanning – but also provides an elegant approach for solving even the most challenging cases of spectral overlap. More reactive gases (e.g., NH_3 or O_2) can now be used in a controlled way, with the certainty the reaction pathways and product ions formed in the cell will not be affected by changes in the sample matrix or by other co-existing analyte ions. This consistency of reaction processes also permits more complex, heavier reaction product ions and clusters to be exploited for ultra-trace element determination. Even for such multi-element reaction product ions, the double mass selection preserves the original isotopic pattern, thereby facilitating spectral interpretation and isotope ratio determination, e.g., in the context of elemental assay via isotope dilution.

The full extent of the capabilities of Triple Quadrupole ICP-MS still needs to be discovered, but it seems clear that with the introduction of the 8800 ICP-QQQ, the analytical community has entered a new era of ICP-MS. This handbook provides an overview of the analytical approaches developed using ICP-QQQ thus far.

Frank Vanhaecke
Department of Analytical Chemistry
Ghent University, Belgium
www.analchem.ugent.be/A&MS

> Search entire document

Table of contents

Title	Elements	Page
Introduction to the Agilent 8800 Triple Quadrupole ICP-MS		5
Semiconductor		
Ultratrace measurement of calcium in ultrapure water	Ca	8
Determination of Ti, V and Cr in 9.8% sulfuric acid	Ti, V, Cr	10
Direct determination of V, Cr, Ge and As in high-purity 20% hydrochloric acid	V, Cr, Ge, As	12
Silicon wafer analysis by ICP-QQQ: Determination of phosphorus and titanium in a high silicon matrix	P, Ti	14
Analysis of sulfur, phosphorus, silicon and chlorine in N-methyl-2-pyrrolidone	S, P, Si, Cl	16
Analysis of silicon, phosphorus and sulfur in 20% methanol	S, P, Si	19
Environmental		
The accurate measurement of selenium in reference materials using online isotope dilution	Se	22
Accurate sulfur quantification in organic solvents using isotope dilution mass spectrometry	S	25
Removal of REE ⁺⁺ interference on arsenic and selenium	As, Se	28
Removal of molybdenum oxide interference on cadmium	Cd	31
Feasibility study of fluorine detection by ICP-QQQ	F	33
ICP-QQQ with oxygen reaction mode for accurate trace-level arsenic analysis in complex samples	As	35
Avoidance of spectral overlaps on reaction product ions with O ₂ cell gas: Comparison of quadrupole ICP-MS and ICP-QQQ	Ti	39
Removal of complex spectral interferences on noble metal isotopes	Ru, Rh Pd, Ag, Os, Ir, Pt and Au	42
Materials		
Arsenic measurement in cobalt matrix using MS/MS mode with oxygen mass-shift	As	48
Determination of sulfur, phosphorus and manganese in high purity iron	Mn, S, P	50
Direct measurement of trace rare earth elements in high purity REE oxides	REEs	53

Title	Elements	Page
The benefits of improved abundance sensitivity with MS/MS for trace elemental analysis of high purity metals	Cu	55
Ultra trace copper analysis in a semiconductor grade organometallic titanium complex	Cu	58
Removal of MH ⁺ interferences in refined REE material analysis	La, Ce, Pr	60
Geology		
Lead isotope analysis: Removal of ²⁰⁴ Hg isobaric interference on ²⁰⁴ Pb using ICP-QQQ MS/MS reaction cell	Pb	64
Nuclear		
Analysis of radioactive iodine-129 using MS/MS with O ₂ reaction mode	I	67
Feasibility study on the analysis of radioisotopes: Sr-90 and Cs-137	Sr, Cs	70
Clinical		
Manganese analysis in whole blood: Expanding the analytical capabilities of ICP-MS	Mn	74
Measurement of titanium in clinical samples: Possible application to monitoring patients with joint replacements	Ti	77
Life science		
Simultaneous quantitation of peptides and phosphopeptides by CapLC-ICP-QQQ	P, S	79
Pharma		
Measurement of selenium in clinical samples in the presence of Gd-based MRI contrasting agents	Se	81
Glossary		83
Index		90

- > Return to table of contents
- > Search entire document

Introduction to the Agilent 8800 Triple Quadrupole ICP-MS

The ground-breaking Agilent 8800 ICP-QQQ is the world's first Triple Quadrupole ICP-MS, representing a major advance in ICP-MS technology and redefining performance for trace element analysis. By taking the proven technology of Agilent's 7700 Series quadrupole ICP-MS and adding the unique power of MS/MS operation, the 8800 ICP-QQQ delivers unmatched analytical capability to handle even the most difficult samples and applications with ease.

Configuration of the Agilent 8800 ICP-QQQ

The IUPAC definition of a triple quadrupole mass spectrometer (term 538 from the 2013 Recommendations) is "*Tandem mass spectrometer comprising two transmission quadrupole mass spectrometers in series, with a (non-selecting) RF-only quadrupole (or other multipole) between them to act as a collision cell.*" The cell containing the ion guide — the Octopole Reaction System (ORS³) in the case of the 8800 ICP-QQQ — can be pressurized with a collision or reaction gas to process the ions selected by the first quadrupole (Q1). The resulting product ions that emerge from the cell are then filtered by the second quadrupole (Q2) before being passed to the detector.

ICP-MS/MS technology

The tandem MS configuration of the Agilent 8800 ICP-QQQ allows unprecedented control over the ions that enter the cell, and therefore ensures consistent reaction processes even when the sample composition varies, a first for ICP-MS. In MS/MS mode, Q1 operates as a unit mass filter, allowing only the ions at a single specific mass to charge ratio (m/z) to enter the collision/reaction cell. All matrix and non-target analyte ions are excluded, and only the analyte ions and any on-mass (same m/z) interferences enter the cell.

This means that the reactions in the ORS³ are simpler and more predictable compared to conventional quadrupole ICP-MS where all ions in the sample enter the cell and react, creating new, sample-dependent interferences and overlaps on the target analytes.

The Agilent 8800 ICP-QQQ unlocks the full potential of reaction cell chemistry to remove spectroscopic interferences, delivering greater accuracy and more consistent results even for previously problematic analytes and in complex and variable matrix samples.




Figure 1. Cutaway diagram of the Agilent 8800 ICP-QQQ

Agilent 8800 ICP-QQQ Application Handbook

- > [Return to table of contents](#)
- > [Search entire document](#)

This page is intentionally left blank.

- > [Return to table of contents](#)
- > [Search entire document](#)

Semiconductor

Ultratrace measurement of calcium in ultrapure water	8
Determination of Ti, V and Cr in 9.8% sulfuric acid	10
Direct determination of V, Cr, Ge and As in high-purity 20% hydrochloric acid	12
Silicon wafer analysis by ICP-QQQ: Determination of phosphorus and titanium in a high silicon matrix	14
Analysis of sulfur, phosphorus, silicon and chlorine in N-methyl-2-pyrrolidone	16
Analysis of silicon, phosphorus and sulfur in 20% methanol	19

Ultratrace measurement of calcium in ultrapure water

Albert Lee, Vincent Yang, Jones Hsu, Eva Wu and Ronan Shin, BASF Taiwan Ltd., Taipei, Taiwan
Katsuo Mizobuchi, Agilent Technologies, Japan

Keywords

semiconductor, process chemicals, ultra pure water, UPW, calcium, method of standard additions, hydrogen on-mass

Introduction

In the semiconductor industry, the control of metal impurities in the process chemicals used in the manufacture of semiconductor devices is critical to achieve the required product performance and yield. As device performance is continually increasing, the required impurity control becomes ever more stringent. For example, metal content of the ultra-pure water (UPW) used in the manufacturing process must be at the sub-ppt level. ICP-MS is the standard technique used for the trace metals analysis of semiconductor chemicals and devices. The most common instrument and measurement technique used in the semiconductor industry is single quadrupole ICP-MS (ICP-QMS) with cool plasma. The cool plasma technique [1], developed in the mid 1990's, enables the quantification of key contaminant elements at the single ppt level. Collision and reaction cell ICP-QMS, developed from 2000 onwards, enabled the direct analysis of more complex semiconductor matrices, but did not improve on the DLs or BECs of cool plasma for low-matrix samples. To achieve measurement at the sub-ppt level, reduction of the BEC is required. As outlined in this paper, the Agilent 8800 ICP-QQQ provides new reaction cell technology that enables a significant reduction in the BEC that can be achieved for Ca, to 100 ppq.

Table 1. Cool plasma operating conditions

Parameter	Unit	Tuning value
RF	W	600
Sampling depth	mm	18
Carrier gas flow	L/min	0.7
Make up gas flow	L/min	1.0
Spray chamber temp.	°C	2

Experimental

Instrumentation: Agilent 8800 #200.

Plasma conditions: For the ultra-trace measurement of Ca, cool plasma operating conditions were used (Table 1). The sample was self-aspirated at a carrier gas flow rate of 0.7 L/min.

Reagents and sample preparation: A Ca standard was prepared in UPW acidified with 0.1% high purity HNO₃. This was used to make 50 ppt and 100 ppt additions to a UPW blank acidified with 0.1% high purity HNO₃.

Results and discussion

Ultra-low BEC for Ca using MS/MS mode

Figure 1 shows the BECs obtained for Ca, measured at its major isotope of ⁴⁰Ca, using the method of standard additions (MSA) under three different operating conditions on the 8800 ICP-QQQ: Single Quad mode with no cell gas, MS/MS mode with no cell gas, and finally MS/MS mode with a H₂ cell gas flow of 1 mL/min. The Single Quad mode uses operating conditions with Q1 acting as an ion guide, to emulate the Agilent 7700 ICP-QMS. The obtained BEC of 6.8 ppt is similar to that routinely achieved with the Agilent 7700 operated in cool plasma mode.

Using MS/MS mode (without cell gas) improved the Ca BEC to 1.4 ppt. MS/MS mode with H₂ at 1 mL/min in the cell further improved the BEC down to 0.041 ppt (41 ppq). The obtained MSA plot is shown in Figure 2. The Agilent 8800 ICP-QQQ in MS/MS mode with H₂ cell gas achieved a BEC for Ca in UPW two orders of magnitude lower than the BEC obtained using conventional ICP-QMS.

Figure 3 shows the spectrum obtained for UPW using cool plasma conditions in Single Quad mode with no cell gas. As can be seen, Ar⁺ (*m/z* 40) is suppressed under the lower temperature plasma conditions, but two intense background peaks are observed at *m/z* = 19 and 30. These are (H₂O)H⁺ and NO⁺ respectively. In Single Quad mode, all ions formed in the plasma,

including these two intense ions, pass through to the cell. Even with no gas added to the cell, a reaction occurs in the cell which causes a new interfering ion at $m/z = 40$. The likely reaction occurring in the cell is: $\text{NO}^+ + \text{Ar} \rightarrow \text{Ar}^+ + \text{NO}$ (charge transfer reaction), which increases the BEC for Ca by several ppt. Although the ionization potential (IP) of NO (IP = 9.26 eV) is lower than that of Ar (IP = 15.7 eV), a metastable ion, NO^+ , exists close to the ionization potential of Ar [2]. So it is reasonable to assume that the charge transfer reaction shown occurs in the cell.

With MS/MS mode on the 8800 ICP-QQQ, Q1 rejects all non-target ions such as NO^+ and $(\text{H}_2\text{O})\text{H}^+$, preventing unwanted reactions from occurring in the cell, which




Figure 1. BECs for Ca obtained using Single Quad mode with no cell gas [6.8 ppt], MS/MS mode with no cell gas [1.4 ppt], and MS/MS mode with an H_2 cell gas flow of 1 mL/min [0.041 ppt].

lowers the Ca BEC. The addition of H_2 in the cell also removes any residual $^{40}\text{Ar}^+$ that is formed even under cool plasma conditions.

References

1. K. Sakata and K. Kawabata, Reduction of fundamental polyatomic ions in inductively coupled plasma mass spectrometry, *Spectrochimica Acta, Part B*, 1994, 49, 1027.
2. R. Marx, Y.M. Yang, G. Mauclaire, M. Heninger, and S. Fenistein, Radioactive lifetimes and reactivity of metastable $\text{NO}^+(\text{a}^3\Sigma^+, \nu)$ and $\text{O}_2^+(\text{a}^4 \text{II}, \nu)$, *J.Chem. Phys.*, Vol. 95, No. 4, 2259-2264, 1991.




Figure 2. MSA calibration plot for Ca using MS/MS mode with H_2 flow of 1 mL/min.




Figure 3. Spectrum of UPW acquired using cool plasma conditions in Single Quad mode with no gas mode.

Determination of Ti, V and Cr in 9.8% sulfuric acid

Junichi Takahashi
Agilent Technologies, Japan

Keywords

semiconductor, process chemicals, sulfuric acid, H_2SO_4 , titanium, vanadium, chromium, ammonia mass-shift, oxygen mass-shift

Introduction

High purity H_2SO_4 is frequently used in the manufacturing of semiconductor devices, in processes such as the removal of organic substances from the surface of silicon wafers. The required metallic impurity level is lower than 100 ppt in the concentrated (usually 98%) acid. ICP-MS is the technique of choice for the measurement of trace metal impurities in semiconductor process chemicals. There are, however, some limitations for the measurement of elements such as Ti, V and Cr in H_2SO_4 . Because of its high viscosity of 27 cP, it is not possible to introduce H_2SO_4 directly into the ICP without dilution. A 10 times dilution in UPW is normally applied, thus the BEC of the calibration curve must be lower than 10 ppt in the 9.8% H_2SO_4 solution measured. In addition, spectral interferences from SO^+ , S_2^+ and ArS^+ originating from H_2SO_4 make it difficult to determine elements such as Ti and Cr at low concentration even by quadrupole ICP-MS (ICP-QMS) equipped with collision/reaction cell (CRC). As outlined in this report, the Agilent 8800 ICP-QQQ with MS/MS mode allows the successful determination of the most problematic elements including Ti, V and Cr in H_2SO_4 .

Table 1. ICP-QQQ operating conditions

		O_2 MS/MS ¹⁾	NH_3 MS/MS ²⁾
RF power	W		1600
Sampling depth	mm		8
CRGS flow rate	L/min		0.8
MUGS flow rate	L/min		0.41
Octopole bias V	V		-20
KED	V		-20
He	mL/min	3	1
O_2	mL/min	0.4	0
NH_3	mL/min	0	3

1) 100% O_2 (purity 99.995%)

2) 10% NH_3 balanced with 90% He (purity 99.995%)

Experimental

Instrumentation: Agilent 8800 #200. Operating parameters are given in Table 1.

Reagents and sample preparation: Highly purified H_2SO_4 , TAMAPURE-AA-100 (98% H_2SO_4) was purchased from Tama Chemicals Co., Ltd. (Kanagawa, Japan). 5 g of H_2SO_4 was diluted by a factor of 10 in a chilled PFA bottle.

Results and discussion

Of the potential polyatomic interferences formed from the H_2SO_4 matrix, the SO^+ ion is very stable and difficult to eliminate because its dissociation energy is as high as 5.44 eV. In addition, its ionization potential is 10.3 eV, which is almost the same as that of S, 10.36 eV. The spectral interferences caused by SO^+ and SOH^+ overlap with ^{48}Ti ($^{32}S^{16}O$), ^{51}V ($^{33}S^{18}O$, $^{34}S^{16}OH$ and $^{32}S^{18}OH$) and ^{52}Cr ($^{34}S^{18}O$). Quadrupole ICP-MS operating in He collision mode provides BECs of 60 ppt for ^{47}Ti (the BEC for the preferred isotope ^{48}Ti is much higher), 3 ppt for V and 8 ppt for Cr in 9.8% H_2SO_4 . The BEC of Ti, in particular, is not acceptable for producers and users of semiconductor grade H_2SO_4 .

Appropriate reaction gases to remove SO^+ successfully in ICP-QMS are difficult to find. NH_3 can reduce SO^+ by two orders of magnitude but the background signal remains too high for this application. Additionally, cluster ions of NH_3 such as N_mH_n produced by the reaction between Ar^+ and the NH_3 cell gas lead to new reaction product ion interferences that increase the background at m/z 51, for example.

The 8800 ICP-QQQ operating in MS/MS mass-shift mode with NH₃ or O₂ reaction gas provides reliable and consistent measurement of Ti as ⁴⁸Ti¹⁴NH(¹⁴NH₃)₃ (Figure 1) and Cr as ⁵²Cr¹⁶O in H₂SO₄. Furthermore, in MS/MS mode, the Ar⁺ ion is removed by Q1, preventing it from reacting with NH₃ to form new product ion interferences in the cell. This reduces the background at *m/z* 51 improving the BEC for V, as shown in Figure 2. The final BECs obtained by ICP-QQQ in 9.8% high purity H₂SO₄ are summarized in Table 2.

Table 2. BECs of Ti, V and Cr in 10x diluted 98% H₂SO₄, measured by ICP-QQQ

Element	Ti	V	Cr
Mode (cell gas)	MS/MS (NH ₃)	MS/MS (NH ₃)	MS/MS (O ₂)
Measured ion	⁴⁸ Ti ¹⁴ NH(¹⁴ NH ₃) ₃ ⁺	⁵¹ V ⁺	⁵² Cr ¹⁶ O ⁺
Mass pair	Q1 = 48, Q2 = 114	Q1 = Q2 = 51	Q1 = 52, Q2 = 68
BEC - ppt	0.72	0.07	3.70




Figure 1. Calibration curve of Ti in 9.8% H₂SO₄

Conclusions

ICP-QQQ operating in MS/MS mode provides a reliable means for manufacturers of high purity H₂SO₄ to guarantee all metallic impurity concentrations at less than 100 ppt in the concentrated acid.

More information

Determination of challenging elements in ultrapure semiconductor grade sulfuric acid by Triple Quadrupole ICP-MS, Agilent application note, 5991-2819EN.




Figure 2. Calibration curve of V in 9.8% H₂SO₄

Direct determination of V, Cr, Ge and As in high-purity 20% hydrochloric acid

Junichi Takahashi
Agilent Technologies, Japan

Keywords

semiconductor, RCA Standard Clean, silicon wafer, hydrochloric acid, vanadium, chromium, germanium, arsenic, ammonia on-mass, ammonia mass-shift, oxygen mass-shift

Introduction

Since the 1970s, the RCA Standard Clean (SC) method has been used extensively in many countries for cleaning silicon wafer surfaces. SC-2 refers to a mixture of HCl and H₂O₂ that is used to remove ionic and metallic contaminants from the surface of silicon wafers. Because cleaning solutions are in direct contact with semiconductor devices, ultra high purity is required for these solutions. The SEMI standard Tier-D protocol for HCl defines the contaminant level to be <10 ppt for each element. Some elements have been very difficult to determine at ppt level by quadrupole ICP-MS (ICP-QMS) due to significant spectral interferences arising from the Cl matrix, even when analyzed by ICP-MS equipped with a collision/reaction cell (CRC). Consequently, some methods for the analysis of high purity HCl by ICP-MS have recommended sample pre-treatment steps to remove the chloride matrix, which can lead to analyte loss and sample contamination. In this study, ICP-QQQ was used to analyze undiluted HCl directly. Using MS/MS mode with mass-shift to remove polyatomic

ions, the most problematic elements, such as V, Cr, Ge and As could be determined in HCl at single-figure ppt detection limits.

Experimental

Instrumentation: Agilent 8800 #200. Operating parameters are given in Table 1.

Reagents: 20% TAMAPURE-AA-100 HCl (metallic impurities are guaranteed to be below 100 ppt) was purchased from Tama Chemicals Co., Ltd. (Kanagawa, Japan). The undiluted HCl was introduced directly into the ICP-QQQ.

Results and discussion

Determination of BECs of V, Cr, Ge and As in high purity HCl

ICP-QMS with a CRC using He collision mode can successfully eliminate some polyatomic ions such as ArCl [1], and the use of NH₃ as a reaction gas also works to remove the ClO⁺ ion for the determination of V. However, ICP-QMS has some serious limitations when highly reactive cell gases (such as NH₃) are used in the CRC. Principal among these limitations is the fact that all ions enter the CRC, so predicted reaction pathways can be disrupted and new reaction product ion overlaps can be formed if the analyte levels in the sample change. ICP-QQQ with MS/MS removes this limitation, as the first quadrupole mass filter (Q1) allows precise selection of the ions that are allowed to enter the cell. This ensures that reaction processes and product ions are strictly controlled, dramatically improving detectability of the analyte ions shown in Table 2.

Table 1. ICP-QQQ operating conditions

		O ₂ MS/MS ¹⁾	NH ₃ MS/MS ²⁾
RF power	W	1600	
Sampling depth	mm	8	
CRGS flow rate	L/min	0.8	
MUGS flow rate	L/min	0.41	
Octopole bias V	V	-20	
KED	V	-20	
He	mL/min	3	1
O ₂	mL/min	0.4	0
NH ₃	mL/min	0	3

1) 100% O₂ (purity 99.995%)

2) 10% NH₃ balanced with 90% He (purity 99.995%)


Table 2. Spectral interferences arising from the Cl matrix on some key elements

Polyatomic interference	m/z	Analyte ion
ClO ⁺	51, 53	⁵¹ V ⁺
ClOH ⁺	52, 54	⁵² Cr ⁺ , (⁵⁴ Fe ⁺) [*]
ClCl ⁺	70, 72, 74	⁷⁰ Ge ⁺ , ⁷² Ge ⁺ , ⁷⁴ Ge ⁺
ArCl ⁺	75, 77	⁷⁵ As ⁺ , (⁷⁷ Se ⁺) [*]

^{*}Alternative isotopes can be chosen to avoid spectral interferences on Fe and Se.

The MS/MS acquisition mode using O₂ or NH₃ as the reaction gas enables the determination of trace ⁵¹V (measured directly as V⁺ using NH₃ cell gas), Cr as ⁵²Cr¹⁶O⁺ (using O₂), Ge as ⁷⁴Ge¹⁴NH₂⁺ (using NH₃) and As as ⁷⁵As¹⁶O⁺ (using O₂). In the case of As, the ⁹¹Zr⁺ ion is removed by Q1 (which is set to the As⁺ precursor ion mass of *m/z* 75), so the potential overlap from Zr on the AsO⁺ product ion at *m/z* 91 is also removed. The complete cut-off of cluster ions by Q1 also eliminates the possibility that ¹⁴NH₂³⁵Cl is created in the cell, so the potential new product ion interference on ⁵¹V is avoided. Representative calibration curves for V and Ge are shown in Figure 1. BECs and DLs determined by the ICP-QQQ for V, Cr, Ge and As are given in Table 3.

51 → 51 V [NH₃,MM]



74 → 90 Ge [NH₃,MM]




Figure 1. Calibration curves of V (NH₃ on-mass mode) and Ge (NH₃ mass-shift mode) in 20% HCl

Investigation of arsenic contamination

As the BEC for arsenic in high purity HCl was relatively high (Table 3), the signal count at *m/z* 91 was investigated further. The signals of the mass-pairs 75/75, 77/77, 75/91 and 77/91 were measured by ICP-QQQ with MS/MS, the mass pair number represents the set mass of Q1 followed by the set mass of Q2, so an MS/MS mode acquisition of mass pair 75/91

Table 3. BECs and DLs for V, Cr, Ge and As in 20% HCl

Element	V	Cr	Ge	As
Mode (cell gas)	MS/MS (NH ₃)	MS/MS (O ₂)	MS/MS (NH ₃)	MS/MS (O ₂)
Measured ion	⁵¹ V ⁺	⁵² Cr ¹⁶ O ⁺	⁷⁴ Ge ¹⁴ NH ₂ ⁺	⁵² Cr ¹⁶ O ⁺
Mass pair	Q1 = Q2 = 51	Q1 = 52, Q2 = 68	Q1 = 74, Q2 = 90	Q1 = 75, Q2 = 91
BEC - ppt	0.3	8.0	1.5	19.7
DL - ppt	0.5	1.1	4.4	3.4

represents a mass-shift mode with Q1 = 75 and Q2 = 91, for example. The four mass pairs were measured in HCl blanks from three different lots, and the results are shown in Table 4. The following observations were made:

1. The ratio of the signal of 75/75 to 77/77 is around four, which is close to the ratio of the abundance of ³⁵Cl to ³⁷Cl, i.e. 3.13.
2. The ratio of the signal of 75/91 to 77/93 is 200–1000, which is far in excess of the ratio of ³⁵Cl to ³⁷Cl.
3. While the signals of 75/75 and 77/77 are similar for the three HCl blanks, those of 75/91 and 77/93 vary.

Table 4. Comparison of background counts (cps) in 3 different lots of 20% HCl*

Mass pair	75->75	77->77	75->91	77->93
Sample A	509.3	133.5	584.4	2.5
Sample B	508.4	126.0	1172.6	1.9
Sample C	612.7	130.0	3175.6	2.6

*All the samples were obtained from the new bottles of high purity HCl

Finding #1 suggests that the remaining signal on 75/75 and 77/77 was mostly from ArCl⁺. This is a reasonable assumption since ArCl⁺ doesn't react with O₂ very efficiently so most ArCl⁺ remains at the original masses of 75 and 77. Finding #2 suggests that the signal of 75/91 is not due to ArCl⁺. Assuming that all counts of 77/93 arise from ⁴⁰Ar³⁷Cl, the contribution of ⁴⁰Ar³⁵Cl to the signal of 75/91 in the HCl blank is estimated to be just 7-8 cps, which is two orders of magnitude lower than the signal that is actually observed. Observation #3, together with #1 and #2, suggests the high count obtained for 75/91 in HCl is due to As impurity in the acid.

Reference

1. Direct analysis of trace metallic impurities in high purity hydrochloric acid by Agilent 7700s ICP-MS, Agilent application note, 5990-7354EN.

Silicon wafer analysis by ICP-QQQ: Determination of phosphorus and titanium in a high silicon matrix

Junichi Takahashi
Agilent Technologies, Japan

Keywords

semiconductor, silicon wafer, phosphorus, titanium, Vapor Phase Decomposition, VPD, oxygen mass-shift

Introduction

The semiconductor industry first used ICP-MS for trace element analysis in the early 1980s. Nowadays the technique is widely used for control of trace impurities in materials and chemicals, particularly by silicon device manufacturers. The major challenge for quadrupole ICP-MS (ICP-QMS) is the presence of spectroscopic interferences on key contaminant elements, although performance has been gradually improved through developments such as cool plasma and collision/reaction cells (CRC), and improved performance has also been provided by high resolution ICP-MS. Consequently metallic impurity control of silicon wafers can be successfully monitored by ICP-MS in the case of low silicon samples such as Vapor Phase Decomposition (VPD) of native silicon wafers. However, difficulties of Si-based spectral interferences, particularly on P and Ti, still affect the analysis of samples that contain high concentrations of Si, such as VPD samples of thermally oxidized wafers and samples relating to bulk silicon wafers. These interferences cannot be reduced adequately by ICP-QMS and have required HR-ICP-MS. In this paper, we evaluate triple quadrupole ICP-MS with MS/MS technology for the determination of ultratrace P and Ti in a high Si matrix.

Table 1. Robust tuning conditions

		O ₂ MS/MS	H ₂ MS/MS
RF power	W	1600	
Sampling depth	mm	8	
CRGS flow rate	L/min	0.6	
MUGS flow rate	L/min	0.6	
He	mL/min	3	0
O ₂	mL/min	0.4	0
H ₂	mL/min	0	10

Experimental

Instrumentation: Agilent 8800 #200 with an inert sample introduction kit including a low flow nebulizer (PFA-20) and a Pt/Ni skimmer cone. The actual sample uptake rate was 36 μ L/min. The sample was self-aspirated from an Agilent I-AS autosampler.

Plasma conditions: Robust tuning conditions were applied as summarized in Table 1.

Ion lens tune: Extract 1 = 0 V was used and other lens voltages were optimized using Auto tune.

Sample preparation: Silicon wafer samples were dissolved in TAMAPURE HF/HNO₃ and the final Si concentration was adjusted to 2000 ppm.

Results and discussion

Phosphorus is monoisotopic at m/z 31, and suffers an interference from ³⁰SiH. While P⁺ can be detected as PO⁺ under cool plasma conditions, it is difficult to maintain cool plasma when the matrix concentration is high. Si sample solutions always contain HF, so Si will form SiF (IP: 7.54 eV) that also interferes with Ti. Table 2 shows the Si-based spectral interferences on P and Ti. Using the 8800 ICP-QQQ operating in MS/MS mode with O₂ mass-shift, P and Ti can be determined as their oxide ions, avoiding the Si-based interferences.

Table 2. Spectral interferences of Si on P and Ti

Polyatomic interference	m/z	Analyte ion
³⁰ SiH ⁺	31	³¹ P ⁺
³⁰ Si ¹⁶ O ⁺	46	⁴⁶ Ti ⁺
²⁸ Si ¹⁹ F ⁺ , ³⁰ Si ¹⁶ OH ⁺	47	⁴⁷ Ti ⁺ , ³¹ P ¹⁶ O ⁺
²⁹ Si ¹⁹ F ⁺ , ³⁰ Si ¹⁸ O ⁺	48	⁴⁸ Ti ⁺
³⁰ Si ¹⁹ F ⁺	49	⁴⁹ Ti ⁺

For Ti analysis, Q1 is set to m/z 48, and so will transmit $^{48}\text{Ti}^+$ and any other interfering ions at mass 48, such as $^{29}\text{Si}^{19}\text{F}^+$ and $^{30}\text{Si}^{18}\text{O}^+$. But only ^{48}Ti reacts with oxygen in the CRC, producing the product ion $^{48}\text{Ti}^{16}\text{O}^+$, which is transmitted by setting Q2 to m/z 64. NH_3 can be used as an alternative reaction gas, as it produces $^{48}\text{Ti}^{14}\text{NH}^+$ that can be detected at m/z 63.

$^{31}\text{P}^+$ reacts readily with O_2 to form $^{31}\text{P}^{16}\text{O}^+$. The selection of ions at m/z 31 by Q1 eliminates the spectral interference of $^{28}\text{Si}^{19}\text{F}$. However, ^{30}SiH passes through Q1 and reacts with O_2 to create $^{30}\text{Si}^{16}\text{OH}$. In order to determine P in a high Si matrix, H_2 mass-shift is a preferred option, despite the relatively low efficiency of production of PH_3^+ or PH_4^+ ions. The MSA calibration curves for P and Ti in a matrix of 2000 ppm Si are shown in Figure 1. The calculated BECs are summarized in Table 3. A long term stability test was carried out by analyzing a spiked sample repeatedly over five hours (Figure 2).

Table 3. BECs of P and Ti in 2000 ppm Si

Element	P	Ti
Mode (cell gas)	MS/MS (H_2)	MS/MS (O_2)
Measured ion	$^{31}\text{PH}_4^+$	$^{48}\text{Ti}^{16}\text{O}^+$
Mass pair	Q1 = 31, Q2 = 35	Q1 = 48, Q2 = 64
BEC - ppt	227	13

Conclusions

The MS/MS mass-shift mode of the ICP-QQQ is effective for the determination of P, Ti and other trace elements in high purity silicon matrices, providing effective removal of the potential Si-based polyatomic interferences.

More information

Improvement of ICP-MS detectability of phosphorus and titanium in high purity silicon samples using the Agilent 8800 Triple Quadrupole ICP-MS, Agilent application note, 5991-2466EN.




Figure 1. MSA curves of P and Ti in 2000 ppm Si




Figure 2. Five-hours test of P and Ti spiked in 2000 ppm Si

Analysis of sulfur, phosphorus, silicon and chlorine in N-methyl-2-pyrrolidone

Naoki Sugiyama
Agilent Technologies, Japan

Keywords

N-methyl-2-pyrrolidone, NMP, semiconductor, process chemicals, sulfur, phosphorus, silicon, chlorine, method of standard additions, oxygen mass-shift

Introduction

N-Methyl-2-Pyrrolidone (NMP), chemical formula: C_5H_9NO , is a stable, water-soluble organic solvent that is widely used in the pharmaceutical, petrochemical, polymer science and especially semiconductor industries. Electronic grade NMP is used by semiconductor manufacturers as a wafer cleaner and photo resist stripper and as such the solvent comes into direct contact with wafer surfaces. This requires NMP with the lowest possible trace metal (and non-metal) contaminant levels. ICP-MS is the technique of choice for the measurement of trace metal impurities in semiconductor process chemicals. It is a challenge, however for ICP-MS to measure non-metallic impurities such as sulfur, phosphorus, silicon, and chlorine in NMP. The low ionization efficiency of these elements greatly reduces analyte signal, while the elevated background signal (measured as background equivalent concentration, BEC) due to N-, O-, and C-based polyatomic ions formed from the NMP matrix makes low-level analysis even more difficult (Table 1).

Table 1. ICP-QMS BECs obtained in no gas mode for selected analytes in NMP

Element	m/z	Ionization potential (eV)	Ionization ratio (%)	BEC without cell (ppm)	Interference
Si	28	8.152	87.9	>100	$^{14}N_2^+$, $^{12}C^{16}O^+$
P	31	10.487	28.8	0.39	$^{14}N^{16}OH^+$, COH_3^+
S	32	10.360	11.5	9.5	$^{16}O_2^+$, NOH_2^+
Cl	35	12.967	0.46	0.26	$^{16}O^{18}OH^+$

Experimental

Instrumentation: Agilent 8800 #200 with narrow injector (id =1.5 mm) torch (G3280-80080) typically used for the analysis of organic solvents. A C-flow 200 PFA nebulizer (G3285-80000) was used in self-aspiration mode. An option gas flow of 20% O_2 in Ar was added to the carrier gas to prevent carbon build up on the interface cones.

Plasma conditions: NMP analysis requires hotter plasma conditions than normal. This was achieved by reducing Make-up Gas (MUGS) by 0.2 L/min. Plasma tuning conditions are summarized in Table 2.

Table 2. Plasma conditions for NMP analysis

Parameter	Unit	Tuning value
RF	W	1550
Sampling depth	mm	8.0
Carrier gas flow	L/min	0.50
Make up gas flow	L/min	0.10
Option gas flow	L/min	0.12 (12% of full scale)
Spray chamber temperature	°C	0

CRC conditions: Table 3 summarizes the cell tuning parameters (gas flow rate and voltages) used.

Reagents and sample preparation: Electronic industry grade NMP was distilled at 120 °C and acidified by adding high purity HNO_3 to a concentration of 1% w/w.

Table 3. CRC operating conditions

Parameter	Unit	O_2 reaction cell		H_2 reaction cell	
		On-mass	Mass-shift	On-mass	Mass-shift
Method	-	On-mass	Mass-shift	On-mass	Mass-shift
Cell gas	-	O_2		H_2	
Gas flow rate	mL/min	0.30		4.0	
Octopole bias	V	-14		-10	
KED	V	-5	-5	0	-5

Results and discussion

NMP was analyzed directly using the method of standard additions (MSA). Three replicate measurements (ten replicates for the blank) were acquired for S, P, Si and Cl using an integration time of 1 s per isotope.

P and S measurement in NMP

The mass-shift method using O₂ worked well for P and S measurement in NMP. The reactions of P and S with O₂ are exothermic, indicated by the negative value for ΔH, as shown below; therefore P⁺ and S⁺ are efficiently converted to their oxide ions, PO⁺ and SO⁺. P and S can be measured as the product ions, avoiding the original spectroscopic interferences on their elemental masses, *m/z* 31 and *m/z* 32.



In MS/MS mode, Q1 rejects ³⁶ArC⁺ before it can enter the cell, preventing it from overlapping SO⁺. This allows ICP-QQQ to control the reaction chemistry pathways and reaction product ions, ensuring that the analyte product ion is measured free from overlap, regardless of the levels of other co-existing analyte (or matrix) elements.





Figure 1. Calibration curve using MS/MS with O₂ mass-shift for P (top) and S (bottom) in NMP

MS/MS mode with the O₂ mass-shift method achieved BECs of 0.55 ppb and 5.5 ppb for P and S respectively in NMP. The low BECs and linear calibration plots achieved in MS/MS mode also prove that the matrix-based interferences do not react with O₂, allowing the analytes to be separated from the interferences.

Si measurement in NMP

H₂ cell gas was applied to the measurement of Si in NMP. The reaction kinetics for Si and its major interferences with H₂ cell gas are shown below. The reaction rate data suggests that Si does not react with H₂ cell gas (endothermic reaction indicated by the positive value for ΔH), and so could be measured in NMP using the direct, on-mass method. While the reaction of Si⁺ with H₂ is endothermic, the reactions of the major interfering ions on Si at mass 28 (N₂⁺ and CO⁺) are exothermic, and these interferences are therefore neutralized or reacted away.



The results obtained are shown in Figure 2 (top). The H₂ on-mass method achieved a BEC of 15.8 ppb for Si in NMP.





Figure 2. Calibration plots for Si in NMP. Top: H₂ on-mass method (Q1=Q2=28). Bottom: O₂ mass-shift method (Q1=28, Q2=60)

Oxygen cell gas was also tested of the measurement of Si in NMP. As shown below, the reaction of Si⁺ with O₂ to form SiO⁺ is endothermic. However, collisional processes in the cell provide additional energy which promotes the reaction, enabling the O₂ mass-shift method to be applied.



Unfortunately a major interference on Si at *m/z* 28 (CO⁺) also reacts with O₂, so the BEC achieved using the O₂ mass-shift method to measure Si as SiO⁺ (Q1 = 28, Q2 = 44) was not satisfactory. Fortunately, another Si reaction product ion (SiO₂⁺) also forms and this can be measured at *m/z* 60 (Q1 = 28, Q2 = 60) giving a BEC of 11.9 ppb for Si in NMP (Figure 2, bottom).

Cl in NMP

Cl⁺ reacts exothermically with H₂ to form HCl⁺ as shown below. HCl⁺ continues to react via a chain reaction to form H₂Cl⁺.




Figure 3 (left) shows calibration plots obtained for Cl in NMP using the H₂ mass-shift method. The plot obtained using the O₂ mass-shift method (Figure 3, right) is also shown for comparison. A slightly better BEC of 34.2 ppb was achieved with much higher sensitivity for Cl in NMP using the H₂ mass-shift method.

More information

Trace level analysis of sulfur, phosphorus, silicon and chlorine in NMP using the Agilent 8800 Triple Quadrupole ICP-MS, 2013, Agilent application note, 5991-2303EN.




Figure 3. Calibration plots for Cl in NMP. Left: H₂ mass-shift method (Q1 = 35, Q2 = 37). Right: O₂ mass-shift method (Q1=35, Q2=51)

Analysis of silicon, phosphorus and sulfur in 20% methanol

Emmett Soffey, Bert Woods and Steve Wilbur
Agilent Technologies, Japan

Keywords

organic solvents, methanol, silicon, phosphorus, sulfur, hydrogen on-mass, oxygen mass-shift

Introduction

Analysis of organic solvents for trace metals presents a number of challenges to ICP-MS, many of which have been overcome to varying degrees on Agilent's 7700 Series quadrupole ICP-MS systems. However, even with these advances, several elements remain challenging in organic solvents, particularly silicon, phosphorus and sulfur. All three elements are subject to intense interferences from polyatomic ions based on carbon, nitrogen and oxygen, which are difficult to completely remove using conventional quadrupole ICP-MS (ICP-QMS). Examples include CO^+ , COH^+ , N_2^+ and NO^+ on silicon 28, 29 and 30; COH^+ , NOH^+ , N_2H^+ , NO^+ and CO^+ on phosphorus 31 and O_2^+ , NO^+ , NOH^+ and NOH_2^+ on sulfur 32 and 34. Additionally, phosphorus and sulfur have high first ionization potentials (IP) of 10.5 eV and 10.4 eV respectively, resulting in relatively poor sensitivity compared to more typical elements whose IPs are in the range of ~6 – 8 eV.

Experimental

Instrumentation: Agilent 8800 #200.

Plasma conditions and ion lens tune: RF power = 1550 W, Sampling depth = 8.0 mm and CRGS flow rate = 1.05 L/min were used with soft extraction tune, Extract 1 = 0 V and Extract 2 = -190 V.

Ultra pure methanol was spiked with silicon (Si), phosphorus (P) and sulfur (S) at 1, 5, 10 and 50 ppb and measured using the ICP-QQQ in several operational modes in order to evaluate the optimum conditions for the simultaneous analysis of all three analytes. Hydrogen and oxygen reaction gases were evaluated, with H_2 cell gas used in both Single Quad (SQ) and MS/MS modes. In addition, helium collision gas was investigated in both SQ and MS/MS mode to determine the effects of using MS/MS with a non-reactive cell gas.

The CRC conditions are outlined in Table 1, which includes the five analysis modes evaluated. Two Single Quad modes were tested, using both He and H_2 in the cell, to simulate the capability of a single quadrupole ICP-MS. In addition, three conditions using MS/MS mode were tested using H_2 , He and O_2 as cell gases.

Table 1. 8800 ICP-QQQ acquisition conditions tested, including five operational modes

Parameter	Unit	He SQ	H_2 SQ	H_2 MS/MS	He MS/MS	O_2 MS/MS
Acquisition mode		SQ	SQ	MS/MS	MS/MS	MS/MS
Cell gas		He	H_2	H_2	He	O_2
Cell gas flow rate	mL/min	5.0	7.0	7.0	7.0	0.40
KED	V	5	0	0	5	-7

Results and discussion

The BECs and DLs results are summarized in Tables 2–4, for silicon, phosphorus and sulfur respectively, for all 5 analysis modes tested. SQ and optimum MS/MS results are in bold type for comparison. A few mass-pairs were measured in each mode as shown. For example, Table 2 shows silicon monitored in MS/MS mode with O_2 cell gas, using a mass-pair of Q1 = 28 and Q2 = 44. With Q1 set to m/z 28, only silicon 28 and any on-mass interferences are allowed to enter the ORS cell. The silicon 28 in the cell reacts with the oxygen cell gas to form SiO^+ , and Q2 is set to measure at Q1 + 16 (m/z = 44), ensuring that only the $\text{M} + {}^{16}\text{O}$ reaction transition is measured.

Table 2. DLs and BECs for silicon. Silicon was not measurable at the spiked concentrations in helium mode

Mode	Mass or mass pair	BEC (ppb)	DL (ppb)
H_2 SQ	Q2=28	25.46	0.12
H_2 MS/MS	Q1=28, Q2=28	2.17	0.03
O_2 MS/MS	Q1=28, Q2=44	85.54	28.21
O_2 MS/MS	Q1=29, Q2=45	N/A	N/A
O_2 MS/MS	Q1=30, Q2=46	99.09	21.26

Table 3. DLs and BECs for phosphorus. Phosphorus was not measurable at the spiked concentrations in H₂ Single Quad mode

Mode	Mass or mass pair	BEC (ppb)	DL (ppb)
He SQ	Q2=31	3.81	0.63
He MS/MS	Q1=31, Q2=31	2.99	0.72
H ₂ MS/MS	Q1=31, Q2=33	0.56	0.07
H ₂ MS/MS	Q1=31, Q2=34	0.58	0.67
O ₂ MS/MS	Q1=31, Q2=47	0.40	0.05
O ₂ MS/MS	Q1=31, Q2=63	0.41	0.02

Table 4. DLs and BECs for sulfur using MS/MS mode with O₂ cell gas. Sulfur was not measurable at the spiked concentrations in helium or hydrogen mode

Mode (cell gas)	Mass or mass pair	BEC (ppb)	DL (ppb)
MS/MS (O ₂)	Q1 = 34, Q2 = 34	51.17	4.37
MS/MS (O ₂)	Q1 = 32, Q2 = 48	3.13	0.10
MS/MS (O ₂)	Q1 = 34, Q2 = 50	3.11	0.20

Sample calibration plots are displayed in Figure 1. They are displayed in pairs showing the results obtained using SQ mode with a typical cell gas (upper calibration), compared to MS/MS mode using the optimum conditions (lower calibration).

Conclusions

It can be seen that in all cases the use of MS/MS mode significantly improves both the BEC and instrument detection limit when compared to Single Quad mode. The most notable improvement was for sulfur which cannot be measured at the spiked concentrations (1, 5, 10, 50 ppb) in SQ He mode due to the intense polyatomic background resulting from the methanol matrix. By contrast, on the 8800 ICP-QQQ using MS/MS mode with O₂ mass-shift, S can be measured with a DL of 0.1 ppb.

More information

Analysis of silicon, phosphorus and sulfur in 20% methanol using the Agilent 8800 Triple Quadrupole ICP-MS, Agilent application note, 5991-0320EN.




Figure 1. Calibration curves for Si, P and S showing SQ results (upper) compared with MS/MS results (lower). ³²S was not measurable at the spiked concentrations in methanol in SQ mode due to the intense ¹⁶O₂⁺ polyatomic interference.

- > Return to table of contents
- > Search entire document

Environmental

The accurate measurement of selenium in reference materials using online isotope dilution	22
Accurate sulfur quantification in organic solvents using isotope dilution mass spectrometry	25
Removal of REE ⁺⁺ interference on arsenic and selenium	28
Removal of molybdenum oxide interference on cadmium	31
Feasibility study of fluorine detection by ICP-QQQ	33
ICP-QQQ with oxygen reaction mode for accurate trace-level arsenic analysis in complex samples	35
Avoidance of spectral overlaps on reaction product ions with O ₂ cell gas: Comparison of quadrupole ICP-MS and ICP-QQQ	39
Removal of complex spectral interferences on noble metal isotopes	42

The accurate measurement of selenium in reference materials using online isotope dilution

Naoki Sugiyama
Agilent Technologies, Japan

Keywords

selenium, environmental, agricultural, human health, online isotope dilution analysis, OIDA, oxygen mass-shift

Introduction

Selenium (Se) is an important element in environmental and agricultural studies and in human health, as it is an essential trace nutrient but is toxic in excess. The role of certain chemical forms of Se is the subject of significant research into cancer prevention. ICP-MS is the analytical method of choice for both total and speciated Se measurements, but Se is a difficult element to quantify accurately at trace levels by ICP-MS for several reasons:

- The signal for Se is low, since it is poorly ionized in the plasma due to its high Ionization Potential (IP) of 9.75 eV.
- Because Se is poorly ionized, it suffers signal suppression in high matrix samples — an issue that is further compounded by the lack of a suitable internal standard element with a similar mass and IP.
- All the analytically useful Se isotopes suffer from multiple spectral interferences, as summarized in Table 1.
- The resolution required to separate all of the spectral interferences is beyond the capabilities of sector-type high resolution (HR-)ICP-MS.

The Agilent 8800 ICP-QQQ in MS/MS mode has a unique ability to remove the complex spectral interferences

from all the Se isotopes shown in Table 1, allowing the use of Isotope Dilution (ID) analysis, which requires at least two interference-free isotopes. ID is the most accurate quantification technique as it is based on direct measurement of isotopic abundances in each sample, rather than a relative measurement of analyte response compared to a standard. As a result, it offers better traceability and improved correction of non-spectroscopic interferences encountered in high matrix sample analysis. This note describes the application of the Agilent 8800 ICP-QQQ using ID for the accurate quantification of Se in a range of certified reference materials (CRMs).

Experimental

Instrumentation: Agilent 8800 #100.

Plasma conditions: Preset plasma/General purpose.

Ion lens tune: Soft extraction tune: Extract 1 = 0 V, Extract 2 = -180 V.

CRC conditions: O₂ gas at 0.4 mL/min plus H₂ gas at 2.0 mL/min, Octopole bias = -18 V and KED = -6 V.

Acquisition parameters: MS/MS O₂ mass-shift method. The reaction of Se⁺ with O₂ to form SeO⁺ is endothermic ($\Delta H_r = 0.71$ eV), but the reaction is efficiently promoted using high collision energy using a low octopole bias voltage setting [1]. Preliminary studies have shown that low BEC for Se isotopes can be achieved via the addition of a small amount of H₂ in MS/MS O₂ mass-shift method.

Table 1. Spectral interferences on Se isotopes

Se isotope		Interference						
m	Abundance %	Isobaric	Argide	Oxides	Hydride	Chloride	Doubly charged	Dimer
77	7.63		³⁹ K ³⁸ Ar ⁺	⁶¹ Ni ¹⁶ O ⁺ , ⁵⁹ Co ¹⁸ O ⁺	⁷⁶ GeH ⁺ , ⁷⁶ SeH ⁺	⁴⁰ Ar ³⁷ Cl ⁺ , ⁴⁰ Ca ³⁷ Cl ⁺	¹⁵⁴ Sm ⁺⁺ , ¹⁵⁴ Gd ⁺⁺	
78	23.77	⁷⁸ Kr ⁺	⁴⁰ Ca ³⁸ Ar ⁺	⁶² Ni ¹⁶ O ⁺	⁷⁷ SeH ⁺	⁴¹ K ³⁷ Cl ⁺	¹⁵⁶ Gd ⁺⁺ , ¹⁵⁶ Dy ⁺⁺	³⁸ Ar ⁴⁰ Ar ⁺ , ³⁹ K ³⁹ K ⁺
80	49.61	⁸⁰ Kr ⁺	⁴⁰ Ca ⁴⁰ Ar ⁺	⁶⁴ Ni ¹⁶ O ⁺ , ⁶⁴ Zn ¹⁶ O ⁺ , ³² S ₂ ¹⁶ O ⁺ , ³² S ¹⁶ O ₃ ⁺	⁷⁹ BrH ⁺	⁴⁵ Sc ³⁵ Cl ⁺	¹⁶⁰ Gd ⁺⁺ , ¹⁶⁰ Dy ⁺⁺ ,	⁴⁰ Ar ⁴⁰ Ar ⁺ , ⁴⁰ Ca ⁴⁰ Ca ⁺
82	8.73	⁸² Kr ⁺	⁴² Ca ⁴⁰ Ar ⁺	⁶⁶ Zn ¹⁶ O ⁺	⁸¹ BrH ⁺	⁴⁵ Sc ³⁷ Cl ⁺	¹⁶⁴ Dy ⁺⁺ , ¹⁶⁴ Er ⁺⁺	

Method: Online isotope dilution analysis (OIDA) [2] was used. OIDA is a useful development of traditional isotope dilution, as it removes the time consuming step of spiking enriched-isotope standards into each individual sample. A ^{82}Se enriched standard purchased from Oak Ridge National Laboratory (USA) was prepared at the appropriate concentration and added via the standard online ISTD mixing kit to the samples. Product ions derived from the ^{16}O -atom addition transition were measured for the three most analytically useful isotopes of Se. On the 8800 ICP-QQQ, this is simply achieved by defining the acquisition method with Q1/Q2 settings: Q1=78/Q2=94, Q1=80/Q2=96 and Q1=82/Q2=98 for the Se isotopes at m/z 78, 80 and 82 respectively.

It should be noted that the use of MS/MS (where Q1 acts as a 1 amu mass filter) is essential for this measurement, as it ensures that only one Se isotope enters the cell for any given mass pair measurement, and only the ^{16}O atom addition is measured because the mass difference between Q1 and Q2 is 16 amu. This ensures that there is no overlap due to the precursor ions from different Se^+ isotopes giving SeO^+ product ions at the same mass, such as the $^{80}\text{Se}^{18}\text{O}^+$ product ion overlap on $^{82}\text{Se}^{16}\text{O}^+$, both of which appear at m/z 98.

ICP-QQQ in MS/MS mode thereby removes one of the critical limitations of reaction chemistry with ICP-QMS, where all the sample ions enter the cell together so no specific reaction transition can be defined. Each Se isotope mass pair was measured with an integration time of 1 s and three replicates.

Sample preparation: The CRMs were microwave digested using a Milestone ETHOS closed vessel microwave digestion system (Milestone, Sorisole, Italy) and following the manufacturer's recommended procedures. The final dilution factor of the samples varied from 250 to 500x.

Results and discussion

Study of cell gases for spectroscopic interference removal

Figure 1 shows the result of a preliminary study of the effects of the choice of cell gas on interference removal. The findings of the study showed that O_2/H_2 mass-shift (Figure 1) enables the measurement of ^{78}Se , ^{80}Se and ^{82}Se relatively free from interferences in the range of synthetic matrices tested.




Figure 1. Preliminary study of the effectiveness of different cell gases for interference removal. Four synthetic matrices likely to give rise to interferences on the Se isotopes were measured using each of the 4 different cell gas modes.

Measurement of Se in CRMs

The concentration of Se was determined in 12 different CRMs using the OIDA method. The CRMs were obtained from NIST (Gaithersburg MD, USA), GSJ Geochemical Reference Samples (Tokyo, Japan), Japan Society for Analytical Chemistry (Tokyo, Japan), and National Institute of Metrology (Beijing, China). The matrices included environmental waters (NIST 1643e and JSAC 0302-3 River Water), rock (JB-3 basalt), sedimentary rock (JSI-1 and NIST 1646a Estuarine Sediment), soil (JSAC0411 Volcanic Ash Soil), biological samples (NIST 1566a Oyster Tissue, NCSZC 81002 Human Hair, NIST 2976 Mussel Tissue), and plant materials (NIST 1575a Pine Needles, NIST 1515 Apple Leaves, NIST1573a Tomato Leaves).

Figure 2 shows the Se results for each CRM expressed as % recovery relative to the certified value. The measured results for Se were in good agreement with the CRM values (90%-112%), using two Se isotope pairs: 78/82 and 80/82. This demonstrates the effectiveness of the Agilent 8800 ICP-QQQ in MS/MS mode for the removal of multiple interferences on ⁷⁸Se, ⁸⁰Se and ⁸²Se.




Figure 2. Result of Se quantification using OIDA in various CRMs

References

1. Agilent 8800 Triple Quadrupole ICP-MS: Understanding oxygen reaction mode in ICP-MS/MS, 2012, Agilent technical overview, 5991-1708EN
2. On-line isotope dilution analysis with the 7700 Series ICP-MS: Analysis of trace elements in high matrix samples, Giuseppe Centineo, Jose Angel Rodriguez Castrillon and Esther Munoz Agudo, 2011, Agilent application note, 5990-9171EN

More information

The accurate measurement of selenium in twelve diverse reference materials using on-line isotope dilution with the 8800 Triple Quadrupole ICP-MS in MS/MS mode, Agilent application note, 5991-0259EN.

Accurate sulfur quantification in organic solvents using isotope dilution mass spectrometry

Glenn Woods
Agilent Technologies (UK) Ltd.

Based upon the published work "Accurate determination of S in organic matrices using isotope dilution ICP-MS/MS" JAAS 2012 DOI: 10.1039/c2ja30265a by:

Lieve Balcaen and Frank Vanhaecke, Ghent University, Department of Analytical Chemistry, Ghent, Belgium

Martin Resano, University of Zaragoza, Department of Analytical Chemistry, Zaragoza, Spain

Glenn Woods, Agilent Technologies UK Ltd., 5500 Lakeside, Cheadle Royal Business Park, SK8 3GR, UK

Keywords

sulfur, ID-MS, biodiesel, environmental, ethanol, NIST SRM 2773, oxygen mass-shift

Introduction

Accurate measurement of sulfur in aqueous and organic media is relatively difficult for ICP-MS due to intense spectral interferences from polyatomic ions formed mainly from oxygen and nitrogen. Sulfur is an important element in environmental terms as it forms SO_x when combusted, contributing to acid rain and photochemical smog. It is also a catalyst poison for some industrial processes and its accurate measurement can be critical.

Experimental

A quadrupole ICP-MS (ICP-QMS) with a collision/reaction cell set up for O_2 mass-shift reaction chemistry can be used to avoid the $^{16}\text{O}_2^+$ overlap on $^{32}\text{S}^+$ by converting the S^+ to SO^+ reaction product ions that are then measured at a new mass (m/z 48) that is free from the O_2^+ overlap. However, in practice, this approach has been of relatively limited use, as ICP-QMS has no way to reject existing ions at the mass of the new analyte product ions, so not all of the interferences are eliminated, particularly when complex or variable matrices are investigated. There has also been some limited success reported by using Xe as a reaction gas to attenuate the O_2 -based interference particularly on

the ^{34}S isotope. Neither of these approaches reduces the backgrounds significantly enough to allow reliable trace level measurement of S, and they do not necessarily preserve the S isotopic abundances. In this investigation, ethanol was used as an example organic solvent and the Agilent 8800 ICP-QQQ was used to determine S by ID-MS in a biodiesel reference material to assess the measurement accuracy of MS/MS mode with O_2 mass-shift for S determination.

Instrumentation: Agilent 8800 #100 with Micromist nebulizer (free aspiration). For organic solvent analysis, a narrow injector torch with id 1.0 mm (G3280-80005) and Pt cones were used. 20% O_2 balanced in Ar was introduced via an option gas flow line to prevent carbon build up.

Plasma conditions: Plasma conditions were optimized manually. (RF power = 1450 W, CRGS flow rate = 0.98 L/min, Option gas flow rate = 0.75 L/min and spray chamber temp. = -5°C).

CRC conditions: O_2 gas at 0.4 mL/min, Octopole bias = -9 V, KED = -8 V.

Sample: Biodiesel certified reference material NIST SRM 2773.

Results and discussion

When using mass-shift mode for sulfur (or any element) it is important to eliminate any potential interferences at the target mass of the reaction product ion, as well as on the primary element mass (the precursor ion). If the target mass suffers from interferences then the measurement would still be compromised. For sulfur, the corresponding isotopes are shifted as follows using M + 16 amu mass-shift:

$^{32}\text{S} \rightarrow \text{SO}$ at 48 amu

$^{33}\text{S} \rightarrow \text{SO}$ at 49 amu

$^{34}\text{S} \rightarrow \text{SO}$ at 50 amu

Unfortunately, the SO^+ product ion masses (m/z 48, 49 and 50) can suffer from multiple interferences including Ca^+ , Cr^+ , V^+ , Ti^+ , ArC^+ and CCl^+ in natural samples. Furthermore the ^{33}S and ^{34}S isotopes can suffer from overlaps due to other combinations of SO^+ product ions, as well as pre-existing ions at the target mass. For example, the $^{34}\text{S}^{16}\text{O}^+$ product ion formed at m/z 50 is overlapped by $^{32}\text{S}^{18}\text{O}^+$ and $^{33}\text{S}^{17}\text{O}^+$, as well as $^{50}\text{Cr}^+$, $^{50}\text{V}^+$, $^{50}\text{Ti}^+$, $^{38}\text{Ar}^{12}\text{C}^+$, and $^{13}\text{C}^{37}\text{Cl}^+$. When operating the 8800 ICP-QQQ in MS/MS mass-shift mode, these overlaps are eliminated and the sulfur isotope pattern is

preserved. Figure 1 provides a graphical representation of the ICP-QQQ setup and the method of interference elimination.

This method would not be useful if the reaction were not quantitative, so to check for linearity, a blank ethanol sample was spiked with sulfur – see Figure 2. Despite the wide variation in absolute sensitivity for the different S isotopes, the BEC was the same for all three isotopes, indicating that the background is due to sulfur in the ethanol.




Figure 1. Mechanism of MS/MS mass-shift for sulfur isotope analysis. The mass difference between Q1 and Q2 is fixed at 16 amu, so only the $+^{16}\text{O}$ -atom transition is observed – the other oxygen isotope transitions are eliminated so the original sulfur isotopic pattern is preserved.




Figure 2. Ethanol with 0, 5, 10, 50 and 100 $\mu\text{g/L}$ sulfur spikes and corresponding calibration curves.

An isotope dilution (ID) method was used to evaluate the accuracy of the 8800 ICP-QQQ MS/MS method, using a biodiesel certified reference material (NIST SRM 2773) and an enriched ³⁴S spike. The biodiesel sample was simply diluted into the ethanol solvent and the appropriate spike added. Reproducibility was tested by analyzing three separate samples of the CRM. The results are presented in Table 1. Repeat measurements were within the expected recovery limits for the material.

Table 1. Isotope dilution analysis of S in diluted biodiesel reference material NIST 2773

Sample	Certified value (µg/g)	Determined conc. (µg/g)
SRM 2773 - 1	7.39 ± 0.39	7.234
SRM 2773 - 2		7.227
SRM 2773 - 3		7.231
Average		7.231
Standard Deviation		0.003
95% confidence interval		7.231 ± 0.015

Conclusions

Until the introduction of ICP-QQQ with MS/MS capability, it was impossible to obtain reliable results for reaction chemistry methods combined with an ID approach, using a quadrupole-based ICP-MS. The novel QQQ configuration of the 8800 ICP-QQQ enables operation in MS/MS mode, which ensures precise control over the reaction chemistry in the cell. This allows the unique isotopic information of the analyte to be retained, while removing the interferences that could affect both precursor and product ions of the target analyte.

Removal of REE⁺⁺ interference on arsenic and selenium

Kazumi Nakano and Yasuyuki Shikamori
Agilent Technologies, Japan

Keywords

Rare Earth Elements, REE, arsenic, selenium, environmental, food, CRMs, oxygen mass-shift

Introduction

Trace analysis of arsenic (As) and selenium (Se) in environmental and food samples is of a great interest, since both elements can be toxic even at quite low levels. It is difficult to quantify As and Se accurately at trace levels in some matrices by quadrupole ICP-MS as all the analytically useful isotopes can suffer from multiple spectral interferences, as summarized in Table 1. This application investigates ICP-QQQ in MS/MS reaction mode to remove interferences on As and Se,

Table 1. Selected spectral interferences on As and Se isotopes

As and Se isotope		Interference			
Element	m	Abundance %	Doubly charged	Chloride	Dimer
As	75	100	¹⁵⁰ Sm ⁺⁺ , ¹⁵⁰ Nd ⁺⁺	⁴⁰ Ar ³⁷ Cl ⁺ , ⁴⁰ Ca ³⁷ Cl ⁺	
Se	77	7.63	¹⁵⁴ Sm ⁺⁺ , ¹⁵⁴ Gd ⁺⁺	⁴⁰ Ar ³⁷ Cl ⁺ , ⁴⁰ Ca ³⁷ Cl ⁺	
	78	23.77	¹⁵⁶ Gd ⁺⁺ , ¹⁵⁶ Dy ⁺⁺	⁴¹ K ³⁷ Cl ⁺	³⁸ Ar ⁴⁰ Ar ⁺ , ³⁹ K ³⁹ K ⁺
	80	49.61	¹⁶⁰ Gd ⁺⁺ , ¹⁶⁰ Gd ⁺⁺	⁴⁵ Sc ³⁵ Cl ⁺	⁴⁰ Ar ⁴⁰ Ar ⁺ , ⁴⁰ Ca ⁴⁰ Ca ⁺
	82	8.73	¹⁶⁴ Dy ⁺⁺ , ¹⁶⁴ Er ⁺⁺	⁴⁵ Sc ³⁷ Cl ⁺	

with an emphasis on the removal of the doubly-charged ions arising from Rare Earth Elements (REE⁺⁺). While the concentration of REEs in environmental and food samples is usually low, some plants will accumulate REEs from the soil, and a high concentration will lead to false positive results for As and Se.

Experimental

Instrumentation: Agilent 8800 #100.

Plasma conditions: Preset plasma/Low matrix.

Ion lens tune: Soft extraction tune: Extract 1 = 0 V, Extract 2 = -180 V.

CRC conditions: O₂ gas flow rate of 0.2 mL/min, Octopole bias = -8 V and KED = -6 V.

Acquisition parameters: MS/MS O₂ mass-shift method to measure As⁺ (as AsO⁺) and Se⁺ (as SeO⁺), as illustrated in Figure 1. Unlike conventional quadrupole ICP-MS, the 8800 ICP-QQQ mass-shift method can be applied to complex matrix samples that may contain Zr and/or Mo. The MS/MS configuration prevents undesired ions such as ⁹¹Zr⁺ and ⁹⁴Mo⁺ from overlapping the MO⁺ product ions, as they are rejected by Q1.




Figure 1. ICP-QQQ MS/MS O₂ mass-shift method for measuring ⁷⁵As (top) and ⁷⁸Se (bottom)

Samples and sample preparation: SPEX XSTC-1 (a mixture of 10 ppm each of Ce, Dy, Er, Eu, Gd, Ho, La, Lu, Nd, Pr, Sm, Sc, Tb, Tm, Yb and Y) purchased from SPEX CertiPrep Ltd. (UK) was used. Four certified reference materials (CRMs): NIST 1515 Apple Leaves, NIST 1573a Tomato Leaves, NIST 1575a Pine Needles and NMIJ 7531a Brown Rice, were used for the method validation. It should be noted that NIST 1515 contains 3 mg/kg Sm and Gd, and 0.2 mg/kg Eu. NIST 1573a contains 0.19 mg/kg Sm, 0.17 mg/kg Gd, 5% Ca and 2.7% K, a combination of matrix elements that might be expected to cause severe interferences on As and Se. All CRMs were microwave-digested in HNO_3 and H_2O_2 , diluted and analyzed.

Results and discussion

Effectiveness of O_2 mass-shift method for removing REE⁺⁺ interferences

To investigate the effectiveness of interference removal modes on the 8800 ICP-QQQ, As and Se were measured in a mixed REE solution containing 1 ppm each of Ce, Dy, Er, Eu, Gd, Ho, La, Lu, Nd, Pr, Sm, Sc, Tb, Tm, Yb and Y.

Three different 8800 ICP-QQQ cell modes were used:

- Single Quad (SQ); no gas
- Single Quad (SQ); reaction mode using hydrogen (H_2) cell gas
- MS/MS; reaction mode using O_2 cell gas with + 16 amu mass-shift

“Single Quad” represents the performance of conventional ICP-QMS while MS/MS mode is unique to ICP-QQQ.

Figure 2 shows the BECs of As and Se in each of the measurement modes. The results in Figure 2 illustrate the excellent interference removal performance of the O_2 mass-shift method for the detection of As and Se in a matrix containing REEs.




Figure 2. BEC of As and Se in 1ppm REE mixed solution with three measurement modes: no gas, H_2 on-mass and O_2 mass-shift mode.

Figure 3 shows the product ion scan spectra obtained using O_2 mass-shift mode for a solution containing 1 ppm REEs without (left) and with (right) a 1 ppb As spike. As illustrated in the schematic, Q1 was fixed at $m/z = 75$ and Q2 was scanned across the selected mass range, to monitor all existing and cell-formed ions derived from precursor ions at $m/z 75$. Figure 3 (left) shows the product ions from $m/z 75$ in the blank REE matrix; the signal at Q2 $m/z = 75$ (mass of As) is due to REE^{2+} . The absence of a signal at $m/z = 91$ (the mass of AsO^+) in the blank REE matrix, indicates that the REEs do not react with O_2 in the cell to give rise to product ions (such as $REEO_2^{2+}$) that overlap AsO^+ at $m/z 91$. Consequently, As can be successfully measured as AsO^+ at $m/z = 91$ as shown in Figure 3 (right).

Method validation with CRMs

The ICP-QQQ method was applied to the measurement of As and Se in four CRMs. Table 2 summarizes the results. The measured concentrations of As and Se in the CRMs were all in good agreement with the certified values.




Figure 3. Product ion scan spectrum of O_2 mass-shift method. (Left) 1 ppm mixed-REE solution and (right) 1 ppm mixed-REE plus 1 ppb As spike.

Table 2. Results of the determination of As and Se in four CRMs using MS/MS O_2 mass-shift mode on the 8800 ICP-QQQ

	As (as AsO^+ at $m/z 91$)			Se (as SeO^+ at $m/z 94$)		
	Certified	Found average	Recovery	Certified	Found average	Recovery
	mg/kg	mg/kg	%	mg/kg	mg/kg	%
NIST1515 Apple Leaves	0.038±0.007	0.037	97	0.050±0.009	0.050	100
NIST1575a Pine Needles	0.039±0.002	0.038	97	0.099±0.004	0.099	100
NIST1573a Tomato Leaves	0.112±0.004	0.113	101	0.054±0.003	0.058	107
NMIJ 7531a Brown Rice	0.280±0.009	0.258	92	NA	0.032	NA

Removal of molybdenum oxide interference on cadmium

Michiko Yamanaka
Agilent Technologies, Japan

Keywords

cadmium, molybdenum oxide, environmental, food, CRMs, hydrogen on-mass

Introduction

Cadmium (Cd) is a well-known toxic element along with As, Hg and Pb. The maximum contamination level of these elements in food, pharmaceuticals, drinking water, wastewater and other matrices is strictly controlled under national and international legislation. Out of the eight natural isotopes of Cd, only ^{111}Cd is free from direct overlap by an atomic isobar (an isotope of a different element at the same mass as the Cd isotope), and even ^{111}Cd is potentially subject to spectroscopic interference by $^{95}\text{MoO}^+$. Fortunately, the concentration of Mo is low in most samples, and quadrupole ICP-MS (ICP-QMS) operating in helium collision mode can remove the interference, allowing the accurate measurement of Cd. However, there are some cases where the Mo concentration is high and a better interference removal technique is required in order to accurately determine Cd. This paper describes the application of MS/MS H_2 reaction mode on the Agilent 8800 ICP-QQQ for the determination of trace Cd in the presence of a high concentration of Mo.




Figure 1. (Left): Signal for m/z 111 with 10 ppm Mo and 10 ppm Mo + 1 ppb Cd, plotted as a function of H_2 flow rate. (Right): Estimated Cd BEC in the presence of 10 ppm Mo as a function of H_2 flow rate.

Experimental

Instrumentation: Agilent 8800 #100. Indium (In) was introduced as the internal standard using the on-line ISTD kit.

Plasma conditions and ion lens tune: RF power = 1550 W; sampling depth = 8.0 mm; carrier gas = 1.01 L/min; make up gas/dilution gas (MUGS/DiGS) = 0.0 L/min; Soft extraction tune: Extract 1 = 0 V, Extract 2 = -165 V, Omega bias = -100 V, Omega = 11.4 V.


CRC conditions: H_2 flow rate 9.0 mL/min, Octopole bias = -22 V, KED = +5 V.

Acquisition parameters: MS/MS H_2 on-mass method i.e. ^{111}Cd was measured at m/z 111 using quadrupole settings of (Q1 = 111, Q2 = 111).

Results and discussion

Optimization of H_2 flow rate

Figure 1 (left) shows the signal at m/z 111 for a 10 ppm Mo solution and a 10 ppm Mo + 1 ppm Cd solution, plotted as a function of H_2 flow rate. Figure 1 (right) shows the calculated BEC of Cd in the presence of 10 ppm Mo. The optimum cell gas flow rate of 9.0 mL/min was used for subsequent experiments.



In order to test the effectiveness of MS/MS mode with H₂ cell gas in comparison to no gas mode, a spike recovery test of 1 ppb Cd in a series of Mo matrix solutions ranging from 0.1 to 100 ppm was conducted. Figure 2 summarizes the results. In no gas mode, the error in quantification of the 1 ppb Cd spike dramatically increases with the concentration of Mo; in contrast, H₂ reaction mode delivers a consistent and accurate result for Cd even in the presence of 100 ppm Mo.

Method validation: Analysis of Cd in CRMs

The concentration of Cd was determined in four different CRMs: NIST 1515 Apple Leaves, NIST 1573a Tomato Leaves, NIST 1575a Pine Needles and NMIJ 7531a

Brown Rice Flour (National Metrology Institute of Japan). Each sample was microwave digested following the manufacturer's recommended procedures, then diluted and analyzed by ICP-QQQ; the final dilution factor was around 100–200. For each CRM, the digested sample was analyzed using the developed method. A second sample of each CRM was prepared and analyzed after the addition of a 10 ppm Mo spike. As summarized in Table 1, good recoveries were obtained for all four references materials both for the unspiked samples and the duplicates with the high added Mo concentration, demonstrating the validity of the method for real sample analysis.




Figure 2. 1 ppb Cd spike recovery in a series of Mo matrix solutions using no gas mode and MS/MS H₂ mode

Table 1. Measurement of Cd in four CRMs using the 8800 ICP-QQQ in MS/MS mode with H₂ reaction gas

CRMs	Certified mg/kg	Without Mo addition		With 10 ppm Mo addition	
		Determined mg/kg	Recovery %	Determined mg/kg	Recovery %
NIST 1515 Apple Leaves	0.014	0.013	93	0.016	115
NIST1573a Tomato Leaves	1.52	1.496	98	1.475	97
NIST1575a Pine Needles	0.223	0.220	99	0.224	100
NMIJ 7531a Brown Rice Flour	0.308	0.298	97	0.293	95

Feasibility study of fluorine detection by ICP-QQQ

Noriyuki Yamada
Agilent Technologies, Japan

Keywords

fluorine-containing polyatomic ions, barium, oxygen on-mass, ammonia mass-shift

Introduction

Fluorine (^{19}F) cannot be directly detected by conventional quadrupole ICP-MS (ICP-QMS) because of severe water-derived interferences at m/z 19 from $^1\text{H}_3^{16}\text{O}^+$ and $^1\text{H}^{18}\text{O}^+$, and extremely low sensitivity due to the fact that it is very difficult to convert fluorine atoms to the positive ions that are measured in ICP-MS. The interference problem can be resolved by high resolution ICP-MS, but the sensitivity issue remains a challenge because almost no F atoms are ionized in an argon plasma due to F having an ionization potential (17.423 eV) that is higher than that of Ar (15.760 eV).

However, fluorine-containing polyatomic ions (XF^+) can be formed in the plasma and they may be used to determine fluorine. Candidate ions are those with a high bond-dissociation energy for the $\text{X}^+\text{-F}$ bond and low ionization potential of X or XF. Since oxygen is present in the plasma (from the water matrix or from air entrainment), the formation of XO^+ or XO often competes against that of XF^+ . Therefore, a low bond-dissociation energy for $\text{X}^+\text{-O}$ and X-O bonds (low affinity of X^+ and X for O) is also desirable for the efficient formation of XF^+ . Barium was selected as "X" for this feasibility study, based on its thermochemical properties (Table 1).

Experimental

Instrumentation: Agilent 8800 #200 with a Micromist nebulizer.

Plasma conditions and ion lens tune: RF power = 1500 W; Sampling depth = 8 mm; CRGS flow rate = 1.00 L/min; sample uptake rate 0.33 mL/min; 100 ppm Ba uptake rate = 0.03 mL/min; MUGS flow rate = 0.32 L/min; Extract 1 = -150 V, Extract 2 = -4 V.

CRC conditions: O_2 gas at 1 mL/min (100%), Octopole bias = -60 V, Energy discrimination = -10 V in O_2 mode; 10% NH_3 /90% He flow rate 8.5 mL/min (85%), Octopole bias = -20 V, Energy discrimination = -10 V in NH_3 mode.

Acquisition parameters: MS/MS O_2 on-mass and MS/MS NH_3 mass-shift. Integration time per mass for BaF and $\text{BaF}(\text{NH}_3)_3$ = 1 sec; integration time per mass for $\text{BaF}(\text{NH}_3)_4$ = 10 sec.

In order to produce BaF^+ in the plasma, Ba solution was mixed online with fluorine standards per a fixed mixing ratio of 1:10. The mixing occurred just before the nebulizer. BaF^+ was efficiently formed under general plasma conditions with the BaO^+/Ba^+ ratio at about 11%. Under hotter plasma conditions, the formation of BaF^+ decreases because it tends to break apart. Under cooler plasma conditions, the formation of BaF^+ also decreases because of the formation of BaO^+ or, possibly, BaO . The signal intensity of BaF^+ was proportional to the concentration of Ba, which was fixed at about 10 ppm (after mixing).

Table 1. Gas phase thermochemical properties of elements having an affinity for fluorine*

Element X	$D_0(\text{X}^+\text{-F})$	IP (X)	$D_0(\text{X-F})$	IP (XF)	$D_0(\text{X}^+\text{-O})$	$D_0(\text{X-O})$
C	7.77	11.27	5.60	9.11	8.35	11.15
Al	3.16	5.99	6.99	9.73	1.81	5.31
Si	7.01	8.15	5.69	7.54	4.99	11.49
Ba	6.39	5.21	5.98	4.70	5.60	5.80
La	6.83	5.61	6.86	5.56	8.73	8.50
Eu	6.05	5.67	5.59	5.90	4.00	5.90

*Unit: eV. $D_0(\text{A-F})$ is the bond-dissociation energy for A-F bond (affinity of A for F) and $\text{IP}(\text{B})$ is the ionization potential of B.

Interference removal using MS/MS mode

$^{138}\text{Ba}^{19}\text{F}^+$ ($m/z=157$) suffers an interference from $^{138}\text{Ba}^{18}\text{O}^+\text{H}^+$. O_2 and NH_3 were tested as reaction gases to reduce the interference. It was found that O_2 reacts with BaOH^+ more efficiently than it reacts with BaF^+ in high energy reaction mode (octopole bias < -50 V). Therefore, using MS/MS mode, a mass pair ($\text{Q1} \rightarrow \text{Q2}$) = ($157 \rightarrow 157$) was selected to detect BaF^+ in O_2 mode. With Q1 set to 157 amu, $^{138}\text{Ba}^+$ was prevented from entering the cell and forming new interferences through unwanted reactions.

NH_3 was found to react with BaF^+ at a high NH_3 flow rate to form $\text{BaF}(\text{NH}_3)_n^+$, where $n = 2, 3, 4$. The most abundant complex ion was $\text{BaF}(\text{NH}_3)_3^+$ at $m/z = 208$, but $\text{BaF}(\text{NH}_3)_4^+$ at $m/z = 225$ was preferable in terms of signal to background ratio or BEC. Mass pairs ($\text{Q1} \rightarrow \text{Q2}$) = ($157 \rightarrow 208$) and ($157 \rightarrow 225$) were selected in NH_3 mode.




Figure 1. Calibration curve for F measured as BaF^+ in O_2 mode

Results and discussion

Figures 1 and 2 show calibration curves up to 2 mg/L (ppm) for fluorine in deionized water. The lowest detection limit (27 ppb) was obtained in O_2 mode. The lowest BEC (87 ppb) was obtained by measuring $\text{BaF}(\text{NH}_3)_4^+$ in NH_3 mode. Table 2 shows the BEC and DL results for F obtained from this study in comparison with the literature values.

Conclusions

Based on this preliminary study, it is clear that the controlled reaction chemistry that is possible with MS/MS mode on the 8800 ICP-QQQ can provide a novel approach to the measurement of F by ICP-MS. In addition to demonstrating detection limits that are comparable with published data measured using conventional quadrupole ICP-MS or high-resolution ICP-MS, the 8800 ICP-QQQ also allows unprecedented flexibility to monitor specific reaction transitions, making it invaluable for method development.




Figure 2. Calibration curve for F measured as $\text{BaF}(\text{NH}_3)_4^+$ in NH_3 mode

Table 2. Analytical performance for fluorine detection by ICP-MS

Analyte ion	Sensitivity [cps/ppm]	BEC [ppm]	DL [ppm]	Technique, reference
F^-	60,000	NA	0.11	Negative ion mode ICP-MS, Appl. Spectrosc, 42, 425 (1988)
F^+	3,000	NA	0.023	He-ICP-MS, Japan analyst 52(4), 275-278, 2003
Al^+ (AlF^{2+} complex)	NA	0.0033	0.0001	IC-ICP-MS (indirect determination), Analyst. 1999 Jan;124(1):27-31
F^+	26	2.05	5.07	HR-ICP-MS, J. Anal. At. Spectrom, 18, 1443, 2003
BaF^+	4,073	0.23	0.027	ICP-QQQ, O_2 mode, this work
$\text{BaF}(\text{NH}_3)_3^+$	929	0.17	0.043	ICP-QQQ, NH_3 mode, this work
$\text{BaF}(\text{NH}_3)_4^+$	67	0.087	0.031	ICP-QQQ, NH_3 mode, this work

ICP-QQQ with oxygen reaction mode for accurate trace-level arsenic analysis in complex samples

Ed McCurdy and Glenn Woods
Agilent Technologies (UK) Ltd.

Keywords

arsenic, zirconium, doubly-charged ion interferences, oxygen mass-shift

Introduction

Arsenic (As), with its high first ionization potential and single isotope at mass 75, is one of the most difficult elements to measure accurately by ICP-MS, particularly in complex matrices. The polyatomic interferences from ArCl^+ and CaCl^+ that overlap As^+ at mass 75 can be removed effectively using quadrupole ICP-MS (ICP-QMS) in helium collision mode, but collision mode cannot resolve the doubly-charged ion interferences from $^{150}\text{Nd}^{++}$ and $^{150}\text{Sm}^{++}$. A quadrupole mass spectrometer separates ions based on their mass to charge ratio (m/z), so doubly-charged ions appear at half their true mass; $^{150}\text{Nd}^{++}$ and $^{150}\text{Sm}^{++}$ therefore give an apparent overlap on As at mass 75.

Oxygen reaction mode (O_2 mode) offers a solution to these doubly-charged ion overlaps, since As can be converted to a reaction product ion $^{75}\text{As}^{16}\text{O}^+$, measured at m/z 91, where it is separated from the doubly charged Nd and Sm, which do not form such product ions. However, the new mass of the AsO^+ product ion is also overlapped by an isotope of zirconium ($^{91}\text{Zr}^+$). The presence of Zr in a sample may therefore cause an error in the results for As measured as AsO^+ using O_2 reaction mode on ICP-QMS.

ICP-QQQ solves this problem, as MS/MS mode allows all masses apart from m/z 75 (including the $^{91}\text{Zr}^+$ ions) to be rejected by the first quadrupole (Q1), ensuring that the AsO^+ product ions can be measured free from overlap. ICP-QQQ with MS/MS therefore allows the accurate determination of As in complex samples that contain any combination of Cl, Ca, Nd, Sm and Zr.

Experimental

Reagents and sample preparation: All of the sample matrices used for this work were prepared using single-element stock solutions (Spex CertiPrep, Claritas grade). The acid matrix and elemental standard concentrations are shown in the caption for each spectrum and are representative of the acid matrix (dilute HNO_3/HCl) and matrix levels commonly found in ICP-MS samples.

The sample matrices investigated were:

- Dilute nitric acid (1% HNO_3)
- Dilute hydrochloric acid (5% HCl)
- Calcium (100 ppm)
- Neodymium and samarium (1 ppm each element)
- Zirconium (0.5 ppm)

Instrumentation: Agilent 8800 #100.

Plasma conditions and ion lens tune: Preset plasma/General purpose, Soft extraction tune: Extract 1 = 0 V, Extract 2 = -170 V.

Acquisition conditions: Four operational modes were used, to investigate the different interference removal performance provided by the different cell modes:

- Single Quad (SQ); no gas
- Single Quad (SQ); collision mode (using helium (He) cell gas at a flow rate of 4 mL/min)
- Single Quad (SQ); reaction mode (using oxygen (O_2) cell gas at a flow rate of 0.2 mL/min).
- MS/MS; reaction mode (using O_2 cell gas at a flow rate of 0.2 mL/min)

KED bias voltage was +5 V in no gas and He mode, and -8 V in O_2 mode.

The three "Single Quad" modes represent the performance available on conventional ICP-QMS operating in collision or reaction mode. MS/MS mode is unique to the tandem mass spectrometer configuration of the 8800 ICP-QQQ.

Results and discussion

Figures 1a, 1b and 1c illustrate how Single Quad mode with He cell gas is effective at removing the common ArCl⁺ and CaCl⁺ polyatomic interferences on As⁺ at *m/z* 75, but is ineffective against the Nd⁺⁺/Sm⁺⁺ interferences.




Figure 1. a) As⁺ (*m/z* 75) in no gas mode, showing polyatomic interferences from ArCl⁺ and CaCl⁺; **b)** ArCl⁺ and CaCl⁺ polyatomics are removed in He collision mode; **c)** He collision mode fails to remove Nd⁺⁺ and Sm⁺⁺ interferences at *m/z* 75.

Figures 2a and 2b show how Single Quad mode with O₂ reaction gas successfully avoids the doubly-charged Nd and Sm interferences by mass-shifting the As to the new AsO⁺ product ion mass at *m/z* 91; but O₂ reaction mode on ICP-QMS cannot remove the ⁹¹Zr⁺ overlap on the AsO⁺ product ion.




Figure 2. a) Nd⁺⁺ and Sm⁺⁺ interferences at *m/z* 75 are avoided in SQ O₂ reaction mode, by measuring As as the AsO⁺ product ion at *m/z* 91; **b)** SQ O₂ reaction mode fails to remove ⁹¹Zr⁺ overlap on the AsO⁺ product ion.

Figure 3 shows that the 8800 ICP-QQQ in MS/MS mode with O₂ reaction gas provides reliable and consistent measurement of As (as AsO⁺) in all matrices. All the original polyatomic and doubly-charged interferences at *m/z* 75 are avoided by mass-shifting the As to *m/z* 91; and in MS/MS mode the ⁹¹Zr⁺ ion is removed by Q1, so the potential overlap on the AsO⁺ product ion at *m/z* 91 is also removed.

Conclusions

With the combination of O₂ reaction mode and MS/MS operation, the 8800 ICP-QQQ provides a reliable approach to the accurate measurement of As in complex samples. All the polyatomic and doubly-charged interferences that affect As measurement at its native mass (*m/z* 75) are avoided by using O₂ mode to mass-shift the As to its AsO⁺ product ion, measured at *m/z* 91. Furthermore, uniquely to the 8800 ICP-QQQ, MS/MS mode also eliminates potential native ion overlaps at *m/z* 91, as they are rejected by Q1 that is set to *m/z* 75 when measuring As.




Figure 3. MS/MS mass-shift with O₂ reaction mode provides consistent, interference-free measurement of As as AsO in all the matrices.

Avoidance of spectral overlaps on reaction product ions with O₂ cell gas: Comparison of quadrupole ICP-MS and ICP-QQQ

Ed McCurdy
Agilent Technologies (UK) Ltd.

Keywords

titanium, reaction chemistry, oxygen mass-shift

Introduction

The 8800 ICP-QQQ opens up many new analytical possibilities and novel methodologies for interference removal based on reaction chemistry. The major benefit provided by the 8800 ICP-QQQ is its unique tandem mass spectrometer configuration, which permits operation in MS/MS mode. In MS/MS, the first quadrupole (Q1) operates as a 1 amu mass filter, providing precise selection of the ions that can enter the reaction cell, and therefore control of the reaction processes that occur. This level of reaction process control is fundamentally different to the operation of conventional quadrupole ICP-MS (ICP-QMS) when using these same reaction chemistries, as ICP-QMS has no way to reject ions before they enter the cell, and so cannot select which ions are involved in the reactions.

This difference is apparent in many reaction chemistries, including both on-mass measurements (where the interfering ions are reactive and are moved away from the analyte ions, which are then measured at the natural mass), and mass-shift methods (where the analyte ions are reactive and are moved to a new product ion

Table 1. Potential native ion overlaps on TiO⁺ product ions in O₂ reaction mode

Precursor ion (Q1)	Product ion (Q2)	Potential overlaps from other analytes		
Ti	TiO	Ni	Cu	Zn
46	62	⁶² Ni	–	–
47	63	–	⁶³ Cu	–
48	64	–	–	⁶⁴ Zn
49	65	–	⁶⁵ Cu	–
50	66	–	–	⁶⁶ Zn

mass that is free from the original overlap). Overlaps on analyte product ions commonly occur in ICP-QMS and can give severe errors in results, especially in cases where the sample matrix or co-existing analyte levels vary from sample to sample.

In this note, we compare the performance of ICP-QMS (the 8800 ICP-QQQ operated in Single Quad mode with Q1 as a bandpass filter) and ICP-QQQ (the 8800 ICP-QQQ operated in MS/MS mode) for the measurement of titanium (Ti) as TiO⁺ product ions, using oxygen reaction mode (O₂ mode).

The native ion overlaps that could affect the measurement of TiO⁺ product ions with oxygen reaction gas are shown in Table 1. It should be noted that these native ion overlaps cannot be rejected by the cell bandpass settings of a conventional quadrupole ICP-MS, because they occur at the same mass as the analyte product ion being measured.

Experimental

For the spectral comparison, scan data were collected for the mass range from *m/z* 60 to 69, covering the TiO⁺ product ions formed from Ti in O₂ reaction mode.

Instrumentation: Agilent 8800 #100.

Plasma conditions and ion lens tune: Preset plasma/General purpose, Soft extraction tune: Extract 1 = 0 V, Extract 2 = -180 V.

CRC conditions: Cell gas = O₂ gas at 0.3 mL/min, Octopole bias = -5 V, KED = -7 V.

Acquisition parameters: Scan range = *m/z* 60 to 69; points per peak = 20; integration time per mass = 1 sec.

Results and discussion

The comparative results for TiO^+ measured in Single Quad (SQ) mode and MS/MS mode are shown in the overlaid spectra in Figures 1 and 2. In both cases, the TiO^+ ions at mass 62, 63, 64, 65 and 66 (from the 5 isotopes of Ti at 46, 47, 48, 49 and 50, respectively) are shown, measured using the same O_2 reaction mode conditions for both modes. The four solutions measured for the overlaid spectra are:

1. 1 ppb Ti in 1% HNO_3
2. 1 ppb Ti + 10 ppb Ni in 1% HNO_3
3. 1 ppb Ti + 10 ppb Cu in 1% HNO_3
4. 1 ppb Ti + 10 ppb Zn in 1% HNO_3

The overlaid spectra in Single Quad mode, shown in Figure 1, show that the peaks for the five TiO^+ isotopes match the theoretical isotopic template in the 1 ppb Ti sample. However, in the other samples containing the elements Ni, Cu and Zn, all of the TiO^+ product

ions suffer significant overlap from the native Ni (m/z 62), Cu (m/z 63 and 65) and Zn (m/z 64 and 66) ions. Unexpected or variable levels of these common elements would lead to an error in the reported results for Ti measured as TiO^+ using quadrupole ICP-MS in O_2 reaction mode.

In contrast, the overlaid spectra for MS/MS mode, shown in Figure 2, demonstrate consistent measurement of all five TiO^+ product ions in all four solutions. The presence of the other elements Ni, Cu and Zn has no impact on the TiO^+ peaks and all five TiO^+ product ion isotopes could be used to give reliable results for Ti in these variable samples. This illustrates how MS/MS mode on the 8800 ICP-QQQ can simplify method development, because consistent cell conditions, acquisition parameters and isotope selection can be used for a range of variable sample types. A further benefit is that interferences are removed from all isotopes under the same cell conditions, so secondary (or qualifier) isotopes become available for data confirmation or isotope analysis.




Figure 1. Overlaid spectra for TiO^+ product ions in variable samples measured using SQ mode (see text for sample composition).




Figure 2. Overlaid spectra for TiO^+ product ions in variable samples measured using MS/MS mode (see text for sample composition)

Conclusions

The comparative spectra presented in this note illustrate the improved accuracy and consistency delivered by ICP-QQQ operating in MS/MS mode, compared to a conventional quadrupole ICP-MS using a reaction cell with bandpass filter. By rejecting non-target native ions that would occur at the same mass as analyte product ions, potential interferences can be eliminated by MS/MS. This allows simpler, more consistent method development, as well as improving accuracy for interfered elements in complex and variable samples.

Removal of complex spectral interferences on noble metal isotopes

Naoki Sugiyama
Agilent Technologies, Japan

Keywords

Platinum Group Elements, gold, silver, ore exploration, geochemical, environmental, catalytic converter, ammonia

Introduction

The precise determination of the noble metals, comprising the Platinum Group Elements (PGEs: Ru, Rh, Pd, Os, Ir and Pt), Au and Ag, is of great interest in areas such as ore exploration and geochemical studies, and these metals are increasingly used for industrial applications including advanced materials and alloys, medical devices, and catalysts for pharmaceutical manufacturing. Environmental monitoring is also required as some of these elements are used in automobile catalytic converters. ICP-MS is widely used for these applications due to its high sensitivity

and multi-element capability. However, the analysis is challenging because the metal concentrations are often low and they are subject to severe spectral overlaps.

Table 1 summarizes the interferences and abundance (%) of each isotope of the elements (the isotopes highlighted in yellow represent the recommended isotope for determination by ICP-MS). Several methods have been developed to resolve the interferences, such as mathematical correction, matrix removal and high-resolution magnetic sector (HR-)ICP-MS. However the mass resolution required to separate some of the interferences is beyond the capability of current commercial HR-ICP-MS. For example separation of $^{103}\text{Rh}^+$ from $^{87}\text{Sr}^{16}\text{O}^+$, $^{105}\text{Pd}^+$ from $^{89}\text{Y}^{16}\text{O}^+$, and $^{109}\text{Ag}^+$ from $^{93}\text{Nb}^{16}\text{O}^+$ requires mass resolution ($M/\Delta M$) of 102900, 27600 and 31500, respectively; commercial HR-ICP-MS instruments

Table 1. Abundance (%) of each isotope of noble metals and the interference on each isotope

m		96	97	98	99	100	101	102	103	104	105	106	107	108	109	110
Element	Ru	5.54		1.87	12.76	12.6	17.06	31.55		8.62						
	Rh								100							
	Pd							1.02		11.14	22.33	27.33		26.46		11.72
	Ag												51.84		48.16	
Interference	Atomic	Zr, Nb, Mo		Mo		Mo						Cd		Cd		Cd
	MH			MoH	MoH		MoH									
	MO, MOH	SeO, BrOH		SeO		SrO	RbO	SrO	SrO, RbO	SrO	YO, SrOH	YOH, ZrO	ZrO	ZrO, MoO	NbO	ZrO, MoO
	Argide			NiAr		NiAr	NiAr	NiAr	CuAr	ZnAr	CuAr	ZnAr	ZnAr	ZnAr		
	others			CuCl	ZnCl	CuCl	ZnCl	CuCl, ZnCl	ZnCl, Pb ⁺⁺	ZnCl	ZnCl					

m		184	185	186	187	188	189	190	191	192	193	194	195	196	197	198
Element	Os	0.02		1.59	1.96	13.24	16.15	26.26		40.78						
	Ir								37.3		62.7					
	Pt							0.014		0.782		32.97	33.83	25.24		7.163
	Au														100	
Interference	Atomic	W		W	Re									Hg		Hg
	MH				WH											
	MO, MOH			YbO	YbO	YbO	YbO	YbO	LuO	YbO, LuD, HfO	HfO	HfO	HfO	HfO	TaO, HfOH	Wo, TaOH
	Argide	NdAr	NdAr	NdAr	SmAr	SmAr, NdAr	SmAr	SmAr, NdAr	EuAr	SmAr	EuAr	SmAr, GdAr	GdAr	GdAr	GdAr	GdAr
	others															

are limited to a maximum resolution of 10,000. To remove the multiple, complex interferences on noble elements, the Agilent 8800 ICP-QQQ was used in MS/MS mode, using ammonia as the reaction gas.

Experimental

Instrumentation: Agilent 8800 #100.

Plasma conditions: Preset plasma/Low matrix.

Ion lens tune: Soft extraction tune: Extract 1 = -3 V, Extract 2 = -200 V.

CRC conditions: NH₃ (10% NH₃ in He) was used as CRC gas in MS/MS mode.

Following a preliminary optimization study, three different NH₃ gas flow rates (low (L), medium (M) and high (H)) were used. Cell conditions are given in Table 2. No gas mode was also applied for comparison purposes.

Table 2. CRC conditions

	No gas	NH ₃ -L	NH ₃ -M	NH ₃ -H
Cell gas	na	NH ₃	NH ₃	NH ₃
Gas flow rate (mL/min)	na	2.0	3.0	5.0
Octopole bias (V)	-8	-5	-10	-12
KED (V)	+5	-8		
Cell exit (V)	-90			
Deflect lens (V)	20	10	6	2
Plate lens (V)	-110			

Method

The BECs of the noble metals were determined in a series of synthetic-matrix samples, using an external calibration method. Indium (In) internal standard (ISTD) was mixed online with the sample via the standard ISTD mixing T-connector. An integration time of 1 s per isotope was used with 3 replicates (7 replicates for blank).

Samples and sample preparation

Standards and matrix samples were prepared from single element stock solutions purchased from Kanto Chemical Co., Inc. (Saitama, Japan) and a REE mixture standard, XSTC-1 purchased from Spex certiPrep. All solutions were diluted into a final acid mix of 1% HNO₃ and 3% HCl.

Results

Matrix interference study

Tables 3 and 4 summarize the results of the spectral interference study obtained by analyzing individual synthetic matrix blank solutions. Table 3 shows the observed interferences, expressed as BEC (ppb), in each matrix blank measured using no gas mode. As expected from Table 1, the synthetic matrices caused significantly elevated BECs (>> 1 ppb) on all the primary and secondary isotopes of all the analytes except for Ru; Rh suffered a relatively minor increase in BEC of ~0.5 ppb in the 10 ppm Pb/1 ppm Hg matrix.

Table 4 shows the results obtained using NH₃ reaction mode. The optimum gas flow rate for NH₃ for each element was investigated and three gas flow rates (Low: 2.0, Medium: 3.0, and High: 5.0 mL/min) were used. The best isotope and method is highlighted in bold in the Table. It can clearly be seen that NH₃ reaction mode effectively removes the interferences on all the analytes, giving BECs of << 0.1 ppb for the preferred isotope/cell mode in all the matrices. The mechanism for the removal of each interference using the MS/MS capability of the 8800 ICP-QQQ is as follows:

- **Ru:** slight interferences from Zn and Mo were resolved using on-mass method with NH₃-M.
- **Rh:** Pb⁺⁺ interference was resolved using on-mass method with NH₃-M.
- **Pd:** significant interferences from SrOH⁺ and YO⁺ were seen on ¹⁰⁵Pd, the only isotope free from atomic isobar. On-mass method with NH₃-H removed the interferences.
- **Ag:** significant ZrO⁺ interference on both ¹⁰⁷Ag and ¹⁰⁹Ag was resolved using on-mass method with NH₃-H.
- **Os:** YbO⁺ interference was observed on both ¹⁸⁸Os⁺ and ¹⁸⁹Os⁺. Since Os⁺ sensitivity in NH₃ mode is low, but Os⁺ forms a product ion of OsNH₃⁺, NH₃-L with mass-shift gave the best result.
- **Ir:** LuO⁺ and HfO⁺ interfere with ¹⁹¹Ir⁺ and ¹⁹³Ir⁺ respectively. NH₃-M with mass-shift method worked for ¹⁹¹Ir⁺ as Ir⁺ forms a product ion of IrNH₃⁺.
- **Pt:** ¹⁹⁵Pt⁺ suffers a significant interference from HfO⁺. While the overlap is less significant on ¹⁹⁸Pt⁺, ¹⁹⁸Pt⁺ suffers an atomic isobar interference from ¹⁹⁸Hg⁺. However Hg⁺ is effectively neutralized by NH₃ so ¹⁹⁸Pt⁺ can be measured free from interference.
- **Au:** significant interferences by TaO⁺ and HfOH⁺ are resolved by mass-shift method with NH₃-M. Au⁺ forms a product ion of Au(NH₃)₂⁺.

Table 3. Summary of spectral interferences in no gas mode, showing analyte BECs (ppb) in each matrix blank. Matrix overlaps that made a significant contribution to the analyte BECs are indicated in red (BEC > 10 ppb) and orange (BEC > 1 ppb).

	Ru		Rh	Pd	Ag		Os			
Isotope	99	101	103	105	107	109	188	189		
NH ₃ flow rate mL/min	NA									
Method	on-mass		on-mass							
Mass pair	99-99	101-101	103-103	105-105	107-107	109-109	188-188		189-189	
10 ppm Cu Zn	0.058	0.041	0.138	0.328	0.064	0.061	0.000		0.000	
10 ppm Sr Rb	0.000	0.034	0.150	4.39	0.005	0.001	0.000		0.000	
10 ppm Ni	0.007	0.019	0.000	0.022	0.012	0.016	0.000		0.000	
10 ppm Mo	0.059	0.018	0.000	0.004	0.000	0.018	0.000		0.000	
10 ppm Pb, 1 ppm Hg	0.000	0.000	0.472	0.002	0.033	0.034	0.000		0.000	
10 ppm Zr Nb	0.000	0.000	0.000	0.022	21.9	1.59	0.000		0.000	
10 ppm REE	0.004	0.000	0.009	165	0.147	0.005	2.78		2.99	
10 ppm Ta	0.008	0.000	0.000	0.004	0.003	0.000	0.000		0.000	
10 ppm Hf	0.000	0.000	0.000	0.004	0.312	0.026	0.000		0.000	
10 ppm W	0.000	0.000	0.000	0.003	0.001	0.001	0.000		0.000	

Table 3. continued

	Ir			Pt			Au	
Isotope	191	193		195	198		197	
NH ₃ flow rate mL/min	NA							
Method	on-mass		on-mass		on-mass		on-mass	
Mass pair	191-191		193-193		195-195		198-198	
							197-197	
10 ppm Cu Zn	0.003		0.002		0.000		0.279	0.001
10 ppm Sr Rb	0.002		0.000		0.001		0.310	0.004
10 ppm Ni	0.009		0.004		0.002		0.444	0.011
10 ppm Mo	0.000		0.000		0.000		0.295	0.000
10 ppm Pb, 1 ppm Hg	0.002		0.002		0.000		1293	0.000
10 ppm Zr Nb	0.002		0.775		1.98		3.17	0.417
10 ppm REE	123		0.712		0.788		2.17	0.138
10 ppm Ta	0.000		0.000		0.244		114	284
10 ppm Hf	0.071		28.1		70.9		2.34	14.1
10 ppm W	0.001		0.000		0.000		19.6	0.002

Table 4. Summary of spectral interferences in MS/MS NH₃ reaction cell mode, showing analyte BECs (ppb) in each matrix blank. Matrix overlaps that made a significant contribution to the analyte BECs are indicated in red (> 10 ppb) and orange (> 1 ppb).

Isotope	Ru		Rh	Pd	Ag		Os			
	99	101	103	105	107	109	188		189	
NH ₃ flow rate mL/min	3.0		3.0	5.0	5.0		2.0			
Method	on-mass	mass-shift	on-mass	mass-shift						
Mass pair	99-99	101-101	103-103	105-105	107-107	109-109	188-188	188-203	189-189	189-204
10 ppm Cu Zn	0.000	0.000	0.000	0.001	0.061	0.057	0.000	0.000	0.001	0.002
10 ppm Sr Rb	0.000	0.005	0.016	0.033	0.000	0.000	0.000	0.000	0.002	0.000
10 ppm Ni	0.000	0.000	0.000	0.000	0.010	0.009	0.000	0.000	0.000	0.000
10 ppm Mo	0.005	0.001	0.000	0.000	0.001	0.000	0.000	0.000	0.001	0.000
10 ppm Pb, 1 ppm Hg	0.000	0.000	0.000	0.001	0.033	0.035	0.000	0.000	0.001	0.000
10 ppm Zr Nb	0.000	0.000	0.000	0.000	0.003	0.000	0.000	0.000	0.002	0.001
10 ppm REE	0.000	0.000	0.000	0.014	0.004	0.004	2.79	0.003	5.85	0.010
10 ppm Ta	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.001	0.003	0.000
10 ppm Hf	0.000	0.000	0.000	0.000	0.000	0.000	0.000	0.001	0.058	0.000
10 ppm W	0.000	0.000	0.000	0.001	0.001	0.001	0.000	0.000	0.000	0.000

Table 4. continued

Isotope	Ir				Pt				Au	
	191		193		195		198		197	
NH ₃ flow rate mL/min	3.0				5.0		3.0		3.0	
Method	on-mass	mass-shift								
Mass pair	191-191	191-206	193-193	193-208	195-195	195-229	198-198	198-232	197-197	197-231
10 ppm Cu Zn	0.004	0.000	0.003	0.004	0.003	0.000	0.002	0.000	0.000	0.000
10 ppm Sr Rb	0.002	0.000	0.001	0.001	0.000	0.000	0.000	0.000	0.000	0.000
10 ppm Ni	0.004	0.000	0.001	0.004	0.000	0.000	0.000	0.002	0.000	0.000
10 ppm Mo	0.000	0.000	0.000	0.000	0.000	0.000	0.003	0.000	0.000	0.000
10 ppm Pb, 1 ppm Hg	0.002	0.001	0.000	0.000	0.000	0.000	0.005	0.001	0.000	0.000
10 ppm Zr Nb	0.017	0.000	0.679	0.066	0.031	0.009	0.003	0.001	0.000	0.000
10 ppm REE	44.3	0.019	1.56	0.019	0.031	0.002	0.000	0.000	0.044	0.003
10 ppm Ta	0.000	0.000	0.000	0.000	0.004	0.000	0.261	0.009	4.11	0.046
10 ppm Hf	0.690	0.095	21.4	2.40	0.904	0.115	0.141	0.070	0.070	0.003
10 ppm W	0.001	0.000	0.000	0.000	0.001	0.000	0.479	0.000	0.001	0.000

Analysis of complex synthetic matrix sample using optimized NH₃ reaction mode

A complex synthetic matrix sample containing 10 ppm each of Cu, Zn, Sr, Rb, Ni, Mo, Pb, Zr, Nb, REEs, Ta, Hf, W and 1 ppm Hg was prepared, and this matrix was spiked with 1 ppb each of Ru, Rh, Pd, Ag, Os, Ir, Pt and Au as analytes. The concentration of the noble metals was determined in two modes: No gas mode and NH₃ reaction cell mode, and the spike recovery results are displayed in Figure 1 for each mode. The results demonstrate that MS/MS mode with NH₃ reaction cell gas successfully removes multiple interferences on all the noble metals, providing accurate results for these analytes even in a complex and challenging matrix.




Figure 1. Result of synthetic matrix sample test. 1 ppb noble metals were measured in a multi-matrix sample containing 10 ppm of each Cu, Zn, Sr, Rb, Ni, Mo, Pb, Zr, Nb, REEs, Ta, Hf, W and 1 ppm Hg

- > Return to table of contents
- > Search entire document

Materials

Arsenic measurement in cobalt matrix using MS/MS mode with oxygen mass-shift	48
Determination of sulfur, phosphorus and manganese in high purity iron	50
Direct measurement of trace rare earth elements in high purity REE oxides	53
The benefits of improved abundance sensitivity with MS/MS for trace elemental analysis of high purity metals	55
Ultra trace copper analysis in a semiconductor grade organometallic titanium complex	58
Removal of MH^+ interferences in refined REE material analysis	60

Arsenic measurement in cobalt matrix using MS/MS mode with oxygen mass-shift

Katsuo Mizobuchi
Agilent Technologies, Japan

Keywords

arsenic, high purity metals, cobalt, zirconium, oxygen mass-shift

Introduction

Measuring the purity of materials such as high purity metals is of interest across advanced technology industries, to support the development of new materials and/or improve the performance of existing products. ICP-MS is widely used for determining elemental impurities in these materials due to its unique features: High sensitivity, low DLs, multi element analysis capability, wide dynamic range, fast analysis and minimal sample preparation requirements.

For many applications, the errors caused by spectral interferences in quadrupole ICP-MS have been adequately addressed by the introduction of CRC technology. However, the analysis of trace contaminants in high purity materials presents a particular challenge due to the high matrix levels and the need to determine impurities at the trace level. For example, the determination of As in Co is difficult for quadrupole ICP-MS due to the signal from CoO^+ that overlaps the only isotope of arsenic at m/z 75. Although only about 0.01% of the Co ions in the plasma are present as CoO^+ ions, the Co concentration in a 1000 ppm solution is 6 or 7 orders of magnitude higher than the trace levels of As that are of interest in this application. Consequently the CoO^+ interference is still very significant relative to the As^+ signal. This note describes the measurement of trace As in a 1000 ppm Co solution using an Agilent 8800 ICP-QQQ in MS/MS mass-shift mode, using oxygen as the reaction gas.

Experimental

Instrumentation: Agilent 8800 #100.

Plasma conditions: Preset plasma/HMI-mid.

CRC conditions: O_2 gas at 0.3 mL/min, Octopole bias = -5 V, KED = -7 V.

Acquisition conditions: Three oxygen (O_2) mass-shift operational modes were compared:

- Single Quad mode A with low mass cut off at $m/z < 59$


- Single Quad mode B with low mass cut off at $m/z \approx 59$
- MS/MS mode with Q1 as a 1 amu mass filter, Q1 = 75 and Q2 = 91

Sample: SPEX CLC02-2Y (SPEX CertiPrep Ltd., UK) was used as 1000 ppm Co solution.

Results and discussion


BEC of As in 1000 ppm Co solution using O_2 mass-shift method

From the equation and reaction enthalpy below, it can be seen that arsenic reacts readily with O_2 cell gas via an O-atom transfer reaction. This creates the reaction product ion AsO^+ at m/z 91, moving the analyte away from the CoO^+ interference on As^+ at m/z 75.



The reaction enthalpy for CoO^+ with the O_2 cell gas is much less favorable, so the overlap from the CoO^+ polyatomic interference is successfully avoided. To evaluate the effectiveness of MS/MS mode for this application, the 8800 ICP-QQQ was operated in three acquisition modes: MS/MS mode and two "Single Quad" modes, in which Q1 functions as a bandpass filter rather than a unit (1 amu) mass filter. In Single Quad mode A, Q1 was set to allow most of the plasma-formed ions to enter the cell; in Single Quad mode B, Q1 was set with a low mass cutoff around m/z 59 to allow only ions with a mass greater than 59 to enter the cell (most $^{59}\text{Co}^+$ ions are rejected); and finally in MS/MS mode Q1 was set to allow only ions at m/z 75 to enter the cell (all $^{59}\text{Co}^+$ rejected).

The BECs for As obtained using the three acquisition modes are shown in Figure 1. MS/MS mode achieved the lowest BEC for As of 330 ppt in 1000 ppm Co. The BEC obtained by the Single Quad modes were orders of magnitude higher, which suggests the occurrence of the following undesired reactions in the cell and indicates the incomplete rejection of Co^+ by Single Quad mode B:



Note that this sequential reaction chemistry leads to a relatively intense signal for CoO_2^+ , because the number




Figure 1. BEC of As in 1000 ppm Co solution with O₂ mass-shift method using three acquisition modes (note log intensity scale).

of precursor ions for the reaction (the Co⁺ ions from the plasma) is so high (10000 times higher intensity than the CoO⁺ signal in the plasma). Consequently, in Single Quad mode, the CoO⁺ overlap cannot be successfully avoided by moving the As⁺ to its AsO⁺ product ion at *m/z* 91 using O₂ cell gas, because CoO₂⁺ (also at *m/z* 91) is formed relatively easily when a large number of Co⁺ ions are allowed to enter the cell.

AsO⁺ in the presence of zirconium

To successfully avoid interferences using the mass-shift method, the mass of the analyte product ion must itself be free from interference. For example, in this application the AsO⁺ product ion is measured at *m/z* 91 where it could be overlapped by an isotope of zirconium (⁹¹Zr⁺). The presence of Zr in a sample may therefore

cause an error in the results for As measured as AsO⁺ using O₂ reaction mode on ICP-QMS. The potential effect of co-existing Zr on AsO⁺ measurement using ICP-QMS was investigated.

Figure 2 is a spectrum of 10 ppm Zr obtained using Single Quad mode A with O₂ mass-shift. Zr reacts with O₂ very efficiently ($\Delta H_r = -3.84$) and is converted to ZrO⁺. However not all the Zr⁺ ions are converted to ZrO⁺ so some Zr⁺ remains, interfering with the measurement of AsO⁺ at *m/z* 91. In contrast, in MS/MS mode the ⁹¹Zr⁺ ion is rejected by Q1, so the potential overlap on the AsO⁺ product ion at *m/z* 91 is removed.

Conclusions

Trace levels of arsenic in a 1000 ppm cobalt matrix can be successfully measured (BEC of 330 ppt) using the 8800 ICP-QQQ operating in MS/MS mass-shift mode, with oxygen as the reaction gas. There are two main advantages of using MS/MS compared to ICP-QMS:

1. In MS/MS mode, Co⁺ is prevented from entering the cell by Q1, which is set to *m/z* 75. In ICP-QMS, CoO₂⁺ is formed in the cell via a chain reaction, and will interfere with AsO⁺ at *m/z* 91.
2. In MS/MS mode, the potential ⁹¹Zr⁺ overlap on the AsO⁺ product ion at *m/z* 91 is eliminated, as ⁹¹Zr⁺ ions (and all other ions apart from *m/z* 75) are rejected by Q1.




Figure 2. Spectrum of 10 ppm Zr obtained using Single Quad mode A with O₂ mass-shift-method

Determination of sulfur, phosphorus and manganese in high purity iron

Yasuyuki Shikamori and Kazumi Nakano
Agilent Technologies, Japan

Keywords

phosphorus, sulfur, manganese, iron, steel, JSS 001-6, JSS 003-6, abundance sensitivity, oxygen mass-shift

Introduction

ICP-MS is the analytical technique of choice for the analysis of trace elements in iron and steel. However, the sensitivity and interference removal performance of quadrupole ICP-MS (ICP-QMS) is not sufficient for the determination of difficult analytes such as phosphorus (P) and sulfur (S) at the low levels required. Furthermore, the determination of manganese (Mn) in an iron matrix is extremely challenging for ICP-QMS due to overlap (or tailing) from the very intense ^{54}Fe and ^{56}Fe peaks that occur either side of the single isotope of manganese at m/z 55.

The Agilent 8800 Triple Quadrupole ICP-MS (ICP-QQQ) provides more effective and reliable removal of polyatomic interferences, such as $^{14}\text{N}^{16}\text{OH}$ on ^{31}P and $^{16}\text{O}_2$ on ^{32}S , using controlled chemical reaction in the CRC. This note describes the performance of the 8800 ICP-QQQ operating in MS/MS mode, for the determination of the trace elements S, P and Mn in two high purity iron CRMs (JSS 001-6 and 003-6).

Experimental

Instrumentation: Agilent 8800 #100.

Plasma condition: Preset plasma/HMI-mid.

Ion lens tune: Soft extraction tune: Extract 1 = 0 V, Extract 2 = -175 V.

CRC and acquisition conditions:

- MS/MS O_2 mass-shift method to remove the $^{14}\text{N}^{16}\text{OH}^+$ and $^{16}\text{O}_2^+$ interferences on $^{31}\text{P}^+$ and $^{32}\text{S}^+$ respectively: O_2 gas at 0.3 mL/min, Octopole bias = -5 V and KED = -7 V.
- MS/MS He on-mass mode to measure $^{55}\text{Mn}^+$: He gas at 5.0 mL/min, Octopole bias = -18 V and KED = 4 V.

All other parameters were optimized by Autotune in the MassHunter software. Figure 1 shows the mechanism used on the 8800 ICP-QQQ to avoid the $^{14}\text{N}^{16}\text{OH}^+$ and $^{16}\text{O}_2^+$ interferences on $^{31}\text{P}^+$ and $^{32}\text{S}^+$ by mass-shift mode ($Q1 \neq Q2$) using O_2 reaction gas.

Sample preparation: Two Steel CRMs, JSS-001 and JSS-003 were purchased from The Japan Iron and Steel Federation (Tokyo, Japan). 0.1 g of each Steel CRM was digested in a mixture of 1 mL HCl and 2 mL HNO_3 and diluted to 100 mL with UPW. No further chemical matrix separation, e.g., solvent extraction, ion exchange, etc. was applied. The digested CRM samples containing 0.1% (1000 ppm) Fe were analyzed directly on the ICP-QQQ using the robust plasma conditions provided by Agilent's HMI aerosol dilution system.





Figure 1. Mechanism of MS/MS mass-shift mode ($Q2 = Q1 + 16$) using O_2 reaction gas for the measurement of P as $^{31}\text{P}^{16}\text{O}^+$ and S as $^{32}\text{S}^{16}\text{O}^+$ at m/z 47 and 48 respectively.

Results and discussion

BEC and DL of P and S

The calibration plots shown in Figure 2 demonstrate that the 8800 ICP-QQQ with MS/MS mass-shift mode can successfully perform the trace level (single ppb) quantitation of P and S in 0.1% Fe solutions. The BEC and DL achieved for P were 0.14 ppb and 0.05 ppb respectively, and the BEC and DL for S were 6.45 ppb and 0.75 ppb respectively.



Trace Mn analysis in Fe matrix

The abundance sensitivity (AS, a measure of peak separation) of ICP-QQQ in MS/MS mode is the product of the Q1 AS x Q2 AS. This means the AS of the 8800 ICP-QQQ is theoretically about 2x that achievable on ICP-QMS, and the ICP-QQQ is therefore able to successfully separate the ^{55}Mn peak from the very intense overlaps from ^{54}Fe and ^{56}Fe in a high iron matrix. This is demonstrated in Figure 3 which shows the spectra of 10 ppb Mn in a 0.1% Fe matrix sample solution measured in Single Quad mode (left) and MS/MS mode on the ICP-QQQ (right). Helium was used as the cell gas in both cases to remove $^{54}\text{FeH}^+$ and ArNH^+ interferences by KED.




Figure 2. Calibration curve for P (left) and S (right) in 0.1% Fe matrix, obtained using O_2 mass-shift mode under robust plasma conditions.




Figure 3. Spectra of 10 ppb Mn in a 0.1% Fe matrix sample solution obtained in Single Quad mode (left) and ICP-QQQ in MS/MS mode (right).

Determination of P, S and Mn in high purity iron CRMs

Trace elements including P, S and Mn were determined by ICP-QQQ in high purity iron CRMs: JSS 001-6 and 003-6, using O₂ mass-shift mode (for P and S) and He mode (for Mn). As summarized in Table 1, excellent agreement was obtained between the measured (found) and certified values for all three elements, indicating the effective interference removal offered by the 8800 ICP-QQQ in MS/MS mode. Excellent spike recovery at the 50 ppb level was also confirmed with JSS 003.

Table 1. Analytical results for P, S and Mn in two high purity iron CRMs

Element	Q1	Q2	ORS	JSS 001-6			JSS 003-6			
				Certified value [mg/kg]	Uncertainty	Found [mg/kg]	Certified value [mg/kg]	Uncertainty	Found [mg/kg]	50 ppb spike recovery %
P	31	47	O ₂	0.5*		0.458	3.5	0.7	3.170	103
S	32	48	O ₂	1.5	0.3	1.512	1.3	0.5	1.287	92
Mn	55	55	He	0.03*		0.036	3.2	0.2	3.432	101

*non-certified value

Direct measurement of trace rare earth elements in high purity REE oxides

Kazumi Nakano
Agilent Technologies, Japan

Keywords

Rare Earth Elements, REE, rare earth oxide, REO, samarium oxide, gadolinium oxide, oxygen mass-shift, ammonia on-mass

Introduction

The rare earth elements (REEs) are widely used in advanced technologies including high-power permanent magnets, lasers, phosphors used in fluorescent lamps, radar screens and plasma displays. REEs are also used in petroleum refining, automobile catalytic converters and batteries, and in high-technology glasses. It is clear from these examples that REEs play a key role in many types of materials used in high-technology industries. However, the presence of other REEs as contaminants in a purified single-element REE material often impacts the functionality of the final product, so impurities in the REE oxide raw material must be carefully controlled.

ICP-MS is the most commonly used atomic spectrometry technique for the measurement of trace REEs due to its simple REE spectra — particularly when compared to emission techniques. The measurement of mid- and high-mass REEs in a low-mass REE matrix is, however, very challenging for ICP-MS because REEs have among the highest metal-oxide (M-O) bond strengths of any element, and the oxide ions of the low mass REE overlap the preferred isotopes of the mid-mass and high-mass REEs. Table 1 shows the interferences observed in the analysis of trace REEs in high-purity samarium (Sm) oxide and gadolinium (Gd) oxide.

Separation of the trace REE analytes from the REE matrix can be performed utilizing a chelating resin, but this technique is time-consuming and customization is needed according to the analyte and matrix element. The direct analysis of trace REEs in a variety of high-purity REE matrices is therefore desired. In this work, an Agilent 8800 Triple Quadrupole ICP-MS was used for the direct analysis of trace REE in two high-purity REE materials: Sm₂O₃ and Gd₂O₃. Operating the ICP-QQQ in MS/MS mode effectively removes the challenging interferences, enabling the determination of REE impurities at trace levels in these two materials.

Experimental

Instrumentation: Agilent 8800 #100.

Plasma conditions: Preset plasma/General purpose.

Ion lens tune: Soft extraction tune: Extract 1 = 0 V, Extract 2 = -180 V.

Acquisition parameters: Three cell modes were used with MS/MS acquisition: No gas, O₂ mass-shift, and NH₃ on-mass mode. In MS/MS O₂ mass-shift mode, the REEs were determined as their oxide ions. REE ions react efficiently with the O₂ cell gas and are converted to the oxide ion REE-O⁺. For example, in the measurement of ¹⁵³Eu⁺, Q1 is set to *m/z* 153 (¹⁵³Eu⁺) and Q2 is set to *m/z* 169 (¹⁵³Eu¹⁶O⁺). Cell tuning parameters are summarized in Table 2.

Table 1. Preferred isotope for ICP-MS analysis of each REE, and the potential interferences caused by Sm₂O₃ and Gd₂O₃ matrices

Element	La	Ce	Pr	Nd	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu
Mass	139	140	141	146	147	153	157	159	163	165	166	169	172	175
Gd ₂ O ₃							N/A	GdH ⁺					GdO ⁺	GdOH ⁺
Sm ₂ O ₃					N/A	SmH ⁺			SmO ⁺	SmO ⁺	SmO ⁺	SmOH ⁺		

Table 2. CRC tuning parameters

Cell mode	Unit	No gas	O ₂	*NH ₃
Scan mode		MS/MS		
Cell gas		N/A	O ₂	NH ₃
Cell gas flow rate	mL/min	N/A	0.35	9.0
Octopole bias	V	-8	-5	-18
KED	V	5	-8	-8
Cell exit	V	-80	-90	-110
Deflect lens	V	20	10	-3
Plate	V	-80	-90	-110

*10% NH₃ balanced in Ar

Results and discussion

Two high purity REE oxide materials Gd₂O₃ (5N) and Sm₂O₃ (4N8) were gently dissolved in semiconductor grade HNO₃ and diluted to a concentration of 1 ppm (as the REE). The other (trace) REEs were measured in each matrix solution using the three cell modes. The results are given in Figure 1 and Figure 2. As expected, analysis of the 1 ppm Gd solution in no gas mode gave positive errors on some elements due to interferences from Gd polyatomic ions: GdH⁺ interferes with ¹⁵⁹Tb⁺, GdO⁺ interferes with ¹⁷²Yb⁺ and GdOH⁺ interferes with ¹⁷⁵Lu⁺.

Preliminary studies showed that NH₃ cell gas reacts with many of the polyatomic ions that interfere with the REE. However, NH₃ also reacts quickly with some of the REE ions, leading to reduced sensitivity of < 1 cps/ppt [1],




Figure 1. Measured concentration of REE impurities in 1 ppm Gd solution. Gd based interferences are observed on Tb, Yb and Lu. Only the elements in red were measured in NH₃ on-mass mode.




Figure 2. Measured concentration of REE impurities in 1 ppm Sm solution. Sm based interferences are observed on Eu, Dy, Ho, Er, Tm and Yb. Only the elements in red were measured in NH₃ on-mass mode.

so this mode is only suitable for the measurement of the less reactive analytes: Pr, Eu, Dy, Ho, Er, Tm and Yb. For these elements, NH₃ on-mass mode gave excellent results, including for Yb in the Gd matrix, where the measured Yb background concentration was reduced by four orders of magnitude (Figure 1) indicating effective removal of the GdO⁺ overlap. Background signals for Dy, Ho, Er and Tm in the Sm matrix were also dramatically improved (Figure 2).

For the REEs that react with NH₃ (La, Ce, Nd, Sm, Gd, Tb and Lu), O₂ mass-shift mode and measurement of the target analyte as its REE-O⁺ ion is the preferred approach. Most REEs are effectively converted to the oxide ion via reaction with O₂ cell gas [1], and this mode was applied to the measurement of Lu in the Gd matrix, avoiding the GdOH⁺ interference on the Lu⁺ isotope and giving a good improvement in the background signal. Compared to no gas mode, O₂ mass-shift mode also gave a good improvement in the background signals for Dy, Ho, Er, Tm and Yb in the Sm matrix, but for all these analytes the backgrounds in NH₃ mode were lower still.

Reference

1. Direct measurement of trace rare earth elements (REEs) in high-purity REE oxide using the Agilent 8800 Triple Quadrupole ICP-MS with MS/MS mode, Agilent application note, 5991-0892EN.

The benefits of improved abundance sensitivity with MS/MS for trace elemental analysis of high purity metals

Fred Fryer
Agilent Technologies, Australia

Keywords

quadrupole, resolution, hyperbolic, abundance sensitivity, copper, high purity zinc

Introduction

The use of a quadrupole mass filter for the separation of compounds in mass spectrometry is well established. Initially used for organic mass spectrometry and residual gas analysis, the quadrupole spectrometer was adopted for the earliest ICP-MS systems, and has remained the default choice throughout the history of ICP-MS. However, the performance characteristics of the quadrupole mass filter do impose several limitations on quadrupole ICP-MS (ICP-QMS).

The resolution (R) of a mass filter (meaning its ability to separate adjacent masses) is defined as $M/\Delta M$, the mass of the target peak/the mass difference to nearest adjacent peak that can be distinguished (separated). However, for practical specifications, the resolution is often simply quoted as the width of the peak at a given peak height. The quadrupole mass filter used in an ICP-QMS instrument is typically operated with a nominal peak width of about 0.75 amu at 10% peak height, illustrated in Figure 1.




Figure 1. Illustration of resolution calculation for a mass spectrometer

For two peaks within the normal signal range of the instrument, this allows the complete, baseline separation of masses 1 amu apart, within the elemental mass range from Li (7 amu) to U (238 amu) and beyond. Higher resolution of 0.4 amu peak width is possible by adjusting the quadrupole voltages, but the signal is reduced (less ion transmission) due to rejection of a higher proportion of the ions that are nominally “on-mass”. Typically the signal loss at higher resolution is around 10-50%, depending on the design and operating characteristics of the quadrupole.

Both the efficiency of transmission of ions at the set-mass (i.e., the sensitivity) and the rejection of ions at other masses (i.e., the resolution of adjacent peaks) are affected by the shape of the field within the quadrupole, and the frequency of the alternating RF fields.

A hyperbolic field (generated by rods with a hyperbolic profile) alternating at high frequency gives more effective filtering of the ion beam than a lower frequency field generated by round quadrupole rods. The practical benefit of hyperbolic rods and high frequency RF voltage is therefore better ion transmission at higher resolution. Some of the many real-world applications where the combination of high sensitivity and good peak separation is required for adjacent low/high concentration elements measured by ICP-MS are shown in Table 1.

Table 1. Elements that require improved separation of adjacent peaks for low level analysis

Low concentration	High signal	Example matrix
^{31}P	$^{16}\text{O}_2, ^{32}\text{S}$	Soil, plants, biological
^{55}Mn	$^{56}\text{Fe}, ^{40}\text{Ar}^{16}\text{O}$	Blood, iron and steel, soil
$^{63}\text{Cu}, ^{65}\text{Cu}$	$^{64}\text{Zn}, ^{66}\text{Zn}$	Metal refining
^{11}B	^{12}C	Soils, solvents, petrochem
^{13}C	^{14}N	Laser imaging of biological

Due to the ion transmission characteristics of a quadrupole, the peak that is generated from the ion signals at each mass forms a non-symmetric Gaussian distribution with a negative skew; i.e., the peak has a longer tail on the leading edge (low mass side) than the trailing edge (high mass side). These "tails" may extend significantly beyond the limits of the nominal 0.75 amu peak width, but since they are at intensities far below 10% of the peak height, they cannot be measured using the simple resolution figure quoted above. The contribution that a peak at mass M makes to its neighbors at M-1 amu and M+1 amu can be quantified, however, and this figure is referred to as the abundance sensitivity (AS) of the quadrupole, illustrated in Figure 2.

For a good quadrupole mass spectrometer in ICP-QMS, the AS would typically be of the order of 10^{-7} , meaning that for an on-mass signal of 10^7 counts, there is a contribution of one count at the adjacent mass (M+/-1 amu).

In applications where the trace analyte must be separated from a very intense matrix peak at the adjacent M+1 mass, such as the examples shown in Table 1, the matrix peak may be at an intensity greater than 10^9 or 10^{10} , and the AS of a quadrupole mass spectrometer is insufficient for accurate trace measurement of an adjacent overlapped analyte at low/sub ppb levels.




Figure 2. Illustration of abundance sensitivity calculation for a mass spectrometer

The problem of adjacent mass overlaps now has an elegant solution in the Agilent 8800 Triple Quadrupole ICP-MS (ICP-QQQ). The 8800 ICP-QQQ uses a tandem mass spectrometer configuration with two quadrupole mass filters (Q1 and Q2) separated by a collision/reaction cell. In MS/MS mode, both quadrupoles are operated as unit mass filters, so the overall AS of the instrument is the product of the Q1 AS x the Q2 AS. With two research-grade, high frequency, hyperbolic quadrupoles, each operating with AS of 10^{-7} , the combined AS of the 8800 ICP-QQQ is theoretically 10^{-14} , although this cannot be verified experimentally as the magnitude of the signal difference exceeds the dynamic range of the detector.

Experimental

Trace copper in high purity zinc

Major uses of Zn include galvanized coating to protect steel, die castings, and solder. Impurities in the metal cause Zn plating to lift, die casts to crack, or solder to 'de-wet', hence high purity zinc (>99.995 %) is a preferred commodity. Common impurities are Cu, Au and Sb, but may also include Cd, Al, Fe, Ag, Bi, As, In, Ni, P, and S.

Instrumentation: Agilent 8800 #100.

Plasma conditions: Preset plasma/Low matrix.

Ion lens tune: Auto tune was used for optimization.

CRC conditions: Helium cell gas at 4.8 mL/min with KED of 4 V.

Results and discussion

Analysis was performed on high purity Zn, dissolved to give a 0.1% (1000 mg/L) Zn solution in a final acid concentration of 2% HNO_3 . The sensitivity of the 8800 ICP-QQQ was reduced to bring the signal for the major Zn isotopes (64, 66) within the detector's upper limit of dynamic range ($\sim 10^{10}$ cps). The intense Zn signals were measured automatically in analog detector mode, while the Cu isotopes were measured in pulse mode. It can clearly be seen in Figure 3 that the intense Zn peaks at m/z 64 and 66 made no contribution to the signal at the two adjacent trace Cu isotopes at m/z 63 and 65. The Cu isotope ratio matched the theoretical abundances ($^{63}\text{Cu}/^{65}\text{Cu}$ natural ratio of 69.17/30.83), at 1 $\mu\text{g/L}$ concentration. If there was a contribution from an adjacent Zn mass then the isotope ratio would be biased.

From the ^{63}Cu calibration (Figure 4), the BEC and DL measured for Cu in the 0.1% Zn matrix were 1.7 ppb and 0.01 ppb respectively, indicating a low and stable background signal. The sensitivity of Cu was 7700 cps/ppb in 1000 $\mu\text{g/L}$ Zn, under the “de-tuned” conditions used to bring the Zn peaks within the detector range. This represents about a 2x reduction in the signal that would be obtained under normal tuning conditions for this type of matrix, if measurement of the matrix element peaks was not required.




Figure 3. 1000 mg/L Zn spiked with 1 $\mu\text{g/L}$ Cu. (1 ppb: 63/total = 0.667, 65/total = 0.333, no mass bias correction performed. If there was some contribution from Zn it would influence m/z 63 differently from m/z 65, because of the different abundance of the Zn isotopes).




Figure 4. Standard addition calibration of ^{63}Cu in 1000 mg/L Zn

Ultra trace copper analysis in a semiconductor grade organometallic titanium complex

Akio Hashizume, Toshiya Shingen, ADEKA Corporation, Japan
Katsuo Mizobuchi, Agilent Technologies, Japan

Keywords

semiconductor, organometallic, copper, titanium, ammonia mass-shift

Introduction

Most quadrupole ICP-MS (ICP-QMS) instruments use CRC technology to resolve spectroscopic interferences. Helium collision mode is widely accepted due to its versatility and ease of use for multi-element analysis of complex and variable samples. While the performance achievable with He mode meets the requirements for most applications, there are some applications, for example impurity analysis of semiconductor materials, that require improved interference removal capability. For these applications, a reactive cell gas (reaction mode) may be used, but the use of highly reactive cell gases in quadrupole ICP-MS is prone to unexpected interferences and overlaps, especially when the matrix is complex, or other analytes are present at varying concentrations. The new Agilent 8800 Triple Quadrupole ICP-MS (ICP-QQQ) eliminates the variability associated with reactive cell gases in ICP-QMS, by using the first quadrupole (Q1) to control the ions that enter the CRC. This ensures that the reactions are predictable and the product ion spectrum is simple and consistent.

This report describes the measurement of trace Cu in a semiconductor grade organometallic Ti complex used in advanced semiconductor processing. It is a challenging application for quadrupole ICP-MS since both isotopes of copper, ^{63}Cu and ^{65}Cu , suffer interference from TiO and TiOH ions, and the use of reactive cell gases to avoid the overlap leads to a very complex product ion spectrum, particularly for organic samples. We demonstrate that the Agilent 8800 ICP-QQQ, operating in MS/MS mass-shift mode using ammonia as a reaction gas, was able to separate Cu^+ from the Ti-based interferences and measure Cu at low ppt levels in a matrix of 500 ppm Ti. Results were also acquired using MS/MS He collision mode, for comparison.

Experimental

Instrumentation: Agilent 8800 #200 with narrow injector (id = 1.5 mm) torch (G3280-80080) used for organic solvent analysis. A low flow PFA nebulizer (G3285-80002) was used in self-aspiration mode. An option gas flow of 20% O_2 balanced in Ar was added to the carrier gas via the standard option-gas line to prevent carbon build up on the interface cones.

Operating conditions: Table 1 summarizes plasma, ion lens and cell tuning conditions.

Acquisition conditions: MS/MS mode was used; cell gas was either NH_3 or He.

Sample and sample preparation: Semiconductor grade organometallic Ti complex (ADEKA Corp., Japan) was diluted with high purity IPA (Tokuyama Corp., Japan) to 500 ppm Ti solution. A spiked standard was prepared from the multi-element standard, xstc-331, purchased from SPEX CertiPrep Ltd. (UK).

Table 1. Experimental conditions

		Units	He collision cell mode	NH_3 reaction cell mode
Cell conditions	Cell gas		He	NH_3 (10% NH_3 in He)
	Cell gas flow rate	mL/min	8.0	6.5
	Octopole bias	V	-18	-18
	KED	V	4	-10
	Cell exit	V	-100	-70
	Deflect	V	-3	-12
	Plate	V	-70	-60
Plasma conditions	RF	W	1600	
	SD	mm	12.0	
	CRGS	L/min	0.70	
	MUGS		0.20	
	Opt gas flow rate		0.20	
Ion lens	Extract 1	V	-60	
	Extract 2	V	-10	

Results and discussion

He collision mode

The He cell gas flow rate was optimized for the lowest BEC of Cu in a 500 ppm Ti solution. As the BEC for ^{63}Cu was lower than the BEC for ^{65}Cu due to the higher abundance of the 63 isotope, and the more significant interference from TiO^+ at m/z 65, Cu was determined on-mass at m/z 63. In MS/MS mode, this is achieved by the acquisition conditions: Q1 = 63; Q2 = 63 (63, 63).

Two solutions were analyzed: 500 ppm Ti solution and 500 ppm Ti + 1 ppb Cu spike. Figure 1 (left) shows the signal at m/z 63 obtained from the analysis of the two solutions, plotted as a function of He flow rate. The BEC calculated from these signals is also given in the figure. It shows that the lowest Cu BEC in He mode was 46 ppt, achieved at a flow rate of 8.0 mL/min He.

NH_3 reaction cell mode

Cu^+ reacts efficiently with NH_3 to form NH_3 cluster ions with the general form $\text{Cu}(\text{NH}_3)_n^+$. TiO^+ does not follow the same reaction pathway as Cu^+ , so the Cu product ion can be measured free from Ti overlap. Based on a

preliminary study, one of the intense product ions, $\text{Cu}(\text{NH}_3)_2^+$, was selected to measure Cu separated from the original TiO^+ interference. A mass pair of Q1 = 63, Q2 = 97 was used with NH_3 as the reaction gas. Figure 1 (right) shows the result. A BEC of 11 ppt for Cu in 500 ppm Ti solution was achieved in NH_3 mode (10% NH_3/He mixed gas), at a flow rate of 6.5 mL/min NH_3 .

Conclusions

Table 2 summarizes the analytical performance achieved by the 8800 ICP-QQQ operating in MS/MS mode with He collision and NH_3 reaction gas. As can be seen, NH_3 reaction mode is more effective than He collision mode for the removal of the TiO^+ interference on Cu. The BEC obtained for Cu in a Ti matrix by NH_3 reaction mode is four times lower than He mode, with seven times higher sensitivity.

Table 2. Summary of Cu measurement in Ti matrix

	Flow rate (mL/min)	BEC (ppt) of Cu in 500 ppm Ti	Sensitivity (cps/ppb)
He collision mode	8.0	45.5	810
NH_3 reaction mode	6.5	10.9	5900




Figure 1. (Left) Cu signal (mass pair 63, 63) vs. He cell gas flow rate, for 500 ppm Ti matrix unspiked and with 1 ppb Cu spike, and calculated BEC.

(Right) Cu signal (mass pair 63, 97) vs. NH_3 cell gas flow rate, for 500 ppm Ti matrix unspiked and with 1 ppb Cu spike, and calculated BEC.

Removal of MH^+ interferences in refined REE material analysis

Naoki Sugiyama

Agilent Technologies, Japan

Keywords

Rare Earth Elements, REE, geochemistry, mining, material science, lanthanum, barium, cerium, method of standard additions, MSA, oxygen mass-shift

Introduction

The measurement of Rare Earth Elements (REEs) is of great importance in geochemistry, mining and material science. Manufacturers of high purity REE materials need to quantify metal impurities, including trace levels of the other REEs, in the refined, single element REE matrix. ICP-MS is the technique of choice for the measurement of REEs, but most of the REE isotopes suffer from interference by polyatomic species (predominantly hydride ions, MH^+ and oxide ions, MO^+) derived from other, lower-mass REE elements. While MH^+ interferences are lower in intensity than MO^+ interferences, they present a more challenging problem for REEs that have no isotope free from interference. For example $^{139}\text{La}^+$ is interfered by $^{138}\text{BaH}^+$ and $^{140}\text{Ce}^+$ by $^{139}\text{LaH}^+$. These interferences are too close in mass to be resolved by high-resolution (HR-)ICP-MS [1]. In this paper, we describe the removal of the MH^+ interferences using an Agilent 8800 ICP-QQQ in MS/MS mass-shift mode, with oxygen as the reaction gas.

Experimental

Instrumentation: Agilent 8800 #100. The standard glass nebulizer was replaced with a C-flow nebulizer (G3285-80000) for optimal washout between the high matrix samples.

Plasma conditions: Preset plasma/General purpose.

Ion lens tune: Soft extraction tune: Extract 1 = 0 V, Extract 2 = -180 V.

CRC conditions: O_2 gas at 0.3 mL/min, Octopole bias = -5 V, KED = -5 V.

Acquisition parameters: MS/MS mode with O_2 mass-shift method.

Figure 1 illustrates the mechanism of MS/MS O_2 mass-shift mode used for measuring Ce in a La matrix sample. The major isotope of Ce at m/z 140 suffers an interference from $^{139}\text{LaH}^+$. Q1 is set to m/z 140, allowing only the analyte ion $^{140}\text{Ce}^+$ and any other ions at m/z 140 to pass through to the cell. All other ions not at m/z 140 are rejected. In the cell, Ce reacts with oxygen to form CeO^+ at m/z 156. Q2 is set to m/z 156, allowing CeO^+ to pass to the detector. Since $^{139}\text{LaH}^+$ does not react with O_2 to form $^{139}\text{LaOH}^+$, it remains as LaH^+ at m/z 140 and is rejected by Q2. The same principle is used for the separation of $^{139}\text{La}^+$ from $^{138}\text{BaH}^+$ in a Ba matrix.




Figure 1. MS/MS mass-shift method with O_2 reaction gas; used for the measurement of Ce, as CeO at m/z 156, in a La matrix.

Results and discussion

Using the Method of Standard Addition (MSA), the BECs and DLs of La in a matrix of 50 ppm Ba, and Ce in a matrix of 50 ppm La were determined. Data was acquired using MS/MS mode with O₂ mass-shift, and also using Single Quad (SQ) mode with O₂ reaction gas to emulate conventional quadrupole ICP-MS (ICP-QMS) for comparison.

As shown in Figures 2A and 2B, SQ mode with O₂ reaction gas suffers from interferences that prevent the measurement of La in the Ba matrix and Ce in the La matrix, respectively. In contrast, the calibration plots shown in Figures 2C and 2D demonstrate that MS/MS mode with O₂ mass-shift can successfully remove the matrix overlaps to permit the trace quantitation of La in a Ba matrix and Ce in a La matrix. The BECs and DLs achieved were 8.5 ppt and 2.5 ppt respectively for La in a 50 ppm Ba solution, and 10.6 ppt and 0.8 ppt respectively for Ce in a 50 ppm La solution.

Investigation of unexpected product ion observed at *m/z* 156 in the 50 ppm La matrix

The background signals that contributed to the poor result obtained for Ce in the La matrix using SQ mode with O₂ reaction gas (Figure 2B) were investigated by carrying out a precursor ion scan for product ion mass 156. The precursor ion scan capability of the 8800 ICP-QQQ provides a uniquely powerful approach to identifying the source of potential polyatomic and reaction product interferences. Oxygen cell gas was introduced into the cell and a precursor ion spectrum was obtained by scanning Q1 from 2 to 260 amu (Figure 3) with Q2 fixed at mass 156. From the spectrum, we can identify which precursor ions react with O₂ to produce product ions at mass 156, overlapping ¹⁴⁰CeO⁺ in SQ mode.




Figure 2. Top: Calibration plots up to 0.2 ppb for La in 50 ppm Ba matrix (A) and Ce in 50 ppm La matrix (B), acquired in SQ mode with oxygen reaction gas (emulating conventional quadrupole ICP-MS). Bottom: Calibration plots up to 0.2 ppb for La in 50 ppm Ba matrix (C) and Ce in 50 ppm La matrix (D) acquired in MS/MS mode with oxygen mass-shift.

Figure 3 shows the precursor ion scan spectrum for product ion mass 156 for the 50 ppm La matrix, with intense peaks at m/z 139 ($^{139}\text{La}^+$) and 156 ($^{139}\text{La}^{16}\text{OH}^+$). In SQ mode, as with conventional ICP-QMS, these ions all enter the cell, and with Q2 set to 156 amu, the $^{139}\text{La}^{16}\text{OH}^+$ polyatomic ions contribute to the signal measured at m/z 156 (^{140}Ce measured as analyte product ion $^{140}\text{CeO}^+$). These unwanted precursor ions cannot be rejected by a CRC operating as a bandpass filter in ICP-QMS, as they are too close in mass to the target analyte precursor ion. Only by using MS/MS mode on the 8800 ICP-QQQ, where Q1 operates as a unit mass filter, can non-target masses (like $^{139}\text{La}^{16}\text{OH}^+$ in this example) be prevented from entering the cell.

Reference

1. Sabine Becker and Hans Joachim Dietze, *Journal of Analytical Atomic Spectrometry*, 1997, vol.12, p881.

More information

Removal of hydride ion interferences (MH^+) on Rare Earth Elements using the Agilent 8800 Triple Quadrupole ICP-MS, Agilent application note, 5991-1481EN.




Figure 3. Precursor ion scan from 2-260 amu for product ion mass 156, in a 50 ppm La matrix. Six peaks are seen at $m/z = 138, 139, 140, 155, 156$ and 157 , with the intense peaks at m/z 139 and m/z 156 being due to $^{139}\text{La}^+$ and $^{139}\text{La}^{16}\text{OH}^+$ respectively.

- > Return to table of contents
- > Search entire document

Geology

Lead isotope analysis: Removal of ^{204}Hg isobaric interference on ^{204}Pb using ICP-QQQ MS/MS reaction cell 64

Nuclear

Analysis of radioactive iodine-129 using MS/MS with O_2 reaction mode 67

Feasibility study on the analysis of radioisotopes: Sr-90 and Cs-137 70

Lead isotope analysis: Removal of ^{204}Hg isobaric interference on ^{204}Pb using ICP-QQQ MS/MS reaction cell

Glenn Woods

Agilent Technologies (UK) Ltd.

Keywords


lead, isotope, ratio, geochronology, dating, mercury, artifacts, precious metals, food, ammonia, on-mass

Introduction

Lead isotope ratio analysis is important as it is used for Pb-Pb dating in geochronology, and to trace the origin of artifacts, precious metals and even foodstuffs. The natural isotopic pattern of lead varies more than any other element in the periodic table, because three of its isotopes are formed from the radioactive decay of uranium ($^{235}\text{U} \rightarrow ^{207}\text{Pb}$; $^{238}\text{U} \rightarrow ^{206}\text{Pb}$) and thorium ($^{232}\text{Th} \rightarrow ^{208}\text{Pb}$). The Pb isotopic pattern can therefore vary depending upon the geology of the rocks and minerals from which the lead was extracted, and the age of the material. In geochronology, the constant rate of U/Th decay allows the Pb/Pb, U/Pb and Th/Pb ratios to be used to date the age of rocks using a so-called geological clock.

When Pb ratios are measured, it is often necessary to correct for the lead naturally present in the sample, and the only non-radiogenic isotope of Pb (^{204}Pb ; natural or common lead), is used for this purpose. For Pb-Pb dating, ^{204}Pb is the reference isotope against which the radiogenic isotopes are compared ($^{206}\text{Pb}/^{204}\text{Pb}$; $^{207}\text{Pb}/^{204}\text{Pb}$). Unfortunately ^{204}Pb is directly overlapped by an isotope of Hg (^{204}Hg), which makes accurate measurement of ^{204}Pb impossible by ICP-MS. Mass resolution of ^{204}Pb from ^{204}Hg is far beyond the capability of any commercial high-resolution (HR-) ICP-MS system, and until recently there has been no reliable chemical means to remove the Hg interference, so mathematical correction has been employed, which introduces error.

Mercury does however undergo a gas-phase charge-transfer reaction with ammonia gas (NH_3), a reaction that can be utilized in the collision/reaction cell of a suitably equipped ICP-MS as follows:



This reaction offers the potential to remove the ^{204}Hg interference from ^{204}Pb , and could be applied to either solution or laser-based ICP-MS analysis.

Experimental

Instrumentation: Agilent 8800 #100.

Plasma conditions: Preset plasma/General purpose.

Ion lens tune: Soft extraction tune: Extract 1 = 0 V, Extract 2 = -170 V.

CRC conditions: NH_3 gas (10% in He) at 1.7 mL/min, Octopole bias = -8 V, KED = -8 V.

Acquisition parameters: Three acquisition modes were compared:

- **No gas:** No reaction cell gas; Single Quad (SQ) mode with Q1 operating as an ion guide
- **NH_3 bandpass:** Ammonia reaction gas; SQ mode with Q1 operating as a bandpass filter
- **NH_3 MS/MS:** Ammonia reaction gas; MS/MS mode with Q1 operating as a mass filter at unit mass resolution

Results and discussion

Removal of $^{204}\text{Hg}^+$ interference on $^{204}\text{Pb}^+$

A preliminary study showed that Pb is almost unreactive with NH_3 cell gas (<0.5% loss of Pb signal) indicating that on-mass sensitivity for Pb should be maintained. On-mass measurement of Pb in NH_3 cell gas mode was therefore investigated in the presence of Hg at 10 ppb. Figure 1 displays the spectra obtained in no gas (left) and NH_3 cell gas (right) modes. The ^{204}Hg interference on ^{204}Pb can be clearly seen in the no gas spectrum, while it has been completely removed under NH_3 reaction mode with MS/MS. A perfect isotopic pattern match was confirmed for Pb in NH_3 mode.

Effectiveness of MS/MS

The NH_3 reaction that removes the ^{204}Hg interference would also work in the reaction cell of a quadrupole ICP-MS (ICP-QMS), but ammonia is a highly reactive gas and can produce many adduct cluster ions, for example from Rare Earth Elements (REEs), see Table 1. The complex matrix composition of many natural samples means that the results obtained with NH_3 cell gas in ICP-QMS are often extremely unreliable. With the 8800 ICP-QQQ, MS/MS mode allows all the co-existing matrix

Table 1. Some possible Rare Earth Element cluster ions that can form in the CRC of an ICP-QMS when using NH_3 reaction gas – the list is by no means exhaustive.

Mass	Potential Cluster Ions of REE
204	$\text{Eu}(\text{NH}_3)_3$; $\text{Yb}(\text{NH}_3)_2$; $\text{Ce}(\text{NH}_3)_4$
205	$\text{Yb}(\text{NH}_3)_2$; $\text{Gd}(\text{NH}_3)_3$
206	$\text{Yb}(\text{NH}_3)_2$; $\text{Lu}(\text{NH}_3)_2$; $\text{La}(\text{NH}_3)_4$; $\text{Ce}(\text{NH}_3)_4$; $\text{Gd}(\text{NH}_3)_3$
207	$\text{La}(\text{NH}_3)_4$; $\text{Yb}(\text{NH}_3)_2$; $\text{Gd}(\text{NH}_3)_3$
208	$\text{Ce}(\text{NH}_3)_4$; $\text{Gd}(\text{NH}_3)_3$; $\text{Tb}(\text{NH}_3)_2$; $\text{Yb}(\text{NH}_3)_2$; $\text{Gd}(\text{NH}_3)_3$

elements to be rejected by Q1, so only the target ions (^{204}Pb and ^{204}Hg) enter the CRC. The NH_3 reactions are therefore controlled and consistent, and no overlapping reaction product ions are formed from other elements in the sample.

To check the formation of cluster ions, the ICP-QQQ was operated with NH_3 cell gas; “Single Quad bandpass” and MS/MS modes were compared for the measurement of a 50 ppb REE mix. Figures 2a and 2b display the spectra obtained using bandpass and MS/MS conditions, respectively.




Figure 1. Standard solution (1 ppb each of Pb and Tl) spiked with 10 ppb Hg without cell gas (left) and with NH_3 (right) using MS/MS mode; Note the ^{204}Hg interference on ^{204}Pb in no gas mode.




Figure 2. Cluster ion formation for 50 ppb REE standard in ammonia mode. Figure 2a (left): REE cluster ion formation using ammonia cell gas in bandpass mode; REE's are allowed into the cell if Q1 is operated as a bandpass filter. The REE cluster ions can be seen at all masses including those for Hg, Tl, Pb and Bi. Figure 2b (right): The identical sample under the same ammonia conditions but this time with Q1 operated at unit mass resolution (MS/MS mode). The REE's are removed from the ion beam before they can enter the cell and form reaction by-products.

$^{204}\text{Pb}/^{208}\text{Pb}$ isotope ratio analysis in presence of Hg

To check the effectiveness of the ^{204}Hg removal, the $^{204}\text{Pb}/^{208}\text{Pb}$ ratio was measured in a 1 ppb lead solution spiked with increasing Hg concentration. Table 2 displays the measured Pb ratio results (without any mass bias correction), showing that the Pb isotope ratio remained constant, regardless of the Hg content.

Table 2. Uncorrected isotopic ratios measured in 1ppb Pb solutions containing mercury at varying concentrations. The lead isotopic ratio 204/208 is not influenced by the presence of Hg.

	^{204}Pb	^{208}Pb	IR (204/208)
Sample	CPS		
Theoretical	NA	NA	0.02671
Pb	3518.5	136124.8	0.02585
Pb Hg 5 ppb	3510.0	139585.9	0.02515
Pb Hg 10 ppb	3439.2	132796.4	0.02590
Pb Hg 20 ppb	3464.8	134417.7	0.02578

Conclusions

With the successful removal of the ^{204}Hg interference on the natural ^{204}Pb isotope, ICP-QQQ displays great promise for Pb/Pb and U/Pb dating and for other applications where accurate measurement of ^{204}Pb is required.

Analysis of radioactive iodine-129 using MS/MS with O₂ reaction mode

Yasuyuki Shikamori, Kazumi Nakano and Naoki Sugiyama
Agilent Technologies, Japan

Keywords

radionuclide, iodine, ¹²⁹I, environmental, nuclear, xenon, NIST 3231 Level I and II, abundance sensitivity, oxygen on-mass

Introduction

Iodine-129 is a long-lived radionuclide (half-life of 15.7 My) which has been released into the environment as a result of human activities such as nuclear weapons testing, accidents at nuclear power plants and especially by emissions from spent nuclear fuel reprocessing plants. The determination of iodine-129 in environmental samples is very difficult by ICP-MS due to the element's relatively low sensitivity, the very low concentrations at which ¹²⁹I must be determined, relative to potentially high levels of ¹²⁷I, the high background caused by ¹²⁹Xe impurities in the argon plasma gas, and possible polyatomic interference from ¹²⁷IH₂⁺. Iodine analysis is further complicated by the fact that it is rapidly volatilized from samples prepared using the acid digestions that are normal for ICP-MS analysis, so an alternative, alkaline sample solubilization and stabilization strategy is required. The isobaric interference from ¹²⁹Xe⁺ can be significantly reduced using ICP-QMS with an Octopole Reaction Cell operated in O₂ reaction mode, resulting in a measured ratio for ¹²⁹I/¹²⁷I of 10⁻⁷ in NIST 3231 SRM Level I^[1]. However, the problem of potential overlap due to tailing from ¹²⁷I and ¹²⁷IH remains, as the relative abundance of the ¹²⁹I to ¹²⁷I will typically exceed 10⁻⁷, which is of the same order as the abundance sensitivity (ability to separate adjacent peaks) of quadrupole ICP-MS (ICP-QMS). In order to overcome these challenges, ICP-QQQ operating in MS/MS mode with O₂ reaction gas was applied to determine ultratrace levels of iodine-129 in aqueous samples.

Experimental

Instrumentation: Agilent 8800 #100.

Plasma conditions: Preset plasma/Low matrix.

Ion lens tune: Soft extraction tune: Extract 1 = 0 V, Extract 2 = -190 V.

CRC conditions: O₂ gas at 0.8 mL/min, Octopole bias = -18 V and KED = -1.5V. MS/MS O₂ on-mass mode was applied to measure iodine-127 and iodine-129 (Q1 = Q2 = 127 for iodine-127, Q1 = Q2 = 129 for iodine-129).

Reference materials and calibration standards:

Calibration standards were prepared by diluting ¹²⁹I isotopic standards NIST SRM 3231 Level I and II (NIST, Gaithersburg MD, USA) with 0.5% TMAH in deionized water. The Level I Certified Value for ¹²⁹I/¹²⁷I = 0.981x10⁻⁶ ± 0.012x10⁻⁶, Level II = 0.982x10⁻⁸ ± 0.012x10⁻⁸. These reference materials were used to check the calibration linearity of the iodine isotopes and to validate the isotopic ratio of iodine-129 and iodine-127.

Results and discussion

Optimization of oxygen cell gas flow

The oxygen gas flow rate was optimized by varying the O₂ flow over the full range of the mass flow controller (0–1.12 mL/min), while monitoring the ^{127,129}I signal and blank intensity, as shown in Figure 1. As the flow rate of O₂ increases, the background signal (due to ¹²⁹Xe) at *m/z* = 129 decreases rapidly, and the iodine signal remains high, dramatically improving the DL for ¹²⁹I.




Figure 1. Profile of $^{127}\text{I}^+$, $^{129}\text{Xe}^+$ signals and estimated ^{129}I BEC. For the BEC calculation, the sensitivity of ^{129}I was assumed to be the same as ^{127}I . Scale of O_2 flow: 100% = 1.12 mL/min.

Abundance sensitivity

Scan spectra over the mass range 127 to 129, covering both ^{127}I and ^{129}I , were acquired for the two SRMs, NIST 3231 Level I and II, using the Agilent 8800 ICP-QQQ in MS/MS on-mass mode with O_2 reaction gas. The overlaid spectra are shown in Figure 2. Excellent abundance sensitivity can be seen, with the sides of the intense ($>10^9$ cps) ^{127}I peak reaching baseline with no tailing of $^{127}\text{I}^+$ or $^{127}\text{IH}^+$ on $^{129}\text{I}^+$.

MS/MS, O_2 reaction ($\text{Q1}=\text{Q2}$)




Figure 2. Iodine spectra showing both ^{127}I and ^{129}I acquired using MS/MS on-mass mode with O_2 cell gas. $^{127}\text{IH}^+$ remains to some extent while $^{127}\text{IH}_2^+$ is completely removed, as noted later.

Calibration curves for ^{127}I and ^{129}I

In order to check the linearity of both iodine isotopes, different concentration solutions of NIST 3231 SRM Level I were prepared in 0.5% TMAH and analyzed as calibration standards, as shown in Figure 3. The BECs for ^{127}I and ^{129}I were 2.9 $\mu\text{g/L}$ and 0.04 ng/L respectively, and the detection limits (3σ , $n=10$) were 0.26 $\mu\text{g/L}$ for ^{127}I and 0.07 ng/L for ^{129}I .





Figure 3. Calibration curve for iodine-127 (top) and iodine-129 (bottom) obtained from multiple dilutions of NIST 3231 SRM.

Analysis of NIST 3231 SRM Level I and Level II

The $^{129}\text{I}/^{127}\text{I}$ ratio in 10x diluted NIST 3231 SRM Levels I ($^{129}\text{I}/^{127}\text{I} = 0.981 \times 10^{-6}$) and II ($^{129}\text{I}/^{127}\text{I} = 0.982 \times 10^{-8}$) was measured using ICP-QQQ in MS/MS on-mass mode with O_2 cell gas. The results are summarized in Table 1. After subtracting the ^{129}I blank, the measured $^{129}\text{I}/^{127}\text{I}$ ratio of NIST 3231 SRM Levels I and II corresponded well with the certified values of 0.981×10^{-6} and 0.982×10^{-8} respectively. The good agreement with the certified ratio indicates that the potential interference of $^{127}\text{IH}_2^+$ on $^{129}\text{I}^+$ is completely removed by O_2 reaction with MS/MS mode.

Reference

1. The ultratrace determination of iodine 129 in aqueous samples using the 7700x ICP-MS with oxygen reaction mode, Agilent application note, 5990-8171EN.

More information

1. The ultratrace determination of iodine 129 using the Agilent 8800 Triple Quadrupole ICP MS in MS/MS mode, Agilent application note, 5991-0321EN.

Table 1. Analytical results for NIST 3231 Level I and Level II

Sample name	Dilution factor	Q1=O2=127 CPS	Q1=O2=129 CPS	$^{129}\text{I}/^{127}\text{I}$	$^{129}\text{I}/^{127}\text{I}$ (average n = 5)	RSD (%)
NIST 3231 10^{-6} ($^{129}\text{I}/^{127}\text{I} = 0.981 \times 10^{-6}$)	10	594,277,896	585.6	0.971×10^{-6}	0.981×10^{-6}	0.8
		592,633,576	597.4	0.994×10^{-6}		
		590,000,723	586.5	0.980×10^{-6}		
		593,387,443	588.5	0.978×10^{-6}		
		592,834,056	588.9	0.979×10^{-6}		
NIST 3231 10^{-8} ($^{129}\text{I}/^{127}\text{I} = 0.982 \times 10^{-8}$)	10	608,737,949	15.1	1.12×10^{-8}	1.02×10^{-8}	7.2
		608,536,242	14.8	1.07×10^{-8}		
		602,626,536	14.2	0.979×10^{-8}		
		603,091,763	13.9	0.929×10^{-8}		
		603,250,003	14.5	1.03×10^{-8}		
NIST Blank	10	600,444,851	8.3	—	—	—

Feasibility study on the analysis of radioisotopes: Sr-90 and Cs-137

Yasuyuki Shikamori and Kazumi Nakano
Agilent Technologies, Japan

Keywords

radioisotopes, radioactive, environmental, nuclear, strontium, ⁹⁰Sr, zirconium, cesium, ¹³⁷Cs, barium, abundance sensitivity, oxygen and hydrogen on-mass, nitrous oxide on-mass

Introduction

ICP-MS can be an effective analytical tool for the analysis of long half-life radioisotopes due to its high sensitivity, speed of analysis, low sample consumption, and ease of sample preparation. The challenge for ICP-MS analysis of radioisotopes arises from interferences; not only by polyatomic ions but also atomic isobar ions that cannot be separated even by high-resolution (HR-) ICP-MS.

Trace analysis of the radionuclide ⁹⁰Sr (half-life = 28.74 years) in environmental samples is of great interest. ⁹⁰Sr is a main fission product that may be present in the environment following accidental releases from nuclear power plants. Geiger-Muller (GM) detectors or Liquid Scintillation Counters (LSC) are used to measure ⁹⁰Sr, though both techniques require complex chemical separation prior to analysis, or long integration times. ICP-MS is also used to measure ⁹⁰Sr, especially when a quick turn-around time is desired. However detection limits of quadrupole ICP-MS are compromised by a spectral overlap from ⁹⁰Zr; in common with all direct isobaric interferences, the ⁹⁰Zr overlap is too close in mass to the ⁹⁰Sr to be resolved using sector field HR-ICP-MS, which is limited to a maximum resolution (M/ΔM) of 10,000. This note describes a method for measuring trace ⁹⁰Sr in the presence of ⁹⁰Zr using ICP-QQQ in MS/MS reaction mode. Since it isn't possible to obtain ⁹⁰Sr, a natural isotope of strontium (⁸⁸Sr) was used to estimate the DL for ⁹⁰Sr. A similar approach was applied to ¹³⁷Cs (half-life = 30.0 years).

Experimental

Instrumentation: Agilent 8800 #100.

Plasma conditions: Preset plasma/Low matrix.

Ion lens tune: Soft extraction tune: Extract 1 = 0 V, Extract 2 = -190 V.

CRC and acquisition conditions: The following conditions were used for the analysis of ⁹⁰Sr and ¹³⁷Cs:

- For ⁹⁰Sr: MS/MS on-mass mode (Q1 = Q2 = 90) with O₂ + H₂ cell gas: 1 mL/min of O₂ and 10 mL/min of H₂, Octopole bias = -5 V and KED = -13 V.
- For ¹³⁷Cs: MS/MS on-mass mode (Q1 = Q2 = 137) with N₂O cell gas: 7 mL/min of N₂O (10% N₂O balanced in He, introduced via the 3rd cell gas flow line), Octopole bias = -5 V and KED = -13 V.

Results and discussion

Radioactive Sr-90 (O₂ + H₂ on-mass mode)

Figure 1 shows spectra of a solution containing Sr and Zr (natural isotopes) acquired on the 8800 ICP-QQQ operated in Single Quad mode (Q1 operated as an ion guide to emulate conventional quadrupole ICP-MS) with no cell gas (left), and in MS/MS mode with O₂ + H₂ cell gas (right). As can be seen in the left hand spectrum, the overlap of ⁹⁰Zr⁺ on ⁹⁰Sr⁺ precludes the low-level determination of ⁹⁰Sr by conventional quadrupole ICP-MS. The spectrum on the right indicates that ⁹⁰Sr⁺ could be measured on-mass at m/z = 90 free from interference by ⁹⁰Zr⁺, since Zr⁺ reacts readily with the O₂ + H₂ gas to form ZrO⁺ and ZrO₂⁺. The signal-to-noise ratio for ⁹⁰Sr was improved by six orders of magnitude using MS/MS O₂ + H₂ reaction cell mode.




Figure 1. Mass spectra of a solution containing 20 ppb Sr + 5 ppm Zr: (left) SQ no gas mode and (right) MS/MS $O_2 + H_2$ reaction mode.

Figure 2 is a spectrum of 100 ppm Sr acquired using MS/MS on-mass mode with $O_2 + H_2$ reaction gas. The excellent abundance sensitivity (peak separation) of MS/MS mode can be confirmed. The peak sides reach the baseline with no tailing from the intense peak of the natural isotope of $^{88}Sr^+$. In addition, no $^{88}SrHH^+$ at $m/z = 90$ is formed in cell, even in a solution containing 100 ppm natural Sr.




Figure 2. Spectrum of 100 ppm Sr solution acquired using MS/MS $O_2 + H_2$ on-mass mode




Figure 3. Mass spectra of a solution containing 20 ppb Cs + 5 ppm Ba: (left) SQ no gas and (right) MS/MS N_2O reaction mode.

Radioactive Cs-137 (N_2O on-mass mode)

Figure 3 shows spectra of a solution containing Cs and Ba (natural isotopes) acquired on the 8800 ICP-QQQ operated in Single Quad mode with no gas mode (left), and in MS/MS mode with N_2O cell gas (right). As can be seen in the left hand spectrum, the $^{137}Ba^+$ overlap on $^{137}Cs^+$ is a problem in conventional quadrupole ICP-MS. As with ^{90}Sr , the right hand spectrum shows that $^{137}Cs^+$ could be measured on mass at $m/z = 137$, free from the $^{137}Ba^+$ interference. Ba^+ reacts readily with N_2O to form BaO^+ and $BaOH^+$ while a part of the Cs^+ analyte ion signal remains at its original mass (as shown by the substantial peak for ^{133}Cs in the right-hand spectrum).

Estimated BEC and DL for Sr-90 and Cs-137

The BEC and DL for two radioisotopes, ^{90}Sr and ^{137}Cs , were estimated from these spectra as summarized in Table 1. This feasibility study demonstrates the potential of ICP-QQQ for the measurement of radioisotopes such as ^{90}Sr and ^{137}Cs .

Table 1. Estimated BEC and DL for ^{90}Sr and ^{137}Cs

Radioisotope	BEC (ng/L)	DL (ng/L)
^{90}Sr	0.08	0.23
^{137}Cs	2.9	15

Agilent 8800 ICP-QQQ
Application Handbook

- > [Return to table of contents](#)
- > [Search entire document](#)

This page is intentionally left blank.

- > [Return to table of contents](#)
- > [Search entire document](#)

Clinical

Manganese analysis in whole blood: Expanding the analytical capabilities of ICP-MS	74
Measurement of titanium in clinical samples: Possible application to monitoring patients with joint replacements	77

Life sciences

Simultaneous quantitation of peptides and phosphopeptides by CapLC-ICP-QQQ	79
--	----

Pharma

Measurement of selenium in clinical samples in the presence of Gd-based MRI contrasting agents	81
--	----

Manganese analysis in whole blood: Expanding the analytical capabilities of ICP-MS

Amir Liba
Agilent Technologies, USA

Keywords

manganese, whole blood, iron, abundance sensitivity, helium MS/MS

Introduction

Analysis of clinical samples is challenging due to their complex matrices. While ICP-MS is an immensely powerful multi-element analytical technique, it does suffer from some well-documented spectral interferences. Achieving low detection limits is limited by background signal from low level impurities and the presence of polyatomic interferences, which require the use of CRC technology for their removal. Although the use of CRC-ICP-MS has alleviated many of these analytical challenges, some spectral interferences remain problematic for quadrupole ICP-MS (ICP-QMS). One such interference is the signal overlap on ^{55}Mn due to peak tailing from both ^{54}Fe and ^{56}Fe . Whole blood contains an average of 500 ppm of Fe, and with the level of Mn in whole blood being roughly 10 ppb, analytical results for Mn tend to bias high due to the significant signal tailing and overlap from the adjacent Fe peaks. In this work, we use the superior abundance sensitivity of the 8800 ICP-QQQ to remove any signal overlap from Fe on Mn in whole blood.

Experimental

Instrumentation: Agilent 8800 #100.

Plasma conditions and ion lens tune: Preset plasma/General purpose with soft extraction tune:
Extract 1 = 0 V.

Method: Samples were analyzed using the 8800 ICP-QQQ in both Single Quad (SQ) mode and MS/MS mode. In this study, the mass range of interest (from m/z 50 to 60) was scanned at twenty points per peak in both no gas and helium (He) modes. For the analysis of Mn in blood, MS/MS mode with on-mass measurement (Q2 set to the same mass as Q1) was used, with helium cell gas (typical flow of 4.3 mL/min) to remove polyatomic ion interferences such as FeH^+ and ArOH^+ .

Sample preparation: A 5 ppb solution of Mn was prepared from a stock of 1000 ppm Mn and either analyzed separately or spiked into “base” whole blood (low level Mn). Whole blood was diluted using an alkali matrix containing ammonium hydroxide, EDTA, Triton X-100, and butanol.

Abundance sensitivity

The abundance sensitivity (AS) of a mass spectrometer is the contribution that the signal at mass M makes to the signals at the adjacent masses ($M\pm 1$), expressed as a ratio ($M-1/M$ on the low-mass side and $M+1/M$ on the high-mass side). Simply put, AS is the measure of the “peak tailing” to adjacent masses, which will contribute to a false positive signal, such as that seen on ^{55}Mn (present at trace levels) from the large contribution from ^{54}Fe and ^{56}Fe (which exists at very high concentration) in whole blood. The abundance sensitivity of the best quadrupole ICP-MS systems is of the order of 10^{-7} .

Results and discussion

Abundance sensitivity study in SQ and MS/MS mode

SQ and MS/MS spectra for a 500 ppm Fe solution acquired in no gas mode are shown in Figure 1. The spectrum on the right illustrates the superior peak-to-peak resolution of the 8800 ICP-QQQ operated in MS/MS mode. Although no interference removal for polyatomic ions was employed, the elimination of the contribution to mass 55 from adjacent peaks is clearly evident in MS/MS mode. The “flat-top” peak shapes are the result of the logarithmic scale.

Abundance sensitivity plays an important role when samples contain a large concentration of Fe. Figure 1 looks at the contribution of “peak tailing” on ^{55}Mn due to high levels of Fe. The high concentration of Fe together with the ArN^+ and ArO^+ contribution in no gas mode resulted in the signals at 54 and 56 being over the range of the detector, and so they were automatically skipped. However, the signal contribution from ^{56}Fe on mass 55 is clearly visible in the SQ mode (indicated by the red box) while it is absent in the MS/MS mode.




Figure 1. Comparison of no gas mode spectra for 500 ppm Fe solution, measured in SQ mode (left) and MS/MS mode (right). The signal colored blue was obtained in pulse counting while the green signal was obtained in analog mode. The dotted lines indicate over-range peaks (automatically skipped to protect the EM detector).

Figure 2 shows three spectra obtained in MS/MS mode with He cell gas. When He cell mode is used for interference removal, precise and accurate analysis is easily achieved. In He MS/MS mode, all interferences (arising from signal overlap from tailing of adjacent peaks and polyatomic ions isobaric interferences) are removed, yielding unbiased analysis and accurate results.




Figure 2. He MS/MS mode spectra: (top to bottom) 5 ppb Mn, 500 ppm Fe, and 5 ppb Mn + 500 ppm Fe

Figure 3 is an overlay of three spectra measured using He MS/MS mode; 1) Blank, 2) 10x whole blood, and 3) 500 ppt Mn spike in 10x whole blood. Table 1 summarizes the results of 10x diluted whole blood analysis and 500 ppt Mn spike recovery test. As shown, very low blank levels were achieved, ⁵⁵Mn was clearly resolved in the spectrum and good spike recoveries were obtained.




Figure 3. Spectra of three samples: blank, 10x diluted blood and 10x diluted blood spiked with 500 ppt Mn

Table 1. 10x diluted whole blood analysis results for Mn

	Blood sample	Blood sample + 500 ppt Mn	Spike recovery
	µg/L		%
Sample A	0.413	0.983	114
Sample B	0.432	0.924	98

Conclusions

Quadrupole ICP-MS has been almost universally accepted for low level analysis of trace analytes in complex matrices. However, many challenging interferences remain unresolved, especially when trace analytes must be measured close to matrix element peaks in complex samples. The Agilent 8800 ICP-QQQ with MS/MS capability has abundance sensitivity better than 10^{-10} , which enables the analysis of trace analytes (such as Mn) in the presence of a high concentration of adjacent elements (such as Fe).

Measurement of titanium in clinical samples: Possible application to monitoring patients with joint replacements

Glenn Woods
Agilent Technologies (UK) Ltd.

Keywords

titanium, biological, serum, urine, joint-replacement, Seronorm, ammonia mass-shift

Introduction

Although titanium (Ti) has little or no direct biological role, it is widely used in dental, artificial/replacement joints and surgical reconstruction applications due to its high strength, light weight and the fact it is biocompatible. It is also used extensively as a pigment and abrasive polishing agent (as TiO_2) and is often found in foods and toothpaste due to its inert nature (this form is usually passed unaltered in faecal matter and is not normally transported or expressed through body fluids).

More recent applications with certain types of metal-on-metal (rather than ceramic or polymer based) joint replacements can lead to the release of wear metal particles or ions within the body of the patient. These can become highly concentrated in the synovial fluid (lubricating fluid of the joint), pass into the bloodstream and be expressed through urine. Unusually “high” concentrations of Ti can indicate a premature failure of a Ti-based joint and such failure can lead to infection or constant pain for the patient. It is therefore important to reliably determine the concentration of Ti within biological fluids at normal endogenous levels in order to obtain a basal concentration. An increase from this concentration could indicate an imminent failure of the joint.

Experimental

The determination of Ti in biological matrices is challenging for conventional ICP-MS, due to its low natural concentration and the presence of spectral interferences on all the Ti isotopes e.g. sulfur (as SO), P (as PO) and Ca. It is possible to use reaction chemistry with NH_3 cell gas in the CRC of a quadrupole ICP-MS (ICP-QMS), to mass-shift the Ti^+ to a higher mass product ion, leaving the interfering species behind.

However, the use of highly reactive cell gases in ICP-QMS is prone to severe errors, as there is no way to control the ions that enter the CRC. This means that the reaction chemistry and the product ions created can change dramatically, with even slight differences in sample matrix or co-existing analyte concentrations. For this application, the 8800 ICP-QQQ was used to provide controlled reaction chemistry with ammonia as the reaction gas and measuring Ti as the $\text{TiNH}_2(\text{NH}_3)_4^+$ cluster ion at the $M + 84$ amu transition.

Instrumentation: Agilent 8800 #100.

Plasma conditions: Preset plasma/General purpose.

Ion lens tune: Soft extraction tune: Extract 1 = 0 V, Extract 2 = -170 V.

CRC conditions: NH_3 gas (10% in He) at 1.7 mL/min, Octopole bias = -8 V, KED = -8 V.

Samples and sample preparation: Certified reference materials of human serum and urine were purchased from Seronorm (Norway). They were prepared in duplicate by 10x dilution into a basic diluent consisting of NH_4OH (0.5%), $\text{H}_4\text{-EDTA}$ (0.01%), BuOH (2%) & Triton X-100 (0.01%) in ultrapure water. No further matrix matching was applied for the standards.

Results and discussion

Selection of product ion for Ti measurement

In order to select the most appropriate Ti cluster ions in NH_3 mode, a product ion scan was performed for the ^{48}Ti isotope by introducing a 10 ppb Ti solution (Figure 1). Q1 was set to m/z 48, allowing only ions at the mass of the target precursor ion to enter the cell; Q2 was scanned over a selected mass range to measure all the product ions formed in the cell by NH_3 reactions with ^{48}Ti . Based upon this scan, the two most abundant cluster ions (Q1 + 84 amu [$\text{TiNH}_2(\text{NH}_3)_4$] and Q1 + 102 amu [$\text{Ti}(\text{NH}_3)_6$]) were selected for further study. For each of the two reaction transitions identified above, neutral gain scans




Figure 1. Product ion scan for $^{48}\text{Ti}^+$ in NH_3 mode




Figure 2. Neutral gain scan for two $\text{Ti} \rightarrow \text{Ti}$ cluster ion transitions: For $\text{TiNH}_2(\text{NH}_3)_4^+$ cluster ions, $Q2 = Q1 + 84$ amu, and for $\text{Ti}(\text{NH}_3)_6^+$ cluster ions, $Q2 = Q1 + 102$ amu. The preservation of the natural Ti isotope abundance pattern ($^{46}\text{Ti}^+$, $^{47}\text{Ti}^+$, $^{48}\text{Ti}^+$, $^{49}\text{Ti}^+$ and $^{50}\text{Ti}^+$) can be seen, confirming that MS/MS mode provides complete control over the complex Ti-NH_3 reaction chemistry.

(where Q1 and Q2 are scanned synchronously, with a set mass difference between them ($Q2 = Q1 + 84$ and $Q2 = Q1 + 102$ in this case)) were performed. These scans are shown in Figure 2 confirming the correct natural isotopic abundances for the different Ti isotopes. Without MS/MS capability, it would be impossible to preserve the isotopic information for this element due to the relatively complex nature of the Ti-ammonia adducts. The instrument cell conditions were optimized using simple HNO_3 acidified Ti standards and applied to the analysis of the CRMs.

Table 1 displays the results for both serum and urine sample types measured against a single calibration. The 8800 ICP-QQQ was operated also under no gas and He mode to provide comparative data, and three Ti isotopes were monitored for the same cluster ion transition, to give confirmation of the results.

Conclusions

Titanium was only certified in two of the four materials measured but the 8800 ICP-QQQ data were all comfortably within the measured ranges when operating under ammonia MS/MS mode, in contrast to no gas and He mode data. Importantly, the three Ti isotopes measured under ammonia MS/MS mode all gave equivalent data; this could indicate applicability of the method to the use of isotope-based analysis such as isotope dilution (ID) or isotope tracer analysis. The use of ammonia combined with MS/MS greatly simplifies the analysis of Ti in biological media for several isotopes. Furthermore, because MS/MS mode provides control over the reaction chemistry, no special attention needs to be paid to specific matrix matching regardless of the fluid investigated.

Table 1. Urine and serum sample recovery ($\mu\text{g/L}$) for Ti in Seronorm CRM using $\text{TiNH}_2(\text{NH}_3)_4^+$ cluster

Sample Name	Target	^{47}Ti [No gas]	^{47}Ti [He]	47 -> 131 Ti [NH_3]	48 -> 132 Ti [NH_3]	49 -> 133 Ti [NH_3]
Urine blank	4.6 (2.2-7.0)	1989.79	41.44	2.80	2.79	2.92
Urine blank	4.6 (2.2-7.0)	2004.91	44.30	3.50	2.93	3.33
Urine trace elements		1789.92	51.41	14.81	15.27	14.42
Urine trace elements		1749.13	52.58	14.99	15.49	15.50
Serum L1	1.28 (0.86-1.80)	144.18	3.79	1.21	1.15	1.14
Serum L1	1.28 (0.86-1.80)	128.97	2.95	1.27	1.18	1.09
Serum L2		100.16	3.95	1.76	1.92	1.61
Serum L2		95.65	3.02	1.82	1.64	1.76

Simultaneous quantitation of peptides and phosphopeptides by CapLC-ICP-QQQ

Silvia Diez Fernández¹, Naoki Sugiyama², Jorge Ruiz Encinar¹ and Alfredo Sanz-Medel¹

¹Department of Physical and Analytical Chemistry, University of Oviedo, Spain, ²Agilent Technologies, Japan

Keywords

proteins, peptides, phosphorus, phosphopeptides, sulfur, S-containing peptides, heteroatom, isotope ratio, pharmaceutical, clinical, drugs, metabolites, environmental, pesticides, nanotechnology, nanoparticles, oxygen mass-shift

Introduction

LC-MS/MS is used for the quantification of target proteins in pharma/biopharma and clinical research. The approach generally relies on the use of synthetic, isotopically-labeled forms of each target protein and peptide, which are used as internal standards for the specific quantitation of the corresponding target compound. In contrast, the high temperature plasma ionization source used in ICP-MS ensures that elemental response is practically independent of the original form of the compound, which enables non-species-specific (or compound-independent) quantitation of compounds by measuring the signal for an element contained in the target compound. In this way, different proteins and peptides containing the heteroatoms S and P can be quantified using a single S- or P-containing compound as a generic standard. Unfortunately, with conventional quadrupole ICP-MS, the DLs for P and S are compromised by their high ionization potential and

by multiple polyatomic interferences. The Agilent 8800 ICP-QQQ can effectively remove those interferences using reaction cell chemistry combined with the unique MS/MS mode, achieving excellent DLs for P and S even in organic solvents. This paper demonstrates the advantage of ICP-QQQ for the determination of proteins and peptides by measurement of P and S heteroatoms.

Experimental

Instrumentation: An Agilent 8800 Triple Quadrupole ICP-MS was used with an Agilent 1260 Series low flow capillary LC system. The standard 2.5 mm internal diameter (id) injector torch was replaced with the narrow injector, 1.5 mm id torch (G3280-80080) used for the analysis of volatile organic solvents. The exit of the LC column was interfaced to the ICP-MS via an Agilent capillary LC interface kit (G3680A) featuring a total consumption nebulizer and micro-volume spray chamber. O₂ gas (20% O₂ in Ar) was supplied to the plasma as an option gas at 0.08 L/min to prevent carbon build-up on the interface cones. Agilent ICP-MS MassHunter chromatographic software was used for integrated control of the LC-ICP-MS system and for data analysis.

CRC conditions: O₂ cell gas flow rate at 0.35 mL/min, Octopole bias = -18 V and KED = -6 V.




Figure 1. 8800 ICP-QQQ MS/MS operation in mass-shift mode to remove interferences on S (a) and P (b)

Acquisition conditions: MS/MS O₂ mass-shift method was applied for P and S measurement as shown in Figure 1.

LC conditions: An Agilent Zorbax SB C18 (5 μm, 150 x 0.3 mm) reverse phase column was used with a flow rate of 5 μL/min. Mobile phases of water (A) and acetonitrile (B) were used for a gradient elution with the following profile: 0-3 min: 1% B; 3-35 min: 1-60% B linear. Both mobile phases contained 0.1% formic acid and 10 ppb Ge as ISTD and for tuning. The injection volume was 2 μL.

Reagents: Bis-4-nitro-phenyl phosphate (BNPP, 99% purity) and methionine (≥ 99% purity) (Sigma-Aldrich, Steinheim, Germany) were used as calibration standards for phosphopeptides and S-containing peptides respectively. Amino acid sequences of the phosphopeptides were LRRRA-pS-LG and KRS-pY-EEHIP, and the S-containing peptides were A-C-TPER-M-AE and VP-M-LK. All peptides were purchased from AnaSpec (Fremont, CA, USA) with purity ≥95%.

Results and discussion

Calibration and DL

Calibration standards containing 25, 50, 100 and 200 ng/mL of both P and S (BNPP and methionine,




Figure 2. Calibration plots for P and S (in BNPP and methionine, respectively) at 25, 50, 100 and 200 ng/mL (as the element)




Figure 3. Chromatogram of phosphopeptides and S-containing peptides. Sample: 45 ng/mL of two phosphopeptides and two S-containing peptides, and 105 ng/mL of BNPP and methionine (conc. as P or S)

respectively) were injected and measured. Excellent linearity and RSD of <4% was obtained (see Figure 2).

The chromatogram for the 50 ng/mL standard was used for signal to noise (S/N) and DL calculation. The DL achieved was 0.10 ng/mL for P and 0.18 ng/mL for S. As the injection volume was 2 μL, the DLs in absolute weight were calculated to be 6.6 fmol and 11 fmol for P and S, respectively.

Measurement of phosphopeptide and S-containing peptides

Finally, a sample containing a mixture of phosphopeptides and S-containing peptides was analyzed. The sample was also spiked with the standards methionine and BNPP for non-species-specific calibration. The chromatogram shown in Figure 3 illustrates the excellent peak shape and S/N obtained, demonstrating the exciting potential of ICP-QQQ for quantitative protein and peptide analysis using measurement of P- and S-heteroatoms.

More information

Simultaneous quantitation of peptides and phosphopeptides by capLC-ICP-MS using the Agilent 8800 Triple Quadrupole ICP-MS, Agilent application note, 5991-1461EN.

Measurement of selenium in clinical samples in the presence of gadolinium-based magnetic resonance imaging contrasting agents

Glenn Woods
Agilent Technologies (UK) Ltd.

Keywords

selenium, enzyme, blood, serum, urine, MRI contrasting agents, gadolinium, molybdenum, zirconium, neutral gain scan, oxygen mass-shift

Introduction

Selenium is an important micronutrient for human, mammalian, bacterial and plant life and is contained within several co-factors and enzyme systems. It is monitored in blood, serum and urine as part of human health, and a deficiency can indicate an illness (particularly if the levels change suddenly), such as cancer, diabetes and tuberculosis (TB).

For cancer patients, determining the location of the tumour often requires the use of Magnetic Resonance Imaging (MRI). However, for some soft tissues such as the brain, a “contrasting agent” is needed to effectively show the location of the tumour or problem area. There are several contrasting agents which are salts or chelates of gadolinium (III) (Gd(III)), trade names are given in brackets:

Gadodiamide (Omniscan), Gadobenate (MultiHance), Gadopentetate (Magnevist), Gadoteridol (ProHance), Gadofosveset (Ablavar, formerly Vasovist), Gadoversetamide (OptiMARK), Gadoxetate (Eovist), Gadobutrol (Gadavist)

Unfortunately Gd has a relatively low second ionization potential (12.09 eV) meaning it can form Gd^{++} ions in the plasma. These Gd^{++} ions appear at half their original mass (as a quadrupole measures ions based on their mass to charge ratio or m/z) and form interferences on all of the main analytical isotopes of Se. This is complicated to a greater extent as Gd has several odd-mass isotopes which form Gd^{++} interferences at half-mass (e.g., $^{155}Gd^{++}$ would appear at m/z 77.5). This makes the spectrum in the mass region of the Se isotopes quite complex when Gd is present in the sample. In a typical patient's sample, the Gd concentration can vary between zero to

several thousand parts per billion ($\mu\text{g/L}$). Because of the variability from patient-to-patient (which is also time-dependant on a sample-to-sample basis due to the contrasting agent's half-life in the body) a simple mathematical correction cannot always be made or a constant “background” be assumed.

Experimental

In order to remove the Gd-based interference, Se^+ can be reacted with oxygen cell gas in the collision/reaction cell to produce SeO^+ as a product ion. The Se-O reaction is slightly endothermic ($\Delta H_r = 0.71$ eV) which means that the reaction yield for SeO^+ would be relatively low. However the bias voltage on the ORS can be adjusted to increase the ion energy improving reaction yield significantly over a more “thermalized” approach. These conditions are referred to as high ORS bias conditions.

Instrumentation: Agilent 8800 #100.

Plasma conditions: Preset plasma/General purpose.

Ion lens tune: Soft extraction tune: Extract 1 = 0 V, Extract 2 = -170 V.

CRC conditions: O_2 gas at 0.3 mL/min, Octopole bias = -15 V, KED = -8 V.

Results and discussion

O_2 mass-shift method

Using O_2 mass-shift, the analyte is measured at $M + 16$ amu (e.g. $^{78}Se^+$ is measured as $^{78}Se^{16}O^+$ at 94 amu). With conventional quadrupole ICP-MS, any ^{94}Mo or ^{94}Zr present in the sample would interfere with the measurement at this mass. However, with MS/MS mode, ^{94}Mo or ^{94}Zr are removed by Q1 as it is set to the mass of the Se^+ precursor ion at 78 amu, and $^{156}Gd^{++}$ is eliminated as Q2 is set to the SeO^+ product ion mass of 94 amu. Even if Gd did form GdO^{++} this would also be eliminated by Q2 as the apparent mass (m/z) of $^{156}Gd^{16}O^{++}$ is 172/2 (86 amu). Figure 1 is a graphical representation of the MS/MS setup.




Figure 1. Representation of ICP-QQQ setup with Q1 set to 78 amu and Q2 set to 94 amu. Mo and Zr-based interferences are eliminated by Q1 and Gd^{++} is eliminated by Q2 allowing the measurement of ^{78}Se as $^{78}\text{Se}^{16}\text{O}^+$

To check for efficient conversion of Se^+ to SeO^+ , a neutral gain scan covering the mass range of all the SeO^+ product ions was performed for a 5 ppb Se solution. Figure 2 displays the isotope pattern of the $+^{16}\text{O}$ -atom transitions for all the Se isotopes, showing a perfect match with the theoretical isotopic fit.

Se measurement in human serum

Instrument cell conditions were optimized using a Se standard in a simple HNO_3 matrix. A pooled human serum sample was prepared by 10x dilution into a basic diluent consisting of NH_4OH (0.5%), H_4 -EDTA (0.01%), BuOH (2%) & Triton X-100 (0.01%) in ultrapure water. The sample was prepared unspiked and also spiked with Gd equivalent to 250, 500 and 1000 $\mu\text{g}/\text{L}$ in the

original sample, and analyzed using the 8800 ICP-QQQ in no gas and O_2 mass-shift modes of operation for comparison. The data is summarized in Table 1. The results show that, under no gas conditions, the apparent Se concentration is influenced by the variable Gd^{++} interference. Recovery based upon the original unspiked sample demonstrates an over-recovery of almost 130% for the no gas data when Gd is at a concentration of 1000 $\mu\text{g}/\text{L}$. In contrast, the Se data measured with MS/MS mass-shift mode remains essentially constant at all levels of Gd matrix. This would indicate that the O_2 mass-shift reaction is independent of the Gd concentration and is highly applicable to this relatively difficult and important application.




Figure 2. Neutral gain scan for 5 ppb Se solution. Note the excellent spectral fit. The six peaks corresponded to SeO^+ product ions formed from the Se isotopes: ^{74}Se , ^{76}Se , ^{77}Se , ^{78}Se , ^{80}Se and ^{82}Se .

Table 1. Serum sample data and recovery for Se with variable Gd concentration. Recovery is calculated based on determined Se concentration in unspiked serum sample. All data is dilution corrected.

	No gas mode		O_2 mass-shift	
	Conc. ppb	Recovery %	Conc. ppb	Recovery %
Serum	93.64	NA	91.42	NA
Serum with 250 $\mu\text{g}/\text{L}$ Gd	99.97	106.7	91.38	100.0
Serum with 500 $\mu\text{g}/\text{L}$ Gd	112.1	120.0	91.70	100.3
Serum with 1000 $\mu\text{g}/\text{L}$ Gd	121.1	129.3	91.78	100.4

- > Return to table of contents
- > Search entire document

Glossary

8800 ICP-QQQ Shorthand for Agilent 8800 Triple Quadrupole ICP-MS.

A

Agilent 8800 Triple Quadrupole ICP-MS	Official name of the Agilent Triple Quadrupole ICP-MS. It is also referred to as the 8800 ICP-QQQ or 8800. The 8800 ICP-QQQ possess two separate quadrupole mass filters, Q1 and Q2, one before and one after the collision/reaction cell. The term "triple quadrupole" or "triple quad" follows the same naming convention as LC-QQQ and GC-QQQ, used in organic mass spectrometry to describe a tandem mass spectrometer with two quadrupole mass analyzers separated by a chamber (or cell) containing an ion guide.
Agilent 8800 #100	The standard configuration instrument: Equipped with an x-lens, Ni cones, a 2.5mm injector quartz torch and a glass concentric nebulizer. It includes 4 cell gas mass flow controllers and 5 plasma gas mass flow controllers.
Agilent 8800 #200	The semiconductor configuration instrument: Equipped with an s-lens, Pt cones, a 2.5mm injector quartz torch and a PFA concentric nebulizer. All other parts of the instrument are the same as the standard configuration (option #100).
Acquisition conditions	Parameters including: Peak profile, mass, integration time, scan number and replicate.
Ammonia, NH ₃	A reaction gas used in the collision/reaction cell. NH ₃ is a very reactive gas, which is used both in on-mass methods and mass-shift methods to remove/avoid interferences.
amu	Atomic mass unit. An obsolete, non-SI unit that is still in common use in its abbreviated form "amu", meaning the same thing as "unified atomic mass unit" (u) or dalton (Da). All are used to indicate the atomic mass of ions, atoms or molecules, based on the carbon 12 standard.
AS, abundance sensitivity	The measure of an analyzer's ability to separate adjacent peaks differing greatly in intensity. The 8800 ICP-QQQ with MS/MS operation delivers unmatched peak separation (abundance sensitivity <math><10^{-10}</math>), as the resolution performance is the product of the abundance sensitivity of the two quadrupoles.
ASX-520	Autosampler suitable for medium to high sample throughput applications, with rack configurations providing up to 360 vial positions (up to 720 with the extended rack XLR-8 version).

B

Bandpass	Mode of operation of a multipole ion guide, where both a low-mass cut-off and high-mass cut-off are applied, rejecting ions below and above a certain m/z . A bandpass filter passes a "window" of masses (typically covering a 20-30 m/z range) through the ion guide, and is therefore distinct from a mass filter, which is capable of unit mass resolution (single m/z mass selection).
Batch	The acquisition batch contains all the information required for a sample analysis or experiment, including peripump program, tuning conditions, acquisition parameters, sample list and data analysis (calibration) parameters. The data analysis (DA) batch contains the measured results for the batch of samples, and includes the calibration plots, internal standard signals and QC flags
BEC, Background equivalent concentration	The magnitude of a signal in a blank, expressed as a concentration.

C

CRGS, Carrier gas	Carrier gas is an Ar gas supply flowing through the nebulizer to convert a liquid sample into a fine aerosol. It is a tuning parameter of the plasma.
CRC conditions (collision/ reaction cell)	Parameters that affect collision/reaction cell performance such as cell gas, cell gas flow rate, octopole bias voltage, KED bias and deflection lens.
Charge transfer	A reaction mechanism that relies on the exchange of charge between ions and cell gas molecules, e.g., $\text{Ar}^+ + \text{NH}_3 \rightarrow \text{Ar} + \text{NH}_3^+$
Collision mode	A cell mechanism to remove interferences either by collisional dissociation or by kinetic energy discrimination (KED). With KED, ions entering the collision/reaction cell collide with the cell gas (such as helium). Since polyatomic ions have a larger ionic cross-section than mono-atomic analyte ions at the same mass, the polyatomic ions undergo more collisions than the analyte ions, and so lose more energy. By the cell exit, the lower energy ions (the polyatomics) can be separated from the higher energy (analyte) ions by applying a bias voltage "step". This is known as kinetic energy discrimination (KED).
Cool plasma	A technique used to reduce interferences. Under low temperature plasma (cool plasma) conditions, the formation of interferences such as Ar^+ , ArO^+ and ArH^+ is suppressed, allowing the detection of Ca^+ , Fe^+ and K^+ at the trace level. Typical RF power for cool plasma is 600-900 W.
CRC	Abbreviation of Collision Reaction Cell.

D

DL	Abbreviation of Detection Limit. Also called LOD (limit of detection). It is the concentration that is equivalent to 3 times the standard deviation (SD) of the background signal.
Dynamic range or analytical working range	The range of linearity of an analytical instrument. The 8800 ICP-QQQ is fitted with an advanced, dual-mode, discrete dynode electron multiplier (DDEM) that provides a full nine orders dynamic range under standard operating conditions.

DiGS, Dilution gas Argon gas flow added to the carrier gas via a dilution gas port located between the torch and the spray chamber. A dilution gas is used for Aerosol Dilution with HMI. The gas supply used for the DiGS can also be switched automatically to add the gas flow to the spray chamber instead (known as make-up gas or MUGS). It is a tuning parameter of the plasma.

E

Enthalpy of reaction, ΔH_r Amount of energy (heat) absorbed or released by a reaction. When ΔH_r is positive ($\Delta H_r > 0$), the reaction is endothermic, meaning energy is required (absorbed) for the reaction to occur. When ΔH_r is negative ($\Delta H_r < 0$), the reaction is exothermic, meaning energy is released by the reaction, which is spontaneous.

G

GC Interface kit Agilent's GC-ICP-MS interface features a fully heated inert transfer line and separately heated inert torch injector that provides reliable separation of volatile compounds.

H

HMI, High Matrix Introduction HMI Aerosol Dilution technology is standard on the 8800 ICP-QQQ, extending the TDS range to % level, while eliminating the added cost, time and potential errors of conventional liquid dilution.

Hard extraction A tuning condition when a negative voltage is applied to the extraction lens. Hard extraction provides higher sensitivity at lower plasma temperature than soft extraction. Cool plasma conditions require hard extraction.

Helium mode, He mode See collision mode.

HR-ICP-MS, high-resolution ICP-MS Also known as magnetic sector, sector field or double focusing. Magnetic sector based ICP-MS instruments are capable of resolution ($M/\Delta M$) of up to 10,000 and are able to resolve most polyatomic species from analytes at the same nominal mass.

I

I-AS, Integrated Auto sampler Integrated, covered auto sampler with pumped rinse station; ideal for ultra-trace analysis and small sample volumes (as low as 0.5 mL). Flexible rack configurations offer a maximum capacity of 89 vials, plus 3 rinse vials.

ICP Inductively coupled plasma, generated by applying a high-power radio frequency (rf) field to a flow of argon gas. The plasma is a high temperature ion source, up to 10,000 K maximum and around 7,500 K in the central channel.

ICP-MS Inductively coupled plasma mass spectrometer or spectrometry.

ICP-QQQ Abbreviation for triple quadrupole ICP-MS.

IDA, ID, ID-MS Isotope Dilution Analysis or Isotope Dilution Mass Spectrometry is a highly accurate method to quantify elements based on the change in isotope ratio that results from the spiking of an unknown sample with a spike enriched in one isotope of the target analyte. Because each sample result is based on the measurement of the change in ratio in that sample, rather than relative to a response in a separate calibration standard, IDMS results are also directly traceable to certified standards, which reduces uncertainty.

Inert Sample Introduction kit	O-ring-free and manufactured from PFA for the lowest contamination levels. Demountable torch with Pt or sapphire injector options. HF resistant, and suitable for high-purity reagents.
Interferences - spectral	Direct overlap from a different element with an isotope at the same nominal mass (isobar), or overlap from a polyatomic ion, or doubly-charged ions resulting from the loss of two electrons instead of just one. Because the quadrupole separates ions based on m/z (mass over charge ratio), a doubly-charged ion (M^{2+}) will appear at mass $M/2$.
Ion guide	Operation of an ion lens where no mass rejection is performed. Applies to simple electrostatic ion lenses, and also to multipole ion guides operated with no low- or high-mass cut-off.
IP, Ionization potential	The first ionization potential of the element is the energy required to remove one electron from a neutral atom and is specific for each element. Most elements are largely converted (>90%) to singly-charged ions in an argon plasma. Elements with a low second IP will also form some doubly-charged ions.
IR, isotope ratio	Ratio of abundance of two isotopes of an element.
Isobar	Refers to isotopes of different elements that appear at the same nominal mass. These overlaps occur when atoms of two different elements (i.e. different number of protons in the nucleus, so different atomic number) each have an isotope with the same atomic weight (same total number of protons plus neutrons in the nucleus, e.g. ^{204}Pb and ^{204}Hg).
Isobaric interferences	Overlaps that occur at the same mass (see isobar). These overlaps/interferences can be resolved by reaction chemistry (e.g. NH_3 is used to separate Pb from the Hg overlap), but cannot be separated by high-resolution ICP-MS; separation of ^{204}Pb from ^{204}Hg would require a resolution of around 500,000 (50x higher than can be achieved by any commercial high-resolution ICP-MS).
Isotope	A specific form (atomic weight) of an element. Many elements have atoms with different atomic weights, such as Pb 204, 206, 207 and 208; these are called isotopes. The different isotopes of Pb all have 82 protons in the nucleus (Pb has atomic number 82) but a different number of neutrons, so the atomic weight is different for each isotope.
ISTD, internal standard	Internal standards are commonly used in ICP-MS, particularly where samples vary in composition from the calibration standards. Changes in sample transport, nebulization efficiency and signal intensity (long-term drift) would all lead to errors, which may be corrected if an ISTD element with similar behavior is used as a reference.

K

KED, Kinetic Energy Discrimination
 KED is used to discriminate the analyte ion of interest from interfering ion(s) by the difference of kinetic energy. Refer to collision mode. KED is also used as a tuning parameter of CRC conditions: $\text{KED} = (\text{Q2 bias voltage}) - (\text{octopole bias voltage})$.

L

LA, Laser ablation
 Method used for the direct analysis of solid samples using a laser to vaporize the sample before introduction to the plasma.

LC Speciation kits
 Sample introduction kits to facilitate LC coupling and provide turn-key methods for common speciation applications. A Capillary-LC connection kit is also available.

M

MUGS, Make-up gas	Make-up gas refers to Ar gas applied to the spray chamber to increase/adjust total injector gas flow rate. It is a tuning parameter of the plasma.
MS/MS mode	Acquisition mode unique to 8800 ICP-QQQ. MS/MS mode operates Q1 as a unit (1 amu window) mass filter and Q2 is also set to the single mass of the target ion or reaction product ion.
Mass pair	MS/MS mode requires a mass setting for Q1 and Q2. The selected mass settings for Q1 and Q2 are known as the mass pair. For example when As is measured in O ₂ mode, Q1 is set to the precursor ion (As ⁺) at <i>m/z</i> 75 and Q2 is set to the product ion (AsO ⁺) at <i>m/z</i> 91. 75 --> 91 is the mass pair for As in O ₂ mode.
Mass filter	Generic term for any mass analyzer capable of unit mass resolution. Note that the ion guide used in the CRC of some quadrupole mass spectrometers appears physically similar to a quadrupole mass filter. However, these ion guides cannot provide unit mass resolution because of the ion scattering effect at the higher pressures present in the CRC.
Mass spectrum	See spectrum.
Mass-shift method	A method where the analyte ion is reactive and is moved to a new mass free from the original interference. Sometimes referred to as "indirect" measurement, e.g. Se ⁺ reacts with O ₂ in the cell and is converted to SeO ⁺ . It can then be detected free from the original interference of ArAr ⁺ .
MH, Mass Hunter software	Software package that provides comprehensive instrument control for the 8800 ICP-QQQ and accessories, and integrated data processing.
MSA, Method of Standard Additions (also known as StdAdd)	A calibration solution is spiked at multiple levels directly into the unknown sample, giving a calibration of response against added concentration. MSA eliminates matrix effects by calibrating in the sample matrix.

N

Neutral Gain Scan	Q1 and Q2 scan together, a fixed mass-shift apart. For example Q2 scans at Q1 + 16 amu for O-atom addition reactions.
-------------------	---

O

O ₂ , oxygen	A reaction gas used with the 8800 ICP-QQQ. A number of elements can be measured in mass-shift method using O ₂ , e.g. Se ⁺ can be measured as SeO ⁺ using O ₂ cell gas. O ₂ is also added to the plasma carrier gas to decompose the carbon matrix when organic solvents are analyzed.
Octopole bias (OctP Bias)	A CRC parameter. It is the bias voltage applied to the octopole ion guide, which determines the collision energy of analyte ions with cell gas molecules.

OIDA, on-line isotope dilution analysis	A very powerful and useful development of traditional isotope dilution, using on-line addition of the isotope spike. Removes the time consuming step of spiking enriched-isotope standards into each separate sample.
On-mass method	A method where reactive interferences are removed to allow an unreactive analyte to be measured at its original mass. Sometimes referred to as "direct" measurement, e.g., the interference of GdO^+ on Yb^+ can be removed by the reaction of GdO^+ with NH_3 .
Organics kit	Contains the sample introduction parts needed to run volatile organic solvents. Includes organics torch, solvent-resistant drain kit and uptake tubing.
ORS	Agilent's CRC design is known as the Octopole Reaction System. It is a temperature-controlled collision/reaction cell with octopole ion guide and four cell gas lines as standard on the 8800 ICP-QQQ. Provides maximum flexibility in collision and reaction modes, and uses a small internal volume cell to ensure rapid cell gas switching and high ion transmission.
ORS ³	Third generation Octopole Reaction System.
Oxygen atom abstraction or oxygen atom transfer	Reaction mechanism associated with the use of oxygen in the collision reaction cell.
P	
Preset method	Preset methods are provided in the ICP-MS MH software. These built-in methods cover a range of predefined operating conditions to suit different applications. Using a preset method, a user can create a new batch with minimum or no customization.
Preset plasma	Preset plasma conditions are a function of ICP-MS MH software. The software provides several predefined plasma conditions that users can select according to the application. This greatly simplifies system optimization by automatically tuning and calibrating the plasma parameters, rather than the user having to set a number of individual plasma tuning parameters. There are three preset plasma conditions that can be selected depending on the sample matrix: Low matrix, general purpose and HMI.
Polyatomic, Polyatomic ion	A molecular ion (an ion composed of more than one atom) that arises in the plasma or during ion extraction, and can appear at the same nominal mass as analyte ions. Polyatomics are usually interferences (such as ArO^+).
Precursor Ion Scan	Q1 scans a user set mass range, while Q2 is set to a single fixed mass, measuring all the reaction product ions at that mass, formed from the different ions entering the cell as Q1 scans the mass range.
Product Ion Scan	Q1 is set to a fixed precursor ion mass, while Q2 scans a user set mass range measuring all reaction product ions formed from that single precursor ion.

Q

Quadrupole bias (QP Bias or Qpole Bias)	Bias voltage applied to the Q2 rods. Used in conjunction with the Octopole bias to provide a bias voltage "step" at the cell exit, usually to reject unwanted low energy ions from the ion beam.
Quantitation or quantification	Quantitative results are produced by comparing signal intensities of elements in the sample to those generated by calibration standards.
Q1	First quadrupole in the configuration of the 8800 ICP-QQQ. Q1 is positioned in front of the ORS ³ , to control the ions that are passed to the cell and enable MS/MS operation.
Q2	Second quadrupole in the configuration of the 8800 ICP-QQQ. Q2 filters the ions that emerge from the cell exit, passing only the target analyte ions to the detector.

R

Resolution	The ability of a mass filter to separate adjacent masses. Defined as $M/\Delta M$; the mass of the target peak/the mass difference to the nearest adjacent peak that can be separated. Sometimes also quoted as the width of the peak at a given peak height (e.g. 0.75 amu at 10% peak height).
------------	---

S

Single Quadrupole MS, ICP-QMS	Conventional ICP-MS containing a single quadrupole mass filter.
Single Quad mode, SQ mode	Q1 operates as a wide band mass filter. SQ mode emulates conventional quadrupole ICP-MS.
Speciation measurement	Combination of chromatographic techniques with ICP-MS as a detector to determine the chemical form of elements in the sample.
Spectrum (Mass Spectrum)	After separation by the final mass filter (Q2), the ions are detected by an electron multiplier. The detector electronics count and store the total signal for each mass (m/z), creating a mass spectrum. The spectrum that is produced provides a simple and accurate qualitative representation of the sample. The magnitude of each peak is directly proportional to the concentration of an element in a sample.
STS, ShieldTorch System	A technique to eliminate capacitive coupling between the RF coil and plasma, keeping the plasma potential low and energy distribution of ions narrow. The technique is crucial for cool plasma and collision mode.

T

TDS, total dissolved solids	The total summed concentration of all non-volatile, dissolved inorganic and organic substances in a liquid. The matrix tolerance of Agilent ICP-MS instruments is 0.2–2 % TDS.
Triple quadrupole ICP-MS	ICP-MS with a tandem MS configuration, featuring a quadrupole mass filter (Q1) in front of the collision/reaction cell (CRC) followed by a second quadrupole mass filter (Q2).

U

UPW	Ultra Pure Water, purified by ion exchange to remove trace contaminants. Used for preparation of standards and for sample dilution for ultra-trace analysis
-----	---

Agilent 8800 ICP-QQQ Application Handbook

> [Return to table of contents](#)
> [Search entire document](#)

Index

Symbols

⁹⁰ Sr	70
¹²⁹ I	67
¹³⁷ Cs	70

A

abundance sensitivity	50, 55, 67, 70, 74
agricultural	22
ammonia	42, 64
ammonia mass-shift	10, 12, 33, 58, 77
ammonia on-mass	12, 53
arsenic	12, 28, 35, 48
artifacts	64

B

barium	33, 60, 70
biodiesel	25
biological	77
blood	81

C

cadmium	31
calcium	8
catalytic converter	42
cerium	60
cesium	70
chlorine	16
chromium	10, 12
clinical	79
cobalt	48
copper	55, 58
CRMs	28, 31

D

dating	64
doubly-charged ion interferences	35
drugs	79

E

environmental	22, 25, 28, 31, 42, 67, 70, 79
enzyme	81
ethanol	25

F

fluorine-containing polyatomic ions	33
food	28, 31, 64

G

gadolinium	81
gadolinium oxide	53
geochemical	42
geochemistry	60
geochronology	64
germanium	12
gold	42

H

helium MS/MS	74
heteroatom	79
high purity metals	48
high purity zinc	55
human health	22
hydrochloric acid	12
hydrogen on-mass	8, 19, 31
hyperbolic	55

I

ID-MS	25
iodine	67
iron	50, 74
isotope	64
isotope ratio	79

J

joint-replacement	77
JSS 001-6	50
JSS 003-6	50

L

lanthanum	60
lead	64

M

manganese	50, 74
material science	60
mercury	64
metabolites	79
methanol	19
method of standard additions	8, 16, 60
mining	60
molybdenum	81
molybdenum oxide	31
MRI contrasting agents	81

N

nanoparticles	79
nanotechnology	79
neutral gain scan	81
NIST 3231 Level I and II	67
NIST SRM 2773	25
nitrous oxide on-mass	70
N-methyl-2-pyrrolidone	16
NMP	16
nuclear	67, 70

O

OIDA	22
online isotope dilution analysis	22
on-mass	64
ore exploration	42
organic solvents	19
organometallic	58
oxygen and hydrogen on-mass	70
oxygen mass-shift	10, 12, 14, 16, 19, 22, 25, 28, 35, 39, 48, 50, 53, 60, 79, 81
oxygen on-mass	33, 67

P

peptides	79
pesticides	79
pharmaceutical	79
phosphopeptides	79
phosphorus	14, 16, 19, 50, 79
Platinum Group Elements	42
precious metals	64
process chemicals	8, 10, 16
proteins	79

Q

quadrupole	55
------------	----

R

radioactive	70
radioisotopes	70
radionuclide	67
Rare Earth Elements	28, 53, 60
rare earth oxide	53
ratio	64
RCA Standard Clean	12
reaction chemistry	39
REE	28, 53, 60
REO	53
resolution	55

S

samarium oxide	53
S-containing peptides	79
selenium	22, 28, 81
semiconductor	8, 10, 12, 14, 16, 58
Seronorm	77
serum	77, 81
silicon	16, 19
silicon wafer	12, 14
silver	42
steel	50
strontium	70
sulfur	16, 19, 25, 50, 79
sulfuric acid	10

T

titanium	10, 14, 39, 58, 77
----------	--------------------

U

ultra pure water	8
UPW	8
urine	77, 81

V

vanadium	10, 12
Vapor Phase Decomposition	14
VPD	14

W

whole blood	74
-------------	----

X

xenon	67
-------	----

Z

zinc	55
zirconium	35, 48, 70, 81

AGILENT TECHNOLOGIES

Agilent 8800 ICP-QQQ Application Handbook

- > [Return to table of contents](#)
- > [Search entire document](#)

For more information

Learn more:

www.agilent.com/chem

Buy online:

www.agilent.com/chem/store

U.S. and Canada

1-800-227-9770

agilent_inquiries@agilent.com

Europe

info_agilent@agilent.com

Asia Pacific

inquiry_lsca@agilent.com

This information is subject to change without notice.

Publication number 5991-2802EN
© Agilent Technologies, Inc. 2013
September 2, 2013



Agilent Technologies