

Analytical Methodology to Monitor the Environmental Fate of Atrazine and Cyromazine

Anthony Macherone, Ronald Honnold, Stephan Baumann
Agilent Technologies, Inc. 5301 Stevens Creek Road Santa Clara, CA 95051 USA

ASMS 2011
WP 329



Introduction

Atrazine is an herbicide used globally. There is compelling evidence that atrazine is linked to deleterious effects in amphibians. One study demonstrated that upon exposure to atrazine, genetically male frogs became hermaphrodites. Male tadpoles lost the ability to develop masculine characteristics and feminization was observed to the point of developing viable ovaries. The catabolism of atrazine by microbial flora yields cyanuric acid which, in the presence of melamine, yields insoluble and highly toxic melamine cyanurate. Cyromazine, a common insecticide, is metabolized by mammals to melamine. It is not unreasonable to consider coeval application of atrazine and cyromazine in farming communities. Herein is presented a sensitive and selective method to monitor these compounds and their metabolites in soil extract.


Atrazine Pathway

Microbial catabolism of atrazine has been found to occur via several pathways. Monitoring the first metabolite in the pathway may offer insight into the particular species involved. We have chosen to monitor the pathway with desethylatrazine as the first metabolite.

Microbes and first metabolites

Compound	Microbes	First Metabolic Product
Atrazine	<i>Pseudomonas</i> sp. ADP, <i>Ralstonia</i> sp. M91-3, <i>Clavibacter</i> sp., <i>Agrobacterium</i> sp. J14a	Hydroxyatrazine
Atrazine	<i>Rhodococcus</i> spp., N186/21, TE1, <i>Pseudomonas</i> spp. 192/194, <i>Streptomyces</i> sp. PS1/5	Desisopropylatrazine
Atrazine	<i>Rhodococcus</i> spp. N186/21, TE1, <i>Nocardia</i> sp., <i>Alcaligenes</i> sp. SG1, <i>Streptomyces</i> sp. PS1/5	Desethylatrazine

Atrazine microbial metabolic pathway



Experimental

Atrazine, desethylatrazine, cyanuric acid, cyromazine, and melamine were prepared at a stock concentration of 1 mg/mL and spiked into extracted soil at 2 pbp through 100 pbp for analysis. Derivatization of the analytes was performed by adding 100 µl 49:49:1 BSTFA/1% TMCS:anhydrous ethyl acetate:anhydrous pyridine to the dried spiked matrix extract and heating at 70 C for 30 minutes. The GC was configured in a sequential, two-column fashion using a purged union to connect them. The purpose of this column configuration is to facilitate backflush of high boilers and heavy matrix and provide clean, reproducible chromatography and spectra.


Experimental

GC-MS/MS Conditions

Oven Program
65 °C for 1 min
then 10 °C/min to 230 °C for 0 min
then 30 °C/min to 300 °C for 2 min
Injection Volume 1 µL, pulsed splitless
Thermal Aux 280 °C
Column flow 1.2 mL/min
Source 280 °C
EI MS/MS Mode: MRM Transitions

Name	Precursor	Product	Collision Energy
Cyanuric Acid	345	330	10
Cyanuric Acid	345	215	8
Cyanuric Acid	345	188	12
Desethylatrazine	259	244	10
Desethylatrazine	259	217	10
Atrazine	215	200	10
Atrazine	200	104	20
Melamine	342	327	10
Melamine	342	285	20
Melamine	327	171	17
Cyromazine	310	309	10
Cyromazine	310	295	10

GC configuration showing back flush flows




Blue arrows = analysis mode

Red arrows = backflush mode

Experimental

Derivatized Analytes




Results and Discussion

Calibration from 2.0 pbp to 100.0 pbp

Excellent signal to noise for all compounds at 2.0 pbp

Analyte	Retention Order	Peak to Peak S/N at 2 pb in Matrix	r ² for cal curve
Cyanuric Acid	1	42:1	0.994
Atrazine	2	70:1	0.996
Desethylatrazine	3	4:1	0.999
Melamine	4	>1000:1	0.995
Cyromazine	5	39:1	0.999



Calibration Curves show for all analytes.

Results and Discussion

In 2009, soils from the Pacific Northwestern part of the United States was extracted for the purposes of pesticide analysis via GC tandem mass spectrometry. The soil was spiked with fifteen common pesticides. MS scan with AMDIS deconvolution and NIST08 library searching elicited the presence of 20 other possible compounds.


Unknowns and Targeted analysis of soil matrix

Azinphos methyl	>90	p,p DDD	>90	Malathion	Imiprothrin
Benzophenone	76	Permethrin	>90	Chlorpyrifos	p,p'-DDT
Chlorthaloniol	>90	Phenanthrene-D10	76	Allethrin	Piperonyl butoxide
Cyhalothrin	>90	Phenol, 2-methyl-	72	o,p'-DDE	Cypermethrin I
Eugenol	70	phorate	>90	Dieldrin	Cypermethrin II
folpet	>90	squalene	83	o,p'-DDD	Cypermethrin III
gamma-tocopherol	55	Terbutcarb	>90	Endosulfan II	Cypermethrin IV
hexachlorobenzene	>90	triphenyl phosphate	>90	o,p'-DDT	Deltamethrin
Metalaxyl	>90	Vitamin E	93		

Since none of the analytes in this study were expected to be present in the extracted soil matrix, the authors used it for this study. The data shown below gives a comparison of an un-spiked soil matrix blank which unexpectedly elicited the presence of cyanuric acid to a 2pp matrix spiked calibrator. Ion ratios were well within expected values and the corrected concentration of cyanuric acid in the un-spiked soil blank is approximately 2 pbp.

Matrix Blank actually contained Cyanuric Acid!

Sample	Type	Method	Cyanuric Acid (345.0 -> 188.0)	(345.0 -> 215.0)
Name	Type	Level Exp. Conc. Units	RT Resp. Final Conc.	Ratio Ratio
Blank 1	Blank		ppb	
Matrix Blank 1	Sample		ppb 12.898 1487	1.25 35.61 45.43
2 pbp 1	Cal	1 2	ppb 12.903 2110	1.92 34.41 51.34



With the exception of cyanuric acid, no interference from the other compounds present in the matrix was observed and all of the monitored analytes gave a linear curve fit coefficient (R²) greater than 0.994. This simple but elegant method demonstrates the feasibility for the detection of toxic atrazine metabolites in the environment. This method further shows the ability to monitor other metabolites in the environment, especially near irrigation and water run off areas that can feed back into the food chain.

Conclusions

The bioavailability of these compounds and their deleterious effects on various species are unknown. However, the propensity of atrazine and cyromazine and their microbial metabolic end products being present in the environment, especially in or near rural watersheds is high. As we know, cyanuric acid and melamine, the metabolites of atrazine and cyromazine, respectively form toxic, insoluble melamine cyanurate under aqueous conditions. The method design presented herein provides a tool for monitoring these compounds and possibly providing answers to the effects of aquatic species exposure to these compounds. The application lends itself to the analysis of fish and amphibian or any water based fauna tissue or bio-fluid. Detection limits of 2 pbp or less were determined for all analytes spiked into the soil matrix. Mid-column backflush also permits excellent reproducibility in terms of retention times and spectral accuracy.

Bibliography

Zeng Y, Sweeney CL, Stephens S, Kotharu P, Turnbull M. Atrazine Pathway Map. University of Minnesota. December 18, 2008. Retrieved May 25, 2011 from http://umbdd.msi.umn.edu/atr/atr_map.html

Hayes T, Haston K, Tsui M, Hoang A, Haeffle C & Vonk A. Atrazine-induced hermaphroditism at 0.1 pbp in American Leopard frogs (*Rana pipiens*): Laboratory and field evidence. *Environmental Health Perspectives* 111(4):April 2003.