

Application Note

Roger Bardsley, Applications Chemist, Teledyne Tekmar

Page | 1

Abstract

US EPA Method 524.2 is used primarily by environmental labs for the analysis of Volatile Organic Compounds (VOCs) in drinking water. While this method is effective at concentrating the trace levels of VOCs sometimes found in finished drinking water, it also tends to transfer significant quantities of water vapor to the Gas Chromatograph/Mass Spectrometer (GC/MS) due to the 4 minute desorb time recommendation.

Teledyne Tekmar's Lumin Purge and Trap (P&T) Concentrator incorporates a Moisture Control System (MCS) specifically designed to reduce the amount of water transferred during desorb, in comparison to current purge and trap technologies. This application note will demonstrate the ability of the Lumin to reduce the amount of water transferred to the Agilent GC/MS system when using this method.



Introduction

Gas Chromatography (GC) inlets, columns and Mass Spectrometers (MS) are sensitive to water. However, to detect toxic VOC compounds, the VOCs must be purged and trapped from water. This places water vapor on the GC/MS system, which contributes to column and inlet issues and frequent MS cleaning to remove the oxidation from the MS source.

The US EPA method that places a large amount of water vapor on the GC/MS system is US EPA method 524.2. This method recommends approximately 4 minutes of desorb time from the analytical trap, which inundates the GC/MS system with large volumes of water, and leads to poor chromatography for early eluting gases.

Traditionally, purge and trap instruments reduced the amount of water transferred to GC/MS instrument through numerous water management techniques. The Teledyne Tekmar Lumin significantly reduces the amount of water transferred to GC/MS system over the current P&T instrumentation through the design of its MCS. Additionally, improvements in analytical trap cooling have further reduced sample cycle times, permitting more samples to be processed within a 12 hour period.

Sample Preparation

Calibration standards were prepared from Restek® standards from Drinking Water VOA MegaMix, Oxygenates, Ketones and 502.2 Mix. The ketones mix compounds were present at 2.5 times the concentration of other compounds in the mix. The oxygenate compound, t-butanol, was present at 5 times the concentration of other compounds in the mix.

A calibration curve was prepared from 0.5 ppb to 50 ppb for all of the compounds except the ketones and t-butanol. The ketones range was from 1.25 ppb to 125 ppb. The t-butanol range was from 2.5 ppb to 250 ppb. The relative Response Factor (RF) using fluorobenzene as the internal standard was calculated for each VOC.

Seven 0.5 ppb standards were prepared to calculate the Method Detection Limit (MDL), accuracy and precision data. All calibration and MDL samples were analyzed with the Lumin and AQUATek 100 conditions in [Table I](#) and the GC/MS conditions in [Table II](#).

The Relative Standard Deviation (RSD) of the Response Factor (RF), MDL, accuracy, and precision data are shown in [Table III](#). A 0.2 ppb standard was analyzed to indicate the lowest detectable peak for the initial six gases unencumbered by excessive water in [Figure 1](#). [Figure 2](#) displays a 5 ppb standard, indicating excellent peak resolution for all of the VOCs, including the first six gases.

Experimental Instrument Conditions

| Table I Lumin and AQUATek 100 Conditions | | | |
|------------------------------------------|-----------------|-------------------------|-----------------|
| Standby | Variable | Bake | Variable |
| Valve Oven Temp | 150 °C | Bake Time | 2.00 min |
| Transfer Line Temp | 150 °C | Bake Temp | 280 °C |
| Sample Mount Temp | 90 °C | MCS Bake Temp | 200 °C |
| Purge Ready Temp | 35 °C | Bake Flow | 200 mL/min |
| MCS Purge Temp | 20 °C | AQUATek 100 | Variable |
| Standby Flow | 10 mL/min | Sample Loop Time | 0.35 min |
| Purge | Variable | Sample Transfer Time | 0.35 min |
| Purge Time | 11.00 min | Rinse Loop Time | 0.30 min |
| Purge Flow | 40 mL/min | Sweep Needle Time | 0.30 min |
| Dry Purge Temp | 20 °C | Presweep Time | 0.25 min |
| Dry Purge Time | 1.00 min | Water Temp | 90 °C |
| Dry Purge Flow | 100 mL/min | Bake Rinse Drain Cycles | 1 |
| Desorb | Variable | Bake Rinse Drain Time | 0.35 min |
| Desorb Preheat Temp | 245 °C | | |
| Desorb Time | 4.00 min | Trap | K |
| Desorb Temp | 250 °C | | |
| Drain Flow | 300 mL/min | | |

| Table II Agilent 7890B GC / 5977A MS | |
|--------------------------------------|-------------------------------------------------------------------|
| Agilent 7890B Conditions | |
| Column | DB-624 UI, 20 m x 0.18 mm, 1 m Film, Helium – 0.8 mL/min |
| Oven Profile | 35 °C, 4 min, 15 °C/min to 70 °C, 25 °C/min to 225 °C, 2 min hold |
| Inlet | 150 °C, 100:1 Split, helium saver 20 mL/min after 4 min |
| Agilent 5977A Conditions | |
| Temp | Transfer Line 225 °C; Source 230 °C; Quad 150 °C |
| Scan | Range 35 m/z to 260 m/z, Delay 0.8 min |
| Gain | Gain Factor 14 |

Results

| Table III Method 524.2 Calibration, Accuracy and Precision Data | | | | | |
|-----------------------------------------------------------------|---------------------|-----------|---------------------------------------|--------------|------------------|
| Compound | Calibration | | Accuracy and precision (n=7, 0.5 ppb) | | |
| | Linearity RF (%RSD) | MDL (ppb) | Avg. Conc. (ppb) | Accuracy (%) | Precision (%RSD) |
| Dichlorodifluoromethane | 8.0 | 0.12 | 0.48 | 95 | 8.0 |
| Chloromethane | 16.2 | 0.17 | 0.69 | 137 | 8.0 |
| Vinyl Chloride | 3.7 | 0.15 | 0.53 | 105 | 9.3 |
| Bromomethane | 11.7 | 0.29 | 0.82 | 165 | 11.2 |
| Chloroethane | 5.4 | 0.14 | 0.56 | 111 | 8.2 |
| Trichlorofluoromethane | 6.4 | 0.11 | 0.50 | 101 | 7.0 |
| Diethyl ether | 3.8 | 0.14 | 0.53 | 106 | 8.5 |
| 1,1-Dichloroethene | 7.6 | 0.16 | 0.54 | 108 | 9.2 |
| Carbon Disulfide | 4.0 | 0.17 | 0.55 | 109 | 10.0 |
| Iodomethane | 12.4 | 0.08 | 0.67 | 133 | 3.6 |
| Allyl Chloride | 4.4 | 0.15 | 0.57 | 115 | 8.1 |
| Methylene Chloride | 8.4 | 0.18 | 0.53 | 107 | 10.7 |
| Acetone (1.25 ppb) | 5.2 | 1.39 | 2.58 | 207 | 17.1 |
| trans-1,2-Dichloroethene | 4.0 | 0.10 | 0.53 | 107 | 6.1 |
| MTBE-d3 (Surrogate, 12.5 ppb) | 2.2 | | 12.69 | 101 | 1.3 |
| MTBE | 3.5 | 0.13 | 0.50 | 100 | 8.0 |
| t-Butyl Alcohol (2.5 ppb) | 3.9 | 0.68 | 2.56 | 102 | 8.4 |
| Diisopropyl Ether (DIPE) | 5.2 | 0.08 | 0.53 | 105 | 5.1 |
| 1,1-Dichloroethane | 3.4 | 0.16 | 0.51 | 102 | 9.9 |
| Acrylonitrile | 7.9 | 0.22 | 0.62 | 124 | 11.4 |
| t-Butyl Ethyl Ether (ETBE) | 2.9 | 0.10 | 0.49 | 99 | 6.3 |
| cis-1,2-Dichloroethene | 4.5 | 0.17 | 0.54 | 108 | 10.1 |
| 2,2-Dichloropropane | 5.8 | 0.08 | 0.49 | 98 | 5.3 |
| Bromochloromethane | 6.0 | 0.23 | 0.55 | 109 | 13.3 |
| Chloroform | 4.7 | 0.03 | 0.51 | 101 | 1.9 |
| Carbon Tetrachloride | 4.6 | 0.13 | 0.47 | 94 | 8.6 |
| Methyl Acrylate | 5.8 | 0.18 | 0.59 | 118 | 9.9 |
| Tetrahydrofuran | 15.3 | 0.44 | 0.62 | 124 | 22.6 |
| 1,1,1-Trichloroethane | 5.3 | 0.11 | 0.52 | 104 | 6.8 |
| 1,1-Dichloropropane | 7.1 | 0.19 | 0.51 | 102 | 12.1 |
| 2-Butanone (MEK, 1.25 ppb) | 8.3 | 0.49 | 1.41 | 113 | 11.1 |
| 1-Chlorobutane | 5.8 | 0.15 | 0.50 | 99 | 9.9 |
| Benzene | 4.4 | 0.16 | 0.49 | 98 | 10.6 |
| Propionitrile (1 ppb to 50 ppb) ¹ | 6.2 | | | | |

| Table III Method 524.2 Calibration, Accuracy and Precision Data | | | | | |
|-----------------------------------------------------------------|---------------------|-----------|---------------------------------------|--------------|------------------|
| Compound | Calibration | | Accuracy and precision (n=7, 0.5 ppb) | | |
| | Linearity RF (%RSD) | MDL (ppb) | Avg. Conc. (ppb) | Accuracy (%) | Precision (%RSD) |
| Methacrylonitrile | 7.0 | 0.40 | 0.61 | 121 | 20.9 |
| t-Amyl Methyl Ether (TAME) | 3.2 | 0.05 | 0.52 | 104 | 3.2 |
| 1,2-Dichloroethane | 2.6 | 0.16 | 0.50 | 100 | 10.5 |
| Fluorobenzene (IS, 12.5 ppb) | 2.8 | | | | 1.1 |
| Trichloroethene | 3.1 | 0.13 | 0.51 | 101 | 7.9 |
| 1,4-Difluorobenzene (Surr, 12.5 ppb) | 2.5 | | 12.51 | 100 | 1.8 |
| t-Amyl Ethyl Ether (TAEE) | 2.9 | 0.10 | 0.49 | 99 | 6.3 |
| Dibromomethane | 8.9 | 0.18 | 0.52 | 104 | 11.2 |
| 1,2-Dichloropropane | 2.7 | 0.14 | 0.55 | 111 | 8.3 |
| Bromodichloromethane | 7.1 | 0.10 | 0.46 | 92 | 6.6 |
| Methyl Methacrylate | 8.6 | 0.12 | 0.51 | 102 | 7.4 |
| cis-1,3-Dichloropropene | 5.7 | 0.14 | 0.46 | 92 | 9.6 |
| Toluene | 4.1 | 0.08 | 0.48 | 96 | 5.4 |
| 2-Nitropropane | 7.6 | 0.32 | 0.47 | 95 | 21.7 |
| Tetrachloroethene | 10.1 | 0.29 | 0.58 | 115 | 16.2 |
| 1,1-Dichloropropanone (63 m/z, 1.25 ppb) ² | 0.9987 | 2.07 | 1.53 | 123 | 43.0 |
| 4-Methyl-2-pentanone (58 m/z, 1.25 ppb) ² | 16.4 | 0.30 | 1.23 | 98 | 7.8 |
| trans-1,3-Dichloropropene | 7.1 | 0.10 | 0.48 | 97 | 6.3 |
| 1,1,2-Trichloroethane | 5.1 | 0.11 | 0.50 | 101 | 6.8 |
| Ethyl Methacrylate | 10.4 | 0.17 | 0.48 | 97 | 11.3 |
| Dibromochloromethane | 7.1 | 0.11 | 0.46 | 92 | 7.4 |
| 1,3-Dichloropropane | 2.5 | 0.16 | 0.50 | 99 | 10.2 |
| 1,2-Dibromoethane (EDB) | 5.3 | 0.15 | 0.51 | 102 | 9.3 |
| 2-Hexanone (1.25 ppb) | 16.6 | 0.51 | 1.27 | 101 | 12.8 |
| Chlorobenzene-d5 (Surr, 12.5 ppb) | 3.4 | | 12.1 | 97 | 2.2 |
| Chlorobenzene | 2.5 | 0.11 | 0.50 | 99 | 6.9 |
| Ethylbenzene | 4.4 | 0.10 | 0.48 | 96 | 6.5 |
| 1,1,1,2-Tetrachloroethane | 5.7 | 0.07 | 0.45 | 90 | 5.0 |
| m-, p-Xylene | 4.6 | 0.09 | 0.94 | 94 | 3.1 |
| o-Xylene | 5.1 | 0.12 | 0.48 | 96 | 8.2 |
| Styrene | 5.9 | 0.12 | 0.50 | 99 | 7.5 |
| Bromoform | 10.1 | 0.06 | 0.46 | 93 | 3.9 |
| Isopropylbenzene | 6.0 | 0.07 | 0.46 | 91 | 5.2 |
| 4-Bromofluorobenzene (Surr, 12.5 ppb) | 4.2 | | 12.2 | 97 | 1.8 |
| Bromobenzene | 3.5 | 0.08 | 0.50 | 100 | 5.2 |

Table III Method 524.2 Calibration, Accuracy and Precision Data

| Compound | Calibration | | Accuracy and precision (n=7, 0.5 ppb) | | |
|---------------------------------------------|---------------------|-----------|---------------------------------------|--------------|------------------|
| | Linearity RF (%RSD) | MDL (ppb) | Avg. Conc. (ppb) | Accuracy (%) | Precision (%RSD) |
| n-Propylbenzene | 5.7 | 0.06 | 0.48 | 96 | 4.3 |
| 1,1,2,2-Tetrachloroethane | 7.9 | 0.17 | 0.52 | 104 | 10.6 |
| 2-Chlorotoluene | 5.1 | 0.09 | 0.48 | 96 | 5.6 |
| 1,3,5-Trimethylbenzene | 6.3 | 0.11 | 0.47 | 93 | 7.5 |
| 1,2,3-Trichloropropane | 4.9 | 0.12 | 0.52 | 103 | 7.4 |
| trans-1,4-Dichloro-2-butene | 8.2 | 0.34 | 0.54 | 108 | 20.2 |
| 4-Chlorotoluene | 7.0 | 0.05 | 0.49 | 97 | 3.5 |
| t-Butylbenzene | 8.3 | 0.15 | 0.47 | 95 | 9.7 |
| Pentachloroethane | 13.2 | 0.36 | 0.37 | 75 | 30.2 |
| 1,2,4-Trimethylbenzene | 4.4 | 0.12 | 0.45 | 89 | 8.8 |
| sec-Butylbenzene | 7.0 | 0.12 | 0.47 | 93 | 8.0 |
| 4-Isopropyltoluene | 5.7 | 0.03 | 0.45 | 90 | 1.8 |
| 1,3-Dichlorobenzene | 6.0 | 0.11 | 0.47 | 93 | 7.5 |
| 1,4-Dichlorobenzene-d4 (Surr, 12.5 ppb) | 4.1 | | 12.3 | 98 | 0.5 |
| 1,4-Dichlorobenzene | 5.1 | 0.08 | 0.51 | 103 | 4.9 |
| n-Butylbenzene | 5.2 | 0.13 | 0.47 | 94 | 8.5 |
| Hexachloroethane | 6.8 | 0.13 | 0.49 | 98 | 8.6 |
| 1,2-Dichlorobenzene-d4 (Surr, 12.5 ppb) | 4.1 | | 12.1 | 96 | 2.4 |
| 1,2-Dichlorobenzene | 4.7 | 0.07 | 0.50 | 100 | 4.4 |
| Nitrobenzene (2 ppb to 50 ppb) ¹ | 1.7 | | | | |
| 1,2-Dibromo-3-chloropropane | 12.1 | 0.42 | 0.67 | 133 | 20.0 |
| Hexachlorobutadiene | 7.1 | 0.17 | 0.52 | 104 | 10.6 |
| 1,2,4-Trichlorobenzene | 9.1 | 0.12 | 0.49 | 97 | 7.8 |
| Napthalene | 13.7 | 0.08 | 0.49 | 98 | 5.3 |
| 1,2,3-Trichlorobenzene | 10.3 | 0.13 | 0.50 | 100 | 8.2 |

1. Compound is a poor purger. Data is presented for calibration curve data and the detectable range. MDL samples were not detectable.
2. Co-eluting peaks. Different ions used due to interference on this column between 1,1-dichloropropanone and 4-methyl-2-pentanone.

Figure 1 Primary Characteristic Ions for the First Six Gases of a 0.2 ppb Standard Indicating Excellent Detection Limits with Minimal Interference from Water.





Figure 2 Total Ion Chromatogram of a 5 ppb VOC Standard with an Inset of the Primary Characteristic Ions for the First Six gases Indicating Consistent Peak Shapes for all Compounds with No Water Interference.



Conclusion

The Teledyne Tekmar Lumin Purge and Trap Concentrator and AQUATek 100 was used to process water samples containing VOCs following US EPA Method 524.2 with detection by an Agilent 7890B GC/5977A MS. The %RSD of the calibration curve passed all method requirements with no interference from excessive water. The MDL, precision and accuracy for seven 0.5 ppb standards also indicated no interference from excessive water.

References

1. Munch, J.W., US EPA Method 524.2 - Measurement of Purgeable Organic Compounds in Water by Capillary Column Gas Chromatography/Mass Spectrometry, Revision 4.1, 1995