**Monoclonal Mouse  
Anti-Human  
Terminal Deoxynucleotidyl Transferase/FITC**  
Clone HT-6  
Code F7139

**ENGLISH**

**Intended use** For in vitro diagnostic use.  
F7139 is intended for use in flow cytometry. Terminal deoxynucleotidyl Transferase (TdT) is accepted as an intracellular marker for immature lymphocytes in the hemopoietic system (1) and is considered essential for the initial evaluation of acute leukaemia's together with a panel of other antibodies (1, 2, 3). Interpretation of results must be made within the context of the patient's clinical history and other diagnostic tests by a certified professional.

**Summary and explanation** TdT is a 58 kDa protein (2) that catalyses the polymerisation of deoxynucleotides on the 3' hydroxyl end of a polydeoxynucleotide without template instructions (3, 4). TdT is present in the nuclei of normal T and B lymphocyte precursors and the neoplastic equivalent (5). Positive cells are present in normal thymus, and in few normal bone marrow cells, corresponding to hemopoietic precursors (2). Most cases of ALL are TdT positive, 37 % of AML cases are TdT positive (3).

**Reagent provided** F7139, is a purified monoclonal mouse antibody conjugated with Fluorescein Isothiocyanate Isomer 1 (FITC). The conjugate is provided in liquid form in buffer containing 1% bovine serum albumin (BSA) and 15 mmol/L NaN<sub>3</sub>, pH 7.2. Each vial contains 50 tests (10 µL of conjugate for up to 10<sup>6</sup> MOLT-4 cells).  
Isotype: IgG1, kappa. Conjugate concentration mg/L: See label on vial.

Antibody Code	Fluorochrome	Negative Control Code
F7139	FITC (Fluorescein Isothiocyanate Isomer 1)	X0927

**Immunogen** Purified TdT isolated from human leukaemia cells (CML-blast crisis).

**Specificity** Cells labelled by the antibody display staining confined to the nucleus.  
Anti-TdT, HT-6, labels MOLT-3 and MOLT-4 cells (both immature T-cell lines derived from ALL, known to be TdT positive) (6).  
In flow cytometry Anti-TdT, HT-6, demonstrates high expression in acute lymphoid and lower expression in myeloid leukaemia's. Anti-TdT, HT-6, together with other antibodies allows for the recognition of minimal leukemic blasts during clinical remission in AML (7).

**Precautions**

1. For professional users.
2. This product contains sodium azide (NaN<sub>3</sub>), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused solution should be disposed of according to local, State and Federal regulations.

**Storage** Store in the dark at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact our Technical Services.

**Staining procedure**

1. Transfer 50 µL (up to 10<sup>6</sup> cells) of the cell suspension to be analysed (whole blood, bone marrow or mononuclear cells) to a test tube.
2. Add 100 µL Dako IntraStain, Reagent A (Fixation), Code K2311. Mix gently by using a vortex mixer to ensure that the cells are in suspension.
3. Incubate at room temperature for 15 minutes.
4. Add 2 mL PBS and mix gently by using a vortex mixer.
5. Centrifuge at 300 x g for 5 minutes, then aspirate the supernatant, leaving approximately 50 µL of fluid.
6. Mix thoroughly by using a vortex mixer to ensure that the cells are in suspension and add 100 µL Dako IntraStain, Reagent B (Permeabilisation), Code K2311. Add 10 µL F7139. Mix gently by using a vortex mixer to ensure that the cells are in suspension.
7. Use a non-reactive monoclonal antibody of the same isotype, and conjugated with the same fluorochrome, as a negative control (see table).
8. Incubate in the dark at room temperature for 15 minutes.
9. Repeat steps 4 and 5.

10. Resuspend the cells in an appropriate fluid for flow cytometry, e.g. 0.3 mL 1% paraformaldehyde (fixative) in 0.01 mol/L PBS, pH 7.4.

11. Analyse on a flow cytometer

It is recommended to include a suitable positive and negative control sample with each run for reagent and preparation control. Note that fluorochrome conjugates are light sensitive, and samples should be protected from light during the staining procedure and until the analysis.

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

1. Bollum FJ. Terminal deoxynucleotidyl transferase as a hematopoietic cell marker [Review]. Blood 1979;54:1203-15.
2. Terminal Deoxynucleotidyl Transferase (TdT). In: Leong AS-Y, Cooper K, Leong FJW-M, editors. Manual of diagnostic antibodies for immunohistology. London: Oxford University Press; 2003. p. 413-5.
3. Erber WN, Mason DY. Immunoalkaline phosphatase labeling of terminal transferase in hematologic samples. Am J Clin Pathol 1987;88:43-50.
4. Bearman RM, Winberg CD, Maslow WC, Racklin B, Carlson F, Nathwani BN, et al. Terminal deoxynucleotidyl transferase activity in neoplastic and nonneoplastic hematopoietic cells. Am J Clin Pathol 1981;75:794-802.
5. Adriaansen HJ, Hooijkaas H, Kappers-Klunne MC, Hählen K, van't Veer MB, van Dongen JJM. Double marker analysis for terminal deoxynucleotidyl transferase and myeloid antigens in acute nonlymphocytic leukemia patients and healthy subjects. Hematol. Blod Transfus. 1990;33:41-9.
6. T-cell leukemia and T-cell lymphoma cell lines. In: Drexler HG, editor. The leukemia-lymphoma cell line factsbook. San Diego, San Francisco, New York, Boston, London, Sydney, Tokyo: Academic Press; 2001. p. 387-8.
7. Paietta E, Meenan B, Heavey C, Thomas D. Detection of Terminal Transferase in Acute Myeloid Leukemia by Flow Cytometry. Cytometry 1994;16:256-61.

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