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LC - Acquisition

File Home Control

Take Release Status Method Single Sample Sequence Ad Hoc Mode Copy Delete Run Queue Activity Log Instrument Status Online Signals Solvent Monitor Spectrum Sequence Creation Template Status Windows

Sequence - Sequence1.sqx

General Properties Run Options Injections Table

| # | Action                                     | Sample type |
|---|--------------------------------------------|-------------|
| 1 | <input checked="" type="checkbox"/> Inject | Blank       |
| 2 | <input checked="" type="checkbox"/> Inject | Cal. Std.   |
| 3 | <input checked="" type="checkbox"/> Inject | Cal. Std.   |
| 4 | <input checked="" type="checkbox"/> Inject | Sample      |
| 5 | <input checked="" type="checkbox"/> Inject | Sample      |
| 6 | <input checked="" type="checkbox"/> Inject | Sample      |
| 7 | <input checked="" type="checkbox"/> Inject | Sample      |
| 8 | <input checked="" type="checkbox"/> Inject | Sample      |

Sample custom parameters Compound

\* Indicates Mandatory Entry

| Parameter | Value |
|-----------|-------|
|           |       |

Compartment 1

**1.00 L**  
Actual Volume

Purging bottle

Consumption (est.): 0.00 L  
Leftover volume (est.): 1.00 L

**No Bottle**  
Actual Volume

Consumption (est.): 0.00 L  
Leftover volume (est.): 0.00 L

**No Bottle**  
Actual Volume

Consumption (est.): 0.00 L  
Leftover volume (est.): 0.00 L

**No Bottle**  
Actual Volume

Consumption (est.): 0.00 L  
Leftover volume (est.): 0.00 L

**0.15 L**  
Actual Volume

Acetonitrile

Consumption (est.): 0.00 L  
Leftover volume (est.): 0.15 L

**1.00 L**  
Actual Volume

Water

Consumption (est.): 0.00 L  
Leftover volume (est.): 1.00 L

Result path: C:\CDSP\Projects\SLS\Results  
Result name:   
Save result as: One result set

Validate Run


Current user: SYSTEM (SYSTEM), Active project: SLS

**0.15 L**  
Actual Volume

Acetonitrile

Consumption (est.): 0.00 L  
Leftover volume (est.): 0.15 L

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
## InfinityLab Assist Troubleshooting

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# Instrument Status

Quat. Pump → **Standby**

Switch On Set Flow Set Composition Purge Seal Wash Prime

A 1.0 % B 99.0 % Flow 0.000 mL/min

C 0.0 % D 0.0 % Pressure [Limit 400] 0.11 bar

Sampler → **Idle**

Switch to Bypass Needle Wash Thermostat Reset Injector Home Arm

Injection Volume 10.00 µL Temperature 27 °C [---]

Sample Location ---

Column Oven → **Idle**

Thermostat Identify Device

Temperature 25.61 °C [25.0 °C] Column Info Eclipse Plus C18 4.6x250mm, 5µm

VWD → **Idle**

Lamp Off Balance Identify Device

Wavelength A 254 nm UV Lamp On Absorbance A 104 mAU

Error

←

Assisted Troubleshooting: Select Cause

Most probable cause

Clogged needle or needle seat

Other causes

Clogged column

Misconfigured method parameters

Clogged heat exchanger

Clogged capillary

Cancel

Notifications

?

×

Binary Pump G7132A:LPBBV00014

less than a minute ago

^

▲

✖

Pressure above high pressure limit (600 bar)

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
InfinityLab Certified ID Vials

Sample Linking Station  
with Sample Linking Software

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Analyzing Iohexol by Compendial Method Produces Excellent Precision with Shallow Gradients Using an Agilent 1260 Infinity III LC

**Abstract**  
The separation of Iohexol and its impurities is a common problem for the analysis of pharmaceuticals. Several methods have been reported for the separation of Iohexol and its impurities, but the results are often inconsistent. This paper describes a new method for the separation of Iohexol and its impurities using a compendial method. The method involves the use of a shallow gradient and a compendial method to achieve excellent precision and resolution. The results show that the method is robust and can be used for the separation of Iohexol and its impurities in a wide range of samples.

|          |           |       |
|----------|-----------|-------|
| Gradient | 0 min     | 1% B  |
|          | 60 min    | 13% B |
|          | 60.01 min | 99% B |


Figure 2. Chromatogram overlay (N = 8) of the separation using the compendial method.

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