
## EVERIAS ROCKER 2





With an unprecedented 5 year limited warranty, Benchmark's Everlast 247 Rocker is built to last. Twenty four hours a day, seven days a week... No problem! The powerful motor and heavy-duty drive system remain reliable and consistent, even under constant usage in a cold room (min. 0°C) or incubator (max. 65°C).

The rocking speed and tilt angle can easily be adjusted to vary between gentle and aggressive agitation. Both the standard platform and the (optional) stacking platform are constructed of stainless steel, preventing against any rusting or corrosion. The Everlast 247 is designed to provide years of uninterrupted service.

## Benchmark

Everlast™ 5 year warranty

Stainless steel surface

16 lbs. load capacity

Adjustable speed & tilt angle

Optional stacking platform

## **Technical Data:**

Speed Range: 3 to 80 rpm

Tilt Angle: 0 to 20° (user adjustable) Platform Size: 16 x 12 in. / 40.5 x 30cm

Load Capacity: 12 lb / 4.45 kg 16 x 13 x 6 in. 40.5 x 33 x 15.25 cm

Operating Temp.: 0°C to +65°C Weight: 20 lb / 9 kg

115V or 230V, 50-60Hz

Warranty: 5 years

## Ordering Information:

BR5000\* Everlast™ Rocker 247 with non-slip rubber mat

**Optional Accessories:** 

BR5000-STACK Optional stacking platform for BR5000

Clearance: 3.5 in.

\* To order a product in 230V / 50Hz, please add -E to item number.

Benchmark

**Benchmark Scientific Inc.** PO Box 709 Edison, NJ 08818 Phone: 908-769-5555

Fax: 908-222-1864

Web: www.BenchmarkScientific.com Email: Info@BenchmarkScientific.com