



Scoops, Spoons & LabWare



SteriWare PharmaScoop®

The Industry Standard - High Quality, Single Use Powder Scoops

The Sampling Systems range of disposable PharmaScoops® have been specially designed for dispensaries and manufacturing. They are the ultimate in quality and affordability.

SteriWare® PharmaScoops® are used by pharmaceutical, food and industrial companies all over the world.

- Moulded & packed in a Cleanroom.
- Individually bagged in ultra strong, easy-to-open packaging
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming materials
- BSE/TSE free
- Full batch traceability (batch number is printed on every pack)
- Available sterile (SAL 10⁻⁶) and non-sterile from stock



The GMP PharmaScoop

It is Good Manufacturing Practice to keep your hands out of the product as much as possible. In accordance with this we have positioned the ergonomic handles high up on the scoop.

The Biggest Choice - Make sure you have the right scoop



Don't leave your scoop behind!
A brightly coloured scoop can be easily seen.

The Correct Material and Colour

We have the widest choice of materials and colours available:

PS - great for most applications. Strong and tough

HDPE and PP - extra levels of chemical resistance

Metal Detectable - class leading PP, with a special metal detectable filler

Coloured scoops are ideal when you are concerned that the scoop could be left in the product.

Wide Range of Sizes and Shapes

We know that not just one shape fits all applications so we have a wide range of styles and sizes so you can have the best scoop for your process.

Long Handle

A long handle can improve the ergonomics for the operator and make the scoop a lot more visible if it is left in the product

Still Not Sure?

Try before you buy. Contact us for a free sample.



White PharmaScoop PS

The original PharmaScoop

Part No.	Nominal Capacity	Sterile	Box Quantity
8045A-30	30ml	No	100
8045A-30S	30ml	Yes	100
8045A-60	60ml	No	100
8045A-60S	60ml	Yes	100
8045A-125	125ml	No	100
8045A-125S	125ml	Yes	100
8045A-250	250ml	No	85
8045A-250S	250ml	Yes	85
8045A-500	500ml	No	50
8045A-500S	500ml	Yes	50
8045A-1000	1000ml	No	25
8045A-1000S	1000ml	Yes	25
8045A-2500	2500ml	No	15
8045A-2500S	2500ml	Yes	15



**Best
Seller**

- ✓ Made from Polystyrene (PS)
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

PharmaScoop HDPE

Extra Chemical Resistance

Part No.	Nominal Capacity	Sterile	Box Quantity
8045H-30	30ml	No	100
8045H-60	60ml	No	100
8045H-60S	60ml	Yes	100
8045H-125	125ml	No	100
8045H-125S	125ml	Yes	100
8045H-250	250ml	No	85
8045H-250S	250ml	Yes	85
8045H-500	500ml	No	50
8045H-500S	500ml	Yes	50
8045H-1000	1000ml	No	25
8045H-1000S	1000ml	Yes	25
8045H-2500S	2500ml	Yes	15



- ✓ Made from HDPE
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

SteriWare

PharmaScoop PP

Great Chemical Resistance

Part No.	Nominal Capacity	Sterile	Box Quantity
8045P-30	30ml	No	100
8045P-30S	30ml	Yes	100
8045P-60	60ml	No	100
8045P-60S	60ml	Yes	100
8045P-125	125ml	No	100
8045P-125S	125ml	Yes	100
8045P-250	250ml	No	85
8045P-250S	250ml	Yes	85
8045P-500	500ml	No	50
8045P-500S	500ml	Yes	50
8045P-1000	1000ml	No	25
8045P-1000S	1000ml	Yes	25
8045P-2500	2500ml	No	15
8045P-2500S	2500ml	Yes	15



- ✓ Made from Polypropylene (PP)
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

Red PharmaScoop PS

Highly Visible

Part No.	Nominal Capacity	Sterile	Box Quantity
8045R-30	30ml	No	100
8045R-30S	30ml	Yes	100
8045R-60	60ml	No	100
8045R-60S	60ml	Yes	100
8045R-125	125ml	No	100
8045R-125S	125ml	Yes	100
8045R-250	250ml	No	85
8045R-250S	250ml	Yes	85
8045R-500	500ml	No	50
8045R-500S	500ml	Yes	50
8045R-1000	1000ml	No	25
8045R-1000S	1000ml	Yes	25
8045R-2500	2500ml	No	15
8045R-2500S	2500ml	Yes	15



- ✓ Made from Polystyrene (PS)
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

Blue PharmaScoop PP

Highly Visible with High Chemical Resistance

Part No.	Nominal Capacity	Sterile	Box Quantity
8045B-60	60ml	No	100
8045B-60S	60ml	Yes	100
8045B-125	125ml	No	100
8045B-125S	125ml	Yes	100
8045B-250	250ml	No	85
8045B-250S	250ml	Yes	85
8045B-500	500ml	No	50
8045B-500S	500ml	Yes	50
8045B-1000	1000ml	No	25
8045B-1000S	1000ml	Yes	25
8045B-2500	2500ml	No	15
8045B-2500S	2500ml	Yes	15



- ✓ Made from Polypropylene (PP)
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

Metal Detectable PharmaScoop PP

Metal Detectable and X-ray Visible

Part No.	Nominal Capacity	Sterile	Pack Quantity
M1000B-30	30ml	No	4
M1000B-60	60ml	No	4
M1000B-125	125ml	No	4
M1000B-250	250ml	No	4
M1000B-500	500ml	No	4
M1000B-1000	1000ml	No	4
M1000B-2500	2500ml	No	4



- ✓ Made from metal detectable Polypropylene (PP)
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

See the full range of Metal Detectable products on pages 48 to 57

GMP Volumetric Scoop PS

Ergonomic Flat Scoop

Part No.	Nominal Capacity	Sterile	Box Quantity
8286A-25	25ml	No	100
8286A-25S	25ml	Yes	100
8286A-50	50ml	No	100
8286A-50S	50ml	Yes	100
8286A-100	100ml	No	100
8286A-100S	100ml	Yes	100
8286A-150	150ml	No	100
8286A-150S	150ml	Yes	100
8286A-250	250ml	No	100
8286A-250S	250ml	Yes	100
8286A-500	500ml	No	40
8286A-500S	500ml	Yes	40
8286A-1000	1000ml	No	30
8286A-1000S	1000ml	Yes	30

New



- ✓ Made from Polystyrene (PS)
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

GMP Sample Scoop With Lid PS

Secures the Sample in the Scoop

Part No.	Nominal Capacity	Sterile	Box Quantity
8916A-25	25ml	No	50*
8916A-25S	25ml	Yes	50*
8916A-50	50ml	No	50*
8916A-50S	50ml	Yes	50*
8916A-100	100ml	No	50*
8916A-100S	100ml	Yes	50*

* Proposed Box Quantity - this may change during packaging trials please check with Sampling Systems

New



- ✓ Made from Polystyrene (PS)
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

SteriWare

SteriWare®

Download Certificates



Batch specific Certificates of Conformity and Certificates of Irradiation* are available for download from:

www.sampling.com/certs

Simply enter your batch number into the search box and download your certificates in PDF format. You will find your batch number printed on the box label and on every bag.

*Certificates of Irradiation are not available for products that are sold as "non-sterile"

SteriWare® Long Handled Scoop

Keep hands away from the product

SteriWare Long Handled Scoops are the ideal tool for reaching into deep containers. The extra long 300mm handle will help operators keep their hands out of the product.

The other advantage of a long 300mm handle is that it makes the scoop more visible if it is left in the powder.

Both non-sterile and sterile scoops are available from stock. They are manufactured in a Cleanroom to ensure cleanliness and reduce particulate contamination.

- Moulded & Packed in a Cleanroom
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming HDPE
- BSE/TSE Free
- Available in uncoloured natural HDPE or bright red HDPE
- Full Batch Traceability
- Available sterile (SAL 10⁻⁶) and non-sterile

The handle length is 300mm. Our range of single use scoops are ideal for use in the laboratory, production and warehouse.



HDPE - Natural

Part No.	Sterile	Nominal Scoop Volume	Box Quantity
8341H-25	No	25ml	50
8341H-25S	Yes	25ml	50
8341H-50	No	50ml	50
8341H-50S	Yes	50ml	50
8341H-100	No	100ml	50
8341H-100S	Yes	100ml	50



HDPE - Red

Part No.	Sterile	Nominal Scoop Volume	Box Quantity
8341R-50S	Yes	50ml	50
8341R-100S	Yes	100ml	50



Metal Detectable PP (Polypropylene) - BLUE

Part No.	Sterile	Nominal Scoop Volume	Pack Quantity
M1040B-100	No	100ml	4



SteriWare® Volumetric Spoons

Single use, accurate measuring spoons

Volumetric Spoons allow small volumes to be accurately measured. The Volumetric Spoons are designed to be used straight from the bag.

- Moulded & Packed in a Cleanroom
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming materials
- BSE/TSE Free
- Full Batch Traceability
- Available sterile (SAL 10⁻⁶) and non-sterile

Note: Volumetric Spoons are not certified measuring devices. User validation is required. If in doubt please contact us



HDPE - Natural

Part No.	Flat Fill Volume	Sterile	Box Quantity
8640H-01	1ml	No	100
8640H-01S	1ml	Yes	100
8640H-0123S	1.23ml	Yes	100
8640H-0246S	2.46ml	Yes	100
8640H-05	5ml	No	100
8640H-05S	5ml	Yes	100
8640H-10	10ml	No	100
8640H-10S	10ml	Yes	100
8640H-20	20ml	No	100
8640H-20S	20ml	Yes	100
8640H-30	30ml	No	100
8640H-30S	30ml	Yes	100

New sizes



Metal Detectable Polypropylene (PP) - BLUE

Part No.	Flat Fill Volume	Sterile	Pack Quantity
M1030B-25	2.5ml	No	4

New



SteriWare

See the full range of Metal Detectable products on pages 48 to 57

SteriWare® Sample Spoons

Ideal for dispensing and sampling small quantities



Our range of low cost single use spoons allow quick and efficient surface sampling of powders and granules. Alternatively they can be also used for dispensing small quantities.

- Moulded & Packed in a Cleanroom
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming polymers
- BSE/TSE Free
- Full Batch Traceability
- Available sterile (SAL 10⁻⁶) and non-sterile

Don't lose your spoon! A brightly coloured spoon can be easily seen.



The Correct Material and Colour

Having the widest choice of materials available.

PS - great for most applications. Strong and tough

HDPE and PP - extra levels of chemical resistance

Metal Detectable - class leading PP

Coloured spoons - ideal when you are concerned that the spoon might be left in the powder

Wide Range of Sizes and Shapes

One shape of spoon will not suit all applications. We offer a huge range of spoon styles, sizes and colours so that you can have the best spoon for your application

Still Not Sure?

Try before you buy. Contact us for a free sample.

SteriWare® Certified BSE/TSE Free



We only use high quality polymers for our SteriWare products. All of the polymer resins that we use are certified BSE/TSE free.

During manufacturing the SteriWare products do not come into contact with any products or processes that may pose a BSE/TSE risk.

A BSE/TSE statement is included in our Certificates of Conformity.

Polystyrene (PS) - WHITE

Part No.	Nominal Spoon Volume	Sterile	Box Quantity
8088A-03	3ml	No	100
8088A-03S	3ml	Yes	100
8088A-10	10ml	No	100
8088A-10S	10ml	Yes	100



Polystyrene (PS) - RED

Part No.	Nominal Spoon Volume	Sterile	Box Quantity
8088R-10S	10ml	Yes	100



Polypropylene (PP) - BLUE

Part No.	Nominal Spoon Volume	Sterile	Box Quantity
8088B-10	10ml	No	100
8088B-10S	10ml	Yes	100



Metal Detectable Polypropylene (PP) - BLUE

Part No.	Nominal Spoon Volume	Sterile	Pack Quantity
M1006B-10	10ml	No	4



SteriWare

See the full range of Metal Detectable products on pages 48 to 57

SteriWare® Mini Spoon PS

Small Economical, Double-ended Spoon



SteriWare Mini Spoons have a 5ml spoon at one end a 2.5ml spoon at the other. They are moulded and packed within a Cleanroom and so can be used straight from the bag.

- Moulded & Packed in a Cleanroom
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming polystyrene
- BSE/TSE Free
- Full Batch Traceability

Part No.	Volume	Sterile	Box Quantity
8642A-01	2.5ml and 5ml	No	130
8642A-01S	2.5ml and 5ml	Yes	130