

# Determination of the essential oil content



KOL

## Compact system for determining the essential oil content

Complete compact system for determining the essential oil content in

- Pharmaceuticals,
- Spices
- Seasoning and
- Herbs

according to ISO 6571.

With base frame, heating devices, brackets, coolant hoses and glass apparatus.

Compact system for determining the essential oil content

Make	Item description	Item no.
KOL	with 500 ml flask	B00217736
KOL 2	with 1,000 ml flask	B00602393

Complete compact system for determining the essential oil content in citrus fruits and their derivatives according to Clevenger (3,000 ml distillation flask).

Make	Item description	Item no.
CLE	Compact system for determining the essential oils in citrus fruits and their derivatives	B00217741