

TECHNICAL DATA

Description	KB ECO 240	KB ECO 400	KB ECO 720	KB ECO 1020
Measures				
Interior volume [L]	247	400	700	1020
Net weight of the unit (empty) [kg]	152	196	272	327
Permitted load [kg]	100	120	150	150
Load per rack [kg]	30	30	45	45
Wall clearance back [mm]	100	100	100	100
Wall clearance sidewise [mm]	100	100	100	100
Internal Dimensions				
Innenraum Breite [mm]	650	650	973	973
Innenraum Höhe [mm]	785	1200	1250	1250
Innenraum Tiefe [mm]	485	390	576	836
Housing dimensions not incl. fittings and connections				
Width net [mm]	925	925	1250	1250
Height net [mm]	1461	1946	1925	1925
Depth net [mm]	796	796	885	1149
Performance Data Temperature				
Temperature range (max. 26 °C under ambient temperature) [°C]	–	–	0...70	0...70
Temperature range (max. 28 °C under ambient temperature) [°C]	0...70	0...70	–	–
Temperature variation at 37 °C [± K]	0.1	0.1	0.1	0.1
Temperature fluctuation at 37 °C [± K]	0.1	0.1	0.1	0.1
Recovery time after 30 seconds door open at 37 °C [min]	3	6	4	4
Electrical data				
Nominal power (version 120 V) [kW]	0.9	1.2	1.3	1.3
Nominal power (version 230 V) [kW]	0.9	1.2	1.3	1.3
Environment-specific data				
Energy consumption at 37 °C [Wh/h]	65	70	105	105
Fixtures				
Number of shelves (std./max.)	2/9	2/15	2/15	2/15

ORDERING INFORMATION [Art. No.]

Further details: Go to www.binder-world.com > Search > Enter Article Number

Voltage	Option model	KB ECO 240	KB ECO 400	KB ECO 720	KB ECO 1020
200...230 V 1~ 50/60 Hz	Standard	9020-0423	9020-0465	9020-0424	9020-0425
100...120 V 1~ 50/60 Hz	Standard	9020-0426	9020-0466	9020-0427	9020-0428