

Detergents and additives

ACIDGLC2 (art. 901117)

ACIDGLASS C2

- liquid phosphate-free neutralizer
- containing citric acid monohydrate
- 5 liters can



ACIDGLP2 (art. 901171)

ACIDGLASS P2

- acid liquid neutralizing agent for laboratory glassware
- containing phosphoric acid, citric acid monohydrate
- 5 liters can



DETERGLASS2 (art. 905448)

DETERGLASS2

- universal alkaline powder
- detergent for automatic washing of laboratory glassware
- contains sodium hydroxide, sodium metasilicate pentahydrate, sodium triphosphate pentabasic, triclosene sodium dihydrate
-
- 10 kg pack



DETERLQ2 (art. 901116)

DETERLIQUID C2

- strong universal liquid alkaline detergent phosphate-free for washing steel up to 95 °C
- also suitable for washing rubber and plastic at 60 °C
- containing sodium silicate, polyacrylate, potassium hydroxide
- 5 liters can



DETERLQI2 (art. 904541)

DETERLIQUID I2

- liquid alkaline detergent
- containing sodium hydroxide, phosphonic-carboxylic acid, sodium hypochlorite
- 5 liters can



DETERLQSP2 (art. 903659)

DETERLIQUID SP2

- liquid alkaline detergent phosphate free for water analysis and pharmaceuticals
- containing potassium hydroxide, sodium hydroxide
- 5 liters can

