

Traceability

KITCON232 (art. 903284)

RS232 CONNECTION KIT

- RS232 – USB converter
- male – female straight connecting cable
- female – female crossover connecting cable



KITSOND (art. 905066)

KIT FOR TEMPERATURE PROBE ACCESS

- pass-through hole for access to the washing chamber, for temperature probes
- normally closed with 3/8" cap, useful diameter 1/2"
- for traceability and validation activities
- placed under the top covering, in the upper part of the machine on the right side

Machines

GW0260-00-00J	GW0260-S0-00J	GW1260-00-000	GW1260-S0-000	GW1260-SC-000
GW1245-S0-100	GW3060BX	GW3060BXA	GW3060BXC	GW3060BXCA
GW3060P	GW3060PA	GW3060S	GW3060SA	GW3060SC
GW3060SCA	GW4090	GW4090A	GW4090BXC	GW4090BXCA
GW4090BXP	GW4090C	GW4090CA	GW4090P	GW4090PA
GW4260-00-000	GW4260-S0-000	GW4260-SC-000	GW4290-00-000	GW4290-SC-000

KITSOND10 (art. 904433)

KIT FOR TEMPERATURE PROBE ACCESS

- pass-through hole for access to the washing chamber, for temperature probes
- normally closed with 3/8" cap, useful diameter 1/2"
- for traceability and validation activities
- placed under the top covering, in the upper part of the machine on the right side

Machines

GW6010L	GW6010LX	GW6010ML	GW6010MLX
---------	----------	----------	-----------

WD-232PRINT (art. 905557)

SERIAL PORT RS232

- RS232 serial port, for external printer (art. WD-PRINTE)
- for machines of the 1260, 4260 product families



Machines

GW4260-00-000	GW4260-S0-000	GW4260-S0-060	GW4260-SC-000	GW1260-00-000
GW1260-S0-000	GW1260-SC-000			

WD-LAN60 (art. 903460)

DATA COMMUNICATION CARD

- Device for LAN ethernet connection to the machine
- Makes monitoring and communication possible from every access point to the LAN net
- To be configured as in the instruction manual

Machines

GW6010L	GW6010LX	GW6010ML	GW6010MLX	GW6010T
GW6010TM	GW6290DS	GW7010	GW7010M	GW7010T
GW7010TM	GW7015	GW7015T		

WD-LANE (art. 903091)

DATA COMMUNICATION CARD

- Device for LAN ethernet connection to the machine
- Makes monitoring and communication possible from every access point to the LAN net
- To be configured as in the instruction manual



Machines

GW3060BX	GW3060BXA	GW3060BXC	GW3060BXCA	GW3060P
GW3060PA	GW3060S	GW3060SA	GW3060SC	GW3060SCA

WD-LANI (art. 903090)

DATA COMMUNICATION CARD

- Device for LAN ethernet connection to the machine
- Makes monitoring and communication possible from every access point to the LAN net
- To be configured as in the instruction manual



Machines

GW4090	GW4090A	GW4090BXC	GW4090BXCA	GW4090BXP
GW4090C	GW4090CA	GW4090P	GW4090PA	

WD-LANI2 (art. 905555)

COMMUNICATION LAN CARD

- data communication card
- for machines of the 0260, 1260, 4260 product families
- equipped with RS232 serial port for connection to an external printer
- equipped with a LAN port for ethernet connection which allows communication and monitoring of the machine from any access point in the local network
- to be set using the instructions in the installation manual



Machines

GW1245-S0-100	GW1260-00-000	GW1260-S0-000	GW1260-SC-000	GW4260-00-000
GW4260-S0-000	GW4260-SC-000	GW4290-00-000	GW4290-SC-000	

WD-PRINT (art. 901702)

PANEL PRINTER

- dot matrix thermal printing
- real time clock
- for panel fastening
- dimensions LxDxH: 120 x 170 x 40 mm
- thermal paper roll diameter: 50 mm - width: 57 mm (cod. WD-PAPER)



Machines

GW4090	GW4090A	GW4090BXC	GW4090BXCA	GW4090BXP
GW4090C	GW4090CA	GW4090P	GW4090PA	

WD-PRINT4290 (art. 905444)

4290 PANEL PRINTER

- dot matrix thermal printing
- real time clock
- for panel fastening
- dimensions LxDxH: 120 x 170 x 40 mm
- thermal paper roll diameter: 50 mm - width: 57 mm (art. WD-PAPER)



WD-PRINT6290 (art. 904555)

PANEL PRINTER

- dot matrix thermal printing
- real time clock
- for panel fastening
- dimensions LxDxH: 120 x 170 x 40 mm
- thermal paper roll diameter: 50 mm - width: 57 mm (cod. WD-PAPER)

Machines

GW6290DS

WD-PRINTE (art. 901783)

EXTERNAL PRINTER

- dot matrix thermal printing
- real time clock
- thermal paper roll diameter: 50 mm - width: 57 mm (cod. WD-PAPER)
- equipped with independent power supply



Machines

GW1245-S0-100	GW3060BX	GW3060BXA	GW3060BXC	GW3060BXCA
GW3060P	GW3060PA	GW3060S	GW3060SA	GW3060SC
GW3060SCA	GW4090	GW4090A	GW4090BXC	GW4090BXCA
GW4090BXP	GW4090C	GW4090CA	GW4090P	GW4090PA

WD-PRINTE2 (art. 905556)

EXTERNAL PRINTER

- external printer (such as WD-PRINTE)
- including RS232 serial port card
- for dot matrix thermal printing
- complete with clock function
- to be combined with a roll of thermal paper diameter 50 mm - width 57 mm (art. WD-PAPER)
- equipped with independent power supply
- NOTE: authorized engineer is required for installation



Machines

GW4260-00-000	GW4260-S0-000	GW4260-S0-060	GW4260-SC-000	GW4290-00-000
GW4290-SC-000				

WD-PRINTEINK (art. 903740)

EXTERNAL IMPACT PRINTER

- impact printer
- paper roll diameter: 50 mm - width: 57 mm (cod. CARTADP24)
- equipped with independent power supply



Machines

GW6010L	GW6010LX	GW6010ML	GW6010MLX	GW6010T
GW6010TM	GW7010	GW7010M	GW7010T	GW7010TM
GW7015	GW7015T			

WD-PRINTP6010 (art. 903462)

CLEAN SIDE PANEL PRINTER

- dot matrix thermal printing
- real time clock
- fitted on clean side panel
- dimensions LxDxH: 120 x 170 x 40 mm
- thermal paper roll diameter: 50 mm - width: 57 mm (cod. WD-PAPER)



Machines

GW6010L	GW6010LX	GW6010T
---------	----------	---------

WD-PRINTP7010 (art. 903463)

CLEAN SIDE PANEL PRINTER

- dot matrix thermal printing
- real time clock
- fitted on clean side panel
- dimensions LxDxH: 120 x 170 x 40 mm
- thermal paper roll diam. 50 mm width 57 mm (cod. WD-PAPER)



Machines

GW7010 GW7010T GW7015

WD-PRINTS6010 (art. 903407)

UNCLEAN SIDE PANEL PRINTER

- dot matrix thermal printing
- real time clock
- fitted on unclean side panel
- dimensions LxDxH: 120 x 170 x 40 mm
- thermal paper roll diameter: 50 mm - width: 57 mm (cod. WD-PAPER)



Machines

GW6010L GW6010LX GW6010ML GW6010MLX GW6010T
 GW6010TM

WD-PRINTS7010 (art. 903464)

UNCLEAN SIDE PANEL PRINTER

- dot matrix thermal printing
- real time clock
- fitted on unclean side panel
- dimensions LxDxH: 120 x 170 x 40 mm
- thermal paper roll diam. 50 mm width 57 mm (cod. WD-PAPER)



Machines

GW7010 GW7010M GW7010T GW7010TM GW7015